

Sysmex Reagents

WDF lysercell and fluoroWDF Feb 2022

Reagent	Discrete mode	
	CBC	CBC+WDF+WPC+RET+PLT-F
Total reagent volume	Approx. 29.9 mL	Approx. 62.8 mL
CELLPACK DCL	Approx. 27.9 mL	Approx. 53.7 mL
SULFOLYSER	Approx. 0.5 mL	Approx. 0.5 mL
Lysercell WNR	Approx. 1.5 mL	Approx. 2.5 mL
Fluorocell WNR	Approx. 20 µL	Approx. 20 µL
Lysercell WDF	-	Approx. 1.5 mL
Fluorocell WDF	-	Approx. 20 µL
Lysercell WPC	-	Approx. 1.5 mL
Fluorocell WPC	-	Approx. 20 µL
CELLPACK DFL	-	Approx. 3.0 mL
Fluorocell RET	-	Approx. 20 µL
Fluorocell PLT	-	Approx. 20 µL

Lysercell WDF + Fluorocell WDF = Differential channel



Lysercell WDF

5L per box

1.5ml per test


3333 sample per box

Store 2 - 35°C

Use 15 – 30°C

On-Board 90 days

37 tests per day minus QC,
primes etc



WDF Lysercell™ WDF

REF ZPPAL337564

5 L

Storage and shelf life after first opening:
1. Once opened, this reagent should be used within 90 days.
2. Do not freeze.

Components:
Organic quaternary ammonium salts 0.07 %
Nonionic surfactant 0.17 %

Mfg. Date 2020-03-23
Exp. Date 2021-03-23
Lot No. AM0006

+2°C +35°C



Sysmex Corporation


1-5-1 Wakinohama-Kaigandori,
Chuo-ku, Kobe 651-0073, Japan

Manufacturing site:

Sysmex Asia Pacific Pte Ltd.
19 Jalan Tukang Singapore 619265

MADE IN SINGAPORE

Sysmex is Registered in U.S. Pat. & Tm. Off.





Lysercell WDF

2 x 42ml cassette per box

20ul per test

2000 tests per cassette,
4000 tests per box

Store 2 - 35°C

Use 15 – 30°C

On-Board 90 days (22 tests per day)