

Title:	Microscope Use in Provider Performed Testing		
Department/Service Line:	Laboratory		
Approver(s):	CLIA Director		
Location/Region/Division:	Baylor Scott & White Health		
Document Number:	BSWH.LAB.PPT.002.R_V1		
Last Review/Revision Date:	See Signatures	Origination Date:	08/2015

SCOPE

This document applies to providers that use a microscope for provider performed testing within Baylor Scott & White Health.

DEFINITIONS

None

METHOD/UTILITY

Microscopes are precision instruments. Achieving an accurate result with microscopic techniques requires an understanding of the operating characteristics and limitations of the equipment used. Therefore, analyzing specimens with the microscope requires training, knowledge of standard precautions, and an understanding of the capabilities, use, and care of the microscope.

PROCEDURE

Parts of the Microscope

Lenses

- The objective lens magnifies the specimen a defined amount. The objective produces the primary image and the eyepiece magnifies it. The total magnification of the image is the product of the magnification of the objective multiplied by the magnification of the eyepiece.
- The oculars (eyepieces) are placed in the top openings of the observation tubes of the microscope and magnify the primary image projected by the objective.

Stand

The stand rests on the base of the microscope and carries the arm and the stage on which the specimen is placed.

Stage

The horizontal platform on which the specimen slide is placed.

Coarse and Fine Adjustment Knobs


Used to bring the slide into focus by moving the stage toward or away from the objectives.

Condenser

The condenser is mounted under the stage to concentrate and focus light from the light source. It can be raised or lowered by means of a condenser knob below the stage. The aperture iris diaphragm can be opened or closed to control the amount of light striking the specimen.

Illumination

The built in light source is usually a tungsten bulb or tungsten-halogen bulb. A switch or dimmer may be used to control the intensity of the light. Lowering the intensity before turning off the microscope lengthens the life of the bulb.



Operation

Low Power Objective (10x)

1. Place a properly labeled specimen slide on the stage.
2. Lower the condenser and turn the lamp to low. Open both the aperture diaphragm and the field diaphragm.
3. With the 10x objective in place and while observing from the *side of the slide and not through the eyepiece*, use the coarse adjustment knob to slowly raise the stage until the slide comes close to the objective.

4. Look through the eyepiece and use the coarse and/or fine adjustment until the image is sharply in focus.
5. Close the field diaphragm almost completely and raise the condenser until the edges of the diaphragm are sharply focused (the condenser is usually at its highest position). Then, open the field diaphragm slowly, stopping just as it disappears from view.
6. Open and close the aperture diaphragm to optimize contrast. Contrast is increased by closing the aperture. If more light is needed, turn up the lamp.

High Power Objective (40x)

1. Focus and center the specimen with the 10x objective as instructed above.
2. Rotate the nosepiece slowly to bring the 40x objective into the light path.
3. Use the FINE adjustment knob to bring the specimen into focus.
4. *Note:* Never raise the stage with the coarse adjustment knob when using the 40x lens. Doing so may cause the lens to hit the slide and break.

Care of the Microscope

- o Microscopes are inspected, cleaned, and checked each day of patient use. For daily cleaning use a high quality lens paper. Document daily maintenance on the Microscope Maintenance Log.
- o Annual microscope maintenance may be performed by Biomed or a contracted third party.
- o Cover the microscope when not in use and leave the 10x objective in position.

ATTACHMENTS

Microscope Maintenance Log for PPM Procedures (BSWH.LAB.PPT.002.A_V1)

RELATED DOCUMENTS

Provider Performed Testing Program (BSWH.LAB.PPT.001.P_V1)

REFERENCES

1. CLSI. *Physician and Nonphysician Provider-Performed Microscopy Testing; Approved Guideline – Second Edition*. CLSI document POCT10-A2. Wayne, PA: Clinical and Laboratory Standards Institute; 2011.

REVISION HISTORY

Version #	Effective Date	Description of Change	Revised By	Removed Date

Attachment Title:	Microscope Maintenance Log for PPM Procedures		
Attachment Number:	BSWH.LAB.PPT.002.A_V1	Last Review/Revision Date:	See Signatures

Location _____
 Year _____

Daily: Inspect microscope and evaluate the need for cleaning or service.

- 10x and 40x objectives must be kept clean at all times.
- If necessary, clean microscope objective lenses, eyepieces, and condenser. Use a high quality lens paper.
- Cleaning by a technical representative is performed as needed (at least annually).

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31
JAN																															
FEB																															
MAR																															
APR																															
MAY																															
JUN																															
JUL																															
AUG																															
SEP																															
OCT																															
NOV																															
DEC																															

Document performance of daily maintenance with initials.
X indicates microscope not in use that day.

Monthly Site Review:

Signature		Date	Signature		Date	Signature		Date	Signature		Date
JAN			APR			JUL			OCT		
FEB			MAY			AUG			NOV		
MAR			JUN			SEP			DEC		