



# HealthPartners/GHI

<b>Blood Culture Collection Procedure</b> (For Bacteria, Fungus and/or Yeast)	<b>Attachments</b> <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No
<b>Key words</b> Phlebotomy, Blood Culture, Blood draw, Sepsis	<b>Number</b> <b>GHP-PC-CLINIC LAB- Procedures- Blood Culture Collection v. 11-2009</b>
<b>Category</b> Provision of Care	<b>Effective Date</b> March 1995
<b>Manual</b> Clinic Laboratory Procedure Manual	<b>Last Review Date</b> November 2012
<b>Issued By</b> Clinic Laboratory Administration	<b>Next Review Date</b> November 2013
<b>Applicable</b> Clinic Laboratory Staff	<b>Origination Date</b> March 1995
	<b>Retired Date</b>
<b>Level of Complexity</b> Waived	<b>Approved Date</b> June 2007
<b>Review Responsibility</b> Laboratory Technical Consultants	<b>Contact</b> Laboratory Technical Consultants
<b>APPROVAL(S)</b> Laboratory Medical Director	

## PURPOSE/PRINCIPLE

The purpose of blood culture collection is to obtain uncontaminated blood for the purpose of diagnostic laboratory testing.

Bacteremia occurs when the rate of multiplication of bacteria in the bloodstream exceeds the capacity of the reticuloendothelial system to remove the bacteria. The clinical pattern of bacteremia may be transient (manipulation of infected tissue, manipulation of contaminated surfaces early systemic and localized infection), intermittent (urinary tract infections, respiratory infections, meningitis or undrained abscesses) or continuous (bacterial endocarditis).

The blood specimen is inoculated into the appropriate blood culture bottles. The bottles are inserted into the BACTEC fluorescent series instrument for incubation and periodic reading. Each bottle contains a chemical sensor which can detect increases in CO<sub>2</sub> produced by the growth of microorganisms. The sensor is monitored every ten minutes for an increase in fluorescence, which is proportional to the CO<sub>2</sub> present. A positive reading indicates the presence of viable organisms in the bottle. Resins have been added in the BACTEC Plus Aerobic and Peds Plus Bottles to enhance recovery of organisms.

## POLICY

Laboratory Staff collecting blood specimens will follow the approved techniques for specimen collection outlined in this procedure.

HealthPartners family of care uses single-use needle and devices for all phlebotomy and blood collection procedures. Should it be necessary to re-stick a patient, a new, single-use needle or device will be used.

**PROCEDURES**

**Specimen for Bacteria:**

The specimen must be collected using sterile technique to reduce the chance of contamination. Routinely, two sets of blood cultures are ordered. When multiple cultures are to be collected follow the Provider’s instructions or collect at least 15 minutes apart. Do not collect more than 3 blood culture sets per day.

Since most bacteremia is intermittent, several techniques may be necessary to isolate the organism responsible for the bacteremia. The following affect the positivity rate:

1. The volume of blood collected - 20 cc is preferred.
2. The number of blood cultures drawn - 3 sets are recommended to rule out bacteremia.
3. The timing of the blood cultures drawn - optimum collection time is before the onset of chills. Subsequent sampling should be spaced at wide intervals.

Recommended collection is before the initiation of antimicrobial therapy. If antimicrobial therapy was initiated before collection, document when therapy was started in the computer.

Blood volumes are as follows:

<b>Criteria</b>	<b>Total Volume</b>	<b>Bottles</b>	<b>Volume Each Bottle</b>
Adult and older children (> 5yrs)  <b>Always collect two sets</b>	16-20 cc <i>Note: If 10 cc or less collected, inoculate the BACTEC Plus Aerobic bottle only.</i>	BACTEC Plus Aerobic  BACTEC Lytic Anaerobic	8-10 cc  8-10 cc
Children under 5 years old	6-10 cc optimal	BACTEC Peds Plus  BACTEC Lytic Anaerobic	3 cc  3-7 cc
Children under 5 years old  <b>If 5 cc or less collected</b>	<i>Note: If 5 cc or less collected, inoculate 3 cc in one BACTEC Peds Plus bottle, and if any remaining inoculate into second BACTEC Peds Plus bottle</i>	BACTEC Peds Plus  BACTEC Peds Plus	3 cc  1-3 cc
Infants or Bone Marrow	1-3 cc	BACTEC Peds Plus	1-3 cc

If using a butterfly draw system, carefully observe the direction of blood flow when starting sample collection. The vacuum in the vial will usually exceed 10 cc, so the user must monitor the volume collected by means of the 5 cc graduation marks on the bottle label. If pediatric bottles (pink) are not available, use the adult BACTEC Plus Aerobic bottle.

<b><u>Materials:</u></b>	<b><u>Stock Number</u></b>
BACTEC Plus Aerobic	#109346
BACTEC Lytic Anaerobic	#109348
BACTEC Peds Plus	#109347



- When final vial or tube is filled, withdraw the needle. Cover the puncture site with a gauze pad and apply pressure.

Syringe/Needle System:

- Insert the needle into the prepared vein and collect 10-20 ml blood in syringe.
- Withdraw needle after collection 10-20 ml blood in syringe.
- Distribute blood equally into aerobic and anaerobic vials.

**NOTES:**


- **Adults and older children: Draw 16-20 cc of blood. Place 8-10 cc of blood in each bottle.**
- **Children under age 5: Draw 6-10 cc of blood. Place approximately 3-5 cc in each bottle.**
- **Pediatric collections of less than 1 cc of blood are placed in the aerobic bottle only.**
- **Adult collections of less than 10 cc are placed in the aerobic bottle only.**
- **Bone marrow: Riverside clinic - See procedure on Bone Marrow.**
- **Refer to Blood Collection Procedure for complete instruction and tips for venipunctures.**

**Specimen Labeling and Transport:**

1. Label both blood culture bottles with the patient's name, chart number, date and clinic. **DO NOT PLACE LABEL OVER BARCODE.** Must use large bar code label on each bottle. See enclosed picture for label placement.

Place large Sunquest barcode label as shown in picture.

Do **NOT** cover BACTEC bottle barcode



2. Mix the contents 2 or 3 times by gentle inversion. **DO NOT INCUBATE SPECIMENS.**
3. Send both bottles to Regions Hospital for testing via the Re-Route system.

**Specimen for Fungus and Yeast**

Fungus and/or yeast cultures are sent to Regions Hospital.

**Supplies:**

2 SPS BD #4960 Sterile Vacutainer Tube (yellow top). Order from Central Lab.

**Specimen:**

Two yellow top SPS Vacutainer tubes are collected from separate sites. Follow instructions for venipuncture procedure for collection on page 2 of this procedure. Disinfect the vacutainer stopper. Draw whole blood in amounts according to these guidelines:

- |    |                           |           |
|----|---------------------------|-----------|
| A. | Infants (up to 18 months) | 0.5 ml    |
| B. | Children under age 5      | 1 - 3 ml  |
| C. | Adolescent & Adults       | 8 - 10 ml |

