

# Training Instructions for Reading Vaginal Wet Mounts for Clue Cells



# Evaluating Wet Mounts for the Presence of Clue Cells

- To clinically diagnose bacterial vaginosis (BV), Amsel's criteria is used.
- One of the four criteria, the presence of  $\geq 20\%$  clue cells per field can be difficult to determine.
- The following tips should be utilized to determine if an epithelial cell is a clue cell or not a clue cell and if there is  $\geq 20\%$  in the field:
  1. Count the number of distinguishable epithelial cells in your field of view. Distinguishable epithelial cells are whole cells that are visible with a nucleus present.

# Evaluating Wet Mounts

2. To determine if any of the epithelial cells are clue cells, it is important to study **ONLY THE BORDERS OF THE CELL**. A cell is a clue cell if the borders of the cell are completely obscured with bacteria and have edges that look “grainy” or “fuzzy”. If any of the border is clear, it is not a clue cell.

**Note:** It is important to realize that the surfaces of epithelial cells can look “grainy” normally due to the cell membrane’s pores and this can be confused as bacteria when it is not.

# Evaluating Wet Mounts

3. To determine the percentage of clue cells in your field:

a. Count the number of clue cells and divide that number by the total number of distinguishable epithelial cells.

For example:


A field has 6 epithelial cells and 2 are clue cells.

$2/6 = 33\%$ . The field is positive for clue cells.

# Evaluating Wet Mounts


- View the following slides to see these tips utilized.

Distinguishable Epithelial Cells (4)



Note: This field is clue cell negative. The majority of the borders of the cells are clear although the surface of the cells appear to be “grainy”.

Clue Cells (2)




Distinguishable Epithelial Cells (3)



Note: This field is positive for clue cells.  
 $2/3 = 66\%$






Distinguishable  
Epithelial Cells (5)  
(Nuclei are  
present)

Clue Cells (3)

Note: This  
field is positive  
for clue cells.  
 $3/5 = 60\%$






Distinguishable  
Epithelial Cells  
(3)


Note: This slide is  
negative for clue  
cells.

The majority of the  
borders are clear on  
each cell, although  
the surfaces appear  
"grainy"



Lactobacillus

Note: This slide is negative for clue cells.  
Notice the large number of neutrophils.




Note: This slide is negative for clue cells.

The borders are clear on each cell, but there are few bacteria on the cells but not enough to cover the cell



Note: This slide is  
negative for clue  
cells.

The borders are clear  
on each cell,  
although the surfaces  
appear "grainy"



Distinguishable Epithelial Cells

Clue Cell

Note: This cell is positive for cl  
cells.

$\frac{1}{4} = 25\%$

This is a cell underneath of  
the clue cell, and since the  
whole cell is not visible, it  
should not be counted

Indistinguishable  
epithelial Cells  
cannot be used for  
this field

NOT  
Cell


Please note:  
this is not a  
cell because  
you can see that  
the visible border  
(the border)  
the cell is  
obscured  
cell on top  
are not covered  
by bacteria

Clue Cells

Please note: These are clue cells because  
the entire borders of both cells are  
completely obscured by bacteria. If the  
borders had any free edge, they would not  
be clue cells.



# Gram stain of clue cells



# Gram stain of normal epithelial cells

