


**BROWN CLINIC  
LABORATORY PROCEDURE MANUAL****PROCEDURE:** Alkaline Phosphatase**PURPOSE:**

The ALP method used on the Dimension EXL200 integrated chemistry system is an *in vitro* diagnostic test intended for the quantitative determination of alkaline phosphatase activity in **serum** and **plasma**. Measurements of alkaline phosphatase or its isoenzymes are used in the diagnosis and treatment of liver, bone, parathyroid and intestinal disorders.

**PRINCIPLE:**

The alkaline phosphatase method is based on a procedure published by Bowers and McComb<sup>1</sup> and more recently reviewed by Rej.<sup>2</sup> This method responds to all ALP isoenzymes in human serum.

Alkaline phosphatase catalyzes the transphosphorylation of p-nitrophenylphosphate (p-NPP) to p-nitrophenol (p-NP) in the presence of the transphosphorylating buffer, 2-amino-2-methyl-1-propanol (AMP). The reaction is enhanced through the use of magnesium and zinc ions. The change in absorbance at 405 nm due to the formation of p-NP is directly proportional to the ALP activity, since other reactants are present in non-rate limiting quantities and is measured using a bichromatic (405, 510 nm) rate technique.

**INSTRUMENT, REAGENTS & SUPPLIES:****Reagents**

Wells <sup>a</sup>	Form	Ingredient	Concentration <sup>b</sup>
1-6	Liquid	2-Amino-2Methyl-1-Propanol(AMP)	3.0 M
(Reagent 1)		Magnesium Acetate	8.0 mmol/L
		Zinc Sulfate	4.0 mmol/L
		HEDTA	8.0 mmol/L
7- 8	Liquid	p-NPP Buffer	101.6 mmol/L

## (Reagent 2)

---

- a. Wells are numbered consecutively from the wide end of the cartridge.
- b. Nominal final value per well in a cartridge.

### **Risk and Safety:**

Irritant. Contains 2-amino-2-methyl-1-propanol, Sulfuric Acid, Zinc, salt, Heptahydrate

R36/38: Irritating to eyes and skin

R52/53: Harmful to aquatic organisms, may cause long-term adverse effects in the aquatic environment.

Safety Data Sheets (SDS/MSDS) can be found at [www.siemens.com/healthcare](http://www.siemens.com/healthcare)

### **Precautions:**

Used cuvettes contain human body fluids; handle with appropriate care to avoid skin contact and ingestion.

For *in vitro* diagnostic use

**Reagent Preparation:** All reagents are liquid and ready to use.

### **STORAGE & STABILITY:**

Store at 2 – 8° C

Protect from light after opening.

**Expiration:** Refer to carton for expiration date of individual unopened reagent cartridges. Sealed wells on the instrument are stable for 30 days.

**Open Well Stability:**           2days for wells 1-6  
  4 days for wells 7-8

### **SPECIMEN REQUIREMENTS:**

Normal procedures for collecting and storing serum and plasma may be used for samples to be analyzed by this method.

- patient preparation: no patient preparation required
- specimen type: serum or heparinized plasma
- handling: stable for 8 hours at room temperature, 7 days at 2-8°C, 6 months when frozen at -20°C or colder.<sup>6</sup>
- Avoid repeated freezing and thawing, thawed specimens which are turbid must be clarified by centrifugation prior to testing.

**CONTROLS:**

At least once daily run solutions at two levels of a quality control material with known concentrations.

For further details, refer to your Dimension® system manual. The result obtained should fall within limits defined by the day-to-day variability of the system as measured in the user's laboratory. If the results fall outside the laboratory's acceptable limits, follow the procedure in the quality control policy.

**PROCEDURE:**

ALPI Flex® reagent cartridge, Cat. No. DF150

**Materials Required But Not Provided**

ALPI CAL, Cat. No. DC150

Enzyme Diluent, Cat. No. 790035901

Quality Control Materials

**Test Steps**

Sampling<sup>c</sup>, reagent delivery, mixing, processing and printing of results are automatically performed by the Dimension® system. For details of this processing, refer to your Dimension® system manual.

c. The sample container (if not a primary tube) must contain sufficient quantity to accommodate the sample volume plus the dead volume. Precise container filling is not required.

**Test Conditions**

- Sample Size: 7 µL
- Reagent 1 Volume: 90 µL
- Reagent 2 Volume: 57 µL
- Diluent Volume: 206 µL
- Test Temperature: 37° C
- Reaction Time: 7.2 minutes
- Wavelength: 405 and 510 nm
- Type of Measurement: Bichromatic rate

**Backup:**

Refer to Brown Clinic Backup Policy

**Calibration:**

Assay Range (@ 37° C): 10 – 1000 U/L[0.17-16.70 µkat/L]<sup>d</sup>

Calibration Material: Dimension® ALPI CAL, Cat. No. DC150

Calibration Scheme:	Three levels (n=3)
Units:	U/L [ $\mu$ kat/L] (U/L x 0.0167)= [ $\mu$ kat/L]
Typical Calibration Levels:	Level 1: 0 U/L [0.00 $\mu$ kat/L] Level 2: 500 U/L [8.35 $\mu$ kat/L] Level 3: 1000 U/L [16.70 $\mu$ kat/L]
Calibration Frequency:	Every 90 days for any one lot
A new calibration is required:	<ul style="list-style-type: none"> <li>• For each new lot of Flex® reagent cartridges</li> <li>• After major maintenance or service, if indicated by quality control results</li> <li>• As indicated in laboratory quality control procedures</li> <li>• When required by government regulations</li> </ul>
Assigned Coefficients:	C <sub>0</sub> 0.000 C <sub>1</sub> 1.000

d. Systeme International d'Unites [SI Units] are in brackets.

### **INTERPRETATION:**

The instrument automatically calculates and prints the activity of alkaline phosphatase in U/L using the calculation scheme illustrated in your Dimension® system manual.

### **Expected Values:**

- 45-116 U/L

The reference interval was calculated non-parametrically and represents the central 95% of the population.

Each laboratory should establish its own reference interval for alkaline phosphatase as performed on the Dimension® system.

### **Reporting:**

- report format: U/L

### **LIMITATIONS:**

Results: >1000 U/L

Manual dilution: Results in excess of 1000 U/L should be repeated after diluting the sample with Enzyme Diluent (Cat. No. 790035901) or equivalent to produce a sample within the assay range. Enter dilution factor. Reassay. Resulting readout is corrected for dilution.

Autodilution (AD): Refer to your Dimension® system literature.

- Samples with results in excess of 1000 U/L[16.70 µkat/L]are reported as “Above Assay Rang” and should be repeated on dilution.
- Samples with results less than 10 U/L[0.17 µkat/L] will be reported as “less than 10U/L”
- A system malfunction may exist if the following 5-test precision is observed:

Activity	SD
150 U/L [2.51 µkat/L]	>7 U/L [0.12 µkat/L]
500 U/L [8.35 µkat/L]	>16 U/L [0.27 µkat/L]

**Analytical Specificity**

For Known Interfering Substances section refer to package insert.

For Known Non-Interfering Substance refer to package insert.

For Additional Technical Information refer to package insert.

**Reference: ALPI Flex® reagent cartridge insert sheet PN 717015.002 Issue Date 2015-01-30**

Origination Date: 9-10-07

Date of Implementation: 11-10-10

\*\*\*\*\*

Written By: \_\_\_\_\_Lori Murray MT(ASCP)\_\_\_\_\_ Date: \_8-8-12

Approved By: \_\_\_\_\_Aaron Shives, MD\_\_\_\_\_ Date: \_\_\_10/21/2017\_\_\_\_\_

Laboratory Director

**REVIEW - REVISION SUMMARY DOCUMENTATION**

Date	By	Revision Summary
07/10	Lori Murray	New format
8-8-12	Lori Murray	chg Instrumentation to EXL200
5-8-15	Sam Legg	updated information to package insert dated 09-04-2012
8-7-15	Heather Hall	Updated risk and safety
7/11/17	LMurray	ref range/ formatting
10/17/17	Heather Hall	Modified Backup Method