

Application Sheet



Laboratory Name
Test Name: Folate III

REF	Σ	SYSTEM
07559992 160	100	MODULAR ANALYTICS E170 cobas e 411 cobas e 601 cobas e 602

For USA: Elecsys Folate III Assay

Intended use

Binding assay for the in vitro quantitative determination of folate in human serum. The binding assay is intended for use on Elecsys and **cobas e** immunoassay analyzers. Folic acid measurements are used in the diagnosis and treatment of anemias.

Summary

Nutritional and macrocytic anemias can be caused by a deficiency of folate. This deficiency can result from diets devoid of raw fruits, vegetables or other foods rich in folic acid, as may be the case with chronic alcoholics, drug addicts, the elderly or persons of low socioeconomic status, etc. In addition, low serum folate during pregnancy has been associated with neural tube defects in the fetus.¹ Dietary deficiency and malabsorption are the major causes of folate deficiency in humans.² Folate is necessary for normal metabolism, DNA synthesis and red blood cell regeneration. Untreated deficiencies may lead to megaloblastic anemia.

Since a deficiency of either vitamin B₁₂ or folate can cause megaloblastic anemia, it is advisable to determine the concentration of both vitamin B₁₂ and folate in order to properly diagnose the etiology of anemia.

Radioassays were first reported for folate in 1973.^{3,4,5,6}

The majority utilize ¹²⁵I-folate radiolabeled tracers and natural binding proteins (milk binding protein, folate binding protein). The various commercial assays differ in their free versus bound separation techniques and choice of specimen pretreatment.

The Elecsys Folate assay employs a competitive test principle using natural folate binding protein (FBP) specific for folate. Folate in the sample competes with the added folate (labeled with biotin) for the binding sites on FBP (labeled with ruthenium complex^a).

a) Tris(2,2'-bipyridyl)ruthenium(II)-complex (Ru(bpy)₃²⁺)

Test principle

Competition principle. Total duration of assay: 27 minutes.

- 1st incubation: By incubating 25 µL of sample with the folate pretreatment reagents 1 and 2, bound folate is released from endogenous folate binding proteins.
- 2nd incubation: By incubating the pretreated sample with the ruthenium labeled folate binding protein, a folate complex is formed, the amount of which is dependent upon the analyte concentration in the sample.
- 3rd incubation: After addition of streptavidin-coated microparticles and folate labeled with biotin, the unbound sites of the ruthenium labeled folate binding protein become occupied, with formation of a ruthenium labeled folate binding protein-folate biotin complex. The entire complex becomes bound to the solid phase via interaction of biotin and streptavidin.
- The reaction mixture is aspirated into the measuring cell where the microparticles are magnetically captured onto the surface of the electrode. Unbound substances are then removed with ProCell/ProCell M. Application of a voltage to the electrode then induces chemiluminescent emission which is measured by a photomultiplier.
- Results are determined via a calibration curve which is instrument-specifically generated by 2-point

Laboratory Name

Test Name: Folate III

calibration and a master curve provided via the reagent barcode.

Reagents - working solutions

The reagent rackpack (M, R1, R2) and the pretreatment reagents (PT1, PT2) are labeled as FOL III.

- PT1 Pretreatment reagent 1 (white cap), 1 bottle, 4 mL: Sodium 2-mercaptoethanesulfonate (MESNA) 40 g/L, pH 5.5.
- PT2 Pretreatment reagent 2 (gray cap), 1 bottle, 5 mL: Sodium hydroxide 25 g/L.
- M Streptavidin-coated microparticles (transparent cap), 1 bottle, 6.5 mL: Streptavidin-coated microparticles 0.72 mg/mL; preservative.
- R1 Folate binding protein~Ru(bpy)₃²⁺ (gray cap), 1 bottle, 9 mL: Ruthenium labeled folate binding protein 75 µg/L; human serum albumin (stabilizer); borate/phosphate/citrate buffer 70 mmol/L, pH 5.5; preservative.
- R2 Folate~biotin (black cap), 1 bottle, 8 mL: Biotinylated folate 17 µg/L; biotin 120 µg/L; human serum albumin (stabilizer); borate buffer 100 mmol/L, pH 9.0; preservative.

Precautions and warnings

For in vitro diagnostic use.

Exercise the normal precautions required for handling all laboratory reagents.

Disposal of all waste material should be in accordance with local guidelines.

Safety data sheet available for professional user on request.

For USA: For prescription use only.

This kit contains components classified as follows in accordance with the Regulation (EC) No. 1272/2008:



Danger

H290 May be corrosive to metals.

H314 Causes severe skin burns and eye damage.

Prevention:

P280 Wear protective gloves/ protective clothing/ eye protection/ face protection.

Response:

P301 + P330
+ P331 IF SWALLOWED: Rinse mouth. Do NOT induce vomiting.

P303 + P361
+ P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower.

P304 + P340
+ P310 IF INHALED: Remove person to fresh air and keep comfortable for breathing. Immediately call a POISON CENTER or doctor/physician.

P305 + P351
+ P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P390 Absorb spillage to prevent material damage.

Product safety labeling primarily follows EU GHS guidance.

Contact phone: 1-800-428-2336

All human material should be considered potentially infectious. All products derived from human blood are prepared exclusively from the blood of donors tested individually and shown to be free from HBsAg and antibodies to HCV and HIV. The testing methods applied were FDA-approved or cleared in compliance with the European Directive 98/79/EC, Annex II, List A.

Laboratory Name

Test Name: Folate III

However, as no testing method can rule out the potential risk of infection with absolute certainty, the material should be handled with the same level of care as a patient specimen. In the event of exposure, the directives of the responsible health authorities should be followed.^{7,8}
Avoid foam formation in all reagents and sample types (specimens, calibrators and controls).

Reagent handling

The reagents in the kit have been assembled into a ready-for-use unit that cannot be separated. All information required for correct operation is read in from the respective reagent barcodes.

Storage and stability

Store at 2-8 °C.

Do not freeze.

Store the Elecsys reagent kit **upright** in order to ensure complete availability of the microparticles during automatic mixing prior to use.

Stability:	
unopened at 2-8 °C	up to the stated expiration date
after opening at 2-8 °C	8 weeks
on the analyzers	2 weeks or 4 weeks when stored alternatively in the refrigerator and on the analyzer, with the total time on-board the analyzer not exceeding 10 x 8 hours

Specimen collection and preparation

Note: Hemolysis may significantly increase folate values due to high concentrations of folate in red blood cells. Therefore, hemolyzed samples are not suitable for use in this assay. Samples for folate determinations should be collected from fasting persons.

Only the specimens listed below were tested and found acceptable.

Serum collected using standard sampling tubes or tubes containing separating gel.

Serum: Stable for 2 hours at 15-25 °C, 2 days at 2-8 °C, 4 weeks at (-15)-(-25) °C. Freeze only once.

Protect from light. Store the samples at 2-8 °C if they cannot be measured immediately.

The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.

Samples should not subsequently be altered with additives (biocides, anti-oxidants or substances possibly changing the pH of the sample) in order to avoid erroneous folate recovery.

Centrifuge samples containing precipitates before performing the assay.

Do not use heat-inactivated samples.

Ensure the samples, calibrators and controls are at 20-25 °C prior to measurement.

Due to possible evaporation effects, samples, calibrators and controls on the analyzers should be analyzed/measured within 2 hours.

Materials provided

See "Reagents – working solutions" section for reagents.

Materials required (but not provided)

-  07560001190, Folate III CalSet, for 4 x 1.0 mL
-  03183971122, Diluent Universal, 2 x 36 mL sample diluent

Laboratory Name

Test Name: Folate III

- General laboratory equipment
- MODULAR ANALYTICS E170 or **cobas e** analyzer

Accessories for MODULAR ANALYTICS E170, **cobas e** 601 and **cobas e** 602 analyzers:

- [REF] 04880340190, ProCell M, 2 x 2 L system buffer
- [REF] 04880293190, CleanCell M, 2 x 2 L measuring cell cleaning solution
- [REF] 03023141001, PC/CC-Cups, 12 cups to prewarm ProCell M and CleanCell M before use
- [REF] 03005712190, ProbeWash M, 12 x 70 mL cleaning solution for run finalization and rinsing during reagent change
- [REF] 03004899190, PreClean M, 5 x 600 mL detection cleaning solution
- [REF] 12102137001, AssayTip/AssayCup Combimagazine M, 48 magazines x 84 reaction vessels or pipette tips, waste bags
- [REF] 03023150001, WasteLiner, waste bags
- [REF] 03027651001, SysClean Adapter M

Accessories for all analyzers:

- [REF] 11298500160, ISE Cleaning Solution/Elecsys SysClean, 5 x 100 mL system cleaning solution

In addition, other suitable control material can be used.

Assay

For optimum performance of the assay follow the directions given in this document for the analyzer concerned. Refer to the appropriate operator's manual for analyzer-specific assay instructions. Resuspension of the microparticles takes place automatically prior to use. Read in the test-specific parameters via the reagent barcode. If in exceptional cases the barcode cannot be read, enter the 15-digit sequence of numbers.

MODULAR ANALYTICS E170, **cobas e** 601 and **cobas e** 602 analyzers: PreClean M solution is necessary.

Bring the cooled reagents to approximately 20 °C and place on the reagent disk (20 °C) of the analyzer. Avoid foam formation. The system automatically regulates the temperature of the reagents and the opening/closing of the bottles.

Calibration

Traceability: This method has been standardized against the WHO International Standard NIBSC code: 03/178.

Every Elecsys reagent set has a barcoded label containing specific information for calibration of the particular reagent lot. The predefined master curve is adapted to the analyzer using the relevant CalSet.

Calibration frequency: Calibration must be performed once per reagent lot using fresh reagent (i.e. not more than 24 hours since the reagent kit was registered on the analyzer). Renewed calibration is recommended as follows:

- after 1 month (28 days) when using the same reagent lot
 - after 7 days (when using the same reagent kit on the analyzer)
 - as required: e.g. quality control findings outside the defined limits
-

Quality control

At least once daily run solutions at two levels of a quality control material with known concentrations.

Refer to Brown Clinic Quality Control Requirements, Rules and Reviews Policy

Refer to Brown Clinic Quality Control Specialty and Subspecialty Policy

Calculation

The analyzer automatically calculates the analyte concentration of each sample (either in nmol/L or ng/mL).

Laboratory Name

Test Name: Folate III

Samples with folate concentrations above the measuring range can be diluted manually with Elecsys Diluent Universal. The recommended dilution is 1:2. The concentration of the diluted sample must be > 8.5 ng/mL or 19.3 nmol/L.

After manual dilution, multiply the result by the dilution factor.

Expected values

Referring to "The American Journal of Clinical Nutrition"⁹ serum folate (folic acid) values were found as follows:

Sex	Age	N	Median		2.5 th -97.5 th percentile	
	years		nmol/L	ng/mL	nmol/L	ng/mL
Both	all	23345	29.5	13.0	10.4-78.9	4.6-34.8
Male	all	11387	27.9	12.3	10.2-73.0	4.5-32.2
Female	all	11958	30.1	13.3	10.9-84.5	4.8-37.3
Both	4-11	3595	39.0	17.2	19.5-85.4	8.6-37.7
Both	12-29	6390	27.4	12.1	11.3-61.6	5.0-27.2
Both	20-59	8689	26.3	11.6	10.0-70.2	4.4-31.0
Both	≥ 60	4671	37.6	16.6	12.7-104	5.6-45.8

These values were obtained in the USA during the National Health and Nutrition Examination Survey (NHANES), 1999-2004.

The values shown below were performed on samples from an apparently healthy population, using the Elecsys Folate III assay on the **cobas e 411** analyzer.

The calculation is based on 214 sera (110 men, 104 women). The age range was between 21 and 59 years. Pregnant or lactating women were excluded. The reference population was selected according to normal homocysteine values.

Country	N	Median		2.5 th -97.5 th percentile	
		nmol/L	ng/mL	nmol/L	ng/mL
USA	214	26.8	11.8	10.9 - 54.9	4.78 - 24.2

Please note: These values should only be used as a guideline.

It should be taken into consideration that differences in the expected values may exist with respect to population and dietary status.

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference ranges.

Specific performance data

For Known Interfering Substances section refer to package insert.

For Known Non-Interfering Substance refer to package insert.

For Additional Technical Information refer to package insert.

References

- 1 Rush D. Folate Supplements Prevent Recurrence of Neural Tube Defects, FDA Dietary Supplement Task Force. Nutrition Reviews 1992;50(1):22-28.
- 2 Herbert V. Drugs effective in megaloblastic anemias. In Goodman LS and Gilman A (eds): The Pharmacological Basis of Therapeutics, 5th Ed, MacMillan Co, 1975;1324-1349.
- 3 Dunn RT, Foster LB. Radioassay of serum Folate. Clin Chem 1973;19:1101-1105.

Laboratory Name
Test Name: Folate III

- 4 Rothenberg SP, DaCosta M, Rosenberg BS. A radioassay for serum Folate: Use of a two phase sequential incubation, ligand-binding system. *New Eng J Med* 1972;285(25):1335-1339.
- 5 Gutcho S, Mansbach L. Simultaneous radioassay of serum Folate and folic acid. *Clin Chem* 1977;23:1609-1614.
- 6 BIO RAD Quantaphase B-12/Folate Radioassay Instruction Manual. March 1995.
- 7 Occupational Safety and Health Standards: bloodborne pathogens. (29 CFR Part 1910.1030). Fed. Register.
- 8 Directive 2000/54/EC of the European Parliament and Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work.
- 9 Pfeiffer CM, Johnson CL, Jain RB, et al. Trends in blood folate and vitamin B-12 concentrations in the United States, 1988-2004. *Am J Clin Nutr* 2007;86:718-727.
- 10 Passing H, Bablok W, Bender R, et al. A general regression procedure for method transformation. *J Clin Chem Clin Biochem* 1988 Nov;26(11):783-790.

Alternative method

Refer to Brown Clinic Back-up Testing Policy

Source document

Reagent Name: Folate III
Method Sheet Version: V1.0 English

Effective date

Effective date for this procedure:

Author

Source documentation compiled by Roche Diagnostics

Revised by: Heather J Hall, MBA, MT(ASCP), CG(ASCP)^{cm} Date: 4/9/2018

Approved by: Aaron Shives MD (Signature on file) Date: 4/11/2018

REVIEW – REVISION SUMMARY DOCUMENTATION

Date: _____ By: _____ Revision Summary: _____