

# STANTON TERRITORIAL HEALTH AUTHORITY

## Yellowknife, Northwest Territories

<b>TITLE:</b> Nova Biomedical Glucose / Ketone Meter Linearity Test	<b>Revision Date:</b> 07-July-2016	<b>Issue Date:</b> 07-July-2014
<b>Document Number:</b> POC20601	<b>Status:</b> <b>Approved</b>	
<b>Distribution:</b> Point of Care Testing Manual	<b>Page:</b> 1 of 5	
<b>Approved by:</b> C. Case, Manager Diagnostic Services	<b>Signed by:</b> <i>Cheryl Case</i>	

### **PURPOSE:**

Linearity testing is performed on a bi-annual basis to ensure that results generated are accurate and precise by testing samples with known analyte concentrations and evaluating the line of best fit when those results are plotted on a line graph.

### **POLICY:**

See Quality Manual Section 6. Process Control

### **REAGENTS and/or MEDIA:**

- Nova Biomedical StatStrip® GLU Test Strips
- Nova Biomedical StatStrip® KET Test Strips
- Nova Biomedical StatStrip® GLU KET Linearity Kit

### **SUPPLIES:**

- Nova Biomedical StatStrip® Glucose / Ketone Hospital Meter
- OR**
- Nova Biomedical StatStrip® Glucose Hospital Meter

### **SPECIAL SAFETY PRECAUTIONS:**

- Handle all patient samples and testing sera using "Routine Practices".
- Please refer to the Northwest Territories Infection Prevention and Control Manual, March 2012.

<b>NOTE:</b> This is a CONTROLLED document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.	
<b>FILENAME:</b> POC20601NovaBiomedicalGlucoseKetoneMeterLinearityPRO.doc	<b>PRINT DATE:</b> 7/7/2014 4:43 PM



<b>TITLE:</b> Nova Biomedical Glucose / Ketone Meter Linearity Test	<b>Revision Date:</b> 07-July-2016	<b>Issue Date:</b> 07-July-2014
<b>Document Number:</b> POC20601	<b>Status:</b> <b>Approved</b>	
<b>Distribution:</b> Point of Care Testing Manual	<b>Page:</b> 2 of 5	

**QUALITY CONTROL:**


See POC20701 Nova Biomedical Glucose / Ketone Meter Quality Control Test

**PROCEDURE INSTRUCTIONS:**



Follow the steps in the table below to perform a linearity test.

Step	Action
1	<p>Turn on by pressing the <b>Sleep Mode</b> button.</p> 
2	Press the <b>OK/Enter</b> button to begin testing.
3	Enter your Operator ID number manually or scan your barcode.
4	Press the <b>OK/Enter</b> button to accept your Operator ID number.
5	<p>Touch the <b>Menu</b> soft key on the instrument screen.</p> 

<b>TITLE:</b> Nova Biomedical Glucose / Ketone Meter Linearity Test	<b>Revision Date:</b> 07-July-2016	<b>Issue Date:</b> 07-July-2014
<b>Document Number:</b> POC20601	<b>Status:</b> <b>Approved</b>	
<b>Distribution:</b> Point of Care Testing Manual	<b>Page:</b> 3 of 5	

6	<p>Touch the <b>Linearity</b> soft key on the instrument screen.</p>  <p>The screenshot shows a menu screen titled 'Menu' with the following options: Patient Test, Manual Entry, QC, Proficiency, Linearity, Correlation, Maintenance, Training, and Admin. A red arrow points to the 'Linearity' button. At the bottom, there are buttons for 'Review', 'Welcome', and 'Logout'.</p>
7	Scan the barcode on the side of the vial for the lot number of the test strips you are using. Press <b>OK/Enter</b> .
8	Linearity reagents are received in a kit from Nova Biomedical. These reagents are ordered directly from the manufacturer. The same linearity kit is used for both whole blood glucose and ketone linearity testing. Scan the barcode of the linearity solution you will be testing by pressing and holding the <b>OK/Enter</b> button and passing the scanner over the barcode on the side of the dropper bottle.
9	Press the <b>OK/Enter</b> button to accept the proficiency test reagent number.

<b>TITLE:</b> Nova Biomedical Glucose / Ketone Meter Linearity Test	<b>Revision Date:</b> 07-July-2016	<b>Issue Date:</b> 07-July-2014
<b>Document Number:</b> POC20601	<b>Status:</b> <b>Approved</b>	
<b>Distribution:</b> Point of Care Testing Manual	<b>Page:</b> 4 of 5	

<b>10</b>	<p>Insert a test strip into the top of the testing meter as shown on the screen.</p> 
<b>11</b>	Gently mix the Linearity Solution before each use
<b>12</b>	<p>With the test strip correctly inserted, the Apply Sample screen will display.</p> 
<b>13</b>	Discard the first drop of control solution from the bottle to avoid contamination.
<b>14</b>	With the meter placed on a flat level surface, place a drop of linearity solution from the bottle at the end of the test strip until the solution is drawn into the well of the test strip. When sample has been drawn into the strip, an audible beep is sounded by the meter.
<b>15</b>	Wipe the applicator tip of the bottle and recap the linearity solution. The Testing Sample screen displays. The screen shows a clock with seconds remaining below the clock.

<b>TITLE:</b> Nova Biomedical Glucose / Ketone Meter Linearity Test	<b>Revision Date:</b> 07-July-2016	<b>Issue Date:</b> 07-July-2014
<b>Document Number: POC20601</b>	<b>Status: <span style="color: red;">Approved</span></b>	
<b>Distribution: Point of Care Testing Manual</b>	<b>Page: 5 of 5</b>	

<b>16</b>	The test result will appear in 6 or 10 seconds depending on which test is being performed. Do not remove the test strip until the countdown is complete. Meter results range from 0.6 – 33.3 mmol/L for glucose and <0.1 to 7.0 mmol/L for Ketones.
<b>17</b>	Run each level in duplicate.
<b>18</b>	Repeat steps 8 to 17 for the all other levels in the linearity kit.
<b>19</b>	Dock the meter into a connected docking station to transmit the result to the NovaNet. Results will be compiled and analyzed by the Technologist II.

**RELATED DOCUMENTS:**

- POC20801 Nova Biomedical Glucose / Ketone Meter Patient Test
- POC20401 Nova Biomedical Glucometer Basic Maintenance

**REFERENCES:**

Nova Biomedical, StatStrip® Glucose Hospital Meter, Instructions for Use Manual, Nova Biomedical Corporation, Waltham, MA, 2011

Nova Biomedical, StatStrip™ Glucose and β-Ketone Control Solution Package Insert, Nova Biomedical UK, Cheshire, UK, 2011

**REVISION HISTORY:**

REVISION	DATE	Description of Change	REQUESTED BY
1.0	23 Dec 2013	Initial Release	JGD Bernier
2.0	15 May 2014	Updated to include reference to: NWT Infection Prevention & Control Manual	JGD Bernier

<b>NOTE: This is a CONTROLLED document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.</b>	
<b>FILENAME: POC20601NovaBiomedicalGlucoseKetoneMeterLinearityPRO.doc</b>	<b>PRINT DATE: 7/7/2014 4:43 PM</b>