

Impact to Results on VITROS® 3600/5600 Systems

Testing was conducted for assays that are sensitive to well wash fluid temperature utilizing the MicroWell Wash module. Samples were tested at nominal temperature (37°C) and at ambient room temperature (~25°C). The data provided below summarizes our assessment.

NOTE: Seventeen assays were tested across multiple concentrations; the values deemed to affect a clinical decision are shown below. Data for other assays will be available upon request.

For assays that have cutoff values, only those samples with results close to the cutoff values may potentially lead to a misclassification.

| Potential Bias if Assays Calibrated at Nominal Temperature (37°C), Samples Processed at Ambient Temperature (~25°C) | | |
|--|---------------------------------------|---------------|
| Assay Name | Concentration Range | Bias Observed |
| NT proBNP | Cutoff value: 300 pg/mL | 19.5% |
| TSH | Lower Reference Range: 0.465 mIU/L | 6% |
| Intact PTH | Lower Reference Range: 7.5pg/mL | 18% |
| Rubella IgG | Cutoff value for Negative: 9.99 IU/mL | -2.8% |
| Anti-HBs | Cutoff value for Positive: 10 mIU/mL | 26% |
| Potential Bias if Assays Calibrated at Ambient Temperature (~25°C), Samples Processed at Nominal Temperature (37°C) | | |
| Assay Name | Concentration Range | Bias Observed |
| NT proBNP | Cutoff value: 450 pg/mL | -16% |
| TSH | Upper Reference Range: 4.68 mIU/L | -1.6% |
| Intact PTH | Upper Reference Range: 53.5 pg/mL | -13% |
| Troponin | 0.034 ng/mL | -8.8% |
| Anti-HBc | Cutoff value for Negative: 0.9 | -10% |
| HBsAg | Cutoff value for Negative: 0.9 | -2.9% |
| Rubella IgM | Cutoff value for Negative: 0.8 | -2.8% |
| Syphilis | Cutoff value for Negative: 0.8 | -6.6% |
| Toxoplasma IgG | Cutoff value for Negative: 3.99 | -4.8% |
| Toxoplasma IgM | Cutoff value for Negative: 0.8 | -5.5% |
| HAV Total | Cutoff for Positive: 0.8 | -6.3% |

NOTE: Testing is ongoing for other MicroWell assays and will be available upon request by contacting our Technical Solutions Center representatives.