

**Stanton Territorial Hospital**P.O. Box 10, 550 Byrne Road  
YELLOWKNIFE NT X1A 2N1

Document Number: TMM80700

Version No: 1.0

Page: 1 of 10

**Distribution:**

Transfusion Medicine Manual

Effective: 03 January, 2018

Date Reviewed: 03 January, 2018

Next Review: 03 January, 2020

**Document Name:**

Unit Retype Testing

**Approved By:**

Jennifer G. Daley Bernier, A/ Manager, Laboratory Services

Status: **APPROVED****PURPOSE:**

To perform and result retype testing on Units of Red Blood Cells received from Canadian Blood Services (CBS).

**POLICY:**

Red Blood Cell units must have the ABO & Rh status confirmed upon receipt from CBS. The ABO and Rh status must be recorded as tested in the Laboratory Information System (LIS).

**SAMPLE INFORMATION:**

<b>Type</b>	Blood
<b>Source</b>	Segment removed from Red Blood Cell Unit
<b>Volume</b>	1 mL or less
<b>Stability</b>	42 days from collection date on unit
<b>Storage Requirements</b>	2-8°C
<b>Criteria for rejection and follow up action</b>	Unit does not pass Visual Inspection

**REAGENTS and/or MEDIA:**

- Anti-A and Anti-B (Ortho Clinical Diagnostics Murine Monoclonal Blend) Reagent
- Anti-D (RH1) IgM + IgG BioClone Reagent
- Rhesus Reagent Control BioClone Reagent
- Isotonic Saline

**SUPPLIES:**

- Serological Centrifuge
- Work block(s)
- Transfer pipettes

**NOTE:** This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

FILENAME: TMM80700.1UnitRetypeTestingPRO.docx

Print Date: 03/01/2018 7:25:00 PM

- Tubes (12x75mm)
- Marker
- Blue Segment Piercer or scissors

**EQUIPMENT CALIBRATION AND MAINTENANCE:**

- See BLB70300 - Equipment Maintenance & Alarm System Checks for Blood Bank
- See BLB70700 - Calibration of Incubators, Pipettes and Thermometers

**SPECIAL SAFETY PRECAUTIONS:**

- Handle all products and samples as potentially infectious.
- Some products may contain hazardous chemicals such as sodium azide (Anti-A, Anti-B) – must be disposed of with flushing of large amounts of water if disposed down a drain.






**QUALITY CONTROL:**

- See TMM\*\*\*\*\* – Reagent QC (to be written)





**PROCEDURE INSTRUCTIONS:**

Follow the steps in the table below to perform and result unit retype testing:



<b>Step</b>	<b>Action</b>	
<b>1</b>	Retype testing is ordered automatically on units that are received into the LIS	
<b>2</b>	<b>If:</b>	<b>Then:</b>
	Units are received using the Batch Delivery Method	Unit Retype worksheets are automatically built for test resulting – Go to Step 12 for Performing testing and Entering results

	<p>Units are received using the Delivery Method</p>	<p>Unit Retype worksheet must be built prior to test resulting – Go to Step 3 below</p>														
<p>3</p>	<p>Go to 4-Results &gt; Unit Test Worksheets</p> 															
<p>4</p>	<p>Click: Build</p> 															
<p>5</p>	<p>Use drop down menu to select testing for the worksheet – Select RETYP</p> 															
<p>6</p>	<p>Click: F12-Accept</p>															
<p>7</p>	<p>List of units with un-resulted Retype testing ordered will appear.</p>  <table border="1" data-bbox="305 1381 1474 1535"> <thead> <tr> <th>Unit#</th> <th>Prod.</th> <th>ABO Rh</th> <th>Order#</th> <th>Test</th> <th>Spec.date,time</th> <th>Priority</th> </tr> </thead> <tbody> <tr> <td>C052117436426</td> <td>RCE</td> <td>O NEG</td> <td>C007001I</td> <td>RETYP</td> <td>17-12-07 08:15</td> <td>routine</td> </tr> </tbody> </table>		Unit#	Prod.	ABO Rh	Order#	Test	Spec.date,time	Priority	C052117436426	RCE	O NEG	C007001I	RETYP	17-12-07 08:15	routine
Unit#	Prod.	ABO Rh	Order#	Test	Spec.date,time	Priority										
C052117436426	RCE	O NEG	C007001I	RETYP	17-12-07 08:15	routine										
<p>8</p>	<p>Select units to retype using the F5-Mark all option or Enter to select individually A check mark will appear next to selected units</p> 															
<p>9</p>	<p>Click: F12-Accept</p>															




**NOTE:** This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.







	If:	Then:
10		Yes
11	Unit retype worksheet is now built.	
12	Go to Enter Results 	
13	If not already selected use the drop down menu to select RETYP in worksheet field	
14	Click F12-Accept	
	If:	Then:
15		Click: Ok
16	Retype Resulting worklist will appear 	
17	Perform ABO and RH retype testing using a segment from the unit.	
18	Label all tubes per Procedure: <b>TMM***** – Labeling Tubes, Gel Cards and Block Set-up</b>	
19	Warm reagents and segment to room temperature: 15-30°C	

**NOTE:** This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

20	<p>Make a 3% red cell suspension from the red cells of the segment and add to the appropriate labelled tubes</p> <ul style="list-style-type: none"> <li>• 1 drop to tube labelled A</li> <li>• 1 drop to tube labelled B</li> <li>• 1 drop to tube labelled D</li> <li>• 1 drop to tube labelled DC</li> </ul>																		
21	Mix the tubes by gentle agitation																		
22	Centrifuge tubes for 15-20 seconds at 3400 ± 200 RPM (high setting on serologic centrifuge)																		
23	Remove tubes from centrifuge.																		
24	<p>Re-suspend the cells and read tubes macroscopically. See <b>Procedure TMM***** - Reading and Recording Hemagglutination Reactions: Tube and Gel (to be written)</b></p> <p>Watch for agglutination or hemolysis</p>																		
25	<p>Enter the result using drop down menu, or type in results into the appropriate testing fields.</p> <p>Note: There is no field to result testing for Rh Control – it is not required to be recorded for retype testing</p>  <table border="1" data-bbox="316 1411 1161 1501"> <thead> <tr> <th>Unit</th> <th>ABO</th> <th>A</th> <th>B</th> <th>D</th> <th>ABO</th> <th>RH</th> </tr> </thead> <tbody> <tr> <td>C052117436426 RCE</td> <td>O-</td> <td>0</td> <td>0</td> <td>0</td> <td>0</td> <td>NEG</td> </tr> </tbody> </table>	Unit	ABO	A	B	D	ABO	RH	C052117436426 RCE	O-	0	0	0	0	NEG				
Unit	ABO	A	B	D	ABO	RH													
C052117436426 RCE	O-	0	0	0	0	NEG													
26	ABO and Rh status will automatically populate based on results entered under the A, B and D testing result fields																		
27	Click: F12-Accept																		
28	<p>Batch results pop up window will appear. Ensure status is C for Complete</p>  <table border="1" data-bbox="316 1789 1474 1858"> <thead> <tr> <th>Test</th> <th>Int. 1</th> <th>Int. 2</th> <th>Stat</th> <th>Cm</th> <th>Except.</th> <th>Ex</th> <th>Add. test 1</th> <th>Add. test 2</th> </tr> </thead> <tbody> <tr> <td>C052117436426 RCE</td> <td>RETYP</td> <td>0</td> <td>NEG</td> <td></td> <td></td> <td></td> <td></td> <td></td> </tr> </tbody> </table>	Test	Int. 1	Int. 2	Stat	Cm	Except.	Ex	Add. test 1	Add. test 2	C052117436426 RCE	RETYP	0	NEG					
Test	Int. 1	Int. 2	Stat	Cm	Except.	Ex	Add. test 1	Add. test 2											
C052117436426 RCE	RETYP	0	NEG																

**NOTE:** This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

29	Click: F12-Accept	
30	<p style="text-align: center;"><b>If:</b></p> 	<p style="text-align: center;"><b>Then:</b></p> <ul style="list-style-type: none"> <li>• Yes – will confirm testing – and change unit status to Available in Inventory</li> <li>• No – test status will remain in process – unit will not be available for use</li> </ul>
31	<p>Option to print Unit Confirmation Label appears – Select printer to print – hit enter or doubleclick on line</p> 	
32	Place printed unit confirmation label on tested units – or place handwritten ABO & Rh confirmation label on tested units if no confirmation label is printed – See examples in Appendix A below.	
33	Click: F12-Accept	
34	<p>The worklist now must be Cleared – Go to Clear</p> 	

<p>35</p>	<p>If retyp work sheet is not already selected in the worksheet field – use the drop down menu or type in RETYP</p> 				
<p>36</p>	<p>F12-Accept</p>				
<p>37</p>	<p>The worklist will appear</p> 				
<p>38</p>	<p>Units to be cleared from the worklist require a checkmark in the box to the left of the Unit field. To select and put a checkmark click either F3-Mark/unmark for single unit or F4-All to select all units on list</p> 				
<p>39</p>	<p>F12-Accept</p>				
<p>40</p>	<table border="1" style="width: 100%;"> <tr> <th style="width: 50%; text-align: center;">If:</th> <th style="width: 50%; text-align: center;">Then:</th> </tr> <tr> <td data-bbox="310 1472 792 1703">  </td> <td data-bbox="878 1472 1461 1719"> <p>Pop up will show how many lines or units will be deleted from the worklist. Click Yes to confirm deletion</p> </td> </tr> </table>	If:	Then:		<p>Pop up will show how many lines or units will be deleted from the worklist. Click Yes to confirm deletion</p>
If:	Then:				
	<p>Pop up will show how many lines or units will be deleted from the worklist. Click Yes to confirm deletion</p>				
<p>41</p>	<p>Worksheet will be cleared of any check marked units.</p>				

<b>Document Name: Unit Retype Testing</b>	<b>Document Number: TMM80700</b>	
	<b>Version No: 1.0</b>	<b>Page: 8 of 10</b>
	<b>Effective: 03 January, 2018</b>	

**EXPECTED RESULTS:**

- Unit is retyped and recorded as such in the LIS.
- Unit status will now be set to Available.
- Unit can now be used in crossmatch testing for patient transfusions.

**RELATED DOCUMENTS:**

- TMM80300 - Receiving Blood & Blood Products from CBS at Stanton Territorial Hospital
- TMM\*\*\*\*\* – Labeling Tubes, Gel Cards and Block Set-up (to be written)
- TMM\*\*\*\*\* – ABO Testing (to be written)
- TMM\*\*\*\*\* – Rh Testing (to be written)
- TMM\*\*\*\*\* - Visual Inspection of Blood, Blood Components and Derivatives.
- TMM80500 - Creating Unit Segments
- TMM\*\*\*\*\* - Reading and Recording Hemagglutination Reactions: Tube and Gel (to be written)

**REFERENCES:**

1. Canadian Society of Transfusion Medicine. (April 2017). *Standards for Hospital Transfusion Services, Version 4.*
2. Canadian Standards Association. (February 2016). *Blood and Blood Components CAN/CSA-Z902-15.*
3. SCC Soft Computer. (n.d.). *SoftBank II v25 Super User Training Manual.* SCC Soft Computer.

**APPENDIX (APPENDICIES):**



- Appendix A: Examples of possible unit confirmation labels



**REVISION HISTORY:**

REVISION	DATE	Description of Change	REQUESTED BY
1.0	03 Jan 2018	Initial Release	A. Richardson

Appendix A: Examples of possible unit confirmation labels

LIS printed Unit Confirmation Label Example	Handwritten Confirmation Label Example
 <p>ABO GROUP CONFIRMED Unit#: C052117436426 RCE <b>O NEG</b> Tested: 17-12-07 14:32by AER TEST2 - Stanton Territorial Hospital Laboratory</p>	 <p>ABO &amp; Rh CONFIRMED <u>7 Dec 17 Test</u> Date/Tech.</p>

**NOTE:** This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.