

Document Name: Gram stain resulting in LIS –  
Bacterial Vaginosis

Approved By:

Status: **DRAFT**

**PURPOSE:** To report the Gram stain results from Bacterial Vaginosis screens in a consistent manner.

**SAMPLE INFORMATION:**

<b>Type</b>	<ul style="list-style-type: none"> <li>• Vagina.</li> <li>• Refer to MIC10230 – Microbiology Specimen Processing.</li> </ul>
-------------	--

**REAGENTS INFORMATION:**

<b>Type</b>	BD™ Gram Crystal Violet, 3.8 L, B4312526 BD™ Gram Iodine (Stabilized), 3.8 L, B4312543 BD™ Gram Decolorizer, 3.8 L, B4312528 BD™ Gram Safranin, 3.8 L, B4312531
<b>Source</b>	Fisher Scientific Canada
<b>Storage</b>	Store at 15° to 30°
<b>Stability</b>	As per expiry date listed on bottle

**SUPPLIES:**

- Frosted end glass microscope slide
- QC slide
- Methanol, absolute
- Immersion oil
- Microscope
- Slide storage tray

**NOTE:** This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

FILENAME:

Print Date:

<b>Document Name: Gram stain resulting in LIS – Bacterial Vaginosis</b>	<b>Document Number: MIC20600</b>	
	<b>Version No: 1.0</b>	Page: 2 of 9
	<b>Effective: DRAFT</b>	

**SPECIAL SAFETY PRECAUTIONS:**

Containment Level 2 facilities, equipment, and operational practices for work involving infectious or potential infectious materials or cultures.

- Lab gown must be worn when performing activities with potential pathogens.
- Gloves must be worn when direct skin contact with infected materials is unavoidable.
- Eye protection must be used when there is a known or potential risk of exposure of splashes.
- All procedures that may produce aerosols, or involve high concentrations or large volumes should be conducted in a biological safety cabinet (BSC).
- The use of needles, syringes and other sharp objects should be strictly limited.


All patient specimens are assumed to be potentially infectious. Universal precautions must be followed. Since viable micro-organisms are used, all cultures must be handled with appropriate precautions. All equipment in contact with cultures should be decontaminated by appropriate methods.

**QUALITY CONTROL:**

- Refer to MIC60060 – Microbiology Stain Quality Control.

<b>NOTE:</b> This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.	
FILENAME:	Print Date:

**PROCEDURE INSTRUCTIONS:**

Step	Action										
1	<p>Pending BV Gram stain orders are found in the LIS Resulting Worklist:  <b>Resulting Worklist → GRM/MISC_1of2</b></p> 										
2	<p>Enter the accession number on the slide and select enter to mark the order.                      Select enter again to open Result Entry or double click on accession number to open.</p>										
3	<p>Under the media resulting log: <b>GMBV</b>, use the <b>GMBV</b> keypad to record the quantity of epithelial cells, white blood cells, clue cells, yeast cells, Trichomonas and bacteria significant to Bacterial Vaginosis: <i>Lactobacillus</i> spp., Gram variable bacilli resembling <i>Gardnerella vaginalis</i> and Gram negative bacilli resembling <i>Mobiluncus</i> spp.</p>										
4	<p>Epithelial cells, white blood cells and red blood cells are quantified as follows under LPF (10X):</p> <table border="1" data-bbox="513 1472 1198 1791"> <tbody> <tr> <td>None seen</td> <td>No cells seen</td> </tr> <tr> <td>1+</td> <td>&lt; 1 cells seen</td> </tr> <tr> <td>2+</td> <td>1 – 9 cells seen</td> </tr> <tr> <td>3+</td> <td>10 – 25 cells seen</td> </tr> <tr> <td>4+</td> <td>&gt; 25 cells seen</td> </tr> </tbody> </table> <p><b>NOTE:</b> Only record “None seen” for WBC.</p>	None seen	No cells seen	1+	< 1 cells seen	2+	1 – 9 cells seen	3+	10 – 25 cells seen	4+	> 25 cells seen
None seen	No cells seen										
1+	< 1 cells seen										
2+	1 – 9 cells seen										
3+	10 – 25 cells seen										
4+	> 25 cells seen										

**NOTE:** This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.



<b>5</b>	Add one drop of immersion oil to the slide. Using the oil immersion lens (100X); examine 20 to 40 fields to observe cell morphology and Gram reaction.																																			
<b>6</b>	<p>Bacterial and yeast cells are quantified as follows under OIF (100x):</p> <table border="1" style="margin-left: auto; margin-right: auto;"> <tr> <td style="text-align: center;">None seen</td> <td style="text-align: center;">No cells seen</td> </tr> <tr> <td style="text-align: center;">1+</td> <td style="text-align: center;">&lt; 1 cells seen</td> </tr> <tr> <td style="text-align: center;">2+</td> <td style="text-align: center;">1 – 9 cells seen</td> </tr> <tr> <td style="text-align: center;">3+</td> <td style="text-align: center;">10 – 25 cells seen</td> </tr> <tr> <td style="text-align: center;">4+</td> <td style="text-align: center;">&gt; 25 cells seen</td> </tr> </table> <p><b>NOTE:</b> If no bacteria are seen, record this result.</p>	None seen	No cells seen	1+	< 1 cells seen	2+	1 – 9 cells seen	3+	10 – 25 cells seen	4+	> 25 cells seen																									
None seen	No cells seen																																			
1+	< 1 cells seen																																			
2+	1 – 9 cells seen																																			
3+	10 – 25 cells seen																																			
4+	> 25 cells seen																																			
<b>7</b>	<p>Calculate the Nugent score of the specimen. The N-score is calculated by assessing the quantity of <i>Lactobacillus</i> spp., Gram-variable bacilli and <i>Mobiluncus</i> spp. Examine 20 to 40 fields and interpret as follows:</p> <table border="1" style="margin-left: auto; margin-right: auto;"> <thead> <tr> <th colspan="6" style="text-align: center;"><b>N-Score Table</b></th> </tr> <tr> <th rowspan="2" style="text-align: center;">Quantitation of bacterial morphotype</th> <th colspan="5" style="text-align: center;">Points scored per morphotype</th> </tr> <tr> <th style="text-align: center;">None</th> <th style="text-align: center;">1+</th> <th style="text-align: center;">2+</th> <th style="text-align: center;">3+</th> <th style="text-align: center;">4+</th> </tr> </thead> <tbody> <tr> <td style="text-align: center;"><i>Lactobacillus</i> spp.: Med/large Gram positive bacilli</td> <td style="text-align: center;">4</td> <td style="text-align: center;">3</td> <td style="text-align: center;">2</td> <td style="text-align: center;">1</td> <td style="text-align: center;">0</td> </tr> <tr> <td style="text-align: center;"><i>Gardnerella vaginalis</i>: Small, Gram variable bacilli</td> <td style="text-align: center;">0</td> <td style="text-align: center;">1</td> <td style="text-align: center;">2</td> <td style="text-align: center;">3</td> <td style="text-align: center;">4</td> </tr> <tr> <td style="text-align: center;"><i>Mobiluncus</i> spp.: Curved, Gram variable bacilli</td> <td style="text-align: center;">0</td> <td style="text-align: center;">1</td> <td style="text-align: center;">1</td> <td style="text-align: center;">2</td> <td style="text-align: center;">2</td> </tr> </tbody> </table>	<b>N-Score Table</b>						Quantitation of bacterial morphotype	Points scored per morphotype					None	1+	2+	3+	4+	<i>Lactobacillus</i> spp.: Med/large Gram positive bacilli	4	3	2	1	0	<i>Gardnerella vaginalis</i> : Small, Gram variable bacilli	0	1	2	3	4	<i>Mobiluncus</i> spp.: Curved, Gram variable bacilli	0	1	1	2	2
<b>N-Score Table</b>																																				
Quantitation of bacterial morphotype	Points scored per morphotype																																			
	None	1+	2+	3+	4+																															
<i>Lactobacillus</i> spp.: Med/large Gram positive bacilli	4	3	2	1	0																															
<i>Gardnerella vaginalis</i> : Small, Gram variable bacilli	0	1	2	3	4																															
<i>Mobiluncus</i> spp.: Curved, Gram variable bacilli	0	1	1	2	2																															
<b>8</b>	Finalize <b>STGM1</b> . Preview instant report and save. Refresh <b>GRM/MISC1of2</b> worklist. If finished reading slides, ensure Gram stains remaining on worklist have been prepared to be read at a later time.																																			
<b>9</b>	Gently blot excess oil from slide using paper towel or gauze and save slides in slide box marked BV.																																			

**NOTE:** This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

FILENAME:

Print Date:



**REPORTING OF RESULTS:**

IF	REPORT	
<b>N-Score is 0 - 3</b>	<ul style="list-style-type: none"> <li>In SCBV test, using the SCBV keypad, select Key 1 to state:  <b>“Smear NEGATIVE for Bacterial Vaginosis &gt;SCBV2”</b>.                             <div style="text-align: center; margin: 10px 0;">  </div> </li> <li>Comment will appear as:  <b>“Smear NEGATIVE for Bacterial Vaginosis”</b> </li> <li>Add comments, if appropriate, from SCBV2 keypad:                             <div style="text-align: center; margin: 10px 0;">  </div> </li> </ul>	
	IF	ADD COMMENT
	Patient is > 55 years Key 1	<b>“Results may not be reliable in post-menopausal women. Correlate with the clinical picture”</b> .
	Yeast cells seen Key 2	<b>“Yeast Seen. Candida species are normal flora in the genital area of 30-40% of women. The presence of yeast must be correlated with clinical picture”</b>
	3+ or 4+ WBC seen Key 3	<b>“Purulence suggests the presence of another infection and/or inflammatory condition. Correlate with the clinical picture. Testing for <i>N. gonorrhoeae</i> and <i>C. trachomatis</i> may be indicated”</b> .
Trichomonas seen Key 4	<b>“Trichomonas vaginalis seen in smear”</b>	

**NOTE:** This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

FILENAME:



Print Date:

IF	REPORT	
<b>N-Score is 4 - 6</b>	<ul style="list-style-type: none"> <li>In SCBV test, using the SCBV keypad, select Key 2 to state: <b>“}ALTF (altered indeterminate) &gt;SCBV2”</b>.</li> </ul>  <ul style="list-style-type: none"> <li>Comment will appear as: <b>“Gram stain shows altered vaginal flora. Results are indeterminate for Bacterial Vaginosis”</b>.</li> <li>Add comments, if appropriate, from SCBV2 keypad:</li> </ul> 	
	IF	ADD COMMENT
	Patient is > 55 years Key 1	<b>“Results may not be reliable in post-menopausal women. Correlate with the clinical picture”</b> .
	Yeast cells seen Key 2	<b>“Yeast Seen. Candida species are normal flora in the genital area of 30-40% of women. The presence of yeast must be correlated with clinical picture”</b>
	3+ or 4+ WBC seen Key 3	<b>“Purulence suggests the presence of another infection and/or inflammatory condition. Correlate with the clinical picture. Testing for <i>N. gonorrhoeae</i> and <i>C. trachomatis</i> may be indicated”</b> .
Trichomonas seen Key 4	<b>“Trichomonas vaginalis seen in smear”</b>	

**NOTE:** This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

FILENAME:


Print Date:

IF	REPORT		
<b>N-Score is 7 - 10</b>	<ul style="list-style-type: none"> <li>In SCBV test, using the SCBV keypad, select Key 3 to state:  <b>“Smear consistent with Bacterial Vaginosis &gt;SCBV2”</b></li> </ul> 		
	<ul style="list-style-type: none"> <li>Comment will appear as:  <b>“Smear consistent with Bacterial Vaginosis”</b></li> </ul>		
	<ul style="list-style-type: none"> <li>Add comments, if appropriate, from SCBV2 keypad:</li> </ul> 		
	IF	ADD COMMENT	
	Patient is > 55 years Key 1	<b>“Results may not be reliable in post-menopausal women. Correlate with the clinical picture”.</b>	
Yeast cells seen Key 2	<b>“Yeast Seen. Candida species are normal flora in the genital area of 30-40% of women. The presence of yeast must be correlated with clinical picture”</b>		
3+ or 4+ WBC seen Key 3	<b>“Purulence suggests the presence of another infection and/or inflammatory condition. Correlate with the clinical picture. Testing for <i>N. gonorrhoeae</i> and <i>C. trachomatis</i> may be indicated”.</b>		
Trichomonas seen Key 4	<b>“Trichomonas vaginalis seen in smear”</b>		

**NOTE:** This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

FILENAME:

Print Date:

IF	REPORT
<p><b>Scant or no cells seen on smear</b></p>	<ul style="list-style-type: none"> <li>In SCBV test, using the SCBV keypad, select Key 4 to state: <b>"}INSF (scant or no cells)"</b></li> </ul>  <ul style="list-style-type: none"> <li>Comment will appear as: <b>"Insufficient sample to assess for vaginitis. Please recollect if clinically indicated"</b></li> <li>Finalize order.</li> </ul>

**LIMITATIONS:**

- The presence or absence of clue cells is not part of the Nugent score and not required for diagnosis.
- For post-menopausal patients, laboratory diagnosis of bacterial vaginosis has not been validated and interpretation of Gram stain results needs to be taken into account. Ensure comment is added.
- For pre-pubescent girls (< 13 years), Bacterial Vaginosis should not be reported. Genital culture should be performed and Gram stain should be reported as per routine specimens. Refer to MIC MIC20200 - Gram stain resulting in LIS – Routine specimens.

**NOTE:** This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

FILENAME:

Print Date:



**REFERENCES:**

- Clinical Microbiology Procedures Handbook, 4<sup>th</sup> edition, ASM Press, 2016.

**REVISION HISTORY:**

REVISION	DATE	Description of Change	REQUESTED BY
1.0		Initial Release	L. Steven

**NOTE:** This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

FILENAME:

Print Date: