



**Stanton Territorial Hospital**  
P.O. Box 10, 550 Byrne Road  
YELLOWKNIFE NT X1A 2N1

**Document Number:** MIC70220

**Version No:** 2.0

Page: 1 of 6

**Distribution:**

**Microbiology Instrumentation Manual**

**Effective:** 22 November, 2017

**Date Reviewed:** 22 November, 2017

**Next Review:** 22 November, 2019

**Document Name:** Dispensette Sterility Test

**Approved By:**

Jennifer G. Daley Bernier, A/ Manager, Laboratory Services

**Status:** **APPROVED**

**PURPOSE:** To check the sterility of the saline used to prepare organisms for Vitek 2 and manual susceptibility testing. To be performed each Monday with weekly Vitek and KB QC testing.

**REAGENTS and/or MEDIA:**

- Thioglycolate broth
- Sterile saline (0.45% and 0.9%)

**SUPPLIES:**

- Dispensette
- Dispensette mounting tool
- Bottle-brush
- Towel
- Steris steam autoclave and supplies

**SPECIAL SAFETY PRECAUTIONS:**

Containment Level 2 facilities, equipment, and operational practices for work involving infectious or potentially infectious materials or cultures.

- Lab gown must be worn when performing activities with potential pathogens.
- Gloves must be worn when direct skin contact with infected materials is unavoidable.
- Eye protection must be used where there is a known or potential risk of exposure to splashes.
- All procedures that may produce aerosols, or involve high concentrations or large volumes should be conducted in a biological safety cabinet (BSC).
- The use of needles, syringes, and other sharp objects should be strictly limited.

All patient specimens are assumed to be potentially infectious. Universal precautions must be followed. Since viable micro-organisms are used, all cultures must be handled with appropriate precautions. All equipment in contact with cultures should be decontaminated by appropriate methods

**NOTE:** This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

FILENAME: MIC70220DispensetteSterilityTestPRO

Print Date: 06/10/2019 2:12:00 PM

**PROCEDURE INSTRUCTIONS:**

Step	Action
<b>Weekly Dispensette Sterility Check</b>	
1	A standing order has been set up in the LIS: <ul style="list-style-type: none"> <li>• Patient last name: <b>TEST.</b></li> <li>• Patient first name: <b>MICRO STERILITY.</b></li> <li>• <b>CXENV</b> will be automatically ordered.</li> <li>• Collect, receive and plate order.</li> <li>• In Micro OE comment, state source of saline (KB Dispensette or Vitek Dispensette).</li> </ul>
2	Label 2 Thioglycolate broth tubes, one for the KB Dispensette order and one for the Vitek Dispensette order.
3	Dispense 3mL of saline into each of the THIO broths.
4	Place broths in the Saline Sterility rack in the O <sub>2</sub> incubator.
5	On MIC70110.1 - Maintenance Record – Vitek 2, check off and initial that saline sterility was performed.
6	Read broths daily for 3 days and observe for signs of growth.

**INTERPRETATION OF RESULTS:**

Step	Action
<b>Interpretation of sterility specimens</b>	
1	<b>If THIO broth remains clear (no growth):</b> <ul style="list-style-type: none"> <li>• Record “No Growth” results daily in LIS.</li> <li>• After 3 days, report “No Growth After 3 Days” and finalize.</li> </ul>
2	<b>If THIO broth turns cloudy (growth):</b> <ul style="list-style-type: none"> <li>• Record “Cloudy-AUTOCLAVE dispenser ^CLAVE”.</li> <li>• Result ^CLAVE with date dispensette was autoclaved and initials.</li> <li>• Report “Cloudy-AUTOCLAVE dispenser” and finalize.</li> <li>• Discard saline bottle and follow instructions below to clean the dispensette.</li> <li>• Assemble new dispensette, which are stored under the Vitek 2.</li> </ul>

**NOTE:** This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.



**1. Piston**




**2. Volume adjustment knob**

**3. Discharge tube**



**4. Discharge tube screw cap**

**5. Telescopic filling tube**

**CLEANING THE DISPENSETTE:**

Step	Action
<b>Cleaning the dispensette</b>	
1	It is important that the dispensette is cleaned before it is autoclaved to prevent damage from occurring.
2	Remove the dispensette from saline bottle and screw onto a bottle of sterile water. Rinse the dispensette several times by completely filling and emptying it.
3	Remove dispensette from sterile water and discharge any remaining liquid.
4	<p>Pull out the telescopic filing tube:</p> 
5	<p>Use mounting tool to unscrew the filing valve:</p> 
6	<p>Loosen piston. Hold the housing securely and unscrew the piston completely by turning it to the right. Carefully pull out the piston. <b>Do not remove the housing!</b></p> 

**NOTE:** This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

7	<p>Clean piston and cylinder with a bottle-brush. If necessary remove deposits at the edge of the glass cylinder taking care not to scratch the interior of the cylinder.</p> 
8	<p>Flush all parts with sterile water and allow parts to dry.</p>
9	<p>Insert the piston completely into the cylinder.</p>
10	<p>Check that the filling valve is securely seated and tightened using the mounting tool.</p> 
11	<p>Reassemble the dispensette.</p>
12	<p>Place dispensette and mounting tool on a towel. Ensure parts do not come into contact with metallic surfaces.</p>
13	<p>Autoclave in the Steris steam autoclave at 250°F for 50 minutes.</p>
14	<p>After cycle is complete, allow dispensette to cool for approximately 2 hours.</p>
15	<p>Inspect all parts for any damage.</p>
16	<p>Once dispensette is cooled and inspected for damage, place back into box and store underneath Vitek 2 instrument.</p>

**NOTE:** This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

**REFERENCES:**

- Brand Dispensette Manual, 2015/3

**REVISION HISTORY:**

REVISION	DATE	Description of Change	REQUESTED BY
1.0	22NOV17	Initial Release	L. Steven
2.0	06OCT19	Updated to reflect new steam autoclave	L. Steven

**NOTE:** This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.