

PROGRAM Standard Operating Procedure – Laboratory Services	
Title: MIC20700-Wet mount reporting in LIS- <i>Trichomonas</i> screens	Policy Number:
Program Name: Laboratory Services	
Applicable Domain: Lab, DI and Pharmacy Services	
Additional Domain(s):	
Effective Date:	Next Review Date:
Issuing Authority: Director, Health Services	Date Approved:
Accreditation Canada Applicable Standard: N/A	

**GUIDING PRINCIPLE:**

Wet mount examination of noninvasively collected specimens can be useful for the rapid, inexpensive detection of *Trichomonas vaginalis*. The sensitivity of the technique varies with the organism and experience of the microscopist but is generally low: about 60 to 70% for *T.vaginalis* in wet mounts of vaginal fluid.

**PURPOSE/RATIONALE:**

This standard operating procedure describes how to report the Wet Mount results of genital specimens in the LIS in a consistent manner.

**SCOPE/APPLICABILITY:**

This procedure applies to Medical Laboratory Technologists (MLTs) performing a wet mount on patient samples.

**SAMPLE INFORMATION:**

<b>Type</b>	<ul style="list-style-type: none"><li>• Vagina</li><li>• Cervix</li><li>• Urethra (male and female)</li><li>• Refer to MIC10100-Microbiology Specimen Processing</li></ul>
-------------	--

**REAGENTS and/or MEDIA:**

- Saline

**SUPPLIES:**

- Glass microscope slides
- Glass coverslips
- Sharps container

**Disclaimer Message:** This is a **CONTROLLED** document for internal use only. Any documents appearing in paper form are not controlled and should be checked against the electronic file version prior to use.

## EQUIPMENT

- Microscope
- 35° ambient air incubator


## SPECIAL SAFETY PRECAUTIONS:

Containment Level 2 facilities, equipment, and operational practices for work involving infectious or potential infectious materials or cultures.

- Ensure that appropriate hand hygiene practices be used.
- Lab gown must be worn when performing activities with potential pathogens.
- Gloves must be worn when direct skin contact with infected materials is unavoidable.
- Eye protection must be used when there is a known or potential risk of exposure of splashes.
- All procedures that may produce aerosols or involve high concentrations or large volumes should be conducted in a biological safety cabinet (BSC).
- The use of needles, syringes and other sharp objects should be strictly limited.

All patient specimens are assumed to be potentially infectious. Routine Practices must be followed. Since viable micro-organisms are used, all cultures must be handled with appropriate precautions. All equipment in contact with cultures should be decontaminated by appropriate methods.




## PROCEDURE INSTRUCTIONS:

Step	Action
<b>Reporting wet mounts in the LIS</b>	
<b>1</b>	<ul style="list-style-type: none"> <li>• Pending wet mounts are found in the LIS under:  <b>Resulting Worklist → WET PREPS</b></li> </ul>  <ul style="list-style-type: none"> <li>• Press enter or double click to open worklist</li> </ul>
<b>2</b>	Remove rack of wet mount samples from the O <sub>2</sub> incubator after a 15 minute incubation.
<b>3</b>	<ul style="list-style-type: none"> <li>• Enter the accession number on the test tube and select enter to mark the order</li> <li>• Select enter again to open Result Entry or double click on accession number to open</li> </ul>
<b>4</b>	<ul style="list-style-type: none"> <li>• Verify that the age of the specimen is not ≥72 hours</li> <li>• If specimen is ≥ 72 hours old, specimen is unsuitable for trichomonas analysis and the test needs to be cancelled</li> <li>• Select <b>Cancel Test</b> and insert cancelation comment <b>IXTRICH</b></li> </ul>

**Disclaimer Message:** This is a **CONTROLLED** document for internal use only. Any documents appearing in paper form are not controlled and should be checked against the electronic file version prior to use.

<b>5</b>	<ul style="list-style-type: none"> <li>Gently press the swab onto the microscope slide and return it to the original transport tube</li> <li>Place coverslip on the slide and allow to settle</li> </ul>
<b>6</b>	Using the highest dry objective (40X), observe for trichomonads, which are flagellated. <b>TIP:</b> Decrease the light of the microscope by closing the diaphragm
<b>7</b>	In order to report trichomonads, the flagella must be seen and motile.
<b>8</b>	Under the test code: <b>SCWPT</b> , use the <b>SCWPT</b> keypad to report the presence or absence of Trichomonas vaginalis.

**REPORTING INSTRUCTIONS:**

IF	REPORT
Trichomonas not seen Sample <24 hours old	1. On the SCWPT keypad, select <b>Key 1</b> to state: <b>"Trichomonas vaginalis – Not Seen"</b>  2. Finalize SCWPT 3. Preview instant report and save
Trichomonas not seen Sample >24 hours old and <72 hours old	1. On the SCWPT keypad, select <b>Key 1</b> to state: <b>"Trichomonas vaginalis – Not Seen"</b> 2. On the SCWPT keypad, select <b>Key 2</b> to state: <b>NOTE: The presence of Trichomonas cannot be ruled out since there was a delay &gt;24hrs in transport and/or processing of this specimen."</b>  3. Finalize SCWPT 4. Preview instant report and save
Trichomonas seen (Must be motile)	1. On the SCWPT keypad, select <b>Key 3</b> to state: <b>"Trichomonas vaginalis – SEEN"</b>  2. Finalize SCWPT 3. Preview instant report and save

**Disclaimer Message:** This is a **CONTROLLED** document for internal use only. Any documents appearing in paper form are not controlled and should be checked against the electronic file version prior to use.

Step	Action
<b>Complete reading of wet mounts</b>	
<b>1</b>	<ul style="list-style-type: none"> <li>Refresh <b>WET PREPS</b> worklist</li> <li>If finished reading wet mounts, ensure no remaining samples are on the worklist</li> </ul>
<b>2</b>	Dispose of slides in sharps container. Place swab in rack designated for the day.

**LIMITATIONS:**

- The wet mount is dependent on the expertise of the microscopist.
- White blood cells are often confused with trichomonads that are no longer moving, which can be avoided by careful examination for the characteristics of the parasite.
- Delays in transit decrease the ability to detect trichomonads, as they can rapidly lose their motility.

**CROSS-REFERENCES:**

- MIC10100-Microbiology Specimen Processing

**REFERENCES:**

- Leber, A. (2016). *Clinical microbiology procedures handbook*. (4<sup>th</sup>ed.) Washington, D.C.: ASM Press

**APPROVAL:**

\_\_\_\_\_  
 Date

**REVISION HISTORY:**

REVISION	DATE	Description of Change	REQUESTED BY
1.0	07 Feb 19	Initial Release	L. Steven
2.0	31 Mar 22	Procedure reviewed and added to NTHSSA policy template	L. Steven

**Disclaimer Message:** This is a **CONTROLLED** document for internal use only. Any documents appearing in paper form are not controlled and should be checked against the electronic file version prior to use.