

PROGRAM Standard Operating Procedure – Laboratory Services	
Title: MIC60110 – BIOFIRE RP2.1 Quality Control	Policy Number:
Program Name: Laboratory Services	
Applicable Domain: Lab, DI and Pharmacy Services	
Additional Domain(s): NA	
Effective Date:	Next Review Date:
Issuing Authority: Director, Laboratory and Diagnostic Imaging Services	Date Approved:
Accreditation Canada Applicable Standard: NA	

Uncontrolled When Printed

GUIDING PRINCIPLE:

Quality control is performed on BIOFIRE RP2.1 pouches to ensure proper function on a monthly basis and to ensure new shipments have not deteriorated during shipping. Quality control should be rotated through the different modules of the instrument.

PURPOSE/RATIONALE:

This standard operating procedure describes the quality control procedure for the BIOFIRE RP2.1 test on the BIOFIRE TORCH.

SCOPE/APPLICABILITY:

This standard operating procedure applies to Medical Laboratory Technologists (MLTs) performing quality control for the BIOFIRE RP2.1 test.

REAGENTS and/or MEDIA:

- MMQCI BIOFIRE RP2.1 controls
- BIOFIRE RP2.1 pouches
- Hydration injection vials (blue)
- Sample buffer ampoules
- Sample injection vials (red)
- Accel TB 1 L bottle
- Accel TB wipes
- 70% isopropyl alcohol

SUPPLIES:

- Personal protective equipment
- Wet waste container
- Dry waste container
- Orange autoclave bags
- Spray bottle
- Sharps container
- Transfer pipettes provided in kit

Disclaimer Message: This is a **CONTROLLED** document for internal use only. Any documents appearing in paper form are not controlled and should be checked against the electronic file version prior to use.

EQUIPMENT:

- BIOFIRE Torch
- Class II biosafety cabinet (BSC)
- Vortex mixer
- Freezer

ENVIRONMENTAL CONTROLS:

- Store the test kit, including reagent pouches and buffers, at room temperature (15°C to 25°C) **DO NOT REFRIGERATE**
- Avoid storage of any materials near heating or cooling vents or in direct sunlight
- All kit components should be stored and used together. Do not use components from one kit with those of another kit
- To prevent amplicon contamination: discard used pouches after the run has completed, avoid excessive handling of pouches after test runs and avoid exposing pouches to sharp edges or anything that might cause a puncture

PROCEDURE INSTRUCTIONS:

Step	Action
Performing BIOFIRE RP2.1 Quality Control	
1	BIOFIRE RP2.1 quality control is performed monthly and upon receipt of new pouches.
2	Perform quality control testing using MMQCI BIOFIRE RP2.1 controls. The kit includes a positive control and a negative control and is stored frozen.
3	Quality control is rotated through the different modules of the instrument using the defined schedule. Refer to appendix below.
4	Prepare the BIOFIRE RP2.1 pouch as per MIC73200-BIOFIRE Respiratory 2.1 (RP2.1) Panel substituting the patient specimen with the control solutions. NOTE: Do not run controls with patient samples in order to prevent possible contamination from the positive control solutions
5	When creating a test run, manually enter the sample ID as Positive Control or Negative Control.
6	If all QC results are acceptable, complete MIC60111-QC Results Record-BIOFIRE RP2.1. Ensure the control lot number, control expiry date, kit lot number, kit expiry date and module numbers are filled in along with the results, the performing technologist's initials, and the reason for performing (monthly, new lot or new shipment). The yearly QC Results Record is kept in the PCR Testing Completed QC Records binder under the MIC60110-QC Results Record-BIOFIRE RP2.1 tab.
7	If QC results are not acceptable, ensure controls were not mixed-up during loading (ex. positive control run as the negative control). After preliminary investigation, repeat testing.
8	If repeat QC testing is still not acceptable: <ul style="list-style-type: none"> Notify the Technical Supervisor-Microbiology for resolution Contact bioMerieux to determine if issues with kits exist Until the problem is resolved, patient testing may not be reported

Disclaimer Message: This is a **CONTROLLED** document for internal use only. Any documents appearing in paper form are not controlled and should be checked against the electronic file version prior to use.

APPENDIX:
BIOFIRE RP2.1 Quality Control Schedule

Run	Module	Control	Run	Module	Control
1	1	Positive	11	6	Positive
2	2	Negative	12	7	Negative
3	2	Positive	13	7	Positive
4	3	Negative	14	8	Negative
5	3	Positive	15	8	Positive
6	4	Negative	16	9	Negative
7	4	Positive	17	9	Positive
8	5	Negative	18	10	Negative
9	5	Positive	19	10	Positive
10	6	Negative	20	1	Negative

CROSS-REFERENCES:

- MIC60111-QC Results Record-BIOFIRE RP2.1
- MIC73200-BIOFIRE Respiratory 2.1 (RP2.1) Panel

REFERENCES:

1. MMQCI. (August 2024). *BIOFIRE RP2.1/RP2.1plus Control Panel* insert
2. BIOFIRE *Respiratory Panel 2.1 (RP2.1)* Instructions for Use. BFR0001-7709-03 September 2023
3. BIOFIRE. *BIOFIRE FilmArray Torch Operator's Manual*, HTFA-PRT-0001-05, March 2020
4. Biosafety Advisory Committee. *STHA Biosafety Program Manual*. January 2016

APPROVAL:

Date

Director, Laboratory and Diagnostic Imaging Services

REVISION HISTORY:

REVISION	DATE	Description of Change	REQUESTED BY
1.0	15 Aug 20	Initial Release	L. Steven
2.0	17 Oct 22	Procedure reviewed and updated to reflect QC run weekly as well as when new kits received	L. Steven
3.0	16 Apr 24	Procedure reviewed and updated to reflect QC run monthly not weekly	L. Steven

Disclaimer Message: This is a **CONTROLLED** document for internal use only. Any documents appearing in paper form are not controlled and should be checked against the electronic file version prior to use.