

 <p>Stanton Territorial Hospital P.O. Box 10, 550 Byrne Road YELLOWKNIFE NT X1A 2N1</p>	<p>Document Name: Acridine Orange Stain Job Aid</p> <p>Distribution: Microbiology Microscopy Manual</p>	<p>Document Number: MIC20110</p> <p>Date Issued: September 19, 2022</p>
---	---	---

Performing the acridine orange stain:

1	<p>SoftMic ordering code is:</p> <ul style="list-style-type: none"> • ACROR in the media resulting plate log if used to verify bacteria seen in gram stain • STAU in test reporting log if used to verify no bacteria seen on positive blood culture smear
2	<ul style="list-style-type: none"> • Turn on the fluorescent component of the microscope and allow the lamp to warm up for at least 15 minutes
3	<ul style="list-style-type: none"> • Prepare a slide of the sample and allow to dry
4	<ul style="list-style-type: none"> • Place slide on the slide warmer in the BSC until dry
5	<ul style="list-style-type: none"> • Fix smear with methanol for 1 minute • Drain off remaining methanol without rinsing • Allow slide to air dry again
6	<ul style="list-style-type: none"> • Flood slide with acridine orange stain for 2 minutes
7	<ul style="list-style-type: none"> • Drain the excess stain and rinse thoroughly with tap water
8	<ul style="list-style-type: none"> • Place slides in a slide holder and allow to drip dry • Stained slides may be placed in the O₂ incubator

NOTE: This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.