

## TRAINING UPDATE

**Lab Location:** SGAH  
**Department:** Blood Bank

**Date Implemented:** 2.4.13  
**Due Date:** 2.25.13

### DESCRIPTION OF PROCEDURE REVISION

|   |
|---|
| <b>Name of procedure:</b>   |
| ABO/Rh Testing for RhIG Candidacy   |
| <b>Description of change(s):</b>  |
| <ul style="list-style-type: none"><li>• ABO/Rh test will be eliminated from Cerner and SMS and replaced with test “ABO/Rh for RhIG Candidacy.”</li><li>• This will come to Sunquest as an ABORH test.</li><li>• When this test is ordered.<ul style="list-style-type: none"><li>○ Run an ABO/Rh.</li><li>○ If the patient is Rh-negative, run an antibody screen and panel (if indicated).</li><li>○ If the patient is <math>\geq 20</math> weeks gestation, perform Kleihauer-Betke testing prior to allocating RhIG.</li></ul></li><li>• It is permissible to issue 1 vial of RhIG pending KBT results.</li></ul> |

Non-Technical SOP

|                    |  |                 |
|--------------------|--|-----------------|
| <b>Title</b>       | <b>ABO/Rh Testing for RhIG Candidacy</b> |                 |
| <b>Prepared by</b> | Stephanie Codina                         | Date: 1.31.2013 |
| <b>Owner</b>       | Stephanie Codina                         | Date: 1.31.2013 |

| <b>Laboratory Approval</b>   |                  |                              |
|--|------------------|------------------------------|
| <b>Print Name and Title</b>  | <b>Signature</b> | <b>Date</b>                  |
| <i>Refer to the electronic signature page for approval and approval dates.</i> |                  |                              |
|  |                  |                              |
|  |                  |                              |
| <b>Local Issue Date:</b>   |                  | <b>Local Effective Date:</b> |

| <b>Review:</b>    |                  |             |
|-------------------|------------------|-------------|
| <b>Print Name</b> | <b>Signature</b> | <b>Date</b> |
|                   |                  |             |
|                   |                  |             |
|                   |                  |             |
|                   |                  |             |
|                   |                  |             |
|                   |                  |             |
|                   |                  |             |

Form revised 3/31/00

**TABLE OF CONTENTS**

1. PURPOSE ..... 2  
 2. SCOPE ..... 2  
 3. RESPONSIBILITY ..... 2  
 4. DEFINITIONS ..... 2  
 5. PROCEDURE ..... 2  
 6. RELATED DOCUMENTS ..... 3  
 7. REFERENCES ..... 3  
 8. REVISION HISTORY ..... 3  
 9. ADDENDA AND APPENDICES ..... 3

**1. PURPOSE**

To define the testing process that will take place if the test “ABO/Rh RhIG Candidate” is ordered.

**2. SCOPE**

This procedure applies to patients who have orders for the “ABO/Rh RhIG Candidate” test. This test is generally only ordered on pregnant females who have experienced a possible sensitization of fetal blood.

**3. RESPONSIBILITY**

All blood bank staff members must understand and adhere to this procedure for performing ABO/Rh RhIG Candidate testing.

**4. DEFINITIONS**

N/A

**5. PROCEDURE**

| Step | Action   |
|------|--|
| 1    | Laboratory testing is necessary to determine whether a patient is a candidate for RhIG injection if the patient has experienced an actual or suspected fetal-maternal hemorrhage (FMH). Examples of events that may result in an FMH include, but are not limited to, placenta previa, amniocentesis, chorionic villus sampling, percutaneous umbilical blood sampling, other obstetrical manipulative procedure or abdominal trauma, actual or threatened pregnancy loss at any stage of gestation, ectopic pregnancy, and hydatiform mole. |
| 2    | When an “ABO/Rh RhIG Candidate” test is received in the blood bank, an ABO/Rh test will be performed and reported per procedure, “ABO/Rh Typing (Manual Tube).”  |

Form revised 3/31/00

| Step | Action  |
|------|---|
| 3    | If the ABO/Rh typing reveals the patient is Rh-negative or weak D positive, an antibody screen will be performed per departmental procedure. Antibody identification procedures will be followed if the antibody screen is positive. RhIG is NOT indicated if the patient has a current or historical record of anti-D.   |
| 4    | Identification and quantification of fetal bleed must be performed if the patient is $\geq 20$ weeks gestation. <ul style="list-style-type: none"> <li>A. If the blood type of the fetus is not known, a Kleihauer-Betke test must be performed.</li> <li>B. If the blood type of the fetus is known (unlikely) and both mom and baby meet acceptability specifications outlined in procedure, "Fetal Bleed Screening," a fetal cell screen may be performed to determine if a fetal bleed has occurred before performing Kleihauer-Betke testing.</li> </ul> |
| 5    | RhIG should be administered to any patient who is Rh-negative and has no current or historical record of anti-D. <ul style="list-style-type: none"> <li>A. 1 vial of RhIG is given to patients who are <math>&lt; 20</math> weeks gestation.</li> <li>B. The amount of RhIG given to patients who are <math>\geq 20</math> weeks gestation will be determined by fetal bleed screen and/or Kleihauer-Betke test results. One vial of RhIG may be given prior to completion of testing.</li> </ul>   |

**6. RELATED DOCUMENTS**

- SOP: ABO/Rh Typing (Manual Tube)
- SOP: Galileo Echo Testing Patient Specimens
- SOP: Antibody Identification
- SOP: Fetal Bleed Screening

**7. REFERENCES**

None

**8. REVISION HISTORY**

| Version | Date | Reason for Revision | Revised By | Approved By |
|---------|------|---------------------|------------|-------------|
|         |      |                     |            |             |
|         |      |                     |            |             |
|         |      |                     |            |             |
|         |      |                     |            |             |

**9. ADDENDA AND APPENDICES**

None

Form revised 3/3/1/00