

## TRAINING UPDATE

**Lab Location:** SGAH & WAH  
**Department:** Micro

**Date Distributed:** 1/1/2014  
**Due Date:** 1/31/2014  
**Implementation:** 2/1/2014

### DESCRIPTION OF PROCEDURE REVISION

<b>Name of procedure:</b>
<b>Measurement of Blood Culture Volume SGAH.M40, WAH.M37 v1</b> <b>Blood Culture Volume Monitor form AG.F28 v2</b>
<b>Description of change(s):</b>
<p>SOP: Section 5.2: Changed # of bottles for SGAH (SOP corrected to match current practice.</p> <p>FORM: Instructions changed to indicate only use aerobic bottle (previously stated to use aerobic and anaerobic)</p> <p><b>This revised SOP and form will be implemented on February 1, 2014</b></p>

Document your compliance with this training update by taking the quiz in the MTS system.

**Approved draft for training all sites (version 1)**

Non-Technical SOP

<b>Title</b>	<b>Measurement of Blood Culture Volume</b>	
<b>Prepared by</b>	Ron Master	Date: 4/26/2013
<b>Owner</b>	Ron Master	Date: 4/26/2013

<b>Laboratory Approval</b>		
<b>Print Name and Title</b>	<b>Signature</b>	<b>Date</b>
<i>Refer to the electronic signature page for approval and approval dates.</i>		
Local Issue Date:		Local Effective Date:

<b>Review:</b>		
<b>Print Name</b>	<b>Signature</b>	<b>Date</b>

**TABLE OF CONTENTS**

1. PURPOSE.....	3
2. SCOPE .....	3
3. RESPONSIBILITY.....	3
4. DEFINITIONS.....	3
5. PROCEDURE.....	3
6. RELATED DOCUMENTS .....	4
7. REFERENCES .....	4
8. REVISION HISTORY.....	4
9. ADDENDA AND APPENDICES .....	4

---

**1. PURPOSE**

Larger volumes of blood increase the yield of true positive cultures. The laboratory monitors collected blood volumes.

**2. SCOPE**

This procedure covers the monitoring of blood culture volumes. It does not include performance of blood cultures.

**3. RESPONSIBILITY**

It is the responsibility of technical laboratory personnel who work in microbiology to learn and understand the process described in this procedure.

**4. DEFINITIONS**

ED - Emergency Department

**5. PROCEDURE**

1. The standard used for measurement is an empty aerobic BACTEC bottle filled with 8 mL of stain such as safrannin or carbol fuchsin.
2. Select 2 **aerobic** bottles from 2 different patients each day. Do not use pediatric bottles.  
SGAH only: **select 2 bottles collected by ED staff and 2 bottles collected by a laboratory phlebotomist.**

Form revised 3/31/00

3. Compare the volume in each of the bottles to be measured to the standard.
4. Record the following information on the Blood Culture Volume Monitor form:
  - a. Number checked, number acceptable, and number unacceptable
  - b. For each unacceptable bottle, record the accession number and collector.
  - c. Record the initials of the person performing the evaluation.
5. At the end of the month, calculate the total number checked, number acceptable, and number unacceptable. Calculate the % acceptable.
6. If the percent acceptable for a month is less than 85%, make a copy of the monitor form for the Field Operations Manager or Supervisor.
7. The monitor form is reviewed weekly by the Field Operations Manager, Supervisor or Phlebotomy Group Lead. Monthly review is performed by the Microbiology Manager or Core Laboratory Supervisor.
8. Communication of findings and Corrective Action
  - a. Results are graphed monthly
  - b. Data is distributed to Field Operations Manager and Supervisor for posting
  - c. Data is presented at Lab / ED meetings and documented in minutes
  - d. Follow up with individual staff is performed by the appropriate supervisor

**6. RELATED DOCUMENTS**

Blood Cultures, with Automated Detection, Microbiology procedure  
 Blood Culture Volume Monitor form (AG.F28)

**7. REFERENCES**

None

**8. REVISION HISTORY**

Version	Date	Reason for Revision	Revised By	Approved By
000	12/17/13	Section 5.2: Changed # of bottles for SGAH Section 6: move form from section 9 Footer: version # leading zero's dropped due to new EDCS in use as of 10/7/13	R. Master	R. Master

**9. ADDENDA AND APPENDICES**

None



## Blood Culture Volume Monitor

- Germantown Emergency Center
- Shady Grove Adventist Hospital
- Washington Adventist Hospital

**Month / Year** \_\_\_\_\_

To be done daily: Two **aerobic blood** culture bottles from two different patients will be monitored daily. Randomly choose two patient aerobic blood culture bottles to compare to the standard bottle containing 8 mL. Acceptable range:  $\geq 85\%$  should contain at least 8 mL of blood. Notify [Field Operations Manager or Supervisor](#) if acceptable range is not met.

Date	Number evaluated	Number acceptable	Number unacceptable	Accession number if unacceptable	Drawn by	Initials
1						
2						
3						
4						
5						
6						
7						
8						
9						
10						
11						
12						
13						
14						
15						
16						
17						
18						
19						
20						
21						
22						
23						
24						
25						
26						
27						
28						
29						
30						
31						
Totals						
% Acceptable						

Weekly review:	Weekly review:	Weekly review:
Weekly review:	Weekly review:	Monthly review:

Document: AG.F28[2] Status: PRERELEASED, Effective: 1/1/2009, Check Version Before Use