

TRAINING UPDATE

Lab Location: SGAH & WAH
Department: Chemistry

Date Distributed: 2/27/2014
Due Date: 3/15/2014
Implementation: 3/1/2014

DESCRIPTION OF PROCEDURE REVISION

Name of procedure:
Ammonia by Dimension Vista® System SGAH.C865 / WAH.C864 v0 Carbon Dioxide by Dimension Vista® System SGAH.C93. WAH.C89v1 Ethyl Alcohol by Dimension Vista® System SGAH.C108. WAH.C104v1
Description of change(s):
New reagent for Ammonia Change in calibrator for CO2 and Ethyl Alcohol - only CO2 sop is attached (after ammonia) because changes are the same These SOPs will be implemented on March 1, 2014

Document your compliance with this training update by taking the quiz in the MTS system.

Approved draft for training all sites (version 0)

Technical SOP

Title	Ammonia by Dimension Vista® System	
Prepared by	Ashkan Chini	Date: 2/5/2014
Owner	Robert SanLuis	Date: 2/5/2014

Laboratory Approval		Local Effective Date:
Print Name and Title	Signature	Date
<i>Refer to the electronic signature page for approval and approval dates.</i>		

Review		
Print Name	Signature	Date

TABLE OF CONTENTS

1. Test Information.....3
 2. Analytical Principle4
 3. Specimen Requirements.....4
 4. Reagents5
 5. Calibrators/Standards6
 6. Quality Control7
 7. Equipment And Supplies9
 8. Procedure9
 9. Calculations.....11
 10. Reporting Results And Repeat Criteria.....11
 11. Expected Values.....11
 12. Clinical Significance12
 13. Procedure Notes12
 14. Limitations Of Method12
 15. Safety13
 16. Related Documents14
 17. References.....14
 18. Revision History14
 19. Addenda14

1. TEST INFORMATION

Assay	Method/Instrument	Local Code
Ammonia	Dimension Vista® System	NH3

Synonyms/Abbreviations
NH3

Department
Chemistry


Form revised 3/31/00

2. ANALYTICAL PRINCIPLE

Ammonia measurements are used in the diagnosis and treatment of severe liver disorders such as cirrhosis, hepatitis and Reye’s syndrome.

The Dimension Vista® Ammonia (AMM) method is an enzymatic method that uses glutamate dehydrogenase (GLDH) and a stabilized NADPH analog. Ammonia reacts with

α-ketoglutarate and reduced cofactor to form L-glutamate and the cofactor. The reaction is catalyzed by glutamate dehydrogenase. The decrease in absorbance due to the oxidation of the reduced cofactor is monitored at 340/700 nm and is proportional to the ammonia concentration.



3. SPECIMEN REQUIREMENTS

3.1 Patient Preparation

Component	Special Notations
Fasting/Special Diets	N/A
Specimen Collection and/or Timing	The tube should be completely filled, stored tightly capped on ice, centrifuged immediately and analyzed within 20 minutes.
Special Collection Procedures	N/A
Other	N/A

3.2 Specimen Type & Handling

Criteria	
Type -Preferred -Other Acceptable	Plasma (Heparin) None
Collection Container	Green Top Tube
Volume - Optimum - Minimum	1.0 mL 0.5 mL
Transport Container and Temperature	Collection container or plastic vial on ice.
Stability & Storage Requirements	Room Temperature: Not stable
	Refrigerated: 2 hours
	Frozen: Not stable
	Instrument on board aliquot stability 2 hours
Timing Considerations	Plasma should be physically separated from cells as soon as possible.

FORM REVISION 3/31/00

Criteria	
Unacceptable Specimens & Actions to Take	Specimens that are unlabeled, improperly labeled, or those that do not meet the stated criteria are unacceptable. Request a recollection and credit the test with the appropriate LIS English text code for “test not performed” message. Examples: Quantity not sufficient-QNS; Wrong collection-UNAC. Document the request for recollection in the LIS.
Compromising Physical Characteristics	Do not use hemolyzed samples. Reject sample and request a recollection. Credit the test with the appropriate LIS English text code explanation of HMT (Specimen markedly hemolyzed)
Other Considerations	Concentrations may more than double in plasma when stored at room temperature for 6 hours.

4. REAGENTS

Refer to the Material Safety Data Sheet (MSDS) supplied with the reagents for complete safety hazards. Refer to the section in this procedure covering “SAFETY” for additional information.

4.1 Reagent Summary

Reagents / Kits	Supplier & Catalog Number
Ammonia	Siemens, Flex® reagent cartridge, Cat. No. K3119

4.2 Reagent Preparation and Storage

NOTES: Each container must be labeled with (1) substance name, (2) lot number, (3) expiration date, (4) any special storage instructions; check for visible signs of degradation. When placed onboard the analyzer, the instrument captures the date / time loaded and calculates and tracks the opened expiration.

Refer to the Material Safety Data Sheet (MSDS) for a complete description of hazards. If a specific hazard is present, it will be noted in this procedure when the hazard is first encountered in a procedural step.

Reagent	Ammonia
Container	Reagent cartridge
Storage	Store at 2-8° C
Stability	<ul style="list-style-type: none"> Reagent is stable until expiration date stamped on the reagent cartridges. Sealed wells on the instrument are stable for 60 days. Once wells 1 - 12 have been entered by the instrument, they are stable for 3 days.
Preparation	All reagents are liquid and ready to use.

5. CALIBRATORS/STANDARDS

5.1 Calibrators/Standards Used

Calibrator	Supplier and Catalog Number
CHEM 3 CAL	Siemens Dimension Vista®, Cat. No. KC130A

5.2 Calibrator Preparation and Storage

NOTE: Date and initial all calibrators upon opening. Each container must be labeled with (1) substance name, (2) lot number, (3) date of preparation, (4) expiration date, (5) any special storage instructions; check for visible signs of degradation. When placed onboard the analyzer, the instrument captures the date / time loaded and calculates and tracks the opened expiration.

Calibrator	CHEM 3 CAL
Preparation	Calibrator is ready for use. No preparation is required.
Storage/Stability	<ul style="list-style-type: none"> Store at 2 – 8 ° C Unopened calibrator is stable until expiration date stamped on the box. Opened Calibrator: once the stopper of the vial is punctured, assigned values are stable for 24 hours when stored on board the Dimension Vista System.

5.3 Calibration Parameter

Criteria	Special Notations
Reference Material	CHEM 3 CAL
Assay Range	10 – 750 µmol/L
Suggested Calibration Level	See Reagent Package Insert for lot specific assigned values in µmol/L
Frequency	<ul style="list-style-type: none"> Every new reagent cartridge lot. Every 60 days for any one lot When major maintenance is performed on the analyzer. When control data indicates a significant shift in assay.
Calibration Scheme	Two levels

5.4 Calibration Procedure

Auto Calibration:

- Place the required calibrator vials in a carrier. Make sure the barcode labels are entirely visible through the slots.
- Place the carrier in the loading area.
- Position the carrier with the labels facing away from the user.
- Press the **Load** button.
- Automatic calibration requires that calibrators be on the instrument. As the time for processing approaches, the instrument reviews onboard inventory for the appropriate calibrators.

Manual Calibration:

1. Verify that calibrators and reagents are in inventory on the instrument.
2. Press **System > Method Summary > Calibration**.
3. Select a method from the sidebar menu. Press the **Order Calibration** button on the screen.
4. Verify that the information on the screen is correct. Verify that the calibrator lot is correct using the drop-down menu.
 - a. When calibrating using Vials press **OK**.
 - b. When calibrating using Cups, check the Use Cups box. This displays the rack and cup position fields. For additional cups use the positions in ascending order. Be sure to use the number of calibration levels and cups as specified in the method IFU. Scan the rack barcode and place calibrator cups in an adapter in position 1 on a rack. Press **OK** and load the rack on the instrument.
5. The status field in the calibration screen changes sequentially to Awaiting Scheduling, Preparing Calibrators and Processing.

5.5 Tolerance Limits

IF.....	THEN.....
If result fall within assay-specific specification, and QC values are within acceptable limits,	proceed with analysis
If result falls outside assay-specific specification, or QC values are out of Acceptable limits,	troubleshoot the assay and/or instrument and repeat calibration

6. QUALITY CONTROL

6.1 Controls Used

Controls	Supplier and Catalog Number
Liquichek™ Ethanol/Ammonia Control Levels 1, 2 & 3	Bio-Rad Laboratories Cat. No. 544, 545 & 546

6.2 Control Preparation and Storage

NOTE: Date and initial all controls upon opening. Each container should be labeled with (1) substance name, (2) lot number, (3) date of preparation, (4) expiration date, (5) initials of tech, and (6) any special storage instructions; check for visible signs of degradation. A barcode label is produced and placed on the vial.

Control	Ethanol/Ammonia
Preparation	Before sampling, allow the vial to reach room temperature (18 - 25°C). Gently swirl the vial several times to ensure homogeneity.
Storage/Stability	Store at 2 – 8 °C Once opened and stored tightly, QC vial will be stable for 20 days.

6.3 Frequency

Analyze all levels of QC material after every calibration and each day of testing (notated on the QC frequency sheets posted on the instruments).

Refer to the Dimension Vista® QC Schedule in the Laboratory policy Quality Control Program and in the Dimension Vista® Quick Reference Guide.

6.4 Tolerance Limits

Step	Action
1	Acceptable ranges for QC are programmed into the Laboratory Information System (LIS), and may be posted near the instrument for use during computer downtime.
2	<p>Run Rejection Criteria</p> <ul style="list-style-type: none"> Anytime the established parameters are exceeded (if one QC result exceeds 2 SD), the run is considered out of control (failed) and patient results must not be reported. The technologist must follow the procedure in the Laboratory QC Program to resolve the problem.
3	<p>Corrective Action:</p> <ul style="list-style-type: none"> All rejected runs must be effectively addressed through corrective action. Steps taken in response to QC failures must be documented. Patient samples in failed analytical runs must be <u>reanalyzed according to the Laboratory QC Program</u>. Supervisors may override rejection of partial or complete runs only with detailed documentation and criteria for overrides that are approved by the Medical Director. Consult corrective action guidelines in Laboratory QC Program. Follow corrective action guidelines in the Laboratory QC Program. Corrective action documentation must follow the Laboratory Quality Control Program.
4	<p>Review of QC</p> <ul style="list-style-type: none"> QC must be reviewed weekly by the Group Lead or designee and monthly by the Supervisor/Manager or designee. If the SD and/or CV are greater than established ranges, investigate the cause for the imprecision and document implementation of corrective actions.

6.5 Review Patient Data

Technologist must review each result with error messages. Refer to the Dimension Vista® system manual “Error messages” section for troubleshooting. Check for unusual patterns, trends, or distributions in patient results (such as an unusually high percentage of abnormal results). Resolve any problems noted before issuing patient reports.

6.6 Documentation

Form revised 3/31/00

- QC tolerance limits are programmed into the instrument and the LIS. The LIS calculates cumulative mean, SD and CV and stores all information for easy retrieval.
- Quality control records are reviewed daily at the bench, weekly by the Group Lead or designee, and monthly by the Supervisor/Manager or designee.
- Refer to complete policies and procedures for QC documentation and for record retention requirements in the Laboratory QC Program.

6.7 Quality Assurance Program

- Each new lot number of reagent or new shipment of the same lot of reagent must be tested with external control materials and previously analyzed samples. Performance of the new lot must be equivalent to the previous lot; utilize published TEA for acceptability criteria.
- Training must be successfully completed and documented prior to performing this test. This procedure must be incorporated into the departmental competency assessment program.
- The laboratory participates in CAP proficiency testing. All proficiency testing materials must be treated in the same manner as patient samples.
- Monthly QC must be presented to the Medical Director or designee for review and signature.
- Monthly QC mean and SD are sent to Bio-Rad Laboratories for peer group comparison.
- Consult the Laboratory QC Program for complete details.

7. EQUIPMENT and SUPPLIES

7.1 Assay Platform Dimension Vista® System

7.2 Equipment

- Refrigerator capable of sustaining 2–8°C.
- Freezer capable of sustaining range not to exceed -20 to -70°C.
- Centrifuge

7.3 Supplies

- Aliquot Plates
- System Fluids
- Assorted calibrated pipettes (MLA or equivalent) and disposable tips

8. PROCEDURE

AMM Flex® reagent cartridge Cat. No. K3119 is required to perform this test.

Ammonia is performed on the Dimension Vista® System after the method is calibrated (see Reference Material in Calibration section) and Quality Controls are acceptable.

NOTE: For all procedures involving specimens, buttoned lab coats, gloves, and face protection are required minimum personal protective equipment. Report all accidents to your supervisor.

The package insert for a new lot of kits must be reviewed for any changes before the kit is used. A current Package Insert is included as a Related Document.

8.1	Sample Processing
1.	A sample rack holding tubes or cups is placed on the rack input lane.
2.	The sample shuttle moves the rack to the barcode reader which identifies the rack and samples to the system.
3.	The rack moves into the sample server and to the rack positioner.
4.	At the same time, aliquot plates move from the aliquot loader into position.
5.	The aliquot probe aspirates the sample from the tubes or cups and dispenses it into the wells of the aliquot plates.
6.	After each aspirate-dispense action, the probe is thoroughly rinsed inside and out to prevent sample carryover.
7.	When sample aspiration is completed, the sample server moves the rack back to the sample shuttle, where it is placed on the output lane and can be removed by the operator.

8.2	Specimen Testing
1.	For QC placement and frequency, refer to the Dimension Vista® QC Schedule in the Laboratory QC Program.
2.	Follow the instructions, outlined in the Dimension Vista® Operator’s Manual
3.	The instrument reporting system contains error messages to warn the user of specific malfunctions. Results followed by such error messages should be held for follow-up. Refer to the Dimension Vista® system manual “Error messages” section for troubleshooting.
4.	Follow protocol in Section 10.5 “Repeat criteria and resulting” for samples with results above or below the Analytical Measurement Range (AMR). Investigate any failed delta result and repeat, if necessary.
5.	Append the appropriate English text code qualifier messages to any samples requiring a comment regarding sample quality and/or any other pertinent factors.

Test Conditions	
Sample Volume:	20 µL
Reagent Volume:	130 µL
Reaction Time:	5.9 minutes
Test Temperature:	37° C
Wavelength:	340 & 700 nm
Type of measurement:	Bichromatic rate

Form revised 3/31/00

9. CALCULATIONS

The instrument automatically calculates the concentration of Ammonia in $\mu\text{mol/L}$.

10. REPORTING RESULTS AND REPEAT CRITERIA

10.1 Interpretation of Data

None required

10.2 Rounding

No rounding is necessary. Instrument reports results as a whole number.

10.3 Units of Measure

$\mu\text{mol/L}$

10.4 Clinically Reportable Range (CRR)

10 – 2,250 $\mu\text{mol/L}$

10.5 Repeat Criteria and Resulting

IF the result is ...	THEN...
< 10 $\mu\text{mol/L}$	Assure there is sufficient sample devoid of bubbles, cellular debris, and/or fibrin clots. Report as: < 10 $\mu\text{mol/L}$
\geq 750 $\mu\text{mol/L}$	On Board Automated Dilution: Results \geq 750 $\mu\text{mol/L}$ will automatically have repeat testing performed into the instrument using dilution factor of 2. No multiplication is necessary.
> 1,500 $\mu\text{mol/L}$	Manual Dilution: Using the primary tube, make the smallest dilution possible to bring the raw data within the AMR. Maximum allowable dilution: x 3 Diluent: reagent grade water Enter dilution factor as a whole number on the “Enter Sample Data” screen.
> 2,250 $\mu\text{mol/L}$	If the recommended dilution does not give results within the clinically reportable range, report as: “> 2,250 $\mu\text{mol/L}$ -REP” Bring to the attention of your supervisor prior to releasing result.

Message	Code
Verified by repeat analysis	Append –REP to the result.

11. EXPECTED VALUES

11.1 Reference Ranges

11 – 32 $\mu\text{mol/L}$

11.2 Critical Values
 ≥ 200 µmol/L

11.3 Priority 3 Limit(s)
 None established

12. CLINICAL SIGNIFICANCE

The major source of circulating ammonia is the gastrointestinal (GI) tract. Under normal conditions, ammonia is metabolized to urea by liver enzymes. Several diseases, both inherited and acquired, cause elevated ammonia (hyperammonemia). The inherited deficiencies of urea cycle enzymes are the major cause of hyperammonemia in infants. Acquired hyperammonemia most often results from liver disease, renal failure, and Reye’s syndrome. Elevated ammonia is toxic to the central nervous system.

13. PROCEDURE NOTES

- **FDA Status:** FDA Approved/cleared
- **Validated Test Modifications:** None

The instrument reporting system contains error messages to warn the operator of specific malfunctions. Any report slip containing such error messages should be held for follow-up. Refer to your Dimension Vista Operator’s Guide.

The expected maximum observed standard deviations for repeatability using n = 5 replicates at the following Lithium concentrations are:

AMM Concentration	Acceptable S.D. Maximum
41 µmol/L	6.9 µmol/L
206 µmol/L	25.9 µmol/L

14. LIMITATIONS OF METHOD

14.1 Analytical Measurement Range (AMR)
 10 – 750 µmol/L

14.2 Precision

Material	Mean µmol/L	Standard Deviation (%CV)	
		Repeatability	Within-Lab
Liquichek Ethanol Ammonia			
Level 1	26	1.6	1.9
Level 2	109	1.4	1.9
Level 3	331	2.1	2.9

14.3 Interfering Substances

Dextran 40 at 1500 mg/dL increases AMM results by 35% at an ammonia concentration of 50 µmol/L and increases AMM results by < 10% at an ammonia concentration of 250 µmol/L.

Form revised 3/31/00

Dextran 40 at 250 mg/dL increases AMM results by 10% at an ammonia concentration of 50 µmol/L.

Immunoglobulin G (IgG) at 5 g/dL increases AMM results by 30% at an ammonia concentration 50 µmol/L and increases AMM results by < 10% at an ammonia concentration of 250 µmol/L.

Triglycerides at 3000 mg/dL tripped a test report message; therefore the magnitude of the interference could not be determined.

HIL Interference:

The AMM method was evaluated for interference according to CLSI/NCCLS EP7-A2.8 Bias is the difference in the results between the control sample (without the interferent) and the test sample (contains the interferent) expressed in percent. Bias exceeding 10% is considered interference.

Substance tested	Substance Concentration	AMM µmol/L	Bias %
Hemoglobin (hemolysate)	75 mg/dL	50	12
	500 mg/dL	250	17
Bilirubin (unconjugated)	80 mg/dL	50, 250	<10
Bilirubin (conjugated)	60 mg/dL	50	-16
	80 mg/dL	250	<10
Lipemia Intralipid®	50 mg/dL	50, 250	13, <10

14.4 Clinical Sensitivity/Specificity/Predictive Values
 None

15. SAFETY

The employee has direct responsibility to avoid injury and illness at work. Nearly all harmful exposures to infectious substances and chemicals, and other injuries, can be avoided with effective training and consistent safe work practices.

Become familiar with the Environmental Health and Safety (EHS) Manual to learn the requirements on working safely and protecting the environment from harm.

Although lab work typically focuses on the hazards of working with specimens and chemicals, we must also control other important hazards.

- Slips, trips, and falls cause many serious injuries. Please ensure that spills are cleaned quickly (to avoid slippery floors) and that you can see and avoid obstacles in your path.
- Ergonomic injuries result from performing tasks with too much repetition, force, or awkward position. Ergonomic injuries include strains and back injuries. Learn about ergonomic hazards and how to prevent this type of injury.
- Scratches, lacerations, and needlesticks can result in serious health consequences. Attempt to find ways to eliminate your risk when working with sharp materials.

Report all accidents and injuries immediately to your supervisor or the business unit Environmental Health and Safety Manager or Specialist.

Form revised 3/31/00

16. RELATED DOCUMENTS

1. Dimension Vista® Clinical Chemistry System Operator's Manual
2. Dimension Vista® Calibration/Verification Procedure
3. Dimension Vista® Cal Accept Guidelines
4. Dimension Vista® Calibration summary
5. Dimension Vista® Sample Processing, Startup and Maintenance procedure
6. Laboratory Quality Control Program
7. QC Schedule for Siemens Dimension Vista®
8. Laboratory Safety Manual
9. Material Safety Data Sheets (MSDS)
10. Dimension Vista® Limits Chart (AG.F200)
11. Quest Diagnostics Records Management Procedure
12. Dimension Vista® System Error Messages Chart
13. Centrifuge Use, Maintenance and Functions Checks (Lab policy)
14. Hemolysis, Icteria and Lipemia Interference (Lab policy)
15. Repeat Testing Requirement (Lab policy)
16. Current Allowable Total Error Specifications at http://questnet1.qdx.com/Business_Groups/Medical/qc/docs/qc_bpt_tea.xls
18. Current package insert AMM Flex® Reagent Cartridge K3119

17. REFERENCES

1. Package Insert, AMM Flex® Reagent Cartridge K3119, Siemens Healthcare Diagnostics Inc., 04/08/2013.
2. Package Insert, CHEM 3 CAL, Siemens Healthcare Diagnostics Inc., 12/2012.
3. Package Insert, Liquichek Ethanol/Ammonia Controls, Bio-Rad Laboratories, 03/2013

18. REVISION HISTORY

Version	Date	Section	Reason	Reviser	Approval

19. ADDENDA

None

Technical SOP

Title	Carbon Dioxide by Dimension Vista® System	
Prepared by	Ashkan Chini	Date: 6/25/2012
Owner	Robert SanLuis	Date: 2/5/2014

Laboratory Approval		Local Effective Date:
Print Name and Title	Signature	Date
<i>Refer to the electronic signature page for approval and approval dates.</i>		

Review		
Print Name	Signature	Date

Form revised 3/02/2007

TABLE OF CONTENTS

1.	Test Information.....	2
2.	Analytical Principle.....	3
3.	Specimen Requirements.....	3
4.	Reagents.....	4
5.	Calibrators/Standards.....	5
6.	Quality Control.....	6
7.	Equipment And Supplies.....	9
8.	Procedure.....	9
9.	Calculations.....	10
10.	Reporting Results And Repeat Criteria.....	10
11.	Expected Values.....	12
12.	Clinical Significance.....	12
13.	Procedure Notes.....	12
14.	Limitations Of Method.....	13
15.	Safety.....	14
16.	Related Documents.....	14
17.	References.....	14
18.	Revision History.....	15
19.	Addenda.....	15

1. TEST INFORMATION

Assay	Method/Instrument	Local Code
Carbon Dioxide	Dimension Vista® System	CO2


Synonyms/Abbreviations
ECO2, CO2, included in Batteries/Packages: BMP, COMP, LYTE and RENP

Department
Chemistry

Form revised 3/02/2007

2. ANALYTICAL PRINCIPLE

The enzymatic carbonate method for the Dimension Vista® System employs a phosphoenolpyruvate carboxylase-malate dehydrogenase coupled enzymatic reaction and a stable analog of the cofactor NADH. The bicarbonate anion reacts with phosphoenolpyruvate in the presence of phosphoenolpyruvate carboxylase (PEPC) and Mg⁺⁺ to form oxaloacetate and inorganic phosphate (Pi). The oxaloacetate is reduced to malate by malate dehydrogenase (MDH) with simultaneous oxidation of the reduced form of an analog (aNADH) of the cofactor, NADH.



3. SPECIMEN REQUIREMENTS

3.1 Patient Preparation

Component	Special Notations
Fasting/Special Diets	N/A
Specimen Collection and/or Timing	Normal procedures for collecting and storing serum and plasma may be used for samples to be analyzed by this method.
Special Collection Procedures	N/A
Other	N/A

3.2 Specimen Type & Handling

Criteria	
Type -Preferred -Other Acceptable	Plasma (Heparin) Serum
Collection Container	Plasma: Green top tube Serum: Red top tube, Serum separator tube (SST)
Volume - Optimum - Minimum	1.0 mL 0.5 mL
Transport Container and Temperature	Collection container or Plastic vial at room temperature
Stability & Storage Requirements	Room Temperature: 8 hours
	Refrigerated: 2 days
	Frozen: 6 months

Form revised 12/02/2010

Criteria	
Timing Considerations	Total Carbon Dioxide concentration may be decreased by 6 mmol/L when uncapped specimens are exposed to the air for one hour. Serum or plasma should be physically separated from cells as soon as possible with a maximum limit of two hours from the time of collection.
Unacceptable Specimens & Actions to Take	Specimens that are unlabeled, improperly labeled, or those that do not meet the stated criteria are unacceptable. Request a recollection and credit the test with the appropriate LIS English text code for "test not performed" message. Examples: Quantity not sufficient-QNS; Wrong collection-UNAC. Document the request for recollection in the LIS.
Compromising Physical Characteristics	Gross hemolysis. Reject sample and request a recollection. Credit the test with the appropriate LIS English text code explanation of HMT (Specimen markedly hemolyzed)
Other Considerations	Allow Red Top or SST to clot completely prior to centrifugation.

4. REAGENTS

Refer to the Material Safety Data Sheet (MSDS) supplied with the reagents for complete safety hazards. Refer to the section in this procedure covering "SAFETY" for additional information.

4.1 Reagent Summary

Reagents	Supplier & Catalog Number
Carbon Dioxide	Siemens, Flex® reagent cartridge, Cat. No. K1137

4.2 Reagent Preparation and Storage

NOTES: Date and initial all reagents upon opening. Each container must be labeled with (1) substance name, (2) lot number, (3) date of preparation, (4) expiration date, (5) initials of tech, (6) any special storage instructions; check for visible signs of degradation. Refer to the Material Safety Data Sheet (MSDS) for a complete description of hazards. If a specific hazard is present, it will be noted in this procedure when the hazard is first encountered in a procedural step.

Reagent	Carbon Dioxide
Container	Reagent cartridge
Storage	Store at 2-8° C
Stability	• Reagent is stable until expiration date stamped on the reagent

Form revised 12/02/2010

	cartridges. <ul style="list-style-type: none"> Sealed wells on the instrument are stable for 30 days. Once wells 1 - 12 have been entered by the instrument, they are stable for 1 day.
Preparation	All reagents are liquid and ready to use.

5. CALIBRATORS/STANDARDS

5.1 Calibrators/Standards Used

Calibrator	Supplier and Catalog Number
CHEM 3 CAL	Siemens Dimension Vista®, Cat. No. KC130A

5.2 Calibrator Preparation and Storage

NOTE: Date and initial all calibrators upon opening. Each container must be labeled with (1) substance name, (2) lot number, (3) date of preparation, (4) expiration date, (5) any special storage instructions; check for visible signs of degradation. When placed onboard the analyzer, the instrument captures the date / time loaded and calculates and tracks the opened expiration.

Calibrator	CHEM 3 CAL
Preparation	Calibrator is ready for use. No preparation is required.
Storage/Stability	<ul style="list-style-type: none"> Store at 2 – 8 ° C Unopened calibrator is stable until expiration date stamped on the box. Opened Calibrator: once the stopper of the vial is punctured, assigned values are stable for 24 hours when stored on board the Dimension Vista System.

5.3 Calibration Parameter

Criteria	Special Notations
Reference Material	CHEM 3 CAL
Assay Range	1 – 45 mmol/L
Suggested Calibration Level	See Reagent Package Insert for lot specific assigned values in mmol/L
Frequency	<ul style="list-style-type: none"> Every new reagent cartridge lot. Every 90 days for any one lot When major maintenance is performed on the analyzer. When control data indicates a significant shift in assay.
Calibration Scheme	2 levels, n = 5

From revised 12/02/2010

5.4 Calibration Procedure

Auto Calibration:

- Place the required calibrator vials in a carrier. Make sure the barcode labels are entirely visible through the slots.
- Place the carrier in the loading area.
- Position the carrier with the labels facing away from the user.
- Press the **Load** button.
- Automatic calibration requires that calibrators be on the instrument. As the time for processing approaches, the instrument reviews onboard inventory for the appropriate calibrators.

Manual Calibration:

- Verify that calibrators and reagents are in inventory on the instrument.
- Press **System > Method Summary > Calibration**.
- Select a method from the sidebar menu. Press the **Order Calibration** button on the screen.
- Verify that the information on the screen is correct. Verify that the calibrator lot is correct using the drop-down menu.
 - When calibrating using Vials press **OK**.
 - When calibrating using Cups, check the Use Cups box. This displays the rack and cup position fields. For additional cups use the positions in ascending order. Be sure to use the number of calibration levels and cups as specified in the method IFU. Scan the rack barcode and place calibrator cups in an adapter in position 1 on a rack. Press **OK** and load the rack on the instrument.
- The status field in the calibration screen changes sequentially to Awaiting Scheduling, Preparing Calibrators and Processing.

5.5 Tolerance Limits

IF.....	THEN.....
If result fall within assay-specific specification, and QC values are within acceptable limits,	proceed with analysis
If result falls outside assay-specific specification, or QC values are out of Acceptable limits,	troubleshoot the assay and/or instrument and repeat calibration

6. QUALITY CONTROL

6.1 Controls Used

Controls	Supplier and Catalog Number
Liquichek™ Unassayed Chemistry Control Levels 1 and 2	Bio-Rad Laboratories Cat. No. 691 and 692

From revised 12/02/2010

6.2 Control Preparation and Storage

NOTE: Date and initial all controls upon opening. Each container should be labeled with (1) substance name, (2) lot number, (3) date of preparation, (4) expiration date, (5) initials of tech, and (6) any special storage instructions; check for visible signs of degradation. A barcode label is produced and placed on the vial.

Control	Liquichek Unassayed Chemistry Controls, Level 1 and 2
Preparation	Allow the frozen control to stand at room temperature (18-25°C) until completely thawed. Swirl the contents gently to ensure homogeneity. (Do not use a mechanical mixer) Use immediately. After each use, promptly replace the stopper and return to 2-8°C storage.
Storage/Stability	Once the control is thawed, all analytes will be stable for 15 days at 2-8°C. Unthawed controls are stable until the expiration date at -20 to -70°C.

6.3 Frequency

Analyze all levels of QC material after every calibration and each day of testing (notated on the QC frequency sheets posted on the instruments).

Refer to the Dimension Vista® QC Schedule in the Laboratory policy Quality Control Program and in the Dimension Vista® Quick Reference Guide.

6.4 Tolerance Limits

Step	Action
1	Acceptable ranges for QC are programmed into the Laboratory Information System (LIS), and may be posted near the instrument for use during computer downtime.
2	Run Rejection Criteria <ul style="list-style-type: none"> Anytime the established parameters are exceeded (if one QC result exceeds 2 SD), the run is considered out of control (failed) and patient results must not be reported. The technologist must follow the procedure in the Laboratory QC Program to resolve the problem.
3	Corrective Action: <ul style="list-style-type: none"> All rejected runs must be effectively addressed through corrective action. Steps taken in response to QC failures must be documented. Patient samples in failed analytical runs must be <u>reanalyzed according to the Laboratory QC Program</u>. Supervisors may override rejection of partial or complete runs only with detailed

From revised 3/03/2007

Step	Action
	documentation and criteria for overrides that are approved by the Medical Director. Consult corrective action guidelines in Laboratory QC Program. Follow corrective action guidelines in the Laboratory QC Program. <ul style="list-style-type: none"> Corrective action documentation must follow the Laboratory Quality Control Program.
4	Review of QC <ul style="list-style-type: none"> QC must be reviewed weekly by the Group Lead or designee and monthly by the Supervisor/Manager or designee. If the SD and/or CV are greater than established ranges, investigate the cause for the imprecision and document implementation of corrective actions.

6.5 Review Patient Data

Technologist must review each result print-out with messages. Refer to the Dimension Vista® system manual “Error messages” section for troubleshooting. Check for unusual patterns, trends, or distributions in patient results (such as an unusually high percentage of abnormal results). Resolve any problems noted before issuing patient reports.

6.6 Documentation

- QC tolerance limits are programmed into the instrument and the LIS. The LIS calculates cumulative mean, SD and CV and stores all information for easy retrieval.
- Quality control records are reviewed daily at the bench, weekly by the Group Lead or designee, and monthly by the Supervisor/Manager or designee.
- Refer to complete policies and procedures for QC documentation and for record retention requirements in the Laboratory QC Program.

6.7 Quality Assurance Program

- Each new lot number of reagent or new shipment of the same lot of reagent must be tested with external control materials and previously analyzed samples. Performance of the new lot must be equivalent to the previous lot; utilize published TEA for acceptability criteria.
- Training must be successfully completed and documented prior to performing this test. This procedure must be incorporated into the departmental competency assessment program.
- The laboratory participates in CAP proficiency testing. All proficiency testing materials must be treated in the same manner as patient samples.
- Monthly QC must be presented to the Medical Director or designee for review and signature.

From revised 3/03/2007

- Monthly QC mean and SD are sent to Bio-Rad Laboratories for peer group comparison.
- Consult the Laboratory QC Program for complete details.

7. EQUIPMENT and SUPPLIES

7.1 Assay Platform

Dimension Vista® System

7.2 Equipment

- Refrigerator capable of sustaining 2–8°C.
- Freezer capable of sustaining range not to exceed -20 to -70°C.
- Centrifuge

7.3 Supplies

- Aliquot Plates
- System Fluids
- Assorted calibrated pipettes (MLA or equivalent) and disposable tips

8. PROCEDURE

CO2 Flex® reagent cartridge Cat. No. K1137 is required to perform this test.

Carbon Dioxide is performed on the Dimension Vista® System after the method is calibrated (see Reference Material in Calibration section) and Quality Controls are acceptable.

NOTE: For all procedures involving specimens, buttoned lab coats, gloves, and face protection are required minimum personal protective equipment. Report all accidents to your supervisor.

The package insert for a new lot of kits must be reviewed for any changes before the kit is used. A current Package Insert is included as a Related Document.

8.1	Sample Processing
1.	A sample rack holding tubes or cups is placed on the rack input lane.
2.	The sample shuttle moves the rack to the barcode reader which identifies the rack and samples to the system.
3.	The rack moves into the sample server and to the rack positioner.
4.	At the same time, aliquot plates move from the aliquot loader into position.

From revised 12/02/2007

8.1	Sample Processing
5.	The aliquot probe aspirates the sample from the tubes or cups and dispenses it into the wells of the aliquot plates.
6.	After each aspirate-dispense action, the probe is thoroughly rinsed inside and out to prevent sample carryover.
7.	When sample aspiration is completed, the sample server moves the rack back to the sample shuttle, where it is placed on the output lane and can be removed by the operator.

8.2	Specimen Testing
1.	For QC placement and frequency, refer to the Dimension Vista® QC Schedule in the Laboratory QC Program.
2.	Follow the instructions, outlined in the Dimension Vista® Operator’s Manual
3.	The instrument reporting system contains error messages to warn the user of specific malfunctions. Results followed by such error messages should be held for follow-up. Refer to the Dimension Vista® system manual “Error messages” section for troubleshooting.
4.	Follow protocol in Section 10.5 “Repeat criteria and resulting” for samples with results above or below the Analytical Measurement Range (AMR). Investigate any failed delta result and repeat, if necessary.
5.	Append the appropriate English text code qualifier messages to any samples requiring a comment regarding sample quality and/or any other pertinent factors.

Test Conditions	
Sample Volume:	1.9 µL
Reagent Volume:	38 µL
Reaction Time:	1.9 minutes
Test Temperature:	37° C
Wavelength:	405 & 700 nm
Type of measurement:	Bichromatic rate

9. CALCULATIONS

The instrument automatically calculates the concentration of carbon dioxide in mmol/L.

10. REPORTING RESULTS AND REPEAT CRITERIA

10.1 Interpretation of Data

None required

From revised 12/02/2007

10.2 Rounding

No rounding is necessary. Instrument reports results as a whole number.

10.3 Units of Measure

mmol/L

10.4 Clinically Reportable Range (CRR)

1 – 90 mmol/L

10.5 Repeat Criteria and Resulting

All repeats must replicate the original result within the total allowable error (TEa) of the assay. Refer to TEa listing for specific information.

Values that fall within the AMR or CRR may be reported without repeat. Values that fall outside these ranges must be repeated.

IF the result is ...	THEN...
< 1 mmol/L	Assure there is sufficient sample devoid of bubbles, cellular debris, and/or fibrin clots. Report as: < 1 mmol/L
> 45 mmol/L	Manual Dilution: Using the primary tube, make the smallest dilution possible to bring the raw data within the AMR. Maximum allowable dilution: x 2 Diluent: Water Enter dilution factor as a whole number. Re-assay. Readout is corrected for dilution.
> 90 mmol/L	If the recommended dilution does not give results within the clinically reportable range, report as: "> 90 mmol/L-REP" Bring to the attention of your supervisor prior to releasing result.

Message	Code
Verified by repeat analysis	Append -REP to the result.

11. EXPECTED VALUES

11.1 Reference Ranges

Age	Female / Male
Adult (>18 years):	21 – 32 mmol/L
Pediatric:	
0 – 6 days	13-21
7 – 30 days	13-22
1 – 5 months	13-23
6 – 12 months	14-23
13 - 23 months	16-25
2 – 18 years	21-32

11.2 Critical Values

< 10 mmol/L

11.3 Priority 3 Limit(s)

None established

12. CLINICAL SIGNIFICANCE

Increase in serum CO₂ content for the most part reflects increase in serum bicarbonate concentration rather than dissolved CO₂ gas (which accounts for only a small fraction of the total). Increased serum bicarbonate is seen in compensated respiratory acidosis and in metabolic alkalosis. Diuretics (thiazides, ethacrynic acid, furosemide, mercurials), corticosteroids (in long term use), and laxatives (when abused) may cause increased bicarbonate.

Decrease in blood CO₂ is seen in metabolic acidosis and compensated respiratory alkalosis. Substances causing metabolic acidosis include ammonium chloride, acetazolamide, ethylene glycol, methanol, paraldehyde, and phenformin. Salicylate poisoning is characterized by early respiratory alkalosis followed by metabolic acidosis with attendant decreased bicarbonate.

13. PROCEDURE NOTES

- **FDA Status:** FDA Approved/cleared
- **Validated Test Modifications:** None

The instrument reporting system contains error messages to warn the operator of specific malfunctions. Any report slip containing such error messages should be held for follow-up. Refer to your Dimension Vista Operator's Guide.

The expected maximum observed standard deviations for repeatability using n = 5 replicates at the following carbon dioxide concentrations are:

CO2 Concentration	Acceptable S.D. Maximum
11 mmol/L	2 mmol/L
23 mmol/L	3 mmol/L

14. LIMITATIONS OF METHOD

14.1 Analytical Measurement Range (AMR)

1 – 45 mmol/L

14.2 Precision

Material	Mean mmol/L	Standard Deviation (%CV)	
		Repeatability	Within-Lab
Multiquel Control			
Level 1	12	0.4 (2.9)	0.6 (4.6)
Level 2	25	0.6 (2.4)	1.1 (4.5)

14.3 Interfering Substances

Turbidity at 3000 mg/dL decreases CO2 results by 35% at CO2 concentration of 28 mmol/L.

HIL Interference:

The CO2 method was evaluated for interference according to CLSI/NCCLS EP7-A2. Bias, defined as the difference between the control sample (does not contain interferent) and the test sample (contains interferent), is shown in the table below. Bias exceeding 10% is considered “interference”.

Substance tested	Substance Concentration	CO2 mmol/L	Bias %
Hemoglobin (hemolysate)	1000 mg/dL	30	<10
Bilirubin (unconjugated)	60 mg/dL	31	<10
Bilirubin (conjugated)	60 mg/dL	31	<10
Lipemia Intralipid®	1000 mg/dL	28	<10
	3000 mg/dL		-35

14.4 Clinical Sensitivity/Specificity/Predictive Values

Not available

From revised 2/02/2007

15. SAFETY

The employee has direct responsibility to avoid injury and illness at work. Nearly all harmful exposures to infectious substances and chemicals, and other injuries, can be avoided with effective training and consistent safe work practices.

Become familiar with the Environmental Health and Safety (EHS) Manual to learn the requirements on working safely and protecting the environment from harm. Although lab work typically focuses on the hazards of working with specimens and chemicals, we must also control other important hazards.

- Slips, trips, and falls cause many serious injuries. Please ensure that spills are cleaned quickly (to avoid slippery floors) and that you can see and avoid obstacles in your path.
- Ergonomic injuries result from performing tasks with too much repetition, force, or awkward position. Ergonomic injuries include strains and back injuries. Learn about ergonomic hazards and how to prevent this type of injury.
- Scratches, lacerations, and needlesticks can result in serious health consequences. Attempt to find ways to eliminate your risk when working with sharp materials.

Report all accidents and injuries immediately to your supervisor or the business unit Environmental Health and Safety Manager or Specialist.

16. RELATED DOCUMENTS

1. Dimension Vista® Clinical Chemistry System Operator’s Manual
2. Dimension Vista® Calibration/Verification Procedure
3. Dimension Vista® Cal Accept Guidelines
4. Dimension Vista® Calibration summary
5. Dimension Vista® Sample Processing, Startup and Maintenance procedure
6. Laboratory Quality Control Program
7. QC Schedule for Siemens Dimension Vista®
8. Laboratory Safety Manual
9. Material Safety Data Sheets (MSDS)
10. Dimension Vista® Limits Chart (AG.F200)
11. Quest Diagnostics Records Management Procedure
12. Dimension Vista® System Error Messages Chart
13. Centrifuge Use, Maintenance and Functions Checks (Lab policy)
14. Hemolysis, Icteria and Lipemia Interference (Lab policy)
15. Repeat Testing Requirement (Lab policy)
16. Current Allowable Total Error Specifications at http://questnet1.qdx.com/Business_Groups/Medical/qc/docs/qc_bpt_tea.xls
17. Current package insert CO2 Flex® Reagent Cartridge K1137

From revised 2/02/2007

17. REFERENCES

1. Ghoshal, Amit K. and Soldin, Steven J., Evaluation of the Dade Behring Dimension® RxL: Integrated chemistry system-pediatric reference ranges. Clinica Chimica Acta 2003; 331:144.
2. Package Insert, CO2 Flex® Reagent Cartridge K1137, Siemens Healthcare Diagnostics Inc., 04/29/2013.
3. Package Insert, CHEM 3 CAL, Siemens Healthcare Diagnostics Inc., 12/2012.
4. Package Insert, Unassayed Liquichek Chemistry Controls, Bio-Rad Laboratories, 01/2013.

18. REVISION HISTORY

Version	Date	Section	Reason	Reviser	Approval
000	2/5/14		Update owner	L Barrett	R SanLuis
000	2/5/14	5	Change in Calibrator, update information	A Chini	R SanLuis
000	2/5/14	16	Update titles	L Barrett	R SanLuis
000	2/5/14	17	Update calibrator package insert	A Chini	R SanLuis
000	2/5/14	Footer	Version # leading zero's dropped due to new EDCS in use as of 10/7/13.	L Barrett	R SanLuis

19. ADDENDA

None

Form revised 2/02/2017