

TRAINING UPDATE

Lab Location: GEC
Department: Core

Date Distributed: 4/4/2014
Due Date: 4/21/2014
Implementation: 4/21/2014

DESCRIPTION OF PROCEDURE REVISION

Name of procedure:
XPAND QC SCHEDULE AG.F210v1
Description of change(s):
Added schedule for 2 nd instrument This revised form will be implemented on April 15, 2014

Document your compliance with this training update by taking the quiz in the MTS system.

XPAND 2 QC SCHEDULE

CHECK SCREEN BEFORE GOING TO NEW SAMPLE, WHEN PROGRAMMING QC

0500	0900	1300	1700	2100	0200
QCBIO1S, QCBIO2S QNBIL2 QCCSF1, QCCSF2 QCCK1, QCCK2, QCCK3 QCTD1S, QCTD2S, QCTD3S QCLIQ1, QCLIQ2, QCLIQ3 QCCRP1, QCCRP2, QCCRP3	QCLYTES2	QCBIO1S	QCLYTES1	QCBIO2S	QCLYTES2
		QCCK1		QCCK3	

BIORAD ETHANOL/AMMONIA CONTROLS (Levels 1, 2 & 3) - 20 day Open Stability (OS)	
QCLIQ1— 3 drops QCLIQ2— 3 drops QCLIQ3— 3 drops	ALC
BIORAD LIQUID UNASSAYED CHEMISTRY CONTROLS (Levels 1 and 2) - 6 days (OS)	
QCBIO1S—8 drops QCBIO2S—8 drops	ALB, ALP, ALT, AMY, AST, BUN, CA, CK, CREA, DBIL, ECO2, GLUC, LIPL, LYLES, MG, TBIL, TP, LA
BIORAD PEDIATRIC CONTROL (Level 2 only) - 5 days (OS)	
QCNBIL2—3 drops	TBIL AND DBIL
BIORAD LIQUICHEK CARDIAC MARKER PLUS CONTROLS (Levels 1, 2 and 3) - 8 days (OS)	
QCCK1—8 drops QCCK2—8 drops QCCK3—8 drops	CTNI & MMB
BIORAD LIQUICHEK IMMUNOASSY PLUS CONTROLS (Levels 1, 2 and 3) - 14 days (OS)	
QCTD1S—5 drops QCTD2S—5 drops QCTD3S—5 drops	ACTM, HCG, SAL & TSH
BIORAD LIQUICHEK SPINAL FLD CONT (RUN AS URINEQC1 & URINEQC2) – 30 days (OS)	
QCCSF1—3 drops QCCSF2—3 drops	UCFP & GLUC
BIORAD IMMUNOLOGY CONTOLS (Levels 1, 2, & 3) - 30 days (OS)	
QCCRP1—3 drops QCCRP2—3 drops QCCRP3—3 drops	CRP

REPEAT QC PROTOCOL

If you have to repeat QC, it must be entered in the following way:

Name: BII01 (or the level you are testing – MAKE SURE TO LEAVE QC FROM THE NAME)

Priority: QC

Sample No: C-BIO1S (or the level you are testing)

Fluid: Serum QC1 (or the level you are testing)

ONLY SELECT THE TESTS NEEDED

MAKE SURE YOU HAVE **NOT** PROGRAMMED THE NAME STARTING WITH “QC” SUCH AS BIO1S

Note: Repeat QC should never be retrieved.

XPAND 1 QC SCHEDULE

CHECK SCREEN BEFORE GOING TO NEW SAMPLE, WHEN PROGRAMMING QC

0900	1300	1700	2100	0200	0500
QCBIO1S, QCBIO2S QNBIL2 QCCSF1, QCCSF2 QCCK1, QCCK2, QCCK3 QCTD1S, QCTD2S, QCTD3S QCLIQ1, QCLIQ2, QCLIQ3 QCCRP1, QCCRP2, QCCRP3	QCLYTES2	QCBIO1S	QCLYTES1	QCBIO2S	QCLYTES2
		QCCK1		QCCK3	

BIORAD ETHANOL/AMMONIA CONTROLS (Levels 1, 2 & 3) - 20 day Open Stability (OS)	
QCLIQ1— 3 drops QCLIQ2— 3 drops QCLIQ3— 3 drops	ALC
BIORAD LIQUID UNASSAYED CHEMISTRY CONTROLS (Levels 1 and 2) - 6 days (OS)	
QCBIO1S—8 drops QCBIO2S—8 drops	ALB, ALP, ALT, AMY, AST, BUN, CA, CK, CREA, DBIL, ECO2, GLUC, LIPL, LYLES, MG, TBIL, TP, LA
BIORAD PEDIATRIC CONTROL (Level 2 only) - 5 days (OS)	
QCNBIL2—3 drops	TBIL AND DBIL
BIORAD LIQUICHEK CARDIAC MARKER PLUS CONTROLS (Levels 1, 2 and 3) - 8 days (OS)	
QCCK1—8 drops QCCK2—8 drops QCCK3—8 drops	CTNI & MMB
BIORAD LIQUICHEK IMMUNOASSY PLUS CONTROLS (Levels 1, 2 and 3) - 14 days (OS)	
QCTD1S—5 drops QCTD2S—5 drops QCTD3S—5 drops	ACTM, HCG, SAL & TSH
BIORAD LIQUICHEK SPINAL FLD CONT (RUN AS URINEQC1 & URINEQC2) – 30 days (OS)	
QCCSF1—3 drops QCCSF2—3 drops	UCFP & GLUC
BIORAD IMMUNOLOGY CONTROLS (Levels 1, 2, & 3) - 30 days (OS)	
QCCRP1—3 drops QCCRP2—3 drops QCCRP3—3 drops	CRP

REPEAT QC PROTOCOL

If you have to repeat QC, it must be entered in the following way:

Name: BII01 (or the level you are testing – MAKE SURE TO LEAVE QC FROM THE NAME)

Priority: QC

Sample No: C-BIO1S (or the level you are testing)

Fluid: Serum QC1 (or the level you are testing)

ONLY SELECT THE TESTS NEEDED

MAKE SURE YOU HAVE **NOT** PROGRAMMED THE NAME STARTING WITH “QC” SUCH AS BIO1S

Note: Repeat QC should never be retrieved.