

TRAINING UPDATE

Lab Location:	SGAH & WAH	Date Distributed:	4/10/2014
Department:	Core Group Leads & TIC staff	Due Date:	5/8/2014
		Implementation:	5/8/2014

DESCRIPTION OF PROCEDURE REVISION

Name of procedure:
Troubleshooting Communication Issues with Devices Connected to a Silex Server SGAH.LIS55, WAH.LIS55 v1
Description of change(s):
<p>Sections 1,2 & 4: replace instrument with device</p> <p>Section 3: add in charge tech</p> <p>Section 5: update background to include GEC and printers, add reference to Silex Server list</p> <p>Section 9: add Silex Server List</p> <p>This revised SOP will be implemented on May 8, 2014</p>

Document your compliance with this training update by taking the quiz in the MTS system.

Approved draft for training all sites (version 1)

Non-Technical SOP

Title	Troubleshooting Communication Issues with Devices Connected to a Silex Server	
Prepared by	Marie Sabonis	Date: 2/16/2011
Owner	Marie Sabonis	Date: 2/16/2011

Laboratory Approval		
Print Name and Title	Signature	Date
<i>Refer to the electronic signature page for approval and approval dates.</i>		
Local Issue Date:		Local Effective Date:

Review:		
Print Name	Signature	Date

TABLE OF CONTENTS

1. PURPOSE.....	3
2. SCOPE	3
3. RESPONSIBILITY.....	3
4. DEFINITIONS.....	3
5. PROCEDURE.....	3
6. RELATED DOCUMENTS	4
7. REFERENCES	4
8. REVISION HISTORY.....	4
9. ADDENDA AND APPENDICES.....	4

1. PURPOSE

This procedure provides guidelines to troubleshoot device interfaces that connect to the LIS via Silex Servers.

2. SCOPE

This procedure applies to devices that are interfaced to the LIS via Silex Servers.

3. RESPONSIBILITY

Group Leads and In-charge Techs perform this procedure

4. DEFINITIONS

Silex Server - Communication server; conduit between the device and the LIS (Sunquest)

5. PROCEDURE

A. Background

Some instruments and printers at Washington Adventist, Shady Grove Adventist and Germantown Emergency Center laboratories are connected to the LIS (Sunquest) via Silex Servers. The Silex servers are connected directly to the device and provide a communication tunnel to the LIS (Sunquest) for data to be sent bidirectionally.

Periodically, we may experience issues with the device not communicating to the LIS (Sunquest). Therefore, orders may not download to the instrument, results may not upload to the LIS (Sunquest) and printers may not print.

This may be caused by communication inadvertently turned off on the instrument or device, the Silex Server, the LIS (Sunquest) or network.

B. Process

1. Refer to Silex Server List to determine if device is connected to a Silex Server
2. Check to see if transmission to LIS (Sunquest) is turned off (see attachment)
3. Check status of Silex Server (see attachment)
4. If transmission is on and Silex Server status is OK, the Group Lead or Tech in Charge is responsible for contacting LIS On-Call staff.

6. RELATED DOCUMENTS

None

7. REFERENCES

Silex Server reference manual

8. REVISION HISTORY

Version	Date	Reason for Revision	Revised By	Approved By
000	3/31/14	Sections 1,2 & 4: replace instrument with device Section 3: add in charge tech Section 5: update background to include GEC and printers, add reference to Silex Server list Section 9: add Silex Server List Footer: version # leading zero's dropped due to new EDCS in use as of 10/7/13.	M. Sabonis	M. Sabonis


9. ADDENDA AND APPENDICES

- A. Job Aid for Silex Server (Troubleshooting Instrument Interface Issues)
- B. Silex Server List

Addendum A

Silex Server (Troubleshooting instrument interface issues)

The Silex Server provides the communication tunnel between the instrument and the LIS(Sunquest)



Function	State	Status	
Power	Color	On	Receiving power
	No color	Off	NOT Receiving power
Network Status			
10	Color	On or Blinking	10base-T network active
	No color	Off	10base-T network NOT active
100	Color	On or Blinking	100base-T network active
	No color	off	100base-T network NOT active

Troubleshooting Instrument Interface (follow steps 1-3)

1. Check instrument to determine if Transmission to LIS(Sunquest) is ON
2. If transmission is **ON** then check lights on Silex Server. The power **AND** the 10 or 100 network light should both be lit. (Refer to chart above)
 - a. If power light is not ON then locate the power supply connection and try the following:
 - Disconnect the power support from the Silex server and then reconnect
 - Unplug the power supply from the electrical outlet and reconnect
 - Reconnect power supply to a different electrical outlet.
 - b. If "10" -or- "100" network light is not ON then locate the network cable connection to the Silex
 - Disconnect the network cable from Silex Server and reconnect making sure that it is securely seated in the jack
 - Disconnect the network cable from the network jack(follow cable from Silex Server to wall jack) and reconnect
3. If both lights are on (power and network) then reset the Silex Server by disconnecting/reconnecting the power supply from the Silex.
 - a. Disconnect the power supply on the Silex Server
 - b. Count to 10
 - c. Reconnect power supply and Wait 5 minutes before try interface again
4. If steps 1-3 do not resolve the instrument interface issue then contact LIS On-Call Support

Marie Sabonis 3/13/2014

Addendum B

Silex Server List

HID	Com type	IP Address	Hardware	SQ ID	Comments
Shady Grove Adventist Hospital					
Phlebotomy / Accessioning					
SGAH	Silex Server	172.16.11.230	Intermec Barcode Printer-Central lab SGAH	323	
SGAH	Silex Server	172.16.11.231	Intermec Barcode Printer- Central Lab SGAH	324	
SGAH	Silex Server	172.16.116.56	Intermec Barcode Printer-CS 3rd floor SGAH	305	
SGAH	Silex Server	172.16.117.181	Intermec Barcode Printer-Phleb SGAH-RT	311	
SGAH	Silex Server	172.16.117.180	Intermec Barcode Printer- Phleb SGAH	312	
SGAH	Silex Server	172.16.11.228	Intermec Barcode Printer-Ref Lab	325	
Instruments					
SGAH	Silex Server	172.16.10.72	Bactec 9000(Microbiology instrument)	BACS1	
SGAH	Silex Server	labtroy02:9100 (172.16.11.50)	Stago (Coag Instrument)	ST1S	
SGAH	Silex Server	labtroy01:9100 (172.16.116.90)	Stago(Coag instrument)	ST2S	
SGAH	Silex Server	labtroyNN:9100(172.16.116.80)	IRIS (Urine instrument)	IRIS1S	
SGAH	Silex Server	172.16.10.63	Coulter LH (Hematology instrument)	LH2S	
SGAH	Silex Server	172.16.10.62	Coulter LH (Hematology instrument)	LH1S	
SGAH	Silex Server	172.16.9.55:10002	Galileo Echo (Blood Bank Instrument)	ECOS1	
Others					
SGAH	Silex Server	172.16.11.229	Lexmark dot matrix printer (Blood Bank unit tag printer)	331	
Germantown Emergency Center					
GEC	Silex Server	10.12.2.62	Coulter LH(Hematology Analyzer)	LHG1	
Washington Adventist Hospital					
Phlebotomy / Accessioning					
WAH	Silex Server	172.16.5.214	Intermec Barcode Printer-CS front desk	105	
Instruments					
WAH	Silex Server	labtroy02:9100 172.16.4.27	Stago	STAW2	
WAH	Silex Server	labtroyw01:9100 172.16.4.26	Stago	STAW1	
WAH	Silex Server	172.16.107.60:10001	Galileo Echo	ECOW1	

revised 3/31/00