

TRAINING UPDATE

Lab Location: SGAH and WAH
Department: Blood Bank

Date Implemented: 5.30.2014
Due Date: 6.8.2014

DESCRIPTION OF PROCEDURE REVISION

Name of procedure:

Reconstituted Whole Blood Worksheet

Description of change(s):

New worksheet to help with the calculations for reconstituted whole blood.



Reconstituted Whole Blood Worksheet

Infant's Information

Mom's Information

Name			
Medical Record Number			
ABO/Rh			
Antibody History	N/A		

Product Information

Date/Time of Exchange		<input type="checkbox"/> Needed STAT
Desired Hematocrit		(A)
Desired Volume		(B)

Red Blood Cell Unit Information

Unit Hematocrit (RBC E0226 Product)	To convert a percent to decimal, move the decimal place 2 places to the left. Example: 70% =0.70	(C)
Unit Volume (mL) (RBC E0226 Product)		(D)

Calculate the volume of plasma needed for reconstitution:

Step 1: Multiply the red cell unit hematocrit (C) times the red cell unit volume (D)

$$(C) \underline{\hspace{2cm}} \times (D) \underline{\hspace{2cm}} \text{ mL} = (E) \underline{\hspace{2cm}} \text{ mL}$$

Step 2: Divide the product (E) obtained in step 1 by the desired hematocrit (A)

$$(E) \underline{\hspace{2cm}} \text{ mL} \div (A) \underline{\hspace{2cm}} = (F) \underline{\hspace{2cm}} \text{ mL}$$

Step 3: Subtract the Quotient (F) obtained in step 2 from the red cell unit volume (D)

$$(D) \underline{\hspace{2cm}} \text{ mL} - (F) \underline{\hspace{2cm}} \text{ mL} = \boxed{\hspace{2cm}} \text{ mL} \rightarrow \text{Volume of plasma to add}$$

Final Hematocrit (RWB E5747 Product)	To convert a percent to decimal, move the decimal place 2 places to the left. Example: 70% =0.70
Final Volume (mL) (RWB E5747 Product)	(G)

Step 4: Subtract the Final Volume (G) from the Desired Volume (B)

$$(B) \underline{\hspace{2cm}} \text{ mL} - (G) \underline{\hspace{2cm}} \text{ mL} = (H) \underline{\hspace{2cm}} \text{ mL}$$

- If H is a positive number, issue the unit as is.
- If H is a negative number, notify the patient care area that the volume is lower than requested and ask if a second unit is needed.

Document: AG.F.288[0] Status: INWORKS_Effective: 6/29/2014, Check Version Before Use