

TRAINING UPDATE

Lab Location: SGAH and WAH **Date Implemented:** 5.30.2014
Department: Blood Bank **Due Date:** 6.8.2014

DESCRIPTION OF PROCEDURE REVISION

Name of procedure:

Reconstituted Whole Blood Preparation for Neonatal Exchange Transfusion

Description of change(s):

- | |
|--|
| <ol style="list-style-type: none">1. Changed name of SOP....This used to be "Neonatal Exchange Transfusion."2. Updated LIS instructions for the Sunquest v6.4 upgrade.3. Removed the need to document the hematocrit of the pRBC in the LIS (we still get the Hct, but only write it on the worksheet). We WILL document the post-reconstitution HCT in the LIS.4. Updated the wording for clarity. |
|--|

Non-Technical SOP

Title	Reconstituted Whole Blood Preparation for Neonatal Exchange Transfusion	
Prepared by	Stephanie Codina	Date: 05.29.2014
Owner	Stephanie Codina	Date: 05.29.2014

Laboratory Approval		
Print Name and Title	Signature	Date
<i>Refer to the electronic signature page for approval and approval dates.</i>		
Local Issue Date:		Local Effective Date:

Review:		
Print Name	Signature	Date

Document: SGAH.BB877[0] Status: INWORKS, Effective: 6/29/2014, Check Version Before Use

Form revised 3/3/00

TABLE OF CONTENTS

1. PURPOSE.....	2
2. SCOPE.....	2
3. RESPONSIBILITY.....	2
4. DEFINITIONS.....	2
5. PROCEDURE.....	3
6. RELATED DOCUMENTS	9
7. REFERENCES	9
8. REVISION HISTORY.....	10
9. ADDENDA AND APPENDICES.....	10

1. PURPOSE

Neonatal exchange transfusion is the treatment of choice for Hemolytic Disease of the Newborn (HDN), hyperbilirubinemia, disseminated intravascular coagulation (DIC), and occasionally the elimination of toxins, drugs, and chemicals in neonates. The procedure consists of replacing one to two whole blood volumes and has several desired effects. Removal of the infant's blood reduces antibody coated red cells, unconjugated bilirubin, and the number of unbound antibody molecules available to bind newly-formed antigen-positive red blood cells. Reconstituted whole blood is used for this procedure. The red cells used for replacement are compatible with the infant and/or maternal specimen and provide increased oxygen-carrying capacity. The plasma restores albumin and coagulation factors.

2. SCOPE

This procedure applies to the preparation of reconstituted whole blood for any neonate who requires an exchange transfusion.

3. RESPONSIBILITY

All blood bank staff members must understand and adhere to this procedure for the preparation of reconstituted whole blood for a neonatal exchange procedure.

4. DEFINITIONS

Neonate: An infant <4 months of age.

Document: SGAH.BB877[0] Status: INWORKS, Effective: 6/29/2014, Check Version Before Use

5. PROCEDURE

This procedure must be carried out in steps to ensure proper preparation in the LIS.

- A. Select a red cell unit.
- B. Select a plasma unit.
- C. Thaw the plasma unit.
- D. Irradiate the red cell unit.
- E. Aliquot the plasma unit.
- F. Combine the red cell and plasma together.

Step	Action
1	<p>The patient care area will order reconstituted whole blood using the order, "TWBNEO." Blood bank staff members will need the following information to process the order.</p> <ul style="list-style-type: none"> A. Infant's name and medical record number (Note: If blood products for an exchange transfusion are requested before an infant has been delivered, order the blood product using the mother's medical record number and type a comment stating the exchange is for an infant that has not been delivered). B. Mother's name and medical record number, if available. Often the baby is transferred to us from another hospital. C. Date and time of exchange procedure D. Specimens for T&S on the mother and TSNEO on the infant (if not already available) E. The desired hematocrit of the final product (usually 50-60%). An exact hematocrit is very difficult to achieve, so a range is generally used.
2	<p>Document the following information on the Reconstituted Whole Blood Worksheet.</p> <ul style="list-style-type: none"> A. Infant's full name and medical record number. B. Mom's full name and medical record number. C. Date and time of expected transfusion (write "STAT" if blood is needed as soon as possible). D. Desired volume and hematocrit of whole blood product.

Form revised 3/31/08

Document: SGAH.BB877[0] Status: INWORKS, Effective: 6/29/2014, Check Version Before Use

Step	Action
3	<p>Order blood products from the supplier (if needed) and emphasize delivery time. Blood products should meet the following criteria:</p> <ul style="list-style-type: none"> A. Group O (Rh-negative if the recipient is Rh-negative) B. Fresh red cells (<7 days old) to avoid high levels of potassium and to maximize red cell survival C. CPDA-1 (AS-3 red cells cannot be used for this procedure) D. CMV-seronegative E. Irradiated (note: irradiation should be performed immediately before reconstitution to minimize the potassium leak) F. Sickle-negative G. Negative for any antigens that correspond to clinically-significant maternal antibodies. H. AHG crossmatch compatible with the mother's plasma. If the mother's plasma is not available, the unit should be AHG crossmatch compatible with the infant's plasma.
4	<p>Obtain the supplies necessary to reconstitute the blood product.</p> <ul style="list-style-type: none"> A. Sterile welding device B. 150 mL transfer bag C. Heat sealer D. Test tubes E. Scale F. Hemostats
5	<p>Perform daily QC of the scale if needed.</p>

Form revised 3/31/00

Document: SGAH.BB877[0] Status: INWORKS, Effective: 6/29/2014, Check Version Before Use

Step	Action
6	<p>Obtain a pre-hematocrit of the original unit:</p> <ul style="list-style-type: none"> A. Mix the red cell unit manually by gently rotating it back-and-forth. B. Sterile dock a 150 ml transfer bag to the unit or use an attached satellite bag if available. C. Allow small sample of blood to flow into the tubing of the empty 150 ml sterile bag. D. Apply hemostat on the tubing toward the original unit once red cells fill up the tubing. E. Seal the tubing of the original unit first using the heat sealer leaving enough tubing for the next component prep. F. Make 3 or more segments, 2-3 inches long, label each segment with the unit number of donor identification number (DIN). Use of the labels found on the back of the original red cell unit is preferred. It is not necessary to label the aliquot bag since the labeled segments are used for unit hematocrit testing. G. Separate the segments from the original blood product unit. H. Pierce the end of 1-2 segments and drip the blood into a clean test tube labeled with the unit number. I. Deliver the sample to hematology and request STAT hematocrit testing. Instruct the hematology tech to test the specimen in duplicate and average the results. J. Convert the hematocrit to a decimal (Ex. 75% = 0.75). The hematocrit should be 0.65-0.80.
7	<p>Document the unit's hematocrit on the worksheet.</p>
8	<p>Select and thaw a unit of plasma per procedure.</p> <ul style="list-style-type: none"> A. Group AB plasma is preferred. B. If group AB plasma is not available, select a unit of plasma that is compatible with the infant and mother's types. C. Plasma must be thawed fresh for this procedure (do not use a unit of thawed plasma on the shelf).
9	<p>Approximate the volume of red cells in the primary red cell container.</p> <ul style="list-style-type: none"> A. Weigh the bag of red cells. B. Subtract 100g (the approximate weight of the bag, anticoagulant, and segments). C. Document the red cell volume on the worksheet.
10	<p>Irradiated the red cell per procedure.</p> <ul style="list-style-type: none"> A. Physically irradiate the red cell. B. Perform the blood component prep function to irradiate the red cell in Sunquest. C. Perform the blood label check in Sunquest.

Form revised 3/31/00

Document: SGAH.BB877[0] Status: INWORKS, Effective: 6/29/2014, Check Version Before Use

Step	Action
11	<p>Calculate the amount of thawed plasma to be added to the packed red cell unit for reconstitution per the following formula. Use the worksheet to aid in this calculation.</p> <p>Desired volume of whole blood unit = $\frac{(\text{Hct of original red cell}) \times (\text{Volume of original red cell in mL})}{(\text{Desired Hct of whole blood unit})}$</p> <p>Volume of plasma to add = (Volume of whole blood unit) – (Volume of original red cell)</p> <p>Key: Original red cell = the packed red blood cell that was pulled out of inventory Whole blood unit = reconstituted whole blood unit; the red cell with added plasma</p> <p>Example: You have 200mL of packed red blood cells with a hematocrit of 75%. The physician is requesting reconstituted whole blood with a final hematocrit of 55%.</p> <p>Volume of whole blood = $(0.75 \times 200) \div 0.55 = 273 \text{ mL}$</p> <p>Volume of plasma to add = $273 \text{ mL} - 200 \text{ mL} = 73 \text{ mL}$</p> <p>Add 73 mL of plasma to the red cell unit to obtain whole blood with a hematocrit of ~55%.</p>
12	<p>Prepare a plasma aliquot per procedure.</p> <ol style="list-style-type: none"> A. The volume of the aliquot will be the “volume of plasma to add” as calculated in the step above. B. The plasma aliquot must be prepared per procedure. <ol style="list-style-type: none"> a. Prepare the plasma aliquot. b. Perform the blood component preparation function to create the aliquot in Sunquest. c. Perform the blood label check in Sunquest.
13	<p>Add the plasma aliquot to the red cell unit.</p> <ol style="list-style-type: none"> A. Use the sterile connecting device to connect the thawed plasma aliquot to the red cell unit. B. Transfer the contents of the aliquotted plasma to the red cell unit. C. Weigh the bag. D. Subtract 100g (the approximate weight of the bag, anticoagulant, and segments) to obtain the new volume of the reconstituted whole blood unit. E. Document the volume on the worksheet.



Form revised 3/31/00

Document: SGAH.BB877[0] Status: INWORKS, Effective: 6/29/2014, Check Version Before Use

Step	Action
14	<p>Prepare segments on the new, reconstituted product.</p> <ul style="list-style-type: none"> A. Gently knead the bag to mix the contents. B. Allow the tubing to refill with the well-mixed reconstituted whole blood. C. Strip the blood in the remaining tubing three times back into the cell bag, mixing well each time. D. Heat seal the line beginning at the end of the tubing and moving towards the bag forming at least 4 segments. Place a double-seal next to the bag.
15	<p>Obtain the hematocrit of the new, reconstituted product.</p> <ul style="list-style-type: none"> A. Label a clean test tube with the unit number. B. Separate 1-2 segments from the original blood product unit. C. Pierce the segments and drain the contents into the clean, labeled test tube. D. Deliver the sample to hematology and request STAT hematocrit testing. E. Instruct the hematology tech to test the specimen in duplicate and average the results. F. Ensure the hematocrit is within 5% of the desired hematocrit. Notify a supervisor if discrepancies exist.
16	<p>Reconstitute the whole blood product in the LIS.</p> <ul style="list-style-type: none"> A. Access Sunquest function, "Blood Component Preparation." B. At the "Value" prompt, type "RE0224" then press the "tab" key. C. Press the tab key to default the current date and time as the reconstitution time. Enter the date and time of reconstitution if prepared at an earlier time (such as during a computer downtime). D. Click the "Continue" button. E. A second "Blood Component Prep" screen will appear. <ul style="list-style-type: none"> a. At the "Unit #" prompt, scan the unit number of the plasma aliquot product. b. At the "Component" prompt, scan the E code of the plasma aliquot. This will autofill both the component and division fields. c. The cursor will return to the "Unit #" prompt. Scan the unit number of the irradiated packed red blood cell (E0224) to be reconstituted. d. The LIS will prompt, "Pick input component." Select the E0224 product and click the "OK" button. e. You will see the QA failure message, "Input unit(s) expire before output unit." Acknowledge the QA failure.

Form revised 3/3/06

Document: SGAH.BB877[0] Status: INWORKS, Effective: 6/29/2014, Check Version Before Use

Step	Action
<p>16 Cont</p>	<p>F. Click on the yellow circle that contains the "O."</p> <p>a. Enter the volume of the whole blood unit. The volume should equal the combined volumes of the red cell unit and plasma aliquot as defined in the computer.</p> <p>b. Enter the expiration date of the reconstituted product. The reconstituted product will expire 24 hours from the time the plasma was aliquotted. The expiration date/time of the reconstituted whole blood unit should equal the expiration date/time of the plasma aliquot.</p>  <p>G. Click the "Save" button.</p> <p>H. A "Preview Output / New Units" screen will appear. Review the information to ensure accuracy, then click on the "finish" button to generate a new label for the reconstituted product.</p>
<p>17</p>	<p>Adhere the new label directly over the label of the red cell unit and perform a blood label check in Sunquest per procedure.</p>
<p>18</p>	<p>Document the post-hematocrit and volume in the LIS.</p> <p>A. Access function "Blood Product Testing."</p> <p>B. At the "Unit Number" prompt, scan the unit number of the original unit.</p> <p>C. At the "Component" prompt, scan the E code of the reconstituted unit. This will autofill both the Component and Division # fields.</p> <p>D. Click on the "Add" button.</p> <p>E. Click on the "Continue" button.</p> <p>F. In the "Test" column, type ";UHCT."</p> <p>G. The prompt, "Confirm adding test UHCT" will appear. Click on the "Yes" button.</p> <p>H. In the result column, type a semicolon followed by the unit's hematocrit. (Ex = ;0.55) then press the "tab" key.</p> <p>I. Click the "Save" button.</p> 

Step	Action
19	Allocate and crossmatch the unit per procedure.
20	Notify the patient care area when the product is available for issue. Instruct the patient care area that the unit has a short expiration and should be transfused as soon as possible.
21	Store the product in the refrigerator (1-6°C) until issue or expiration.
22	If the volume of the reconstituted product is: A. GREATER than the volume requested. Issue the entire product. The patient care area will pull of the desired volume. B. LESS THAN the volume requested. Notify the patient care area that you will prepare an additional reconstituted unit if needed.

6. RELATED DOCUMENTS

- Form: Reconstituted Whole Blood Worksheet (AG.F288)
- SOP: Procurement of Blood Products and Desired Inventory Levels
- SOP: Scale Quality Control
- SOP: Plasma for Transfusion
- SOP: Blood Component Irradiation
- SOP: Plasma Aliquot Preparation
- SOP: Blood Label Check
- SOP: Crossmatch

7. REFERENCES

1. Roback, J.D., Grossman, B.J., Harris, T., and Hillyer, C.D. 2011. Technical Manual of the AABB, 17th ed. AABB Publishing, Bethesda, Maryland.
2. Standards for Blood Banks and Transfusion Services, 2014. AABB, 29th ed. AABB Publishing, Bethesda, Maryland.
3. Circular of Information for the Use of Human Blood and Blood Components Prepared by: AABB, the American Red Cross, America's Blood Center and the Armed Services Blood Program. 2009. Bethesda, MD.
4. Code of Federal regulations, 21 CFR, Parts 200 and 600. Washington DC: US Government Printing Office, Current edition.

Document: SGAH.BB877[0] Status: INWORKS, Effective: 6/29/2014, Check Version Before Use

8. REVISION HISTORY

Version	Date	Reason for Revision	Revised By	Approved By

9. ADDENDA AND APPENDICES
N/A

Form revised 3/31/00