

## TRAINING UPDATE

**Lab Location:** SGAH and WAH      **Date Implemented:** 6.3.2014  
**Department:** Blood Bank      **Due Date:** 6.8.2014

### DESCRIPTION OF PROCEDURE REVISION

#### **Name of procedure:**

Plasma Aliquot Preparation

#### **Description of change(s):**

1. When plasma is aliquotted, the ORIGINAL PARENT product will become part "AO" and the first aliquot will be "BO."  
Subsequent aliquots will be Aa, Ab, Ac, etc.
2. Following aliquot, the labels for BOTH the parent and aliquot will print automatically.
3. Ensure the new labels are adhered to BOTH the parent and aliquot BEFORE separating them from each other.
4. You will have to handwrite the volume of anticoagulant on the new plasma label for apheresis plasma.
5. Label check MUST be performed for BOTH units.
6. The aliquots will have an "A" product code instead of an "E" product code. This is so aliquots will automatically bill in Sunquest (if we use the E code, we will have to manually bill).

Non-Technical SOP

<b>Title</b>	<b>Plasma Aliquot Preparation</b>	
<b>Prepared by</b>	Stephanie Codina	Date: 12.21.2012
<b>Owner</b>	Stephanie Codina	Date: 12.21.2012

**Laboratory Approval**

Print Name and Title	Signature	Date
<i>Refer to the electronic signature page for approval and approval dates.</i>		

Local Issue Date: \_\_\_\_\_ Local Effective Date: \_\_\_\_\_

**Review:**

Print Name	Signature	Date

Document: SGAH.BB146[2] Status: INWORKS, Effective: 6/29/2014, Check Version Before Use

Form revised 3/31/00

**TABLE OF CONTENTS**

1. PURPOSE ..... 2  
 2. SCOPE ..... 2  
 3. RESPONSIBILITY ..... 2  
 4. DEFINITIONS ..... 2  
 5. PROCEDURE ..... 2  
 6. RELATED DOCUMENTS ..... 6  
 7. REFERENCES ..... 7  
 8. REVISION HISTORY ..... 7  
 9. ADDENDA AND APPENDICES ..... 7

**1. PURPOSE**

To describe the procedure for making small-volume plasma aliquots from plasma units. This procedure allows small amounts of a plasma unit to be transfused over the lifespan of the original blood product. This process limits donor exposures and decreases donor-related risks in the recipient while minimizing overall blood product wastage.

**2. SCOPE**

This procedure applies to small-volume plasma transfusions requested for a neonate or small child. SGAH staff members perform all plasma aliquot procedures for WAH and SGAH.

**3. RESPONSIBILITY**

All blood bank staff members must be trained and competent in plasma aliquot preparation to ensure the purity, potency, and safety of the aliquotted plasma product.

**4. DEFINITIONS**

N/A

**5. PROCEDURE**

Step	Action
1	The patient care area will order plasma aliquots using test "TPLANE." Review the order, special instructions, and volume.
2	Obtain the supplies necessary to aliquot the plasma. A. 150 mL transfer bag B. Sterile welding device C. Heat sealer D. Scale E. Hemostats

Document: SGAH.BB146[2] Status: INWORKS, Effective: 6/29/2014, Check Version Before Use

Form revised 3/31/00

Document: SGAH.BB146[2] Status: INWORKS, Effective: 6/29/2014, Check Version Before Use

Step	Action
3	Perform daily QC of the scale if needed.
4	<p>Select plasma that meets the recipient's transfusion specifications.</p> <p>A. Only group AB plasma is transfused to neonates.</p> <p>a. If group AB plasma is not available, obtain pathologist approval to transfuse group specific plasma to a neonate.</p> <p>b. We must have 2 independent blood types on file (ABO retype) prior to transfusing group-specific plasma.</p> <p>B. Pediatric patients may receive group AB or group specific plasma products.</p> <p>All plasma transfused to neonates will be transfused within 24 hours of the time thawed. <b>If a previously-thawed unit is used, ensure that the transfusion will take place before the 24-hour limit.</b></p>
5	If the plasma is frozen, thaw per procedure "Plasma for Transfusion." Thaw the plasma physically and BCP the thaw in the computer.
6	<p>Tighten all connections. The hub connection nearest the syringe has disconnected on rare occasions.</p> <p><b>Use aseptic technique for this procedure!</b></p>
7	<p>Document the following on the "Product Modification Log."</p> <p>A. Tech identification</p> <p>B. Date of modification</p> <p>C. Unit number</p> <p>D. E code of original and new units (or A code if applicable)</p> <p>E. Lot number of bag</p> <p>F. Wafer lot number</p>
8	Gently mix the primary bag to ensure the plasma is completely thawed and no frozen bits remain.
9	Connect the transfer bag to the primary plasma per procedure, "Sterile Tubing Welder."
10	Tare the scale using an empty 150 mL transfer bag.
11	Allow the required amount of plasma to flow into the transfer bag via gravity. Include an extra 5 mL of plasma to compensate for the volume that will be lost in the tubing.


Form revised 3/3/00

Document: SGAH.BB146[2] Status: INWORKS, Effective: 6/29/2014, Check Version Before Use

Step	Action
12	<p>Clamp the line when an appropriate volume of plasma has been transferred. Seal the line between the hemostat clamps using a tube sealer at least twice. <b>Do not separate the aliquot from the parent unit at this time.</b></p> <ul style="list-style-type: none"> <li>A. Always make sure the hemostat is clamped between the parent unit and the location in which the tubing will be sealed.</li> <li>B. This will protect the sterility of the unit should the heat seal fail.</li> </ul>
13	<p>Access Sunquest function "Blood Component Preparation."</p> <p>Note: Do not branch to blood component preparation from blood order processing.</p>
14	<p>At the "Value" prompt, type the aliquot function that corresponds to the plasma unit to be aliquotted then press the "Tab" key. The aliquot function is A + the E code of the platelet product. Refer to appendix A for additional information.</p>
15	<p>Press the tab key to default the current date and time as the aliquot time. Enter the date and time on which the aliquot was prepared if prepared at an earlier time (such as during computer downtime).</p>
16	<p>Click the "continue" button.</p>
17	<p>A second "Blood Component Prep" screen will appear.</p> <ul style="list-style-type: none"> <li>A. At the "Unit #" prompt, scan the unit number DIN of the parent plasma to be aliquotted.</li> <li>B. At the "Component" prompt, scan the product code of the parent plasma to be aliquotted. This will autofill both the product code and division fields.</li> </ul>
18	<p>A pop-up screen will appear asking the user to indicate the number of units that will be prepared (ie the number of aliquots being prepared at one time).</p> <ul style="list-style-type: none"> <li>A. Enter 1 in the field.</li> <li>B. Click the "OK" button.</li> </ul> <div data-bbox="703 1478 1195 1873" data-label="Image"> <p>The screenshot shows a software interface for blood component preparation. A main window titled 'Blood Component Prep' is visible, with various fields and buttons. A smaller pop-up dialog box is overlaid on top, titled 'New Entry for Blood Component Prep 1'. This dialog box has a field labeled 'Number of units' with the value '1' entered, and 'OK' and 'Cancel' buttons at the bottom.</p> </div>

Form revised 3/3/00

Document: SGAH\_BB146[Z] Status: INWORKS, Effective: 6/29/2014, Check Version Before Use

Step	Action
19	<p>On the next screen, click on the yellow circle containing the “N” (for new product). Enter the volume of the plasma aliquot being prepared, then press the “Tab” key.</p> 
20	<p>Verify the new expiration dates/times. Document the new expiration date and time on the log.</p> <ul style="list-style-type: none"> <li>A. The expiration date of the parent unit will not change if a closed system is used.</li> <li>B. The expiration date of the parent unit will change to 24 hours from the time of aliquot if an open system is used (ie the sterile connection failed). If an open system is used, the output blood product codes in appendix A do not apply. See a supervisor for guidance.</li> <li>C. The expiration date of the aliquot will always be 24 hours from the time of thaw (not aliquot) regardless of whether an open or closed system is used.</li> </ul>
21	Click the “save” button.
22	<p>A “Preview Output / New Units” screen will appear. Review the information to ensure accuracy, then click on the “finish” button to generate new labels for the parent and aliquotted products.</p> <ul style="list-style-type: none"> <li>A. The first time an aliquot is prepared, the system will convert the parent unit to division “AO” and the aliquot to division “BO.”</li> <li>B. All subsequent divisions will assign a division code to the aliquot using the division labeling convention of “Aa, Ab, Ac, Ad, .....Az.”</li> </ul>
23	Adhere the new labels to BOTH the parent unit and the aliquot. <b>Ensure you adhere the “AO” label to the parent unit.</b>
24	After labeling, disconnect the aliquot from the parent unit.

Form revised 3/3/100

Document: SGAH.BB146[2] Status: INWORKS, Effective: 6/29/2014, Check Version Before Use

Step	Action
25	<p>When the aliquot was made from an <b>apheresis</b> plasma product, calculate the amount of anticoagulant in the aliquoted and parent plasma products using the following formulas. Then, document the volume in the designated space on each product label.</p> $\text{Ratio} = \frac{\text{Amount of Anticoagulant in Plasma}}{\text{Volume of Plasma}}$ <p>Volume of anticoagulant in the aliquot = ratio x aliquot volume</p> <p>Volume of anticoagulant in the parent product = ratio x remaining volume</p> <p>For example,</p> <p>If the plasma label indicated an anticoagulant volume of 5 mL and a total volume of 200 mL. Therefore, the ratio would be:</p> $\text{Ratio} = 5 \text{ mL} \div 200 \text{ mL} = 0.025$ <p>Assume the new aliquot is 15 mL. The amount of anticoagulant in the aliquot is:</p> $0.025 \times 15 \text{ mL} = 0.4 \text{ mL of anticoagulant in the aliquot.}$
26	<p>Document the following on the "Product Modification Log."</p> <ul style="list-style-type: none"> <li>A. Division of new product</li> <li>B. Documentation of the weld inspection</li> </ul>
27	<p>Perform a label check in Sunquest per procedure. Label check must be performed on BOTH the parent unit and the aliquot.</p>
28	<p>Allocate the plasma aliquot per procedure.</p>
29	<p>Store both the parent unit and aliquot in the refrigerator (1-6°C) until issue. Aliquots should be transfused as soon as possible following preparation.</p>

- 6. RELATED DOCUMENTS**  
 SOP: Plasma for Transfusion  
 SOP: Sterile Tubing Welder  
 SOP: Scale Quality Control  
 SOP: Blood Label Check  
 Form: Product Modification Log (AG.F01)

Form revised 3.31/00

7. **REFERENCES**

1. Roback, J.D., Grossman, B.J., Harris, T., and Hillyer, C.D. 2011. Technical Manual of the AABB, 17<sup>th</sup> ed., AABB Publishing, Bethesda, Maryland.
2. Standards for Blood Banks and Transfusion Services, 2014. AABB, 19<sup>th</sup> ed., AABB Publishing, Bethesda, Maryland.

8. **REVISION HISTORY**

Version	Date	Reason for Revision	Revised By	Approved By
		Supersedes SGAH.B408.01		
000	5.9.13	Section 5: Added ISBT-128 labeling information Section 9: Added appendix B	SCodina	NCacciabeve
001	5.27.14	Section 5: Removed references to codabar-labeled units. Updated LIS instructions to include the Sunquest v6.4 upgrade. Moved appendix for blood component preparation in the LIS to the procedure. Section 9: Re-numbered appendix Footer: version # leading zero's dropped due to new EDCS in use as of 10/7/13.	SCodina	NCacciabeve

9. **ADDENDA AND APPENDICES**

- A. Plasma Aliquot Blood Component Prep Functions



**Appendix A**  
**Plasma Aliquot Blood Component Prep Functions**

**Plasma Aliquot Products**

<b>Original Product</b>	<b>Component Prep Function</b>	<b>Final Product</b>
E2121	AE2121	A2121
E2684	AE2684	A2684
E2702	AE2702	A2702
E2720	AE2720	A2720
E5548	AE5548	A5548
E5549	AE5549	A5549
E5550	AE5550	A5550

Document: SGAH.BB146[2] Status: INWORKS, Effective: 6/29/2014, Check Version Before Use

Form revised 3/3/00