

TRAINING UPDATE

Lab Location: GEC
Department: Core

Date Distributed: 7/30/2014
Due Date: 8/20/2014
Implementation: 8/20/2014

DESCRIPTION OF PROCEDURE REVISION

Name of Documents:

Enzymatic Carbonate (ECO₂) by Dimension® Xpand Chemistry Analyzer GEC.C24v1

Ethyl Alcohol by Dimension® Xpand Chemistry Analyzer GEC.C41v1

Description of change(s):

CO₂ and Alcohol SOPs:

Section	Reason
5, 17	Revised to reflect new CHEM III calibrator

CO₂ only

Section	Reason
3.2	Add timing for uncapped specimens exposed to air
6.7	Add use of TEA for lot to lot runs
10.2	Correct rounding to whole number
10.5	Remove code QNSR
11.3	Remove priority 3 reporting for SGAH

These revised SOPs will be implemented on August 20, 2014.

Document your compliance with this training update by taking the quiz in the MTS system.

Technical SOP

Title	Enzymatic Carbonate (ECO₂) by Dimension® Xpand Chemistry Analyzer	
Prepared by	Ashkan Chini	Date: 4/12/2011
Owner	Robert SanLuis	Date: 4/12/2011

Laboratory Approval	Local Effective Date:	
Print Name and Title	Signature	Date
<i>Refer to the electronic signature page for approval and approval dates.</i>		

Review		
Print Name	Signature	Date

TABLE OF CONTENTS

1. Test Information.....2
 2. Analytical Principle3
 3. Specimen Requirements.....3
 4. Reagents4
 5. Calibrators/Standards5
 6. Quality Control6
 7. Equipment And Supplies9
 8. Procedure9
 9. Calculations.....10
 10. Reporting Results And Repeat Criteria.....10
 11. Expected Values.....11
 12. Clinical Significance12
 13. Procedure Notes12
 14. Limitations Of Method12
 15. Safety13
 16. Related Documents14
 17. References.....14
 18. Revision History15
 19. Addenda15

1. TEST INFORMATION

Assay	Method/Instrument	Local Code
Enzymatic Carbonate	Dimension® Xpand Chemistry Analyzer	CO ₂


Synonyms/Abbreviations
CO ₂ , ECO ₂

Department
Chemistry

Form revised 2/02/2007

2. ANALYTICAL PRINCIPLE

The enzymatic carbonate (ECO₂) method for the Dimension® system employs a phosphoenolpyruvate carboxylase-malate dehydrogenase coupled enzymatic reaction and a stable analog of the cofactor NADH. The bicarbonate anion reacts with phosphoenolpyruvate in the presence of phosphoenolpyruvate carboxylase (PEPC) and Mg⁺⁺ to form oxaloacetate and inorganic phosphate (P_i). The oxaloacetate is reduced to malate-by-malate dehydrogenase (MDH) with simultaneous oxidation of the reduced form of an analog (aNADH) of the cofactor, NADH.



3. SPECIMEN REQUIREMENTS

3.1 Patient Preparation

Component	Special Notations
Fasting/Special Diets	N/A
Specimen Collection and/or Timing	Normal procedures for collecting and storing serum and plasma may be used for samples to be analyzed by this method.
Special Collection Procedures	None
Other	N/A

3.2 Specimen Type & Handling

Criteria	
Type -Preferred -Other Acceptable	Plasma (Heparin) Serum
Collection Container	Plasma: Green top tube serum separator Serum: Red top tube, Serum separator tube (SST)
Volume - Optimum - Minimum	1.0 mL 0.5 mL
Transport Container and Temperature	Collection container or Plastic vial at room temperature
Stability & Storage Requirements	Room Temperature: 8 hours
	Refrigerated: (2-8°C) 2 days

Form revised 2/02/2007

Criteria	
	Frozen: (-20°C or colder) 6 months
Timing Considerations	Total CO ₂ concentration may be lowered by as much as 6 mmol/L when uncapped specimens are exposed to air for one hour.
Unacceptable Specimens & Actions to Take	Specimens that are unlabeled, improperly labeled, or those that do not meet the stated criteria are unacceptable. Request a recollection and credit the test with the appropriate LIS English text code for “test not performed” message. Examples: Quantity not sufficient-QNS; Wrong collection-UNAC. Document the request for recollection in the LIS.
Compromising Physical Characteristics	Gross hemolysis. Reject sample and request a recollection. Credit the test with the appropriate LIS English text code explanation of HMT (Specimen markedly hemolyzed)
Other Considerations	Allow to clot completely prior to centrifugation.

4. REAGENTS

Refer to the Material Safety Data Sheet (MSDS) supplied with the reagents for complete safety hazards. Refer to the section in this procedure covering “SAFETY” for additional information.

4.1 Reagent Summary

Reagents	Supplier & Catalog Number
ECO ₂	Siemens, Flex® reagent cartridge, Cat. No. DF137

4.2 Reagent Preparation and Storage

NOTES: Date and initial all reagents upon opening. Each container must be labeled with (1) substance name, (2) lot number, (3) date of preparation, (4) expiration date, (5) initials of tech, (6) any special storage instructions; check for visible signs of degradation.

Refer to the Material Safety Data Sheet (MSDS) for a complete description of hazards. If a specific hazard is present, it will be noted in this procedure when the hazard is first encountered in a procedural step.

Irritant. Contains sodium azide as a preservative. Sodium Azide can react with copper or lead pipes in drain lines to form explosive compounds. May cause sensitization by skin contact.

Reagent	Enzymatic Carbonate
Container	Reagent cartridge
Storage	Store at 2-8°C

Stability	<ul style="list-style-type: none"> • Reagent is stable until expiration date stamped on the reagent cartridges. • Sealed or unhydrated cartridge wells on the instrument are stable for 30 days. • Once wells 1 – 6 have been entered by the instrument, they are stable for 2 days.
Preparation	Reagents are supplied ready for use. No additional preparation is required.

5. CALIBRATORS/STANDARDS

5.1 Calibrators/Standards Used

Calibrator	Supplier and Catalog Number
CHEM III Calibrator	Siemens Dimension®, Cat. No. DC130

5.2 Calibrator Preparation and Storage

NOTE: Date and initial all calibrators upon opening. Each container must be labeled with (1) substance name, (2) lot number, (3) date of preparation, (4) expiration date, (5) initials of tech (6) any special storage instructions; check for visible signs of degradation.

Calibrator	CHEM III Calibrator
Preparation	Calibrator is ready for use. No preparation is required.
Storage/Stability	<ul style="list-style-type: none"> • Store at 2 – 8 °C • Unopened calibrators are stable until the expiration date printed on the label. • Once cap is removed, assigned values are stable for 30 days when recapped immediately after use and stored at 2 – 8 °C

5.3 Calibration Parameter

Criteria	Special Notation
Reference Material	CHEM III Calibrator
Assay Range	5 – 45 mmol/L
Suggested calibration level	See Reagent Package Insert for lot specific assigned values in mmol/L
Frequency	<ul style="list-style-type: none"> • Every new reagent cartridge lot. • Every 90 days for any one lot. • When major maintenance is performed on the analyzer. • When control data indicates a significant shift in assay.
Calibration Scheme	Three levels in triplicate.

Assigned Coefficients	C ₀ 0.000
	C ₁ 1.000

5.4 Calibration Procedure

1. From Operating Menu press F5:Process Control press F1: Calibration Enter Password press F2: SETUP and RUN
2. Select the test method to be calibrated - if lot number is incorrect Press F1: Other Lot
3. Enter all information on screen
4. Press F8: QC yes/no to change to yes
5. Press F4: Assign cups If additional methods need to be calibrated, select the method.
6. Press F7: Load/run
7. Load cups into assigned position
8. Press F4: RUN

5.5 Tolerance Limits

IF.....	THEN.....
If result fall within assay-specific specification, and QC values are within acceptable limits,	proceed with analysis
If result falls outside assay-specific specification, or QC values are out of Acceptable limits,	troubleshoot the assay and/or instrument and repeat calibration

6. QUALITY CONTROL

6.1 Controls Used

Controls	Supplier and Catalog Number
Liquichek Unassayed Chemistry Controls Levels 1 & 2	Bio-Rad Laboratories Cat. No. 691 & 692

6.2 Control Preparation and Storage

NOTE: Date and initial all controls upon opening. Each container should be labeled with (1) substance name, (2) lot number, (3) date of preparation, (4)

expiration date, (5) initials of tech, and (6) any special storage instructions; check for visible signs of degradation.

Control	Liquichek Unassayed Chemistry Controls Levels 1 & 2
Preparation	Allow the frozen control to stand at room temperature (18-25°C) until completely thawed. Swirl the contents gently to ensure homogeneity. (Do not use a mechanical mixer) Use immediately. After each use, promptly replace the stopper and return to 2 to 8°C storage.
Storage/Stability	Open controls are stable for 15 days at 2-8°C. Unopened controls are stable until the expiration date at -20 to -70°C.

6.3 Frequency

Analyze all levels of QC material after every calibration and each day of testing.

Refer to the **Dimension Xpand® QC** Schedule in the Laboratory policy Quality Control Program and in the **Dimension Xpand® Quick Reference** Guide.

6.4 Tolerance Limits

Step	Action
1	Acceptable ranges for QC are programmed into the Laboratory Information System (LIS), and may be posted near the instrument for use during computer downtime.
2	Run Rejection Criteria <ul style="list-style-type: none"> Anytime the established parameters are exceeded (if one QC result exceeds 2 SD), the run is considered out of control (failed) and patient results must not be reported. The technologist must follow the procedure in the Laboratory QC Program to resolve the problem.
3	Corrective Action: <ul style="list-style-type: none"> All rejected runs must be effectively addressed through corrective action. Steps taken in response to QC failures must be documented. Patient samples in failed analytical runs must be <u>reanalyzed according to the Laboratory QC Program</u>. Supervisors may override rejection of partial or complete runs only with detailed documentation and criteria for overrides that are approved by the Medical Director. Consult corrective action guidelines in Laboratory QC Program. Follow corrective action guidelines in the Laboratory QC Program. Corrective action documentation must follow the Laboratory Quality Control Program.

Form revised 2/02/2007

Step	Action
4	Review of QC <ul style="list-style-type: none">• QC must be reviewed weekly by the Group Lead or designee and monthly by the Supervisor/Manager or designee.• If the SD and/or CV are greater than established ranges, investigate the cause for the imprecision and document implementation of corrective actions.

6.5 Review Patient Data

Technologist must review each result with error messages. Refer to the Dimension Xpand[®] system manual “Error messages” section for troubleshooting. Check for unusual patterns, trends, or distributions in patient results (such as an unusually high percentage of abnormal results). Resolve any problems noted before issuing patient reports.

6.6 Documentation

- QC tolerance limits are programmed into the instrument and the LIS. The LIS calculates cumulative mean, SD and CV and stores all information for easy retrieval.
- Quality control records are reviewed daily at the bench, weekly by the Group Lead or designee, and monthly by the Supervisor/Manager or designee.
- Refer to complete policies and procedures for QC documentation and for record retention requirements in the Laboratory QC Program.

6.7 Quality Assurance Program

- Each new lot number of reagent or new shipment of the same lot of reagent must be tested with external control materials and previously analyzed samples. Performance of the new lot must be equivalent to the previous lot; **utilize published TEA for acceptability criteria.**
- Training must be successfully completed and documented prior to performing this test. This procedure must be incorporated into the departmental competency assessment program.
- The laboratory participates in CAP proficiency testing. All proficiency testing materials must be treated in the same manner as patient samples.
- Monthly QC must be presented to the Medical Director **or designee** for review and signature.
- Monthly QC mean and SD are sent to Bio-Rad Laboratories for peer group comparison.
- Consult the Laboratory QC Program for complete details.

7. EQUIPMENT and SUPPLIES

7.1 Assay Platform

Dimension Xpand® System

7.2 Equipment

- Refrigerator capable of sustaining 2–8°C.
- Freezer capable of sustaining range not to exceed -20 to -70°C.
- Centrifuge

7.3 Supplies

- Plastic serum tubes and serum cups
- **Reagent Grade** water (Millipore® or equivalent)
- Calibrated pipettes and disposable tips

8. PROCEDURE

ECO₂ Flex® reagent cartridge Cat. No. DF137 is required to perform this test.

Enzymatic Carbonate is performed on the **Dimension Xpand®** clinical chemistry system after the method is calibrated (see Reference Material in Calibration section) and Quality Controls are acceptable.

NOTE: For all procedures involving specimens, buttoned lab coats, gloves, and face protection are required minimum personal protective equipment. Report all accidents to your supervisor.

The package insert for a new lot of kits must be reviewed for any changes before the kit is used. A current Package Insert is included as a Related Document.

8.1	Instrument Set-Up Protocol
1.	For instrument set up and operation: Refer to Startup and Maintenance, Siemens Dimension Xpand® procedure.
2.	Check reagent inventory
3.	Sampling, reagent delivery, mixing and processing of results are automatically performed by the Dimension Xpand® system. For details of the automated parameters, see below under “Test conditions.”

8.2	Specimen/Reagent Preparation
1.	Centrifuge the specimens.
2.	Specimens are placed in Dimension® Xpand segments for analysis by the instrument. Refer to the Sample Processing, Siemens Dimension® Xpand procedure. The sample container (if not a primary tube) must contain sufficient quantity to accommodate the

8.2	Specimen/Reagent Preparation
	sample volume plus 50 µL of dead volume. Precise container filling is not required.
8.3	Specimen Testing
1.	For QC placement and frequency, refer to the Dimension [®] Xpand QC Schedule in the Laboratory QC Program.
2.	Follow the instructions, outlined in the Dimension [®] Xpand Operators Manual
3.	The instrument reporting system contains error messages to warn the user of specific malfunctions. Results followed by such error messages should be held for follow-up. Refer to the Dimension Xpand [®] system manual “Error messages” section for troubleshooting.
4.	Follow protocol in Section 10.5 “Repeat criteria and resulting” for samples with results above or below the Analytical Measurement Range (AMR). Investigate any failed delta result and repeat, if necessary.
5.	Append the appropriate English text code qualifier messages to any samples requiring a comment regarding sample quality and/or any other pertinent factors.

Test Conditions	
Sample Size:	5 µL
Reagent 1 Volume:	100 µL
Temperature:	37° C
Wavelength:	405 and 700 nm
Type of Measurement:	Bichromatic rate

9. CALCULATIONS

The instrument automatically calculates and prints the concentration of Enzymatic Carbonate in mmol/L.

10. REPORTING RESULTS AND REPEAT CRITERIA

10.1 Interpretation of Data

None required

10.2 Rounding

No rounding is necessary. **Instrument reports results as a whole number.**

10.3 Units of Measure

mmol/L

10.4 Clinically Reportable Range (CRR)

5 - 90 mmol/L

10.5 Repeat Criteria and Resulting

All repeats must replicate the original result within the total allowable error (TEa) of the assay. Refer to TEa policy for specific information.

Values that fall within the AMR or CRR may be reported without repeat. Values that fall outside these ranges must be repeated

IF the result is ...	THEN...
≤5 mmol/L	Assure there is sufficient sample devoid of bubbles, cellular debris, and/or fibrin clots. Report as: <5 mmol/L
>45 mmol/L	Manual Dilution: Using the primary tube, make the smallest dilution possible to bring the raw data within the AMR. Maximum allowable dilution: x 2 Diluent: Reagent Grade Water Enter dilution factor as a whole number on the “Enter Sample Data” screen.
>90 mmol/L	If the recommended dilution does not give results within the clinically reportable range, report as: “>90 mmol/L-REP” Bring to the attention of your supervisor prior to releasing result.

Message	Code
Verified by repeat analysis	Append -REP to the result.

11. EXPECTED VALUES

11.1 Reference Ranges

Age	Female / Male
Adult (>18 years):	21 – 32 mmol/L
Pediatric:	
0 – 6 days	13-21
7 – 30 days	13-22
1 – 5 months	13-23
6 – 12 months	14-23
13 - 23 months	16-25
2 – 18 years	21-32

11.2 Critical Values

< 10 mmol/L

11.3 Priority 3 Limit(s)

None established

12. CLINICAL SIGNIFICANCE

Increase in serum CO₂ content for the most part reflects increase in serum bicarbonate concentration rather than dissolved CO₂ gas (which accounts for only a small fraction of the total). Increased serum bicarbonate is seen in compensated respiratory acidosis and in metabolic alkalosis. Diuretics (thiazides, ethacrynic acid, furosemide, mercurials), corticosteroids (in long term use), and laxatives (when abused) may cause increased bicarbonate.

Decrease in blood CO₂ is seen in metabolic acidosis and compensated respiratory alkalosis. Substances causing metabolic acidosis include ammonium chloride, acetazolamide, ethylene glycol, methanol, paraldehyde, and phenformin. Salicylate poisoning is characterized by early respiratory alkalosis followed by metabolic acidosis with attendant decreased bicarbonate.

13. PROCEDURE NOTES

- **FDA Status:** FDA Approved/cleared
- **Validated Test Modifications:** None

The instrument reporting system contains error messages to warn the operator of specific malfunctions. Any report slip containing such error messages should be held for follow-up. Refer to your Dimension Xpand Operator's Guide.

A system malfunction may exist if the following 5-test precision is observed:

Concentration	S.D.
25 mmol/L	> 1.2 mmol/L
50 mmol/L	> 2.0 mmol/L

14. LIMITATIONS OF METHOD

14.1 Analytical Measurement Range (AMR)

5 – 45 mmol/L

14.2 Precision

Material	Mean mmol/L	Standard Deviation (%CV)	
		Within-run	Total
Dade® Moni-trol® TOTAL Control			
Level 1	13.1	0.5	0.7
Level 2	30.5	0.6	0.9
Plasma Pool	24.5	0.7	1.1
Serum Pool	24.4	0.7	1.1

14.3 Interfering Substances

In rooms with poor ventilation, an open Flex® reagent cartridge well can absorb CO₂ which may cause results to be elevated by up to 30%.

Hemoglobin (hemolysate) of 1000 mg/dL decreases an ECO₂ result of 13 mmol/L by 21%.

Lipemia (Intralipid®) of 3000 mg/dL decreases an ECO₂ result of 13mmol/L by 16%.

HIL Interference:

The ECO₂ method was evaluated for interference from hemolysis, icterus and lipemia according to CLSI/NCCLS EP7-P. Bias, defined as the difference between the control sample (does not contain interferent) and the test sample (contains interferent), is shown in the table below. Bias exceeding 10% is considered “interference”.

Substance tested	Test Concentration SI Units	ECO ₂ Concentration mmol/L	Bias %
Hemoglobin (hemolysate)	500 mg/dL (monomer)	13	<10
Bilirubin (unconjugated)	80 mg/dL	14	<10
Lipemia (Intralipid®)	1000 mg/dL	13	<10

14.4 Clinical Sensitivity/Specificity/Predictive Values

Not available.

15. SAFETY

The employee has direct responsibility to avoid injury and illness at work. Nearly all harmful exposures to infectious substances and chemicals, and other injuries, can be avoided with effective training and consistent safe work practices.

Become familiar with the Environmental Health and Safety (EHS) Manual to learn the requirements on working safely and protecting the environment from harm. Although lab

work typically focuses on the hazards of working with specimens and chemicals, we must also control other important hazards.

- Slips, trips, and falls cause many serious injuries. Please ensure that spills are cleaned quickly (to avoid slippery floors) and that you can see and avoid obstacles in your path.
- Ergonomic injuries result from performing tasks with too much repetition, force, or awkward position. Ergonomic injuries include strains and back injuries. Learn about ergonomic hazards and how to prevent this type of injury.
- Scratches, lacerations, and needlesticks can result in serious health consequences. Attempt to find ways to eliminate your risk when working with sharp materials.

Report all accidents and injuries immediately to your supervisor or the business unit Environmental Health and Safety Manager or Specialist.

16. RELATED DOCUMENTS

1. Dimension Xpand[®] Clinical Chemistry System Operator's Manual
2. Calibration / Verification Siemens Dimension[®] Xpand procedure
3. Dimension Xpand[®] Cal Accept Guidelines
4. Dimension Xpand[®] Calibration summary
5. Sample Processing, Siemens Dimension[®] Xpand procedure
6. Start up and Maintenance, Siemens Dimension[®] Xpand procedure
7. Laboratory Quality Control Program
8. QC Schedule for Siemens Dimension Xpand[®]
9. Laboratory Safety Manual
10. Material Safety Data Sheets (MSDS)
11. Siemens Dimension Xpand[®] Limits Chart (AG.F143)
12. Quest Diagnostics Records Management Procedure
13. Dimension Xpand[®] System Error Messages Chart
14. Centrifuge Use, Maintenance and Functions Checks (Lab policy)
15. Hemolysis, Icteria and Lipemia Interference (Lab policy)
16. Repeat Testing Requirements (Lab policy)
17. Critical Values (Lab policy)
18. Current Allowable Total Error Specifications at
http://questnet1.qdx.com/Business_Groups/Medical/qc/docs/qc_bpt_tea.xls
19. Current package insert ECO₂ Flex[®] Reagent Cartridge DF137

17. REFERENCES

1. Ghoshal, Amit K. and Soldin, Steven J., Evaluation of the Dade Behring Dimension[®] RxL: Integrated chemistry system-pediatric reference ranges. Clinica Chimica Acta 2003; 331:144
2. Package Insert, ECO₂ Flex[®] Reagent Cartridge DF137, Siemens Healthcare Diagnostics Inc., 06/05/2013.
3. Package insert, CHEM III Calibrator DC130, Siemens Healthcare Diagnostics Inc., 04/2013.

4. Package insert, Liquichek Unassayed Serum Chemistry Controls, Bio-Rad Laboratories, 01/2013.

18. REVISION HISTORY

Version	Date	Section	Reason	Reviser	Approval
			Supersedes SOP C062.000		
000	7/15/14	1, 7.1	Add analyzer name	L Barrett	R SanLuis
000	7/15/14	3.2	Add timing for uncapped specimens exposed to air	A Chini	R SanLuis
000	7/15/14	5, 17	Revised to reflect new CHEM III calibrator	A Chini	R SanLuis
000	7/15/14	6.7	Add use of TEA for lot to lot runs	L Barrett	R SanLuis
000	7/15/14	10.2	Correct rounding to whole number	A Chini	R SanLuis
000	7/15/14	10.5	Remove code QNSR	L Barrett	R SanLuis
000	7/15/14	11.3	Remove priority 3 reporting for SGAH	L Barrett	R SanLuis
000	7/15/14	15	Update to standard wording	L Barrett	R SanLuis
000	7/15/14	16	Update document titles	L Barrett	R SanLuis
000	7/15/14	Footer	Version # leading zero's dropped due to new EDCS in use as of 10/7/13	L Barrett	R SanLuis

19. ADDENDA

None

Technical SOP

Title	Ethyl Alcohol by Dimension® Xpand Chemistry Analyzer	
Prepared by	Ashkan Chini	Date: 8/14/2013
Owner	Robert SanLuis	Date: 8/14/2013

Laboratory Approval		Local Effective Date:
Print Name and Title	Signature	Date
<i>Refer to the electronic signature page for approval and approval dates.</i>		

Review		
Print Name	Signature	Date

Form revised 2/02/2007

TABLE OF CONTENTS

1.	Test Information.....	2
2.	Analytical Principle.....	3
3.	Specimen Requirements.....	3
4.	Reagents.....	4
5.	Calibrators/Standards.....	5
6.	Quality Control.....	6
7.	Equipment And Supplies.....	8
8.	Procedure.....	9
9.	Calculations.....	10
10.	Reporting Results And Repeat Criteria.....	10
11.	Expected Values.....	11
12.	Clinical Significance.....	12
13.	Procedure Notes.....	12
14.	Limitations Of Method.....	12
15.	Safety.....	13
16.	Related Documents.....	14
17.	References.....	14
18.	Revision History.....	14
19.	Addenda.....	14

1. TEST INFORMATION

Assay	Method/Instrument	Local Code
Ethyl Alcohol	Dimension® Xpand Chemistry Analyzer	ALCO


Synonyms/Abbreviations
Ethanol, ETOH

Department
Chemistry

Form revised 2/02/2007

2. ANALYTICAL PRINCIPLE

The Ethyl Alcohol Assay is based on an enzymatic reaction. Reagent 1 contains the buffering system. Reagent 2 contains alcohol dehydrogenase (ADH), the coenzyme nicotinamide adenine dinucleotide (NAD), buffer, preservatives, and stabilizers. The ADH catalyzes the oxidation of ethyl alcohol to acetaldehyde. During this reaction, NAD is reduced to NADH with a concomitant increase in absorbance at 340 nm proportional to the concentration of alcohol in the specimen.



3. SPECIMEN REQUIREMENTS

3.1 Patient Preparation

Component	Special Notations
Fasting/Special Diets	N/A
Specimen Collection and/or Timing	Normal procedures for collecting and storing serum and plasma may be used for samples to be analyzed by this method.
Special Collection Procedures	Use non-alcohol germicidal solution to cleanse the skin.
Other	N/A

3.2 Specimen Type & Handling

Criteria	
Type -Preferred -Other Acceptable	Serum Plasma (Heparin)
Collection Container	Serum: Red top tube Plasma: Green top tube
Volume - Optimum - Minimum	1.0 mL 0.5 mL
Transport Container and Temperature	Collection container or Plastic vial at room temperature
Stability & Storage Requirements	Room Temperature: (18 – 28°C) 2 days
	Refrigerated: (2 – 8 °C) 2 weeks
	Frozen: (-20 °C or colder) stable indefinitely
Timing Considerations	Tubes that have been open for any great length of time are unacceptable. Open and process samples in STAT mode.

FORM 2012C REVISED 01/2012

Criteria	
Unacceptable Specimens & Actions to Take	Specimens that are unlabeled, improperly labeled, or those that do not meet the stated criteria are unacceptable. Request a recollection and credit the test with the appropriate LIS English text code for “test not performed” message. Examples: Quantity not sufficient-QNS; Wrong collection-UNAC. Document the request for recollection in the LIS.
Compromising Physical Characteristics	Gross hemolysis. Reject sample and request a recollection. Credit the test with the appropriate LIS English text code explanation of HMT (Specimen markedly hemolyzed)
Other Considerations	Allow Red Top tube to clot completely prior to centrifugation.

4. REAGENTS

Refer to the Material Safety Data Sheet (MSDS) supplied with the reagents for complete safety hazards. Refer to the section in this procedure covering “SAFETY” for additional information.

4.1 Reagent Summary

Reagents	Supplier & Catalog Number
Ethyl Alcohol	Siemens, Flex® reagent cartridge, Cat. No. DF22

4.2 Reagent Preparation and Storage

NOTES: Date and initial all reagents upon opening. Each container must be labeled with (1) substance name, (2) lot number, (3) date of preparation, (4) expiration date, (5) initials of tech, (6) any special storage instructions; check for visible signs of degradation.

Refer to the Material Safety Data Sheet (MSDS) for a complete description of hazards. If a specific hazard is present, it will be noted in this procedure when the hazard is first encountered in a procedural step.

Reagent	Ethyl Alcohol
Container	Reagent cartridge
Storage	Store at 2-8° C
Stability	<ul style="list-style-type: none"> Reagent is stable until expiration date stamped on the reagent cartridges. Sealed or unhydrated cartridge wells on the instrument are stable for 30 days. Open well stability: 5 days for all wells
Preparation	All reagents are liquid and ready to use.

FORM 2012C REVISED 01/2012

5. CALIBRATORS/STANDARDS

5.1 Calibrators/Standards Used

Calibrator	Supplier and Catalog Number
CHEM III Calibrator	Siemens Dimension®, Cat. No. DC130

5.2 Calibrator Preparation and Storage

NOTE: Date and initial all calibrators upon opening. Each container must be labeled with (1) substance name, (2) lot number, (3) date of preparation, (4) expiration date, (5) initials of tech (6) any special storage instructions; check for visible signs of degradation.

Calibrator	CHEM III Calibrator
Preparation	Calibrator is ready for use. No preparation is required.
Storage/Stability	<ul style="list-style-type: none"> Store at 2 – 8 °C Unopened calibrators are stable until the expiration date printed on the label. Once cap is removed, assigned values are stable for 30 days when recapped immediately after use and stored at 2 – 8 °C

5.3 Calibration Parameter

Criteria	Special Notations
Reference Material	CHEM III Calibrator
Assay Range	3 – 300 mg/dL
Suggested Calibration Level	See Reagent Package Insert for lot specific assigned values in mg/dL
Frequency	<ul style="list-style-type: none"> Every new reagent cartridge lot. Every 90 days for any one lot When major maintenance is performed on the analyzer. When control data indicates a significant shift in assay.
Calibration Scheme	Three levels in triplicate
Assigned Coefficients	C ₀ 0.10 C ₁ 3.3

5.4 Calibration Procedure

1. From Operating Menu press F5:Process Control press F1: Calibration

Form revised 12/02/2007

Enter Password press F2: SETUP and RUN
2. Select the test method to be calibrated - if lot number is incorrect Press F1: Other Lot
3. Enter all information on screen
4. Press F8: QC yes/no to change to yes
5. Press F4: Assign cups If additional methods need to be calibrated, select the method.
6. Press F7: Load/run
7. Load cups into assigned position
8. Press F4: RUN

5.5 Tolerance Limits

IF.....	THEN.....
If result fall within assay-specific specification, and QC values are within acceptable limits,	proceed with analysis
If result falls outside assay-specific specification, or QC values are out of Acceptable limits,	troubleshoot the assay and/or instrument and repeat calibration

6. QUALITY CONTROL

6.1 Controls Used

Controls	Supplier and Catalog Number
Liquichek™ Ethanol/Ammonia Control Levels 1, 2 and 3	Bio-Rad Laboratories Cat. No. 544, 545 and 546

6.2 Control Preparation and Storage

NOTE: Date and initial all controls upon opening. Each container should be labeled with (1) substance name, (2) lot number, (3) date of preparation, (4) expiration date, (5) initials of tech, and (6) any special storage instructions; check for visible signs of degradation.

Control	Liquichek Ethanol/Ammonia Controls, Levels 1, 2 and 3
Preparation	Before sampling, allow this product to reach room temperature (18 - 25°C). Gently swirl the vial several times to ensure homogeneity. After each use, promptly replace the stopper and return to 2 - 8°C storage.

Form revised 12/02/2007

Storage/Stability	Once the control is opened, all analytes will be stable for 20 days at 2-8°C. Unopened controls are stable until the expiration date at 2-8°C.
--------------------------	---

6.3 Frequency

Analyze all levels of QC material after every calibration and each day of testing.

Refer to the Dimension Xpand® QC Schedule in the Laboratory policy Quality Control Program and in the Dimension Xpand® Quick Reference Guide.

6.4 Tolerance Limits

Step	Action
1	Acceptable ranges for QC are programmed into the Laboratory Information System (LIS), and may be posted near the instrument for use during computer downtime.
2	Run Rejection Criteria <ul style="list-style-type: none"> Anytime the established parameters are exceeded (if one QC result exceeds 2 SD), the run is considered out of control (failed) and patient results must not be reported. The technologist must follow the procedure in the Laboratory QC Program to resolve the problem.
3	Corrective Action: <ul style="list-style-type: none"> All rejected runs must be effectively addressed through corrective action. Steps taken in response to QC failures must be documented. Patient samples in failed analytical runs must be <u>reanalyzed according to the Laboratory QC Program</u>. Supervisors may override rejection of partial or complete runs only with detailed documentation and criteria for overrides that are approved by the Medical Director. Consult corrective action guidelines in Laboratory QC Program. Follow corrective action guidelines in the Laboratory QC Program. Corrective action documentation must follow the Laboratory Quality Control Program.
4	Review of QC <ul style="list-style-type: none"> QC must be reviewed weekly by the Group Lead or designee and monthly by the Supervisor/Manager or designee. If the SD and/or CV are greater than established ranges, investigate the cause for the imprecision and document implementation of corrective actions.

From revised 12/02/2007

6.5 Review Patient Data

Technologist must review each result with error messages. Refer to the Dimension Xpand® system manual “Error messages” section for troubleshooting. Check for unusual patterns, trends, or distributions in patient results (such as an unusually high percentage of abnormal results). Resolve any problems noted before issuing patient reports.

6.6 Documentation

- QC tolerance limits are programmed into the instrument and the LIS. The LIS calculates cumulative mean, SD and CV and stores all information for easy retrieval.
- Quality control records are reviewed daily at the bench, weekly by the Group Lead or designee, and monthly by the Supervisor/Manager or designee.
- Refer to complete policies and procedures for QC documentation and for record retention requirements in the Laboratory QC Program.

6.7 Quality Assurance Program

- Each new lot number of reagent or new shipment of the same lot of reagent must be tested with external control materials and previously analyzed samples. Performance of the new lot must be equivalent to the previous lot; utilize published TEA for acceptability criteria.
- Training must be successfully completed and documented prior to performing this test. This procedure must be incorporated into the departmental competency assessment program.
- The laboratory participates in CAP proficiency testing. All proficiency testing materials must be treated in the same manner as patient samples.
- Monthly QC must be presented to the Medical Director or designee for review and signature.
- Monthly QC mean and SD are sent to Bio-Rad Laboratories for peer group comparison.
- Consult the Laboratory QC Program for complete details.

7. EQUIPMENT and SUPPLIES

7.1 Assay Platform

Dimension Xpand® System

7.2 Equipment

- Refrigerator capable of sustaining 2–8°C.
- Freezer capable of sustaining range not to exceed -20 to -70°C.
- Centrifuge

From revised 12/02/2007

7.3 Supplies

- Plastic serum tubes and serum cups
- Reagent Grade water (Millipore® or equivalent)
- Calibrated pipettes and disposable tips

8. PROCEDURE

ETOH Flex® reagent cartridge Cat. No. DF22 is required to perform this test.

Ethyl Alcohol is performed on the Dimension Xpand® clinical chemistry system after the method is calibrated (see Reference Material in Calibration section) and Quality Controls are acceptable.

NOTE: For all procedures involving specimens, buttoned lab coats, gloves, and face protection are required minimum personal protective equipment. Report all accidents to your supervisor.

The package insert for a new lot of kits must be reviewed for any changes before the kit is used. A current Package Insert is included as a Related Document.

8.1	Instrument Set-Up Protocol
1.	For instrument set up and operation: Refer to Startup and Maintenance, Siemens Dimension Xpand® procedure.
2.	Check reagent inventory
3.	Sampling, reagent delivery, mixing and processing of results are automatically performed by the Dimension Xpand® system. For details of the automated parameters, see below under "Test conditions."

8.2	Specimen/Reagent Preparation
1.	Centrifuge the specimens.
2.	Specimens are placed in Dimension® Xpand segments for analysis by the instrument. Refer to the Sample Processing, Siemens Dimension® Xpand procedure. The sample container (if not a primary tube) must contain sufficient quantity to accommodate the sample volume plus 50 µL of dead volume. Precise container filling is not required.

8.3	Specimen Testing
1.	For QC placement and frequency, refer to the Dimension® Xpand QC Schedule in the Laboratory QC Program.
2.	Follow the instructions, outlined in the Dimension® Xpand Operators Manual

8.3	Specimen Testing
3.	The instrument reporting system contains error messages to warn the user of specific malfunctions. Results followed by such error messages should be held for follow-up. Refer to the Dimension Xpand® system manual "Error messages" section for troubleshooting.
4.	Follow protocol in Section 10.5 "Repeat criteria and resulting" for samples with results above or below the Analytical Measurement Range (AMR). Investigate any failed delta result and repeat, if necessary.
5.	Append the appropriate English text code qualifier messages to any samples requiring a comment regarding sample quality and/or any other pertinent factors.

Test Conditions	
Sample Volume:	9 µL
Buffer Volume:	225 µL
Enzyme Reagent Volume:	121 µL
Temperature:	37° C
Wavelength:	340 and 383 nm
Type of measurement:	Bichromatic rate

9. CALCULATIONS

The instrument automatically calculates the concentration of Ethyl Alcohol in mg/dL.

10. REPORTING RESULTS AND REPEAT CRITERIA

10.1 Interpretation of Data

None required

10.2 Rounding

No rounding is necessary. Instrument reports results as a whole number.

10.3 Units of Measure

mg/dL

10.4 Clinically Reportable Range (CRR)

3 - 900 mg/dL

10.5 Repeat Criteria and Resulting

All repeats must replicate the original result within the total allowable error (TEa) of the assay. Refer to TEa listing for specific information.

Values that fall within the AMR or CRR may be reported without repeat. Values that fall outside these ranges must be repeated.

IF the result is ...	THEN...
< 3 mg/dL	Assure there is sufficient sample devoid of bubbles, cellular debris, and/or fibrin clots. Report as: < 3 mg/dL
≥ 300 mg/dL	On Board Automated Dilution: Results ≥ 300 mg/dL will automatically have repeat testing performed into the instrument using dilution factor of 1.5. No multiplication is necessary.
> 450 mg/dL	Manual Dilution: Using the primary tube, make the smallest dilution possible to bring the raw data within the AMR. Maximum allowable dilution: x 3 Diluent: Reagent Grade Water Enter dilution factor as a whole number on the "Enter Sample Data" screen.
> 900 mg/dL	If the recommended dilution does not give results within the clinically reportable range, report as: "> 900 mg/dL-REP" Bring to the attention of your supervisor prior to releasing result.

Message	Code
Verified by repeat analysis	Append -REP to the result.

11. EXPECTED VALUES

11.1 Reference Ranges

< 5 mg/dL

11.2 Critical Values

> 400 mg/dL

11.3 Priority 3 Limit(s)

None established

From revised 2/02/2007

12. CLINICAL SIGNIFICANCE

Alcohol (ethyl alcohol, ethanol) is the most frequently performed medico legal test, and is the most common toxic substance encountered. In addition to beverages, products containing alcohol in significant amounts include mouthwashes, colognes, and medicinal preparations. Measurements of alcohol levels are used to determine legal impairment, for forensic purposes, in the diagnosis and treatment of alcohol dependency and in emergency settings to detect alcohol poisoning.

Alcohol's deleterious effects are well documented. It has been linked with birth defects (fetal alcohol syndrome), cardiac conditions, high blood pressure, liver disease, and mental deterioration. It is by far the leading cause of death from hepatic failure. Additionally, alcohol-induced behavior is a contributing factor in the majority of accidents and murders.

Within approximately one hour of ingestion, alcohol will have permeated all tissues of the body in proportion to water content. Some alcohol is absorbed while in the stomach, but the principal site of absorption is the upper portion of the small intestine. Rate of absorption is dependent upon emptying time of the stomach, which is subject to various influences. Since alcohol distributes evenly throughout the body water, its concentration in blood following a known dose may be estimated indirectly by measuring concentrations in urine, serum, or plasma.

About 95% of the elimination of alcohol from the body is accomplished by metabolism in the liver. The remainder is excreted unchanged by the lungs, kidneys, and in the feces. Alcohol is rapidly metabolized so that a moderate dose will clear from the blood in approximately one hour.

13. PROCEDURE NOTES

- **FDA Status:** FDA Approved/cleared
- **Validated Test Modifications:** None

The instrument reporting system contains error messages to warn the operator of specific malfunctions. Any report slip containing such error messages should be held for follow-up. Refer to your Dimension Xpand Operator's Guide.

A system malfunction may exist if the following 5-test precision is observed:

Activity	S.D.
100 mg/dL	> 3 mg/dL
300 mg/dL	> 6 mg/dL

14. LIMITATIONS OF METHOD

14.1 Analytical Measurement Range (AMR)

3 – 300 mg/dL

From revised 2/02/2007

14.2 Precision

Material	Mean mg/dL	Standard Deviation (%CV)	
		Repeatability	Within-Lab
Bio-Rad Ethanol/Ammonia			
Level 1	39	0.3 (0.8)	1.0 (2.6)
Level 2	98	0.7 (0.8)	2.4 (2.5)
Level 3	255	1.5 (0.6)	5.4 (2.1)
Plasma Pool	246	1.4 (0.6)	3.0 (1.2)
Serum Pool	102	1.2 (1.2)	1.4 (1.4)

14.3 Interfering Substances

HIL Interference:

The ETOH method was evaluated for interference according to CLSI/NCCLS EP7-A2. Bias, defined as the difference between the control sample (does not contain interferent) and the test sample (contains interferent), is shown in the table below. Bias exceeding 10% is considered “interference”.

Substance tested	Substance Concentration SI Units	ETOH mg/dL	Bias %
Hemoglobin (hemolysate)	1000 mg/dL	100	<10
Bilirubin (unconjugated)	80 mg/dL	96	<10
Bilirubin (conjugated)	80 mg/dL	97	<10
Lipemia Intralipid®	3000 mg/dL	102	<10

14.4 Clinical Sensitivity/Specificity/Predictive Values

Not available

15. SAFETY

The employee has direct responsibility to avoid injury and illness at work. Nearly all harmful exposures to infectious substances and chemicals, and other injuries, can be avoided with effective training and consistent safe work practices.

Become familiar with the Environmental Health and Safety (EHS) Manual to learn the requirements on working safely and protecting the environment from harm. Although lab work typically focuses on the hazards of working with specimens and chemicals, we must also control other important hazards.

- Slips, trips, and falls cause many serious injuries. Please ensure that spills are cleaned quickly (to avoid slippery floors) and that you can see and avoid obstacles in your path.
- Ergonomic injuries result from performing tasks with too much repetition, force, or awkward position. Ergonomic injuries include strains and back injuries. Learn about ergonomic hazards and how to prevent this type of injury.

- Scratches, lacerations, and needlesticks can result in serious health consequences. Attempt to find ways to eliminate your risk when working with sharp materials.

Report all accidents and injuries **immediately** to your supervisor or the business unit Environmental Health and Safety Manager or Specialist.

16. RELATED DOCUMENTS

1. Dimension Xpand® Clinical Chemistry System Operator’s Manual
2. Calibration / Verification Siemens Dimension® Xpand procedure
3. Dimension Xpand® Cal Accept Guidelines
4. Dimension Xpand® Calibration summary
5. Sample Processing, Siemens Dimension® Xpand procedure
6. Start up and Maintenance, Siemens Dimension® Xpand procedure
7. Laboratory Quality Control Program
8. QC Schedule for Siemens Dimension Xpand®
9. Laboratory Safety Manual
10. Material Safety Data Sheets (MSDS)
11. Siemens Dimension Xpand® Limits Chart (AG.F143)
12. Quest Diagnostics Records Management Procedure
13. Dimension Xpand® System Error Messages Chart
14. Centrifuge Use, Maintenance and Functions Checks (Lab policy)
15. Hemolysis, Icteria and Lipemia Interference (Lab policy)
16. Repeat Testing Requirements (Lab policy)
17. Critical Values (Lab policy)
18. Current Allowable Total Error Specifications at http://questnet1.qdx.com/Business_Groups/Medical/qc/docs/qc_bpt_tea.xls
19. Current package insert ETOH Flex® Reagent Cartridge DF22

17. REFERENCES

1. Package Insert, ETOH Flex® Reagent Cartridge DF22, Siemens Healthcare Diagnostics Inc., 01/31/2009.
2. Package Insert, CHEM III Calibrator, Siemens Healthcare Diagnostics Inc., 04/2013.
3. Package Insert, Liquichek Ethanol/Ammonia Control, Bio-Rad Laboratories, 03/2013.

18. REVISION HISTORY

Version	Date	Section	Reason	Reviser	Approval
000	7/15/14	5, 17	Revised to reflect new CHEM III calibrator	A Chini	R SanLuis
000	7/15/14	Footer	Version # leading zero’s dropped due to new EDCS in use as of 10/7/13	L Barrett	R SanLuis

19. ADDENDA

None