

## TRAINING UPDATE

**Lab Location:** SGAH and WAH      **Date Implemented:** 12.22.2014  
**Department:** Blood Bank      **Due Date:** 01.15.2015

### DESCRIPTION OF PROCEDURE REVISION

#### **Name of procedure:**

Cold Agglutinin Screen

#### **Description of change(s):**

1. Note this procedure replaces the previous procedure, "Short Cold Panel."
2. We now have a test that will be used to enter this directly into the LIS (CAGS).
  - a. The test will be reported as positive, negative, or inconclusive.
  - b. The test will automatically generate billing.
3. When we perform this test, we will no longer include A1 and B cells or cord cells. We will simply run screen cells and an auto control. The antibody ID will be reported as "non-specific cold antibody" or "cold autoagglutinin" only.

**Electronic Document Control System**



**Document No.:** WAHBB893[0]

**Title:** Cold Agglutinin Screen

**Owner:** LESLIE.X.BARRETT LESLIE BARRETT

**Status** INWORKS

**Effective Date:** 18-Jan-2015

**Next Review Date:**

Quest Diagnostics

Site: Washington Adventist Hospital

Title: Cold Agglutinin Screen

Technical SOP

<b>Title</b>	<b>Cold Agglutinin Screen</b>		
<b>Prepared by</b>	Stephanie Codina	<b>Date:</b>	12/18/2014
<b>Owner</b>	Stephanie Codina	<b>Date:</b>	12/18/2014

<b>Laboratory Approval</b>		<b>Local Effective Date:</b>	
<b>Print Name and Title</b>	<b>Signature</b>	<b>Date</b>	
<i>Refer to the electronic signature page for approval and approval dates.</i>			

<b>Review</b>		
<b>Print Name</b>	<b>Signature</b>	<b>Date</b>

Document: WAH.BB893[0] Status: INWORKS, Effective: 1/18/2015, Check Version Before Use

Form revised 10/31/02

Document: WAH.BB893[0] Status: INWORKS, Effective: 1/18/2015, Check Version Before Use

**TABLE OF CONTENTS**

1. Test Information.....2

2. Analytical Principle .....3

3. Specimen Requirements.....3

4. Reagents .....4

5. Calibrators/Standards .....4

6. Quality Control .....4

7. Equipment And Supplies .....4

8. Procedure .....5

9. Calculations.....6

10. Reporting Results And Repeat Criteria.....6

11. Expected Values.....7

12. Clinical Significance.....7

13. Procedure Notes .....7

14. Limitations Of Method .....7

15. Safety .....7

16. Related Documents .....8

17. References.....8

18. Revision History .....8

19. Addenda .....8

**1. TEST INFORMATION**

Assay	Method/Instrument	Local Code
Cold Agglutinin Screen	Tube test	N/A

Synonyms/Abbreviations
Short Cold Panel

Department
Blood Bank

Document: WAH-BB893[0] Status: INWORKS, Effective: 1/18/2015, Check Version Before Use

**2. ANALYTICAL PRINCIPLE**

Cold agglutinins are those antibodies that optimally react at temperatures between 4-25°C. Cold agglutinins can be auto or allo in nature; they are differentiated by running an autologous control. These antibodies rarely cause destruction of red blood cells, because body temperature is closer to 37°C. However, the antibodies do interfere with serologic testing in the blood bank and may mask the reactions of other, clinically-significant antibodies. Identifying the specificity of cold antibodies is generally not necessary as this information tends to be more academic than clinically-necessary.

**3. SPECIMEN REQUIREMENTS**

Refer to procedure 'Sample Specifications for Blood Bank Testing' for labeling requirements.

**3.1 Patient Preparation**

N/A

**3.2 Specimen Type & Handling**

Criteria	
<b>Type</b> -Preferred -Other Acceptable	Preferred: EDTA Other acceptable: ACD, CPD, CPDA-1, CP2D, oxalate, or clotted blood
<b>Collection Container</b>	Vacutainer
<b>Volume</b> - Optimum - Minimum	1 ml 1 ml
<b>Transport Container and Temperature</b>	Transport vacutainer at room temperature
<b>Stability &amp; Storage Requirements</b>	Room Temperature: within 8 hours
	Refrigerated: 1 to 10C for 48 hours
	Frozen: Plasma or serum can be stored frozen indefinitely
<b>Timing Considerations</b>	EDTA samples must be tested within 48 hours of collection
<b>Unacceptable Specimens &amp; Actions to Take</b>	Heparin, sodium citrate, or vacutainers with gel separators are not acceptable and must be recollected.
<b>Compromising Physical Characteristics</b>	Specimens must be aseptically collected
<b>Other Considerations</b>	Not applicable

Document: WAH.BB893[0] Status:INWORKS,Effective:1/18/2015, Check Version Before Use

**4. REAGENTS**

Refer to the Material Safety Data Sheet (MSDS) supplied with the reagents for complete safety hazards. Refer to the section in this procedure covering "SAFETY" for additional information.

**4.1 Reagent Summary**

Reagents / Kits	Supplier & Catalog Number
Screening cells	Immucor, Cat.# 2381 or equivalent

**4.2 Reagent Preparation and Storage**

**NOTES:** Date and initial all reagents upon opening. Each container must be labeled with (1) substance name, (2) lot number, (3) date of preparation, (4) expiration date, (5) initials of tech, (6) any special storage instructions; check for visible signs of degradation.

Refer to the Material Safety Data Sheet (MSDS) for a complete description of hazards. If a specific hazard is present, it will be noted in this procedure when the hazard is first encountered in a procedural step.

Screening Cells	
Container	10ml each
Storage/Stability	1-10C / Stable until manufacturer's expiration date.
Preparation	Resuspend red cells before use by gently inverting each vial several times.

**5. CALIBRATORS/STANDARDS**

N/A

**6. QUALITY CONTROL**

N/A

**7. EQUIPMENT and SUPPLIES**

**7.1 Assay Platform**

N/A

**7.2 Equipment**

Serological centrifuge  
Refrigerator

**7.3 Supplies**

Test tubes, (10 x 75 mm and/or 12 x 75 mm)  
Pipettes

**8. PROCEDURE**

**NOTE: For all procedures involving specimens, buttoned lab coats, gloves, and face protection are required minimum personal protective equipment. Report all accidents to your supervisor.**

Step	Action
1	Confirm the specimen acceptability per procedure, "Sample Specifications for Blood Bank Testing."
2	Label 5 tubes with the patient identifiers. A. At a minimum, tubes will contain the first 3 digits of the patient's last name or the patient's last and first initial. B. Additional identifiers will be used if duplicate initials or numbers are present in a testing batch.
3	Prepare a cell suspension from the patient/test cells in one of the tubes per procedure.
4	Label each of the remaining tubes with one of the following: A. Label one tube "I" or "1." B. Label one tube "II" or "2." C. Label one tube "III" or "3." D. Label one tube "AC" or "AUTO."
5	Add 2 drops of patient's plasma/serum to each tube.
6	Add 1 drop of cells to each of the appropriate tubes. A. Add 1 drop of screen cell I to the tube labeled "I" or "1." B. Add 1 drop of screen cell II to the tube labeled "II" or "2." C. Add 1 drop of screen cell III to the tube labeled "III" or "3." D. Add 1 drop of the patient/test cell suspension to the tube labeled "AC" or "AUTO."
7	Mix well and incubate at 4°C for at least 15 minutes.
8	Serofuge immediately after incubation for the time listed on the serofuge.

Document: WAH.BB893[0] Status: INWORKS, Effective: 1/18/2015, Check Version Before Use

Step	Action
9	Without delay, resuspend the cells and observe macroscopically for agglutination with the use of an agglutination viewer. Any delay in reading the reactions may cause the agglutination to disperse.
10	Record the reactions on an antigram that corresponds to the correct lot of screen cells used.

**9. CALCULATIONS**

N/A

**10. REPORTING RESULTS AND REPEAT CRITERIA**

**10.1 Interpretation**

Positive: Agglutination or hemolysis is present in any screening cell used for testing.  
 Negative: No agglutination or hemolysis is presented in any screening cell used for testing.

The autologous control will determine the specificity of the antibody in a positive cold agglutinin screen.

A. Negative autocontrol: Interpret a positive as a "Non-specific cold antibody."

B. Positive autocontrol: Interpret a positive as a "Cold autoantibody."

A positive autocontrol in the presence of a negative cold agglutinin screen must be investigated further. Interpret as "Inconclusive results."

**10.2 Rounding**

N/A

**10.3 Units of Measure**

N/A

**10.4 Clinically Reportable Range (CRR)**

N/A

**10.5 Repeat Criteria and Resulting**

N/A



Document: WAH.BB893[0] Status: INWORKS, Effective: 1/18/2015, Check Version Before Use

**11. EXPECTED VALUES**

**11.1 Reference Ranges**

N/A

**11.2 Critical Values**

N/A

**11.3 Priority 3 Limit(s)**

N/A

**12. CLINICAL SIGNIFICANCE**

None

**13. PROCEDURE NOTES**

- **FDA Status:** LDT without message
- **Validated Test Modifications:** None

Even when a cold agglutinin is identified, it is possible that there is an underlying alloantibody present. If all significant antibodies cannot be ruled out on a panel when a cold reacting autoantibody is present, further antibody work-up is required. A LISS or prewarm screen should be considered.

**14. LIMITATIONS OF METHOD**

N/A

**15. SAFETY**

The employee has direct responsibility to avoid injury and illness at work. Nearly all harmful exposures to infectious substances and chemicals, and other injuries, can be avoided with effective training and consistent safe work practices.

Become familiar with the Environmental Health and Safety (EHS) Manual to learn the requirements on working safely and protecting the environment from harm. Although lab work typically focuses on the hazards of working with specimens and chemicals, we must also control other important hazards.

- Slips, trips, and falls cause many serious injuries. Please ensure that spills are cleaned quickly (to avoid slippery floors) and that you can see and avoid obstacles in your path.
- Ergonomic injuries result from performing tasks with too much repetition, force, or awkward position. Ergonomic injuries include strains and back injuries. Learn about ergonomic hazards and how to prevent this type of injury.

- Scratches, lacerations, and needlesticks can result in serious health consequences. Attempt to find ways to eliminate your risk when working with sharp materials.

Report all accidents and injuries immediately to your supervisor or the business unit Environmental Health and Safety Manager or Specialist.

**16. RELATED DOCUMENTS**

- SOP: Sample Specifications for Blood Bank Testing
- SOP: Preparing a 2-4% Cell Suspension for Testing
- SOP: Prewarmed Antiglobulin Technique
- SOP: Antibody Screen, LISS Tube Method

**17. REFERENCES**

- A. Bryant N, An Introduction to Immunohematology, WB Saunders Co., 1976, p 196-199.
- B. Johns, G., Gockel-Blessing, E, Zundel, W, and Denesiuk, L. 2015. Pearson Education, Upper Saddle River, NJ.

**18. REVISION HISTORY**

Version	Date	Section	Reason	Reviser	Approval

**19. ADDENDA**

Appendix A: LIS entry of the Cold Agglutinin Screen

Document: WAH.BB893[0] Status:INWORKS,Effective:1/18/2015, Check Version Before Use

Form revised 10/31/02

**Appendix A: LIS entry of the Cold Agglutinin Screen**

Step	Action
1	Access TS order using Sunquest function, "Blood Order Processing."
2	In the "Add Spec Test" field, type ";CAGS" then press the "Tab" key or click on the ellipses and search for the cold agglutinin screen. This will add the cold agglutinin screen to the TS specimen.
3	Result the CAGS test with one of the following: A. Enter "P" for positive. B. Enter "N" for negative. C. Enter ")" for inconclusive result.  Note: Antibody identification should be entered in the ABI field per procedure.
4	Click the "Save" button.
5	Billing for this test will automatically generate when the test is ordered in Sunquest.

Document: WAH.BB893[0] Status:INWORKS,Effective:1/18/2015, Check Version Before Use

Form revised 10/31/02