

TRAINING UPDATE


Lab Location: WAH
Department: Core

Date Distributed: 1/8/2016
Due Date: 1/31/2016
Implementation: 2/1/2016

DESCRIPTION OF PROCEDURE REVISION

Name of procedure:
Core Laboratory Biological Safety Cabinet Preventive Maintenance Chart AG.F84.1
Description of change(s):
<p>Removed old gauge certification, incorrect range and UV light instruction.</p> <p>Added -</p> <ul style="list-style-type: none">• space to record gauge certification date (<i>done every year, usually in October</i>)• space to record gauge acceptable value and instruction (<i>range is greater than or equal to most recent certification value</i>) See photos on next page <p>This revised FORM will be implemented on February 1, 2016</p>

Document your compliance with this training update by taking the quiz in the MTS system.



October 2016 gauge reading was 0.60 therefore current acceptable range is ≥ 0.60

Core Laboratory Biological Safety Cabinet Preventive Maintenance Chart

Equipment Description: Biosafety Hood

Equipment ID: 133886110609
Gauge Certification Date: _____

Gauge Acceptable Performance Range: \geq _____

Gauge acceptability is defined as greater than or equal to the most recent gauge certification value (certification performed annually)

Frequency: Performed by each shift that the cabinet is used. Clean BSC with disinfectant and rinse with dH2O before and after use. Only supplies needed for test being performed should be inside hood. Do not store reagent kits in Cabinet. Blower should run continuously. **DO NOT TURN OFF. Turn on UV for 15 minutes per day on the Night Shift.**

MONTH / YEAR: _____

Date	Day Shift					Evening Shift				Night Shift			
	Temp 22-28C	Visual Airflow	Gauge Reading	Clean & Rinse	Tech	Visual Airflow	Gauge Reading	Clean & Rinse	Tech	Visual Airflow	Gauge Reading	Clean & Rinse	Tech
1													
2													
3													
4													
5													
6													
7													
8													
9													
10													
11													
12													
13													
14													
15													
16													
17													
18													
19													
20													
21													
22													
23													
24													
25													
26													
27													
28													
29													
30													
31													
Weekly review:				Weekly review:					Weekly review:				
Weekly review:				Weekly review:					Monthly review:				