

## TRAINING UPDATE

**Lab Location:**

GEC, SGMC & WAH

**Date Distributed:**

3/13/2016

**Department:**

Mgmt, IT & QA

**Due Date:**

3/31/2016

**Implementation:**

**3/31/2016**

### DESCRIPTION OF PROCEDURE REVISION

**Name of procedure:**

**Verification of Patient Results SGAH/WAH.IT24 v2**

**Medical Director Patient Report Review AG.F300.1**

**Description of change(s):**

Current SOP titled “Review of Cumulative Reports” (QA45) will be retired.

Process for medical director review of patient reports will include printed LIS copies and HIS screen shots. It will be performed in tandem with the process to review result transmission from origin to LIS / HIS posting as described in the SOP titled “Verification of Patient Results.”

Form F300 has been re-titled (replaced ‘Cumulative’ with ‘Patient’)

**This revised SOP & form will be implemented on March 31, 2016**

**Document your compliance with this training update by taking the quiz in the MTS system.**

**Approved draft for training (version 2)**

Non-Technical SOP

<b>Title</b>	<b>Verification of Patient Results</b>	
<b>Prepared by</b>	Marie Sabonis	Date: 12/15/2008
<b>Owner</b>	Marie Sabonis	Date: 12/15/2008

<b>Laboratory Approval</b>		
<b>Print Name and Title</b>	<b>Signature</b>	<b>Date</b>
<i>Refer to the electronic signature page for approval and approval dates.</i>		
Local Issue Date:		Local Effective Date:

<b>Review:</b>		
<b>Print Name</b>	<b>Signature</b>	<b>Date</b>

## **TABLE OF CONTENTS**

1. PURPOSE.....	2
2. SCOPE .....	2
3. RESPONSIBILITY.....	2
4. DEFINITIONS.....	2
5. PROCEDURE.....	2
6. RELATED DOCUMENTS .....	3
7. REFERENCES .....	3
8. REVISION HISTORY.....	3
9. ADDENDA AND APPENDICES.....	4

---

### **1. PURPOSE**

To describe the system in place to verify that patient results are accurately transmitted from the point of data entry (interfaced and manual input) to all types of patient reports.

### **2. SCOPE**

LIS staff verifies that results from the starting point, i.e., the instrument, post correctly in the LIS and then finally posts correctly in HIS. The documentation is provided to the Medical Directory for review.

### **3. RESPONSIBILITY**

LIS staff and Medical Director perform this procedure.

### **4. DEFINITIONS**

N/A

### **5. PROCEDURE**

- a. LIS staff will obtain result examples for Hematology, Coagulation, Manual Entry, Chemistry, Urinalysis, Reference lab (Micro and Non-Micro), Blood Bank, Micro In house and Point of Care. The following is obtained, if applicable
  - Instrument printouts or manual log sheets
  - Instrument LIS trace file, if applicable
  - LIS Interim and Cumulative reports
  - HIS Print screen
  - Copy of Reference Lab report (Micro and Non-Micro)

- b. LIS staff reviews from point of entry to last point of output (HIS). Review is documented on the chart titled HIS and LIS Mini Validation (Addenda A).
- c. All documents, HIS and LIS Mini Validation chart and **Medical Director Patient Report Review (see Related Documents) are submitted for review and sign off by the Medical Director.**
- d. This process will be completed biennially, and in the following cases:

**NEW INSTRUMENTS**

New instruments are tested in TEST AREA before bringing them online in Production. In the test area, all test codes and calculations (if defined) are exercised, the instrument printout verified for expected results and then accepted in the LIS. Using Inquiry in the LIS, the results are checked against the printout from the instrument for accuracy; this must include High (H), Low (L) and C (Critical) flags. A patient sample should also be run on an instrument and checked for accuracy on the HIS side.

**PATCH, SOFTWARE VERSION UPDATES**

All online instruments are tested in the TEST area for accuracy by comparing data sent by instrument to LIS and HIS TEST.

- Instrument printout
- Inquiry
- HIS Result Display

**IMPLEMENTATION OF AN INTERFACE**

- e. Documentation is kept in LIS offices at respective sites.

**6. RELATED DOCUMENTS**

Medical Director Patient Report Review (AG.F300)

**7. REFERENCES**

College of American Pathologists, Laboratory Accreditation Manual, Laboratory General Checklist, current version

**8. REVISION HISTORY**

Version	Date	Reason for Revision	Revised By	Approved By
		Supersedes SOP IT025.001		
000	12/30/13	Section 2 & 5: Replace SMS with HIS Section 9: Revise title and checklist by removing specific systems (SMS & Sunquest) Footer: version # leading zero's dropped due to new EDCS in use as of 10/7/13.	M Sabonis	M Sabonis

Form revised 3/31/00

1	2/22/16	Section 2: add Medical Director review Section 5: add documents for review, Medical Director sign-off and interface implementation Section 6: add form Section 7: remove General CAP Checklist number Section 9: update HIS and LIS Mini Validation checklist to include Medical Director sign-off.	M Sabonis	M Sabonis

**9. ADDENDA AND APPENDICES**

- A. HIS and LIS Mini Validation checklist (see Attachments on the document Profile)



**MEDICAL DIRECTOR PATIENT REPORT REVIEW**

Year:

\_\_\_\_\_  
Medical Director\_\_\_\_\_  
DateNo Issues noted: \_\_\_\_\_  

---

Issues Noted:	Action taken to resolve:	Date Resolved:	By:
1.			
2.			
3.			
4.			
5.			
6.			
7.			