

TRAINING UPDATE

Lab Location: SGMC & WAH
Department: Core

Date Distributed: 7/22/2016
Due Date: 8/16/2016
Implementation: 8/16/2016

DESCRIPTION OF PROCEDURE REVISION

Name of procedure:

Urine Amphetamine/Methamphetamine Screen by Vista® System	SGAH.C120 v1
Urine Barbiturates Screen by Dimension Vista® System	SGAH.C121 v1
Urine Benzodiazepines Screen by Dimension Vista® System	SGAH.C122 v1
Urine Cannabinoids Screen by Dimension Vista® System	SGAH.C123 v1
Urine Cocaine Metabolite Screen by Dimension Vista® System	SGAH.C124 v1
Urine Opiates Screen by Dimension Vista® System	SGAH.C125 v1
Urine Phencyclidine Screen by Dimension Vista® System	SGAH.C126 v1

Note: all have been converted to system SOPs

Description of change(s):

All SOPs have the following changes -

Section	Reason
Header	Add WAH
3.1, 3.2	Add urine collection kit
4.2	Add safety instructions
5.2	Remove uncapped calibrator storage
6.4, 6.6	Replace LIS with Unity Real Time
11.3	Add report comments
16	Update document titles
17	Update PI revision dates

The Amphetamine/Methamphetamine SOP is attached in its entirety.

Only the page that shows report comments is attached for the other drugs –

- the comment when a drug screen is resulted covers all the drugs
- if the drug is ordered individually, then the comment is specific for that drug

These revised SOPs will be implemented on August 16, 2016

Document your compliance with this training update by taking the quiz in the MTS system.

Technical SOP

Title	Urine Amphetamine/Methamphetamine Screen by Dimension Vista® System	
Prepared by	Ashkan Chini	Date: 6/25/2012
Owner	Robert SanLuis, Jean Buss	Date: 7/6/2016

Laboratory Approval		Local Effective Date:
Print Name and Title	Signature	Date
<i>Refer to the electronic signature page for approval and approval dates.</i>		

Review		
Print Name	Signature	Date

TABLE OF CONTENTS

1. Test Information.....2
 2. Analytical Principle3
 3. Specimen Requirements.....3
 4. Reagents4
 5. Calibrators/Standards5
 6. Quality Control7
 7. Equipment And Supplies9
 8. Procedure9
 9. Calculations.....11
 10. Reporting Results And Repeat Criteria.....11
 11. Expected Values.....11
 12. Clinical Significance12
 13. Procedure Notes13
 14. Limitations Of Method13
 15. Safety14
 16. Related Documents14
 17. References.....14
 18. Revision History15
 19. Addenda15

1. TEST INFORMATION

Assay	Method/Instrument	Local Code
Urine Amphetamine, Qualitative	Dimension Vista® System	UAMPT


Synonyms/Abbreviations
“Speed”/AMP Included in Batteries/Packages: Urine Drug Screen: UDRGS (SGAH) & UDRGW (WAH)

Department
Chemistry

Form revised 2/02/2007

2. ANALYTICAL PRINCIPLE

The principle of this test is based on the competition for antibody binding sites between drug in the sample and drug labeled with the enzyme glucose-6-phosphate dehydrogenase (G6PDH). Matched lots of monoclonal antibody reactive to d-amphetamine and d-methamphetamine and d-amphetamine and d-methamphetamine labeled with glucose-6-phosphate dehydrogenase are used in this Syva® Emit® II Plus methodology.



Where: Ab =antibody reactive to d-amphetamine and d-methamphetamine
 AMPH =amphetamines and methamphetamines
 AMPH-G6PDH = d-amphetamine and d-methamphetamine glucose-6-phosphate dehydrogenase conjugates

The concentration of drug in the sample determines the amount of AMPH-glucose-6-phosphate dehydrogenase (AMPH-G6PDH) conjugate that is bound to the antibody. The unbound conjugate catalyzes the oxidation of glucose-6-phosphate, with the simultaneous reduction of NAD⁺ to NADH, more rapidly than does the bound conjugate. The rate of increasing absorbance at 340 nm due to the increase in NADH is related to the concentration of drug in the sample.

3. SPECIMEN REQUIREMENTS

3.1 Patient Preparation

Component	Special Notations
Fasting/Special Diets	N/A
Specimen Collection and/or Timing	Freshly voided urine specimens should be used for testing.
Special Collection Procedures	No additives or preservatives are needed. Adulteration of the urine specimen may cause erroneous results. If adulteration is suspected, obtain a fresh specimen. Urine specimens should be handled and treated as if they are potentially infected. Preferred method is the Urine Collection Kit with specimen transferred to Urine Chemistry Collection Tube (yellow top).
Other	If Urine Collection Kit is not used, submit to Laboratory within 2 hours of collection.

Form revised 2/02/2007

3.2 Specimen Type & Handling

Criteria	
Type -Preferred -Other Acceptable	Urine None
Collection Container	Urine Collection Kit or sterile container
Volume - Optimum - Minimum	15 mL 2 mL
Transport Container and Temperature	Urine Chemistry Collection Tube (yellow top) or container at room temperature.
Stability & Storage Requirements	Room Temperature: Processing/Testing should take place immediately.
	Refrigerated: 24 hours
	Frozen: If storage longer than 24 hours is required.
	Instrument on board aliquot stability 2 hours
Timing Considerations	Deliver specimens to laboratory immediately.
Unacceptable Specimens & Actions to Take	Specimens that are unlabeled, improperly labeled, or those that do not meet the stated criteria are unacceptable. Samples in Urine Analysis Preservative Tubes are NOT acceptable. Request a recollection and credit the test with the appropriate LIS English text code for “test not performed” message. Examples: Quantity not sufficient-QNS; Wrong collection-UNAC. Document the request for recollection in the LIS.
Compromising Physical Characteristics	Turbidity: Centrifuge turbid samples before analysis. Temperature: Specimens should be at a temperature of 20 – 25° C before testing.
Other Considerations	Boric acid should not be used as a preservative. Plastic transfer pipettes should NOT be used for delivering patient specimen to sample cup.

4. REAGENTS

Refer to the Safety Data Sheet (SDS) supplied with the reagents for complete safety hazards. Refer to the section in this procedure covering “SAFETY” for additional information.

4.1 Reagent Summary

Reagents	Supplier & Catalog Number
AMPH	Siemens, Flex® reagent cartridge, Cat. No. K5091

4.2 Reagent Preparation and Storage

NOTES: Each container must be labeled with (1) substance name, (2) lot number, (3) expiration date, (4) any special storage instructions; check for visible signs of degradation. When placed onboard the analyzer, the instrument captures the date / time loaded and calculates and tracks the opened expiration. Refer to the Safety Data Sheet (SDS) for a complete description of hazards. If a specific hazard is present, it will be noted in this procedure when the hazard is first encountered in a procedural step.

**Contains 2-methyl-4-isothiazolin-3-one.
 Harmful to aquatic life with long lasting effects.
 Wear protective clothing, gloves and eye/face protection.**

Reagent	AMPH
Container	Reagent cartridge
Storage	Store at 2-8° C
Stability	<ul style="list-style-type: none"> • Reagent is stable until expiration date stamped on the reagent cartridges. • Sealed wells on the instrument are stable for 30 days. • Once wells 1 - 12 have been entered by the instrument, they are stable for 2 days.
Preparation	All reagents are liquid and ready for use.

5. CALIBRATORS/STANDARDS

5.1 Calibrators/Standards Used

Calibrator	Supplier and Catalog Number
UDAT CAL	Siemens Dimension Vista®, Cat. No. KC510

5.2 Calibrator Preparation and Storage

NOTE: Date and initial all calibrators upon opening. Each container must be labeled with (1) substance name, (2) lot number, (3) date of preparation, (4) expiration date, (5) any special storage instructions; check for visible signs of degradation. When placed onboard the analyzer, the instrument captures the date / time loaded and calculates and tracks the opened expiration.

Calibrator	UDAT CAL
Preparation	Calibrator is ready for use. No preparation is required.
Storage/Stability	<ul style="list-style-type: none"> • Store at 2-8° C • Unopened calibrator is stable until expiration date stamped on the box.

	<ul style="list-style-type: none"> • Opened Calibrator: once the stopper of the vial is punctured, assigned values are stable for 15 days when stored on board the Dimension Vista System.
--	--

5.3 Calibration Parameter

Criteria	Special Notations
Reference Material	UDAT CAL
Assay Range	0 – 1000 ng/mL
Suggested Calibration Level	1000 ng/mL cutoff. Validate the calibration by assaying a positive and negative control.
Frequency	<ul style="list-style-type: none"> • Every new reagent cartridge lot. • Every 30 days for any one lot • When major maintenance is performed on the analyzer. • When control data indicates a significant shift in assay.
Calibration Scheme	4 levels, n = 5

5.4 Calibration Procedure

Auto Calibration:

1. Place the required calibrator vials in a carrier. Make sure the barcode labels are entirely visible through the slots.
2. Place the carrier in the loading area.
3. Position the carrier with the labels facing away from the user.
4. Press the **Load** button.
5. Automatic calibration requires that calibrators be on the instrument. As the time for processing approaches, the instrument reviews onboard inventory for the appropriate calibrators.

Manual Calibration:

1. Verify that calibrators and reagents are in inventory on the instrument.
2. Press **System > Method Summary > Calibration**.
3. Select a method from the sidebar menu. Press the **Order Calibration** button on the screen.
4. Verify that the information on the screen is correct. Verify that the calibrator lot is correct using the drop-down menu.
 - a. When calibrating using Vials press **OK**.
 - b. When calibrating using Cups, check the Use Cups box. This displays the rack and cup position fields. For additional cups use the positions in ascending order. Be sure to use the number of calibration levels and cups as specified in the method IFU. Scan the rack barcode and place calibrator cups in an adapter in position 1 on a rack. Press **OK** and load the rack on the instrument.
5. The status field in the calibration screen changes sequentially to Awaiting Scheduling, Preparing Calibrators and Processing.

5.5 Tolerance Limits

IF.....	THEN.....
If result fall within assay-specific specification, and QC values are within acceptable limits,	proceed with analysis
If result falls outside assay-specific specification, or QC values are out of Acceptable limits,	troubleshoot the assay and/or instrument and repeat calibration

6. QUALITY CONTROL

6.1 Controls Used

Controls	Supplier and Catalog Number
Liquichek™ Urine Toxicology Control Levels S1E and S2E Low Opiate	Bio-Rad Laboratories Cat. No. 423 and 424

6.2 Control Preparation and Storage

NOTE: Date and initial all controls upon opening. Each container should be labeled with (1) substance name, (2) lot number, (3) date of preparation, (4) expiration date, (5) initials of tech, and (6) any special storage instructions; check for visible signs of degradation. A barcode label is produced and placed on the vial.

Control	Liquichek Urine Toxicology Controls, Levels S1E and S2E
Preparation	Before sampling allow the control to reach room temperature (18-25°C) and swirl gently to ensure homogeneity. Use immediately. After each use, promptly replace the stopper and return to 2-8°C storage.
Storage/Stability	Once the control is opened, all analytes will be stable for 30 days at 2-8°C. Unopened controls are stable until the expiration date at 2-8°C.

6.3 Frequency

Analyze all levels of QC material after every calibration and each day of testing (notated on the QC frequency sheets posted on the instruments).

Refer to the Dimension Vista® QC Schedule in the Laboratory policy Quality Control Program and in the Dimension Vista® Quick Reference Guide.

Form revised 2/02/2007

6.4 Tolerance Limits

Step	Action
1	Acceptable ranges for QC are programmed into the instrument's Quality Control software system and Unity Real Time , and may be posted near the instrument for use during computer downtime.
2	Run Rejection Criteria <ul style="list-style-type: none"> Anytime the established parameters are exceeded, the run is considered out of control (failed) and patient results must not be reported. The technologist must follow the procedure in the Laboratory QC Program to resolve the problem.
3	Corrective Action: <ul style="list-style-type: none"> All rejected runs must be effectively addressed through corrective action. Steps taken in response to QC failures must be documented. Patient samples in failed analytical runs must be <u>reanalyzed according to the Laboratory QC Program</u>. Supervisors may override rejection of partial or complete runs only with detailed documentation and criteria for overrides that are approved by the Medical Director. Consult corrective action guidelines in Laboratory QC Program. Follow corrective action guidelines in the Laboratory QC Program. Corrective action documentation must follow the Laboratory Quality Control Program.
4	Review of QC <ul style="list-style-type: none"> QC must be reviewed weekly by the Group Lead or designee and monthly by the Supervisor/Manager or designee. If the SD and/or CV are greater than established ranges, investigate the cause for the imprecision and document implementation of corrective actions.

6.5 Review Patient Data

Each result is reviewed for error messages. Refer to the Dimension Vista system manual "Error messages" section for troubleshooting. Resolve any problems noted before issuing patient reports.

6.6 Documentation

- QC tolerance limits are programmed into the instrument [and Unity Real Time](#); it calculates cumulative mean, SD and CV and stores all information for easy retrieval.

- Quality control records are reviewed daily at the bench, weekly by the Group Lead or designee, and monthly by the Supervisor/Manager or designee.
- Refer to complete policies and procedures for QC documentation and for record retention requirements in the Laboratory QC Program.

6.7 Quality Assurance Program

- Each new lot number of reagent or new shipment of the same lot of reagent must be tested with external control materials and previously analyzed samples. Performance of the new lot must be equivalent to the previous lot; utilize published TEA for acceptability criteria.
- Training must be successfully completed and documented prior to performing this test. This procedure must be incorporated into the departmental competency assessment program.
- The laboratory participates in CAP proficiency testing. All proficiency testing materials must be treated in the same manner as patient samples.
- Monthly QC must be presented to the Medical Director or designee for review and signature.
- Monthly QC mean and SD are sent to Bio-Rad Laboratories for peer group comparison.
- Consult the Laboratory QC Program for complete details.

7. EQUIPMENT and SUPPLIES

7.1 Assay Platform

Dimension Vista® System

7.2 Equipment

- Refrigerator capable of sustaining 2–8°C.
- Freezer capable of sustaining range not to exceed -20 to -70°C.
- Centrifuge

7.3 Supplies

- Aliquot Plates
- System Fluids
- Assorted calibrated pipettes (MLA or equivalent) and disposable tips

8. PROCEDURE

AMPH Flex® reagent cartridge Cat. No. K5091 is required to perform this test.

Urine Amphetamine/Methamphetamine Screen is performed on the Dimension Vista® System after the method is calibrated (see Reference Material in Calibration section) and Quality Controls are acceptable.

NOTE: For all procedures involving specimens, buttoned lab coats, gloves, and face protection are required minimum personal protective equipment. Report all accidents to your supervisor.

The package insert for a new lot of kits must be reviewed for any changes before the kit is used. A current Package Insert is included as a Related Document.

8.1	Sample Processing
1.	A sample rack holding tubes or cups is placed on the rack input lane.
2.	The sample shuttle moves the rack to the barcode reader which identifies the rack and samples to the system.
3.	The rack moves into the sample server and to the rack positioner.
4.	At the same time, aliquot plates move from the aliquot loader into position.
5.	The aliquot probe aspirates the sample from the tubes or cups and dispenses it into the wells of the aliquot plates.
6.	After each aspirate-dispense action, the probe is thoroughly rinsed inside and out to prevent sample carryover.
7.	When sample aspiration is completed, the sample server moves the rack back to the sample shuttle, where it is placed on the output lane and can be removed by the operator.

8.2	Specimen Testing
1.	For QC placement and frequency, refer to the Dimension Vista® QC Schedule in the Laboratory QC Program.
2.	Follow the instructions, outlined in the Dimension Vista® Operator's Manual
3.	The instrument reporting system contains error messages to warn the user of specific malfunctions. Results followed by such error messages should be held for follow-up. Refer to the Dimension Vista® system manual "Error messages" section for troubleshooting.
4.	Follow protocol in Section 10.5 "Repeat criteria and resulting" for samples with results above or below the Analytical Measurement Range (AMR). Investigate any failed delta result and repeat, if necessary.
5.	Append the appropriate English text code qualifier messages to any samples requiring a comment regarding sample quality and/or any other pertinent factors.

Test Conditions	
Sample Volume:	1.2 µL
Reagent 1 Volume:	98 µL
Reagent 2 Volume:	42 µL
Reaction Time:	5.3 minutes
Test Temperature:	37° C
Wavelength:	340 & 600 nm
Type of measurement:	Rate

9. CALCULATIONS

None

10. REPORTING RESULTS AND REPEAT CRITERIA

10.1 Interpretation of Data

None required

10.2 Rounding

N/A

10.3 Units of Measure

N/A

10.4 Clinically Reportable Range (CRR)

N/A

10.5 Repeat Criteria and Resulting

Specimens that give an “Abnormal Reaction” message must be repeated.

Message	Code
Verified by repeat analysis	Append –REP to the result.

11. EXPECTED VALUES

11.1 Reference Ranges

None detected

11.2 Critical Values

None established

11.3 Standard Required Messages

The following comment is automatically added to the report by the LIS when a urine amphetamine test is ordered:

This is a screening assay. Amphetamines are detected in concentrations at or above 1000 ng/mL. The ingestion of certain herbal or plant products containing Ephedra or its metabolites may cause false positive amphetamine/methamphetamine results. A more specific testing method GCMS is available from the lab.

The following comment is automatically added to the report by the LIS when a urine drug screen is ordered:

The drug of abuse panel is a screening assay. It detects the following drugs of abuse in concentrations at or above the concentrations listed below.

Phencyclidine	25 ng/mL
Benzodiazepines	200 ng/mL
Cocaine	300 ng/mL
Amphetamines	1000 ng/mL
THC	50 ng/mL
Opiates	300 ng/mL
Barbiturates	200 ng/mL

The ingestion of certain herbal or plant products containing Ephedra or its metabolites may cause false positive amphetamine/metamphetamine results.

This test is for medical screening purposes ONLY. For confirmation a separate order for Gas Chromatography by Mass Spectrophotometry (GCMS) is required.

12. CLINICAL SIGNIFICANCE

Amphetamines are central nervous system stimulants that produce wakefulness, alertness, increased energy, reduced hunger, and an overall feeling of well-being. Amphetamines can be taken orally, intravenously, by smoking, or by snorting.

Amphetamines are readily absorbed from the gastrointestinal tract and are then either deactivated by the liver or excreted unchanged in the urine. The relative importance of these elimination modes depends on urinary pH. Amphetamine is metabolized to deaminated (hippuric and benzoic acids) and hydroxylated metabolites. Methamphetamine is partially metabolized to amphetamine, its major active metabolite.

Amphetamines appear in the urine within three hours after any type of administration¹, and can be detected by this Emit® assay for as long as 24–48 hours after the last dose.

13. PROCEDURE NOTES

- **FDA Status:** FDA Approved/Modified
- **Validated Test Modifications:** Removed pH testing per communication from manufacturer.

The instrument reporting system contains error messages to warn the operator of specific malfunctions. Any report slip containing such error messages should be held for follow-up. Refer to your Dimension Vista Operator's Guide.

Concentration	S.D.
300 ng/mL	47 ng/mL
500 ng/mL	86 ng/mL
1000 ng/mL	208 ng/mL

14. LIMITATIONS OF METHOD

- A positive result suggests the likely presence of amphetamines but does not indicate or measure intoxication.
- The presence of amphetamines in urine is only an indication of recent exposure to or use of amphetamines.
- The psychological and physiological effects of amphetamines do not necessarily correlate with urinary concentration.
- A positive AMPH result suggests the likely presence of drug and its metabolites. The AMPH method cannot fully quantitate the concentration of individual components.
- Interpretation of results must take into account that urine concentrations vary extensively with fluid intake, and other biological variables.
- There is a possibility that other substances and/or factors not listed above may interfere with the test and cause false results, e.g., technical or procedural errors.
- **NEGATIVE** results for specimens with concentrations below the assay range may be accompanied by an "assay range" or by a "below assay range" message. These results should be reported as **NEGATIVE**.
- **POSITIVE** results for specimens with concentrations above the assay range may be accompanied by an "assay range" or by an "above assay range" message. These results should be reported as **POSITIVE**.

14.1 Analytical Measurement Range (AMR)

Qualitative Assay: 125 – 1800 ng/mL (for 1000 ng/mL cutoff)

14.2 Precision

Material	Mean ng/mL	Standard Deviation (%CV)	
		Repeatability	Within-Lab
Calibrator Control ng/mL			
225	248	19.1 (8)	20.1 (8)
300	320	11.3 (4)	18.5 (6)
375	391	21.4 (5)	28.0 (7)
500	545	20.5 (4)	26.2 (5)
625	684	40.3 (6)	50.2 (7)

14.3 Interfering Substances

None

14.4 Clinical Sensitivity/Specificity/Predictive Values

Not available

15. SAFETY

The employee has direct responsibility to avoid injury and illness at work. Nearly all harmful exposures to infectious substances and chemicals, and other injuries, can be avoided with effective training and consistent safe work practices.

Become familiar with the Environmental Health and Safety (EHS) Manual to learn the requirements on working safely and protecting the environment from harm. Although lab work typically focuses on the hazards of working with specimens and chemicals, we must also control other important hazards.

- Slips, trips, and falls cause many serious injuries. Please ensure that spills are cleaned quickly (to avoid slippery floors) and that you can see and avoid obstacles in your path.
- Ergonomic injuries result from performing tasks with too much repetition, force, or awkward position. Ergonomic injuries include strains and back injuries. Learn about ergonomic hazards and how to prevent this type of injury.
- Scratches, lacerations, and needlesticks can result in serious health consequences. Attempt to find ways to eliminate your risk when working with sharp materials.

Report all accidents and injuries immediately to your supervisor or the business unit Environmental Health and Safety Manager or Specialist.

16. RELATED DOCUMENTS

1. Dimension Vista® Clinical Chemistry System Operator's Manual
2. Dimension Vista® Calibration/Verification Procedure
3. Dimension Vista® Cal Accept Guidelines

4. Dimension Vista® Calibration summary
5. Dimension Vista® Sample Processing, Startup and Maintenance procedure
6. Laboratory Quality Control Program
7. QC Schedule for Siemens Dimension Vista®
8. Laboratory Safety Manual
9. Safety Data Sheets (SDS)
10. Dimension Vista® Limits Chart (AG.F200)
11. Quest Diagnostics Records Management Procedure
12. Dimension Vista® System Error Messages Chart
13. Centrifuge Use, Maintenance and Functions Checks (Lab policy)
14. Hemolysis, Icteria and Lipemia Interference (Lab policy)
15. Repeat Testing Requirement (Lab policy)
16. Current Allowable Total Error Specifications
at http://questnet1.qdx.com/Business_Groups/Medical/qc/docs/qc_bpt_tea.xls
17. Current package insert AMPH Flex® Reagent Cartridge K5091

17. REFERENCES

1. Package Insert, AMPH Flex® Reagent Cartridge K5091, Siemens Healthcare Diagnostics Inc., 05/05/2015.
2. Package Insert, UDAT CAL, Siemens Healthcare Diagnostics Inc., 05/2011.
3. Package Insert, Liquichek Urine Toxicology Controls, Bio-Rad Laboratories, 01/2015.

18. REVISION HISTORY

Version	Date	Section	Reason	Reviser	Approval
000	7/6/16		Update owner	L Barrett	R SanLuis
000	7/6/16	3.1, 3.2	Add urine collection kit	L Barrett	R SanLuis
000	7/6/16	4.2	Add safety instructions	A Chini	R SanLuis
000	7/6/16	5.2	Remove uncapped calibrator storage	A Chini	R SanLuis
000	7/6/16	6.4, 6.6	Replace LIS with Unity Real Time	L Barrett	R SanLuis
000	7/6/16	11.3	Add report comments	A Chini	R SanLuis
000	7/6/16	16	Update titles	L Barrett	R SanLuis
000	7/6/16	17	Update PI dates	A Chini	R SanLuis
000	7/6/16	Footer	Version # leading zero's dropped due to new EDCS in use as of 10/7/13.	L Barrett	R SanLuis

19. ADDENDA

Validated Test Modification statement from manufacturer

From: Brodbeck, Beate (H USA)
Sent: Thursday, June 24, 2010 3:54 PM
To: SanLuis, Robert; Mcmillan, Wendell R
Subject: pH for urine DAU samples

Hi Robert and Wendell,

Below is the response regarding your inquiry of pH testing for drugs of abuse urines.

Beate Brodbeck

Chemistry Instrument Specialist - Western Maryland

Siemens Healthcare Diagnostics
<blocked::http://www.siemens.com/diagnostics>

C: 410-370-4382 | VM: 800-948-3234 x-2684

beate.brodbeck@siemens.com

From Kevin Mulrooney:

pH correction of urine samples prior to running the Drugs of Abuse assays on Dimension is not an absolute requirement. The Dimension IFU's all say that the acceptable pH range is 5-8. The Syva Emit IFU's for the same tests all say that the acceptable pH range is 3-11, except for THC (pH range 4.5-8). The assays all work at pH 3-11. The vast majority of urines will fall in this range. THC is an exception in that at acid pH <4.5, THC recovery is decreased, and at basic pH >8, THC recovery is increased. When the Dimension IFU's were written, I suppose the decision was made to standardize the pH acceptable range to the most narrow (THC). Dimension customers can run the DAT's without checking pH, but there is a slight chance of inaccuracy with THC. Urine pH outside the 5-8 range is not common, either.

Hope this helps. We don't have this in a formal document, but you can share this information with your customer.

Regards,

Customer Care.
With you every step of the way.

Kevin Mulrooney
Kevin Mulrooney
Staff Product Support Specialist
Global Product Support
Siemens Healthcare Diagnostics
700 GBC Drive M/S 707
Glasgow, DE 19714
(302) 631-8854
Fax (302) 631-7487
kevin.l.mulrooney@siemens.com

11.2 Critical Values

None established

11.3 Standard Required Messages

The following comment is automatically added to the report by the LIS when a urine barbiturates test is ordered:

This is a screening assay. Barbiturates are detected in concentrations at or above 200 ng/mL.

The following comment is automatically added to the report by the LIS when a urine drug screen is ordered:

The drug of abuse panel is a screening assay. It detects the following drugs of abuse in concentrations at or above the concentrations listed below.

Phencyclidine	25 ng/mL
Benzodiazepines	200 ng/mL
Cocaine	300 ng/mL
Amphetamines	1000 ng/mL
THC	50 ng/mL
Opiates	300 ng/mL
Barbiturates	200 ng/mL

The ingestion of certain herbal or plant products containing Ephedra or its metabolites may cause false positive amphetamine/metamphetamine results.

This test is for medical screening purposes ONLY. For confirmation a separate order for Gas Chromatography by Mass Spectrophotometry (GCMS) is required.

12. CLINICAL SIGNIFICANCE

Barbiturates, a class of nervous system depressants, are usually taken orally, but are sometimes injected intravenously or intramuscularly. They are absorbed rapidly; 30–40% is bound to plasma protein, and the rest is distributed to muscle, fat, and to the liver (where they are ultimately inactivated). They are classified based on their duration of action, ranging from very short acting (approximately 15 minutes) to long acting (a day or more). Some of the most commonly abused barbiturates are the short-acting ones, including pentobarbital and secobarbital. An example of a longacting barbiturate is phenobarbital. The ratio of unchanged drug to metabolites varies depending upon duration of action. Short-acting barbiturates will generally be excreted in urine as metabolites, while the long-acting barbiturates will primarily appear unchanged.

11.2 Critical Values

None established

11.3 Standard Required Messages

The following comment is automatically added to the report by the LIS when a urine benzodiazepines test is ordered:

This is a screening assay. Benzodiazepines are detected in concentrations at or above 200 ng/mL.

The following comment is automatically added to the report by the LIS when a urine drug screen is ordered:

The drug of abuse panel is a screening assay. It detects the following drugs of abuse in concentrations at or above the concentrations listed below.

Phencyclidine	25 ng/mL
Benzodiazepines	200 ng/mL
Cocaine	300 ng/mL
Amphetamines	1000 ng/mL
THC	50 ng/mL
Opiates	300 ng/mL
Barbiturates	200 ng/mL

The ingestion of certain herbal or plant products containing Ephedra or its metabolites may cause false positive amphetamine/metamphetamine results.

This test is for medical screening purposes ONLY. For confirmation a separate order for Gas Chromatography by Mass Spectrophotometry (GCMS) is required.

12. CLINICAL SIGNIFICANCE

Benzodiazepines are sedative-hypnotic drugs that are structurally similar and include widely used drugs such as chlordiazepoxide, diazepam, and oxazepam. The different benzodiazepines are absorbed at different rates, and the timing of their psychoactive effects varies with the absorption rate. Benzodiazepines are usually taken orally and are metabolized in the liver. Some benzodiazepine metabolites are pharmacologically active. Benzodiazepines potentiate the effect of other central nervous system depressants, such as ethyl alcohol.

13. PROCEDURE NOTES

- **FDA Status:** FDA Approved/Modified

11.2 Critical Values

None established

11.3 Standard Required Messages

The following comment is automatically added to the report by the LIS when a urine cannabinoids test is ordered:

This is a screening assay. THC is detected in concentrations at or above 50 ng/mL.

The following comment is automatically added to the report by the LIS when a urine drug screen is ordered:

The drug of abuse panel is a screening assay. It detects the following drugs of abuse in concentrations at or above the concentrations listed below.

Phencyclidine	25 ng/mL
Benzodiazepines	200 ng/mL
Cocaine	300 ng/mL
Amphetamines	1000 ng/mL
THC	50 ng/mL
Opiates	300 ng/mL
Barbiturates	200 ng/mL

The ingestion of certain herbal or plant products containing Ephedra or its metabolites may cause false positive amphetamine/metamphetamine results.

This test is for medical screening purposes ONLY. For confirmation a separate order for Gas Chromatography by Mass Spectrophotometry (GCMS) is required.

12. CLINICAL SIGNIFICANCE

Marijuana is a mixture of dried leaves and flowering tops of the plant *Cannabis sativa L.* The agents that produce the hallucinogenic and other biological effects of marijuana are called cannabinoids. The cannabinoid $\Delta 9$ -tetrahydrocannabinol ($\Delta 9$ -THC) is the principal psychoactive ingredient in marijuana and hashish. The compound $\Delta 9$ -THC is quickly and effectively absorbed by inhalation or from the gastrointestinal tract, and is almost completely metabolized by liver enzymes. Peak plasma levels of $\Delta 9$ -THC occur within 10 minutes of inhalation and approximately 1 hour after ingestion. Excretion of urinary metabolites and excretion by way of the feces begins within 72 hours after exposure. Concentration depends on the total amount of THC absorbed, frequency of abuse, rate of release from fatty tissue, and time of specimen collection with respect to use. In chronic users, THC may accumulate in fatty tissue faster than it can be eliminated. This accumulation leads to longer detection times in urinalysis for chronic users than for occasional users.

11.2 Critical Values

None established

11.3 Standard Required Messages

The following comment is automatically added to the report by the LIS when a urine cocaine test is ordered:

This is a screening assay. Benzoylcegonine (Cocaine Metabolite) is detected in concentrations at or above 300 ng/mL.

The following comment is automatically added to the report by the LIS when a urine drug screen is ordered:

The drug of abuse panel is a screening assay. It detects the following drugs of abuse in concentrations at or above the concentrations listed below.

Phencyclidine	25 ng/mL
Benzodiazepines	200 ng/mL
Cocaine	300 ng/mL
Amphetamines	1000 ng/mL
THC	50 ng/mL
Opiates	300 ng/mL
Barbiturates	200 ng/mL

The ingestion of certain herbal or plant products containing Ephedra or its metabolites may cause false positive amphetamine/metamphetamine results.

This test is for medical screening purposes ONLY. For confirmation a separate order for Gas Chromatography by Mass Spectrophotometry (GCMS) is required.

12. CLINICAL SIGNIFICANCE

Cocaine is a central nervous system stimulant that is extracted from the coca plant. As a drug of abuse, it is self-administered in a variety of ways, including inhalation and intravenous injection. Cocaine base can be smoked in a form that is commonly known as “crack.” Cocaine is rapidly absorbed, especially when smoked. While all forms are potentially addicting, “crack” is especially likely to lead to dependence because of its more rapid and heightened effect on the abuser. Excretion rate patterns vary with the mode of administration and from individual to individual. Cocaine is almost completely metabolized, primarily in the liver, with only about one percent excreted in the urine unchanged. Most cocaine is eliminated as benzoylcegonine, the major metabolite of cocaine. Cocaine is also excreted in relatively lesser amounts as ecgonine methyl ester and ecgonine. Cocaine metabolites may be detected in urine for up to a couple of days after cocaine is used.

11.2 Critical Values

None established

11.3 Standard Required Messages

The following comment is automatically added to the report by the LIS when a urine opiates test is ordered:

This is a screening assay. Opiates are detected in concentration at or above 300 ng/mL.

The following comment is automatically added to the report by the LIS when a urine drug screen is ordered:

The drug of abuse panel is a screening assay. It detects the following drugs of abuse in concentrations at or above the concentrations listed below.

Phencyclidine	25 ng/mL
Benzodiazepines	200 ng/mL
Cocaine	300 ng/mL
Amphetamines	1000 ng/mL
THC	50 ng/mL
Opiates	300 ng/mL
Barbiturates	200 ng/mL

The ingestion of certain herbal or plant products containing Ephedra or its metabolites may cause false positive amphetamine/metamphetamine results.

This test is for medical screening purposes ONLY. For confirmation a separate order for Gas Chromatography by Mass Spectrophotometry (GCMS) is required.

12. CLINICAL SIGNIFICANCE

Opiates are a class of compounds that includes morphine, codeine, and heroin. Morphine and codeine are naturally occurring alkaloids that are found in opium, a substance exuded from the unripe seed pod of the opium poppy *Papaver somniferum*. Heroin is a semisynthetic derivative of morphine. Morphine is a potent analgesic. Codeine is used in analgesic preparations and as a cough suppressant. Heroin is an even more potent analgesic than morphine. Both morphine and codeine are legitimate drugs. Heroin is a drug of abuse that may be snorted, smoked, or dissolved and injected subcutaneously or intravenously. Opiates are absorbed rapidly. Heroin is converted almost immediately to morphine, which is excreted in urine both unchanged and as a glucuronidated metabolite. Excretion takes place over a period of a couple of days. Codeine is excreted in urine as a glucuronidated conjugate,

11.2 Critical Values

None established

11.3 Standard Required Messages

The following comment is automatically added to the report by the LIS when a urine phencyclidine test is ordered:

This is a screening assay. Phencyclidine is detected in concentrations at or above 25 ng/mL.

The following comment is automatically added to the report by the LIS when a urine drug screen is ordered:

The drug of abuse panel is a screening assay. It detects the following drugs of abuse in concentrations at or above the concentrations listed below.

Phencyclidine	25 ng/mL
Benzodiazepines	200 ng/mL
Cocaine	300 ng/mL
Amphetamines	1000 ng/mL
THC	50 ng/mL
Opiates	300 ng/mL
Barbiturates	200 ng/mL

The ingestion of certain herbal or plant products containing Ephedra or its metabolites may cause false positive amphetamine/metamphetamine results.

This test is for medical screening purposes ONLY. For confirmation a separate order for Gas Chromatography by Mass Spectrophotometry (GCMS) is required.

12. CLINICAL SIGNIFICANCE

Phencyclidine, also known as PCP and “angel dust,” is a synthetic drug that was originally developed for its anesthetic properties but is now a drug of abuse used solely for its potent hallucinogenic effects. It may be self-administered in a variety of ways, including ingestion, inhalation, and intravenous injection. Phencyclidine is absorbed well and quickly, and concentrated in the brain and fatty tissues. Excretion patterns vary widely, ranging from several hours to a couple of weeks. Phencyclidine is excreted in the urine unchanged, as conjugated metabolites, and primarily as unidentified compounds.

13. PROCEDURE NOTES

- **FDA Status:** FDA Approved/Modified