

TRAINING UPDATE

Lab Location: GEC
Department: Core

Date Distributed: 5/19/2017
Due Date: 6/5/2017
Implementation: 5/19/2017

DESCRIPTION OF PROCEDURE REVISION

Name of procedure:
Xpand QC Schedule AG.F210.2
Description of change(s):
Update Xpand 2 to reflect limited testing. Update Cardiac and Pediatric control Open Stability (OS) Note – Because the instrument changes went into effect on 5/17, this revised FORM was implemented on May 19, 2017

Document your compliance with this training update by taking the quiz in the MTS system.

XPAND 2 QC SCHEDULE

CHECK SCREEN BEFORE GOING TO NEW SAMPLE, WHEN PROGRAMMING QC

0500	0900	1300	1700	2100	0200
QC BIO1S, QC BIO2S QC CK1, QC CK2, QC CK3	QCLYTES2	QC BIO1S	QCLYTES1	QC BIO2S	QCLYTES2
		QC CK1		QC CK3	

BIORAD LIQUID UNASSAYED CHEMISTRY CONTROLS (Levels 1 and 2) - 6 days (OS)
QC BIO1S—8 drops ALB, ALP, ALT, AST, BUN, CA, CK, CREA, ECO2, GLUC, LYLES, TBIL, TP, LA QC BIO2S—8 drops
BIORAD LIQUICHEK CARDIAC MARKER PLUS CONTROLS (Levels 1, 2 and 3) - 10 days (OS)
QC CK1—8 drops QC CK2—8 drops CTNI & MMB QC CK3—8 drops

REPEAT QC PROTOCOL

If you have to repeat QC, it must be entered in the following way:

Name: BII01 (or the level you are testing – MAKE SURE TO LEAVE QC FROM THE NAME)

Priority: QC

Sample No: C-BIO1S (or the level you are testing)

Fluid: Serum QC1 (or the level you are testing)

ONLY SELECT THE TESTS NEEDED

MAKE SURE YOU HAVE *NOT* PROGRAMMED THE NAME STARTING WITH “QC” SUCH AS BIO1S

Note: Repeat QC should never be retrieved.

Document: AG.F210[2] Status: PRERELEASED, Effective: 5/19/2017, Check Version Before Use

XPAND 1 QC SCHEDULE

CHECK SCREEN BEFORE GOING TO NEW SAMPLE, WHEN PROGRAMMING QC

Document: AG.F210[2] Status: PRERELEASED, Effective: 5/19/2017, Check Version Before Use

0900	1300	1700	2100	0200	0500
QC BIO1S, QC BIO2S QNBIL2 QCCSF1, QCCSF2 QCCK1, QCCK2, QCCK3 QCTD1S, QCTD2S, QCTD3S QCLIQ1, QCLIQ2, QCLIQ3 QCCRP1, QCCRP2, QCCRP3	QCLYTES2	QC BIO1S	QCLYTES1	QC BIO2S	QCLYTES2
		QCCK1		QCCK3	

BIORAD ETHANOL/AMMONIA CONTROLS (Levels 1, 2 & 3) - 20 day Open Stability (OS)	
QCLIQ1— 3 drops	ALC
QCLIQ2— 3 drops	
QCLIQ3— 3 drops	
BIORAD LIQUID UNASSAYED CHEMISTRY CONTROLS (Levels 1 and 2) - 6 days (OS)	
QC BIO1S—8 drops	ALB, ALP, ALT, AMY, AST, BUN, CA, CK, CREA, DBIL, ECO2, GLUC, LIPL, LYLES, MG, TBIL, TP, LA
QC BIO2S—8 drops	
BIORAD PEDIATRIC CONTROL (Level 2 only) - 14 days (OS)	
QCNBIL2—3 drops	TBIL AND DBIL
BIORAD LIQUICHEK CARDIAC MARKER PLUS CONTROLS (Levels 1, 2 and 3) - 10 days (OS)	
QCCK1—8 drops	CTNI & MMB
QCCK2—8 drops	
QCCK3—8 drops	
BIORAD LIQUICHEK IMMUNOASSY PLUS CONTROLS (Levels 1, 2 and 3) - 14 days (OS)	
QCTD1S—5 drops	ACTM, HCG, SAL & TSH
QCTD2S—5 drops	
QCTD3S—5 drops	
BIORAD LIQUICHEK SPINAL FLD CONT (RUN AS URINEQC1 & URINEQC2) – 30 days (OS)	
QCCSF1—3 drops	UCFP & GLUC
QCCSF2—3 drops	
BIORAD IMMUNOLOGY CONTROLS (Levels 1, 2, & 3) - 30 days (OS)	
QCCRP1—3 drops	CRP
QCCRP2—3 drops	
QCCRP3—3 drops	

REPEAT QC PROTOCOL

If you have to repeat QC, it must be entered in the following way:

Name: BII01 (or the level you are testing – MAKE SURE TO LEAVE QC FROM THE NAME)

Priority: QC

Sample No: C-BIO1S (or the level you are testing)

Fluid: Serum QC1 (or the level you are testing)

ONLY SELECT THE TESTS NEEDED

MAKE SURE YOU HAVE **NOT** PROGRAMMED THE NAME STARTING WITH “QC” SUCH AS BIO1S

Note: Repeat QC should never be retrieved.