

TRAINING UPDATE

Lab Location: SGMC & WAH
Department: Core Lab

Date Distributed: 6/6/2017
Due Date: 6/27/2017
Implementation: 6/27/2017

DESCRIPTION OF PROCEDURE REVISION

| Name of procedure: | |
|--|---|
| Cholesterol, Total by Dimension Vista® System SGAH.C78 v4 Triglycerides by Dimension Vista® System SGAH.C111 v2 | |
| Note: these have been converted to system SOPs | |
| Description of change(s): | |
| <i>(note QC change is already in effect)</i> | |
| Section | Reason |
| Header | Add WAH |
| 3.2 | Remove specimen onboard stability |
| 4,5,6 | Remove individual section labeling instructions and add general one |
| 6.1, 6.2 | Update QC material and storage |
| 7.2 | Change freezer limits to match products |
| 10.5 | Move patient review from section 6 |
| 10.6 | Add requirement to repeat on different analyzer if below CRR |
| 15 | Update to new standard wording |
| 17 | Update QC product |

These revised SOPs will be implemented on June 27, 2017

Document your compliance with this training update by taking the quiz in the MTS system.

Technical SOP

| | | |
|--------------------|--|------------------|
| Title | Cholesterol, Total by Dimension Vista® System | |
| Prepared by | Ashkan Chini | Date: 6/22/2012 |
| Owner | Robert SanLuis | Date: 10/21/2013 |

| Laboratory Approval | Local Effective Date: | |
|--|------------------------------|------|
| Print Name and Title | Signature | Date |
| <i>Refer to the electronic signature page for approval and approval dates.</i> | | |
| | | |
| | | |
| | | |
| | | |
| | | |

| Review | | |
|---------------|-----------|------|
| Print Name | Signature | Date |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |

TABLE OF CONTENTS

1. Test Information.....2
 2. Analytical Principle3
 3. Specimen Requirements.....3
 4. Reagents4
 5. Calibrators/Standards5
 6. Quality Control6
 7. Equipment And Supplies9
 8. Procedure9
 9. Calculations.....10
 10. Reporting Results And Repeat Criteria.....11
 11. Expected Values.....12
 12. Clinical Significance13
 13. Procedure Notes13
 14. Limitations Of Method13
 15. Safety14
 16. Related Documents14
 17. References.....14
 18. Revision History15
 19. Addenda16

1. TEST INFORMATION

| Assay | Method/Instrument | Local Code |
|--------------------|--------------------------|-------------------|
| Cholesterol, Total | Dimension Vista® System | CHOL |


| Synonyms/Abbreviations |
|--|
| CHOL, Included in Batteries/Packages: LIPD |

| Department |
|-------------------|
| Chemistry |

Form revised 2/02/2007

2. ANALYTICAL PRINCIPLE

Cholesterol esterase (CE) catalyzes the hydrolysis of cholesterol esters to produce free cholesterol which, along with preexisting free cholesterol, is oxidized in a reaction catalyzed by cholesterol oxidase (CO) to form cholest-4-ene-3-one and hydrogen peroxide. In the presence of horseradish peroxidase (HPO), the hydrogen peroxide thus formed is used to oxidize N, N-diethylaniline- HCl/4-aminoantipyrine (DEA-HCl/AAP) to produce a chromophore that absorbs at 540 nm. The absorbance due to oxidized DEA-HCl/AAP is directly proportional to the total cholesterol concentration and is measured using a polychromatic (540, 452, 700 nm) endpoint technique.



3. SPECIMEN REQUIREMENTS

3.1 Patient Preparation

| Component | Special Notations |
|-----------------------------------|--|
| Fasting/Special Diets | N/A |
| Specimen Collection and/or Timing | Normal procedures for collecting and storing serum and plasma may be used for samples to be analyzed by this method. |
| Special Collection Procedures | N/A |
| Other | N/A |

3.2 Specimen Type & Handling

| Criteria | |
|--------------------------------------|--|
| Type -Preferred -Other Acceptable | Plasma (Lithium Heparin) Serum |
| Collection Container | Plasma: Mint green top tube (PST) Serum: Red top tube, Serum separator tube (SST) |
| Volume - Optimum - Minimum | 1.0 mL 0.5 mL |
| Transport Container and Temperature | Collection container or Plastic vial at room temperature |

| Criteria | |
|---|---|
| Stability & Storage Requirements | Room Temperature: 8 hours |
| | Refrigerated: 5 - 7 days |
| | Frozen: 3 months |
| Timing Considerations | Serum or plasma should be physically separated from cells as soon as possible with a maximum limit of two hours from the time of collection. |
| Unacceptable Specimens & Actions to Take | Specimens that are unlabeled, improperly labeled, or those that do not meet the stated criteria are unacceptable. Request a recollection and credit the test with the appropriate LIS English text code for “test not performed” message. Examples: Quantity not sufficient-QNS; Wrong collection-UNAC. Document the request for recollection in the LIS. |
| Compromising Physical Characteristics | Lipemic Samples: Ultra-centrifugation removes lipemia. Thus, if lipid testing (CHOL, TRIG, HDL, or LDL) is requested, testing for lipids must be performed prior to ultra-centrifugation. Note: Saved aliquot must be clearly marked as ultra-centrifuged. Gross hemolysis. Reject sample and request a recollection. Credit the test with the appropriate LIS English text code explanation of HMT (Specimen markedly hemolyzed) |
| Other Considerations | Allow Red Top or SST to clot completely prior to centrifugation. |

NOTE: Labeling requirements for all reagents, calibrators and controls include: (1) Open date, (2) Substance name, (3) Lot number, (4) Date of preparation, (5) Expiration date, (6) Initials of tech, and (7) Any special storage instructions. Check all for visible signs of degradation. When placed onboard the analyzer, the instrument captures the date / time loaded and calculates and tracks the opened expiration.

4. REAGENTS

The package insert for a new lot of kits must be reviewed for any changes before the kit is used. A current Package Insert is included as a Related Document.

4.1 Reagent Summary

| Reagents | Supplier & Catalog Number |
|-------------|--|
| Cholesterol | Siemens, Flex® reagent cartridge, Cat. No. K1027 |

4.2 Reagent Preparation and Storage

| | |
|--------------------|--|
| Reagent | Cholesterol |
| Container | Reagent cartridge |
| Storage | Store at 2-8°C |
| Stability | <ul style="list-style-type: none"> • Reagent is stable until expiration date stamped on the reagent cartridges. • Sealed wells on the instrument are stable for 30 days. • Once wells 1 - 12 have been entered by the instrument, they are stable for 7 days. |
| Preparation | Hydrating, mixing and diluting are automatically performed by the instrument. |

5. CALIBRATORS/STANDARDS

5.1 Calibrators/Standards Used

| Calibrator | Supplier and Catalog Number |
|------------|---|
| CHEM 1 CAL | Siemens Dimension Vista®, Cat. No. KC110B |

5.2 Calibrator Preparation and Storage

| | |
|--------------------------|---|
| Calibrator | CHEM 1 CAL |
| Preparation | Allow CHEM 1 Calibrator to thaw and equilibrate to room temperature (22 – 28° C) for 1 hour. Before use, gently invert the calibrator vials at least 10 times to ensure that the contents are thoroughly mixed. Do not vortex. |
| Storage/Stability | <ul style="list-style-type: none"> • Store at -25 to - 15°C • Unopened calibrator is stable until expiration date stamped on the box. • Opened Calibrator: once the stopper of the vial is punctured, assigned values are stable for 7 days when stored on board the Dimension Vista System. |

5.3 Calibration Parameter

| Criteria | Special Notations |
|------------------------------------|--|
| Reference Material | CHEM 1 CAL |
| Assay Range | 50 – 600 mg/dL |
| Suggested Calibration Level | See Reagent Package Insert for lot specific assigned values in mg/dL |

| | |
|---------------------------|--|
| Frequency | <ul style="list-style-type: none"> • Every new reagent cartridge lot. • Every 90 days for any one lot • When major maintenance is performed on the analyzer. • When control data indicates a significant shift in assay. |
| Calibration Scheme | 2 levels, n = 5 |

5.4 Calibration Procedure

Auto Calibration:

1. Place the required calibrator vials in a carrier. Make sure the barcode labels are entirely visible through the slots.
2. Place the carrier in the loading area.
3. Position the carrier with the labels facing away from the user.
4. Press the **Load** button.
5. Automatic calibration requires that calibrators be on the instrument. As the time for processing approaches, the instrument reviews onboard inventory for the appropriate calibrators.

Manual Calibration:

1. Verify that calibrators and reagents are in inventory on the instrument.
2. Press **System > Method Summary > Calibration**.
3. Select a method from the sidebar menu. Press the **Order Calibration** button on the screen.
4. Verify that the information on the screen is correct. Verify that the calibrator lot is correct using the drop-down menu.
 - a. When calibrating using Vials press **OK**.
 - b. When calibrating using Cups, check the Use Cups box. This displays the rack and cup position fields. For additional cups use the positions in ascending order. Be sure to use the number of calibration levels and cups as specified in the method IFU. Scan the rack barcode and place calibrator cups in an adapter in position 1 on a rack. Press **OK** and load the rack on the instrument.
5. The status field in the calibration screen changes sequentially to Awaiting Scheduling, Preparing Calibrators and Processing.

5.5 Tolerance Limits

| IF..... | THEN..... |
|--|---|
| If result fall within assay-specific specification, and QC values are within acceptable limits, | proceed with analysis |
| If result falls outside assay-specific specification, or QC values are out of Acceptable limits, | troubleshoot the assay and/or instrument and repeat calibration |

6. QUALITY CONTROL

6.1 Controls Used

| Controls | Supplier and Catalog Number |
|--|--|
| Liquid Assayed Multiquel® Levels 1 and 3 | Bio-Rad Laboratories Cat. No. 337 and 339 |

6.2 Control Preparation and Storage

| | |
|--------------------------|---|
| Control | Liquid Assayed Multiquel® Levels 1 and 3 |
| Preparation | Allow the frozen control to stand at room temperature (18-25°C) for 30 minutes or until completely thawed. Swirl the contents gently to ensure homogeneity. (Do not use a mechanical mixer) Use immediately. After each use, promptly replace the stopper and return to 2-8°C storage. |
| Storage/Stability | <p>Frozen controls are stable until the expiration date at -20 to -50°C.</p> <p>Thawed and Unopened: When this product is stored at 2-8°C and the stopper is not punctured, it will be stable for 30 days. This product can be used for 7 days when stored on-board the Siemens Dimension Vista at 2- 8°C.</p> <p>Thawed and Opened: Once the product stopper is punctured, all analytes will be stable for 5 days when stored at 2- 8°C.</p> <p>Store away from light.</p> |

6.3 Frequency

Analyze all levels of QC material after every calibration and each day of testing (notated on the QC frequency sheets posted on the instruments).

Refer to the Dimension Vista® QC Schedule in the Laboratory policy Quality Control Program and in the Dimension Vista® Quick Reference Guide.

6.4 Tolerance Limits and Criteria for Acceptable QC

| Step | Action |
|------|--|
| 1 | Acceptable ranges for QC are programmed into the instrument’s Quality Control software system and Unity Real Time, and may be posted near the instrument for use during computer downtime. |
| 2 | <p>Run Rejection Criteria</p> <ul style="list-style-type: none"> Anytime the established parameters are exceeded (if one QC result exceeds 2 SD), the run is considered out of control (failed) and patient results must not be reported. The technologist must follow the procedure in the Laboratory QC |

| Step | Action |
|------|--|
| | Program to resolve the problem. |
| 3 | <p>Corrective Action:</p> <ul style="list-style-type: none"> • All rejected runs must be effectively addressed through corrective action. Steps taken in response to QC failures must be documented. Patient samples in failed analytical runs must be <u>reanalyzed according to the Laboratory QC Program</u>. Supervisors may override rejection of partial or complete runs only with detailed documentation and criteria for overrides that are approved by the Medical Director. Consult corrective action guidelines in Laboratory QC Program. Follow corrective action guidelines in the Laboratory QC Program. • Corrective action documentation must follow the Laboratory Quality Control Program. |
| 4 | <p>Review of QC</p> <ul style="list-style-type: none"> • QC must be reviewed weekly by the Group Lead or designee and monthly by the Supervisor/Manager or designee. • If the SD and/or CV are greater than established ranges, investigate the cause for the imprecision and document implementation of corrective actions. |

6.5 Documentation

- QC tolerance limits are programmed into the instrument and Unity Real Time; it calculates cumulative mean, SD and CV and stores all information for easy retrieval.
- Quality control records are reviewed daily at the bench, weekly by the Group Lead or designee, and monthly by the Supervisor/Manager or designee.
- Refer to complete policies and procedures for QC documentation and for record retention requirements in the Laboratory QC Program.

6.6 Quality Assurance Program

- Each new lot number of reagent or new shipment of the same lot of reagent must be tested with external control materials and previously analyzed samples. Performance of the new lot must be equivalent to the previous lot; utilize published TEA for acceptability criteria.
- Training must be successfully completed and documented prior to performing this test. This procedure must be incorporated into the departmental competency assessment program.
- The laboratory participates in CAP proficiency testing. All proficiency testing materials must be treated in the same manner as patient samples.

- Monthly QC must be presented to the Medical Director or designee for review and signature.
- Monthly QC mean and SD are sent to Bio-Rad Laboratories for peer group comparison.
- Consult the Laboratory QC Program for complete details.

7. EQUIPMENT and SUPPLIES

7.1 Assay Platform

Dimension Vista® System

7.2 Equipment

- Refrigerator capable of sustaining 2–8°C.
- Freezer capable of sustaining range not to exceed -15 to -25°C for calibrator
- Freezer capable of sustaining range not to exceed -20 to -50°C for QC product
- Centrifuge

7.3 Supplies

- Aliquot Plates
- System Fluids
- Assorted calibrated pipettes (MLA or equivalent) and disposable tips

8. PROCEDURE

CHOL Flex® reagent cartridge Cat. No. K1027 is required to perform this test.

Cholesterol is performed on the Dimension Vista® System after the method is calibrated (see Reference Material in Calibration section) and Quality Controls are acceptable.

NOTE: For all procedures involving specimens, buttoned lab coats, gloves, and face protection are required minimum personal protective equipment. Report all accidents to your supervisor.

| 8.1 | Sample Processing |
|-----|--|
| 1. | A sample rack holding tubes or cups is placed on the rack input lane. |
| 2. | The sample shuttle moves the rack to the barcode reader which identifies the rack and samples to the system. |
| 3. | The rack moves into the sample server and to the rack positioner. |
| 4. | At the same time, aliquot plates move from the aliquot loader into position. |

| 8.1 | Sample Processing |
|------------|---|
| 5. | The aliquot probe aspirates the sample from the tubes or cups and dispenses it into the wells of the aliquot plates. |
| 6. | After each aspirate-dispense action, the probe is thoroughly rinsed inside and out to prevent sample carryover. |
| 7. | When sample aspiration is completed, the sample server moves the rack back to the sample shuttle, where it is placed on the output lane and can be removed by the operator. |

| 8.2 | Specimen Testing |
|------------|--|
| 1. | For QC placement and frequency, refer to the Dimension Vista® QC Schedule in the Laboratory QC Program. |
| 2. | Follow the instructions, outlined in the Dimension Vista® Operator's Manual |
| 3. | The instrument reporting system contains error messages to warn the user of specific malfunctions. Results followed by such error messages should be held for follow-up. Refer to the Dimension Vista® system manual "Error messages" section for troubleshooting. |
| 4. | Follow protocol in Section 10.5 "Repeat criteria and resulting" for samples with results above or below the Analytical Measurement Range (AMR). Investigate any failed delta result and repeat, if necessary. |
| 5. | Append the appropriate English text code qualifier messages to any samples requiring a comment regarding sample quality and/or any other pertinent factors. |

| Test Conditions | |
|------------------------|------------------------|
| Sample Volume: | 1.25 µL |
| Reagent 1 Volume: | 36.5 µL |
| Reagent 2 Volume: | 10.8 µL |
| Reaction Time: | 5.6 minutes |
| Test Temperature: | 37° C |
| Wavelength: | 540, 542 & 700 nm |
| Type of measurement: | Polychromatic endpoint |

NOTE: In the event that the test system becomes inoperable, notify supervision or designee for further direction. Patient specimens must be stored in a manner that maintains the integrity of the specimen.

9. CALCULATIONS

Total Cholesterol (TC) = High Density Lipoprotein (HDL) + Low Density Lipoprotein (LDL) + Very Low Density Lipoprotein (VLDL)

The following calculations are performed by the LIS (Sunquest) when a Lipid Panel is performed:

$$TC = HDL + LDL + VLDL$$

$$\text{Calculated LDL} = TC - HDL - VLDL$$

$$VLDL = \text{Triglycerides} / 5$$

$$\text{Total Cholesterol / HDL ratio} = TC / HDL$$

10. REPORTING RESULTS AND REPEAT CRITERIA

10.1 Interpretation of Data

None required

10.2 Rounding

No rounding is necessary. Instrument reports results as a whole number.

10.3 Units of Measure

mg/dL

10.4 Clinically Reportable Range (CRR)

50 – 3,000 mg/dL

10.5 Review Patient Data

Each result is reviewed for error messages. Refer to the Dimension Vista system manual “Error messages” section for troubleshooting. Resolve any problems noted before issuing patient reports.

10.6 Repeat Criteria and Resulting

All repeats must replicate the original result within the total allowable error (TEa) of the assay. Refer to TEa listing for specific information.

Values that fall within the AMR or CRR may be reported without repeat. Values that fall outside these ranges must be repeated.

| IF the result is ... | THEN... |
|----------------------|--|
| < 50 mg/dL | Assure there is sufficient sample devoid of bubbles, cellular debris, and/or fibrin clots and re-run on a different analyzer. If both results are <50 mg/dL, add comment REP and release result. |

| IF the result is ... | THEN... |
|----------------------|--|
| ≥ 600 mg/dL | On Board Automated Dilution: Results ≥ 600 mg/dL will automatically have repeat testing performed into the instrument using dilution factor of 4. No multiplication is necessary. |
| > 2,400 mg/dL | Manual Dilution: Using the primary tube, make the smallest dilution possible to bring the raw data within the AMR. Maximum allowable dilution: x 5 DILUENT: WATER Enter dilution factor as a whole number. Re-assay. Readout is corrected for dilution. |
| > 3,000 mg/dL | If the recommended dilution does not give results within the clinically reportable range, report as: “> 3,000 mg/dL-REP” Bring to the attention of your supervisor prior to releasing result. |

| Message | Code |
|-----------------------------|----------------------------|
| Verified by repeat analysis | Append -REP to the result. |

11. EXPECTED VALUES

11.1 Reference Ranges

| Age | Female | Male |
|------------------------------|-------------|-------------|
| Adult (>18 years): | < 200 mg/dL | < 200 mg/dL |
| Pediatric: | | |
| 16 – 18 years | 101 - 200 | 105 - 200 |
| 14 – 15 years | 125 - 200 | 101 - 200 |
| 12 – 13 years | 120 - 200 | 122 - 200 |
| 10 – 11 years | 122 - 200 | 120 - 200 |
| 7 – 9 years | 107 - 200 | 107 - 200 |
| 4 – 6 years | 103 - 184 | 103 - 184 |
| 1 – 3 years | 37 - 178 | 37 - 178 |
| 7 – 11 months | 68 - 200 | 83 - 200 |
| 2 – 6 months | 59 - 200 | 53 - 194 |
| 0 – 1 month | 56 - 195 | 37 - 174 |

Calculated LDL: < 130 mg/dL
 VLDL: 8 – 32 mg/dL
 Chol / HDL Ratio: < 5.0

11.2 Critical Values

None established

11.3 Standard Report Messages

The following comment is automatically added to the report by the LIS when a lipid panel is ordered:

Lipid Interpretation

| RISK OF CORONARY HEART DISEASE | | |
|--------------------------------|------|-------|
| TOTAL CHOL. / HDL-CHOL. RATIO | | |
| | MEN | WOMEN |
| ½ average risk | 3.4 | 3.4 |
| average risk | 5.0 | 4.4 |
| 2 times average risk | 9.6 | 7.1 |
| 3 times average risk | 23.4 | 11.0 |

12. CLINICAL SIGNIFICANCE

Lipids and lipoproteins in circulation have been strongly associated with coronary heart disease (CHD), associated lipid metabolism disorders, and atherosclerosis, a cause of CHD.

13. PROCEDURE NOTES

- **FDA Status:** FDA Approved/cleared
- **Validated Test Modifications:** None

The instrument reporting system contains error messages to warn the operator of specific malfunctions. Any report slip containing such error messages should be held for follow-up. Refer to your Dimension Vista Operator's Guide.

The expected maximum observed standard deviations for repeatability using n = 5 replicates at the following cholesterol concentrations are:

| CHOL Concentration | Acceptable S.D. Maximum |
|--------------------|-------------------------|
| 180 mg/dL | 17 mg/dL |
| 450 mg/dL | 30 mg/dL |

14. LIMITATIONS OF METHOD

14.1 Analytical Measurement Range (AMR)

50 – 600 mg/dL

14.2 Precision

| Material | Mean mg/dL | Standard Deviation (%CV) | |
|-------------------|---------------|--------------------------|------------|
| | | Repeatability | Within-Lab |
| Multiquel Control | | | |
| Level 1 | 176 | 4 (2) | 5 (3) |
| Level 2 | 278 | 4 (2) | 7 (2) |

14.3 Interfering Substances

Potassium Oxalate/Sodium Fluoride can decrease cholesterol results an average of 12% and should not be used.

HIL Interference:

The CHOL method was evaluated for interference according to CLSI/NCCLS EP7-A2. Bias, defined as the difference between the control sample (does not contain interferent) and the test sample (contains interferent), is shown in the table below. Bias exceeding 10% is considered “interference”.

| Substance tested | Substance Concentration | CHOL mg/dL | Bias % |
|--------------------------|-------------------------|------------|--------|
| Hemoglobin (hemolysate) | 1000 mg/dL | 202 | <10 |
| Bilirubin (unconjugated) | 5 mg/dL | 202 | <10 |
| | 10 mg/dL | | -11 |
| | 20 mg/dL | | -13 |
| | 40 mg/dL | | -26 |
| Bilirubin (conjugated) | 5 mg/dL | 202 | <10 |
| | 10 mg/dL | | -12 |
| | 20 mg/dL | | -13 |
| | 40 mg/dL | | -32 |
| Lipemia Intralipid® | 1000 mg/dL | 202 | <10 |
| | 3000 mg/dL | | -11.3 |

14.4 Clinical Sensitivity/Specificity/Predictive Values

Not available

15. SAFETY

Refer to your local and corporate safety manuals and Safety Data Sheet (SDS) for detailed information on safety practices and procedures and a complete description of hazards.

16. RELATED DOCUMENTS

1. Dimension Vista® Clinical Chemistry System Operator's Manual
2. Dimension Vista® Calibration/Verification Procedure

3. Dimension Vista® Cal Accept Guidelines
4. Dimension Vista® Calibration summary
5. Dimension Vista® Sample Processing, Startup and Maintenance procedure
6. Laboratory Quality Control Program
7. QC Schedule for Siemens Dimension Vista®
8. Laboratory Safety Manual
9. Safety Data Sheets (SDS)
10. Dimension Vista® Limits Chart (AG.F200)
11. Quest Diagnostics Records Management Procedure
12. Dimension Vista® System Error Messages Chart
13. Centrifuge Use, Maintenance and Functions Checks (Lab policy)
14. Hemolysis, Icteria and Lipemia Interference (Lab policy)
15. Repeat Testing Requirement (Lab policy)
16. Current Allowable Total Error Specifications
at http://questnet1.qdx.com/Business_Groups/Medical/qc/docs/qc_bpt_tea.xls
17. Current package insert CHOL Flex® Reagent Cartridge K1027

17. REFERENCES

1. Ghoshal, Amit K. and Soldin, Steven J., Evaluation of the Dade Behring Dimension® RxL: Integrated chemistry system-pediatric reference ranges. Clinica Chimica Acta 2003; 331:144.
2. Package Insert, CHOL Flex® Reagent Cartridge K1027, Siemens Healthcare Diagnostics Inc., 07/28/2014.
3. Package Insert, CHEM I CAL, Siemens Healthcare Diagnostics Inc., 03/2015.
4. Package Insert, Liquid Assayed Multiqual® Chemistry Controls, Bio-Rad Laboratories, 09/2015.

18. REVISION HISTORY

| Version | Date | Section | Reason | Reviser | Approval |
|---------|----------|----------|--|----------------------|-----------|
| 000 | 10/21/13 | | Update owner | L Barrett | R SanLuis |
| 000 | 10/21/13 | 10, 14 | Change values to whole numbers | L Barrett | R SanLuis |
| 000 | 10/21/13 | 11.1 | Change adult range | L Barrett | R SanLuis |
| 000 | 10/21/13 | 16 | Update titles | L Barrett | R SanLuis |
| 000 | 10/21/13 | Footer | Version # leading zero's dropped due to new EDCS in use as of 10/7/13. | L Barrett | R SanLuis |
| 1 | 6/12/14 | 5.2 | Updated open calibrator stability | A Chini | R SanLuis |
| 2 | 2/17/16 | 1 | Add battery code | A Chini | R SanLuis |
| 2 | 2/17/16 | 3.2 | Specify anticoagulant, add instructions for lipemia | L Barrett A Chini | R SanLuis |
| 2 | 2/17/16 | 6.4, 6.6 | Replace LIS with Unity Real Time | L Barrett | R SanLuis |

| Version | Date | Section | Reason | Reviser | Approval |
|---------|---------|----------|---|-----------|-----------|
| 2 | 2/17/16 | 9 | Correct to state performed by LIS | A Chini | R SanLuis |
| 2 | 2/17/16 | 11.1 | Add ranges for calculated values | A Chini | R SanLuis |
| 2 | 2/17/16 | 11.3 | Add report comment for lipid panel | A Chini | R SanLuis |
| 2 | 2/17/16 | 17 | Update package insert information | A Chini | R SanLuis |
| 3 | 5/30/17 | Header | Add WAH | L Barrett | R SanLuis |
| 3 | 5/30/17 | 3.2 | Remove specimen onboard stability | L Barrett | R SanLuis |
| 3 | 5/30/17 | 4,5,6 | Remove individual section labeling instructions and add general one | L Barrett | R SanLuis |
| 3 | 5/30/17 | 6.1, 6.2 | Update QC material and storage | L Barrett | R SanLuis |
| 3 | 5/30/17 | 7.2 | Specify freezer ranges for products | L Barrett | R SanLuis |
| 3 | 5/30/17 | 10.5 | Move patient review from section 6 | L Barrett | R SanLuis |
| 3 | 5/30/17 | 10.6 | Add requirement to repeat on different analyzer if below CRR | L Barrett | R SanLuis |
| 3 | 5/30/17 | 15 | Update to new standard wording | L Barrett | R SanLuis |
| 3 | 5/30/17 | 17 | Update QC product | L Barrett | R SanLuis |

19. ADDENDA

None

Technical SOP

| | | |
|--------------------|---|-----------------|
| Title | Triglycerides by Dimension Vista® System | |
| Prepared by | Ashkan Chini | Date: 6/25/2012 |
| Owner | Robert SanLuis | Date: 2/17/2016 |

| Laboratory Approval | Local Effective Date: | |
|--|------------------------------|------|
| Print Name and Title | Signature | Date |
| <i>Refer to the electronic signature page for approval and approval dates.</i> | | |
| | | |
| | | |
| | | |
| | | |
| | | |

| Review | | |
|---------------|-----------|------|
| Print Name | Signature | Date |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |

TABLE OF CONTENTS

1. Test Information.....2

2. Analytical Principle3

3. Specimen Requirements.....3

4. Reagents4

5. Calibrators/Standards5

6. Quality Control6

7. Equipment And Supplies9

8. Procedure9

9. Calculations.....10

10. Reporting Results And Repeat Criteria.....10

11. Expected Values.....12

12. Clinical Significance13

13. Procedure Notes13

14. Limitations Of Method13

15. Safety14

16. Related Documents14

17. References.....14

18. Revision History15

19. Addenda16

1. TEST INFORMATION

| Assay | Method/Instrument | Local Code |
|---------------|--------------------------|-------------------|
| Triglycerides | Dimension Vista® System | TRIG |


| Synonyms/Abbreviations |
|---|
| TGL, TRIG, Included in Batteries/Packages: LPNL |

| Department |
|-------------------|
| Chemistry |

Form revised 2/02/2007

2. ANALYTICAL PRINCIPLE

The triglycerides method is based on an enzymatic procedure in which combinations of enzymes are employed for the measurement of serum or plasma triglycerides. The sample is incubated with lipoprotein lipase (LPL) enzyme reagent that converts triglycerides into free glycerol and fatty acids. Glycerol kinase (GK) catalyzes the phosphorylation of glycerol by adenosine-5- triphosphate (ATP) to glycerol-3-phosphate. Glycerol-3-phosphate-oxidase (GPO) oxidizes glycerol-3-phosphate to dihydroxyacetone phosphate and hydrogen peroxide (H2O2). The catalytic action of peroxidase (POD) forms quinoneimine from H2O2, aminoantipyrine and 4-chlorophenol. The change in absorbance due to the formation of quinoneimine is directly proportional to the total amount of glycerol and its precursors in the sample and is measured using a bichromatic (510, 700 nm) endpoint technique.



3. SPECIMEN REQUIREMENTS

3.1 Patient Preparation

| Component | Special Notations |
|--|--|
| Fasting/Special Diets | Fasting specimen preferred. Patient should be fasting for 12 hours before collection. |
| Specimen Collection and/or Timing | Normal procedures for collecting and storing serum and plasma may be used for samples to be analyzed by this method. |
| Special Collection Procedures | None |
| Other | Avoid blood collection tubes containing glycerol lubricated stoppers which will falsely elevate results. |

3.2 Specimen Type & Handling

| Criteria | |
|--------------------------|--------------------------|
| Type | Plasma (Lithium Heparin) |
| -Preferred | |
| -Other Acceptable | Serum |

| Criteria | |
|---|---|
| Collection Container | Plasma: Mint green top tube (PST) Serum: Red top tube, Serum separator tube (SST) |
| Volume - Optimum - Minimum | 1.0 mL 0.5 mL |
| Transport Container and Temperature | Collection container or Plastic vial at room temperature |
| Stability & Storage Requirements | Room Temperature: 8 hours |
| | Refrigerated: 7 days |
| | Frozen: 3 months |
| Timing Considerations | Serum or plasma should be physically separated from cells as soon as possible with a maximum limit of two hours from the time of collection. |
| Unacceptable Specimens & Actions to Take | Specimens that are unlabeled, improperly labeled, or those that do not meet the stated criteria are unacceptable. Request a recollection and credit the test with the appropriate LIS English text code for “test not performed” message. Examples: Quantity not sufficient-QNS; Wrong collection-UNAC. Document the request for recollection in the LIS. |
| Compromising Physical Characteristics | Lipemic Samples: Ultra-centrifugation removes lipemia. Thus, if lipid testing (CHOL, TRIG, HDL, or LDL) is requested, testing for lipids must be performed prior to ultra-centrifugation. Note: Saved aliquot must be clearly marked as ultra-centrifuged. Gross hemolysis. Reject sample and request a recollection. Credit the test with the appropriate LIS English text code explanation of HMT (Specimen markedly hemolyzed) |
| Other Considerations | Allow Red Top or SST to clot completely prior to centrifugation. |

NOTE: Labeling requirements for all reagents, calibrators and controls include: (1) Open date, (2) Substance name, (3) Lot number, (4) Date of preparation, (5) Expiration date, (6) Initials of tech, and (7) Any special storage instructions. Check all for visible signs of degradation. When placed onboard the analyzer, the instrument captures the date / time loaded and calculates and tracks the opened expiration.

4. REAGENTS

The package insert for a new lot of kits must be reviewed for any changes before the kit is used. A current Package Insert is included as a Related Document.

4.1 Reagent Summary

| Reagents | Supplier & Catalog Number |
|---------------|--|
| Triglycerides | Siemens, Flex® reagent cartridge, Cat. No. K2069 |

4.2 Reagent Preparation and Storage

| Reagent | Triglycerides |
|-------------|---|
| Container | Reagent cartridge |
| Storage | Store at 2-8°C |
| Stability | <ul style="list-style-type: none"> Stable until expiration date stamped on the cartridges. Sealed wells on the instrument are stable for 30 days. Once wells 1 - 12 have been entered by the instrument, they are stable for 7 days. |
| Preparation | All reagents are liquid and ready to use. |

5. CALIBRATORS/STANDARDS**5.1 Calibrators/Standards Used**

| Calibrator | Supplier and Catalog Number |
|------------|--|
| CHEM 2 CAL | Siemens Dimension Vista®, Cat. No. KC120 |

5.2 Calibrator Preparation and Storage

| | |
|-------------------|---|
| Calibrator | CHEM 2 CAL |
| Preparation | CHEM 2 CAL is ready for use. |
| Storage/Stability | <ul style="list-style-type: none"> Store at 2-8°C Unopened calibrator is stable until expiration date stamped on the box. Opened Calibrator: once the stopper of the vial is punctured, assigned values are stable for 24 hours when stored on board the Dimension Vista System. |

5.3 Calibration Parameter

| Criteria | Special Notations |
|-----------------------------|--|
| Reference Material | CHEM 2 CAL |
| Assay Range | 2 – 1000 mg/dL |
| Suggested Calibration Level | See Reagent Package Insert for lot specific assigned values in mg/dL |

| | |
|---------------------------|--|
| Frequency | <ul style="list-style-type: none"> • Every new reagent cartridge lot. • Every 90 days for any one lot • When major maintenance is performed on the analyzer. • When control data indicates a significant shift in assay. |
| Calibration Scheme | 2 levels, n = 5 |

5.4 Calibration Procedure

Auto Calibration:

1. Place the required calibrator vials in a carrier. Make sure the barcode labels are entirely visible through the slots.
2. Place the carrier in the loading area.
3. Position the carrier with the labels facing away from the user.
4. Press the **Load** button.
5. Automatic calibration requires that calibrators be on the instrument. As the time for processing approaches, the instrument reviews onboard inventory for the appropriate calibrators.

Manual Calibration:

1. Verify that calibrators and reagents are in inventory on the instrument.
2. Press **System > Method Summary > Calibration**.
3. Select a method from the sidebar menu. Press the **Order Calibration** button on the screen.
4. Verify that the information on the screen is correct. Verify that the calibrator lot is correct using the drop-down menu.
 - a. When calibrating using Vials press **OK**.
 - b. When calibrating using Cups, check the Use Cups box. This displays the rack and cup position fields. For additional cups use the positions in ascending order. Be sure to use the number of calibration levels and cups as specified in the method IFU. Scan the rack barcode and place calibrator cups in an adapter in position 1 on a rack. Press **OK** and load the rack on the instrument.
5. The status field in the calibration screen changes sequentially to Awaiting Scheduling, Preparing Calibrators and Processing.

5.5 Tolerance Limits

| IF..... | THEN..... |
|--|---|
| If result fall within assay-specific specification, and QC values are within acceptable limits, | proceed with analysis |
| If result falls outside assay-specific specification, or QC values are out of Acceptable limits, | troubleshoot the assay and/or instrument and repeat calibration |

6. QUALITY CONTROL

6.1 Controls Used

| Controls | Supplier and Catalog Number |
|--|--|
| Liquid Assayed Multiquel® Levels 1 and 3 | Bio-Rad Laboratories Cat. No. 337 and 339 |

6.2 Control Preparation and Storage

| Control | Liquid Assayed Multiquel® Levels 1 and 3 |
|--------------------------|---|
| Preparation | Allow the frozen control to stand at room temperature (18-25°C) for 30 minutes or until completely thawed. Swirl the contents gently to ensure homogeneity. (Do not use a mechanical mixer) Use immediately. After each use, promptly replace the stopper and return to 2-8°C storage. |
| Storage/Stability | <p>Frozen controls are stable until the expiration date at -20 to -50°C.</p> <p>Thawed and Unopened: When this product is stored at 2-8°C and the stopper is not punctured, it will be stable for 7 days for triglycerides.</p> <p>This product can be used for 7 days when stored on-board the Siemens Dimension Vista at 2- 8°C.</p> <p>Thawed and Opened: Once the product stopper is punctured, all analytes will be stable for 5 days when stored at 2- 8°C.</p> <p>Store away from light.</p> |

6.3 Frequency

Analyze all levels of QC material after every calibration and each day of testing (notated on the QC frequency sheets posted on the instruments).

Refer to the Dimension Vista® QC Schedule in the Laboratory policy Quality Control Program and in the Dimension Vista® Quick Reference Guide.

6.4 Tolerance Limits and Criteria for Acceptable QC

| Step | Action |
|------|---|
| 1 | Acceptable ranges for QC are programmed into the instrument’s Quality Control software system and Unity Real Time, and may be posted near the instrument for use during computer downtime. |
| 2 | <p>Run Rejection Criteria</p> <ul style="list-style-type: none"> Anytime the established parameters are exceeded (if one QC result exceeds 2 SD), the run is considered out of control (failed) and patient results must not be reported. |

| Step | Action |
|------|--|
| | <ul style="list-style-type: none"> The technologist must follow the procedure in the Laboratory QC Program to resolve the problem. |
| 3 | <p>Corrective Action:</p> <ul style="list-style-type: none"> All rejected runs must be effectively addressed through corrective action. Steps taken in response to QC failures must be documented. Patient samples in failed analytical runs must be <u>reanalyzed according to the Laboratory QC Program</u>. Supervisors may override rejection of partial or complete runs only with detailed documentation and criteria for overrides that are approved by the Medical Director. Consult corrective action guidelines in Laboratory QC Program. Follow corrective action guidelines in the Laboratory QC Program. Corrective action documentation must follow the Laboratory Quality Control Program. |
| 4 | <p>Review of QC</p> <ul style="list-style-type: none"> QC must be reviewed weekly by the Group Lead or designee and monthly by the Supervisor/Manager or designee. If the SD and/or CV are greater than established ranges, investigate the cause for the imprecision and document implementation of corrective actions. |

6.5 Documentation

- QC tolerance limits are programmed into the instrument and Unity Real Time; it calculates cumulative mean, SD and CV and stores all information for easy retrieval.
- Quality control records are reviewed daily at the bench, weekly by the Group Lead or designee, and monthly by the Supervisor/Manager or designee.
- Refer to complete policies and procedures for QC documentation and for record retention requirements in the Laboratory QC Program.

6.6 Quality Assurance Program

- Each new lot number of reagent or new shipment of the same lot of reagent must be tested with external control materials and previously analyzed samples. Performance of the new lot must be equivalent to the previous lot; utilize published TEA for acceptability criteria.
- Training must be successfully completed and documented prior to performing this test. This procedure must be incorporated into the departmental competency assessment program.
- The laboratory participates in CAP proficiency testing. All proficiency testing materials must be treated in the same manner as patient samples.

- Monthly QC must be presented to the Medical Director or designee for review and signature.
- Monthly QC mean and SD are sent to Bio-Rad Laboratories for peer group comparison.
- Consult the Laboratory QC Program for complete details.

7. EQUIPMENT and SUPPLIES

7.1 Assay Platform

Dimension Vista® System

7.2 Equipment

- Refrigerator capable of sustaining 2–8°C.
- Freezer capable of sustaining range not to exceed -20 to **-50°C.**
- Centrifuge

7.3 Supplies

- Aliquot Plates
- System Fluids
- Assorted calibrated pipettes (MLA or equivalent) and disposable tips

8. PROCEDURE

TRIG Flex® reagent cartridge Cat. No. K2069 is required to perform this test.

Triglycerides is performed on the Dimension Vista® System after the method is calibrated (see Reference Material in Calibration section) and Quality Controls are acceptable.

NOTE: For all procedures involving specimens, buttoned lab coats, gloves, and face protection are required minimum personal protective equipment. Report all accidents to your supervisor.

| 8.1 | Sample Processing |
|------------|--|
| 1. | A sample rack holding tubes or cups is placed on the rack input lane. |
| 2. | The sample shuttle moves the rack to the barcode reader which identifies the rack and samples to the system. |
| 3. | The rack moves into the sample server and to the rack positioner. |
| 4. | At the same time, aliquot plates move from the aliquot loader into position. |
| 5. | The aliquot probe aspirates the sample from the tubes or cups and dispenses it into the wells of the aliquot plates. |

| 8.1 | Sample Processing |
|------------|---|
| 6. | After each aspirate-dispense action, the probe is thoroughly rinsed inside and out to prevent sample carryover. |
| 7. | When sample aspiration is completed, the sample server moves the rack back to the sample shuttle, where it is placed on the output lane and can be removed by the operator. |

| 8.2 | Specimen Testing |
|------------|--|
| 1. | For QC placement and frequency, refer to the Dimension Vista® QC Schedule in the Laboratory QC Program. |
| 2. | Follow the instructions, outlined in the Dimension Vista® Operator's Manual |
| 3. | The instrument reporting system contains error messages to warn the user of specific malfunctions. Results followed by such error messages should be held for follow-up. Refer to the Dimension Vista® system manual "Error messages" section for troubleshooting. |
| 4. | Follow protocol in Section 10.5 "Repeat criteria and resulting" for samples with results above or below the Analytical Measurement Range (AMR). Investigate any failed delta result and repeat, if necessary. |
| 5. | Append the appropriate English text code qualifier messages to any samples requiring a comment regarding sample quality and/or any other pertinent factors. |

| Test Conditions | |
|------------------------|----------------------|
| Sample Volume: | 1.6 µL |
| Reagent Volume: | 55 µL |
| Reaction Time: | 5.6 minutes |
| Test Temperature: | 37° C |
| Wavelength: | 510 & 700 nm |
| Type of measurement: | Bichromatic endpoint |

NOTE: In the event that the test system becomes inoperable, notify supervision or designee for further direction. Patient specimens must be stored in a manner that maintains the integrity of the specimen.

9. CALCULATIONS

The instrument automatically calculates the concentration of Triglycerides in mg/dL.

10. REPORTING RESULTS AND REPEAT CRITERIA

10.1 Interpretation of Data

None required

10.2 Rounding

No rounding is necessary. Instrument reports results as a whole number.

10.3 Units of Measure

mg/dL

10.4 Clinically Reportable Range (CRR)

2 – 5000 mg/dL

10.5 Review Patient Data

Each result is reviewed for error messages. Refer to the Dimension Vista system manual “Error messages” section for troubleshooting. Resolve any problems noted before issuing patient reports.

10.6 Repeat Criteria and Resulting

All repeats must replicate the original result within the total allowable error (TEa) of the assay. Refer to TEa listing for specific information.

Values that fall within the AMR or CRR may be reported without repeat. Values that fall outside these ranges must be repeated.

| IF the result is ... | THEN... |
|----------------------|--|
| < 2 mg/dL | Assure there is sufficient sample devoid of bubbles, cellular debris, and/or fibrin clots and re-run on a different analyzer. If both results are <2 mg/dL, add comment REP and release result. |
| ≥ 1000 mg/dL | On Board Automated Dilution: Results ≥ 1000 mg/dL will automatically have repeat testing performed into the instrument using dilution factor of 4. No multiplication is necessary. |
| > 4,000 mg/dL | Manual Dilution: Using the primary tube, make the smallest dilution possible to bring the raw data within the AMR. Maximum allowable dilution: x 5 DILUENT: Water Enter dilution factor as a whole number. Re-assay. Readout is corrected for dilution. |
| > 5,000 mg/dL | If the recommended dilution does not give results within the clinically reportable range, report as: “> 5,000 mg/dL-REP” Bring to the attention of your supervisor prior to releasing result. |

Form revised 2/02/2007

| Message | Code |
|-----------------------------|----------------------------|
| Verified by repeat analysis | Append –REP to the result. |

11. EXPECTED VALUES

11.1 Reference Ranges

| Age | Female | Male |
|------------------------------|---------------|---------------|
| Adult (>18 years): | 0 – 149 mg/dL | 0 – 149 mg/dL |
| Pediatric: | | |
| 16 – 18 years | 35 - 134 | 32 - 134 |
| 14 – 15 years | 36 - 129 | 32 - 158 |
| 12 – 13 years | 35 - 124 | 22 - 138 |
| 10 – 11 years | 37 - 134 | 22 - 131 |
| 7 – 9 years | 26 - 123 | 26 - 123 |
| 4 – 6 years | 30 - 110 | 30 - 110 |
| 1 – 3 years | 25 - 119 | 25 - 119 |
| 1 – 11 months | 34 - 340 | 42 - 279 |
| 8 – 30 days | 33 - 270 | 37 - 279 |
| 0 – 7 days | 26 - 159 | 19 - 174 |

11.2 Critical Values

None established

11.3 Standard Required Messages

The following comment is automatically added to the report by the LIS when a lipid panel is ordered:

Lipid Interpretation

| RISK OF CORONARY HEART DISEASE | | |
|--------------------------------|------|-------|
| TOTAL CHOL. / HDL-CHOL. RATIO | | |
| | MEN | WOMEN |
| ½ average risk | 3.4 | 3.4 |
| average risk | 5.0 | 4.4 |
| 2 times average risk | 9.6 | 7.1 |
| 3 times average risk | 23.4 | 11.0 |

12. CLINICAL SIGNIFICANCE

Triglycerides are water-insoluble lipids consisting of three fatty acids linked to one glycerol molecule. Triglycerides are transported in the blood as core constituents of all lipoproteins, but the greatest concentration of these molecules is carried in the triglycerides-rich chylomicrons and very low density lipoproteins (VLDL). Through the action of lipases and bile acids, triglycerides are hydrolyzed into glycerol and fatty acids which are absorbed by adipose tissue for storage or by other tissues requiring a source of energy. A peak concentration of chylomicron-associated triglycerides occurs within 3–6 hours after ingestion of a fat-rich meal; however, the rate of absorption of fats is highly variable, depending on the individual and dietary composition of the fat. After absorption, triglycerides are resynthesized in the epithelial cells and combined with cholesterol and a number of apolipoproteins to form chylomicrons.

Triglycerides measurements obtained are used in the diagnosis and treatment of patients with diabetes mellitus, nephrosis, liver obstruction, other diseases involving lipid metabolism, or various endocrine disorders.

13. PROCEDURE NOTES

- **FDA Status:** FDA Approved/cleared
- **Validated Test Modifications:** None

The instrument reporting system contains error messages to warn the operator of specific malfunctions. Any report slip containing such error messages should be held for follow-up. Refer to your Dimension Vista Operator's Guide.

The expected maximum observed standard deviations for repeatability using n = 5 replicates at the following Triglycerides concentrations are:

| TRIG Concentration | Acceptable S.D. Maximum |
|---------------------------|--------------------------------|
| 70 mg/dL | 9 mg/dL |
| 375 mg/dL | 26 mg/dL |

14. LIMITATIONS OF METHOD**14.1 Analytical Measurement Range (AMR)**

2 – 1000 mg/dL

14.2 Precision

| Material | Mean mg/dL | Standard Deviation (%CV) | |
|-------------------|-----------------------|---------------------------------|-------------------|
| | | Repeatability | Within-Lab |
| Multiquel Control | | | |
| Level 1 | 68 | 2 (3) | 3 (4) |
| Level 2 | 384 | 6 (2) | 9 (2) |

14.3 Interfering Substances

Ascorbic acid at a concentration of 5 mg/dL decreases triglycerides results by 11.8% at triglyceride concentration of 180 mg/dL.

HIL Interference:

The TRIG method was evaluated for interference according to CLSI/NCCLS EP7-A2. Bias, defined as the difference between the control sample (does not contain interferent) and the test sample (contains interferent), is shown in the table below. Bias exceeding 10% is considered “interference”.

| Substance tested | Substance Concentration | TRIG mg/dL | Bias % |
|--------------------------|-------------------------|------------|--------|
| Hemoglobin (hemolysate) | 1000 mg/dL | 180 | <10 |
| Bilirubin (unconjugated) | 5 mg/dL | 180 | <10 |
| | 10 mg/dL | | 11 |
| | 20 mg/dL | | 20 |
| | 60 mg/dL | | 24 |
| Bilirubin (conjugated) | 60 mg/dL | 180 | <10 |

14.4 Clinical Sensitivity/Specificity/Predictive Values

Not available

15. SAFETY

Refer to your local and corporate safety manuals and Safety Data Sheet (SDS) for detailed information on safety practices and procedures and a complete description of hazards.

16. RELATED DOCUMENTS

1. Dimension Vista® Clinical Chemistry System Operator’s Manual
2. Dimension Vista® Calibration/Verification Procedure
3. Dimension Vista® Cal Accept Guidelines
4. Dimension Vista® Calibration summary
5. Dimension Vista® Sample Processing, Startup and Maintenance procedure
6. Laboratory Quality Control Program
7. QC Schedule for Siemens Dimension Vista®
8. Laboratory Safety Manual
9. Safety Data Sheets (SDS)
10. Dimension Vista® Limits Chart (AG.F200)
11. Quest Diagnostics Records Management Procedure
12. Dimension Vista® System Error Messages Chart
13. Centrifuge Use, Maintenance and Functions Checks (Lab policy)
14. Hemolysis, Icteria and Lipemia Interference (Lab policy)
15. Repeat Testing Requirement (Lab policy)
16. Current Allowable Total Error Specifications
at http://questnet1.qdx.com/Business_Groups/Medical/qc/docs/qc_bpt_tea.xls
17. Current package insert TRIG Flex® Reagent Cartridge K2069

17. REFERENCES

1. Ghoshal, Amit K. and Soldin, Steven J., Evaluation of the Dade Behring Dimension® RxL: Integrated chemistry system-pediatric reference ranges. Clinica Chimica Acta 2003; 331:144.
2. Package Insert, TRIG Flex® Reagent Cartridge K2069, Siemens Healthcare Diagnostics Inc., 10/20/2016.
3. Package Insert, CHEM 2 CAL, Siemens Healthcare Diagnostics Inc., 03/18/2015.
4. Package Insert, Liquid Assayed Multiqual® Chemistry Controls, Bio-Rad Laboratories, 09/2015.

18. REVISION HISTORY

| Version | Date | Section | Reason | Reviser | Approval |
|---------|---------|----------|--|----------------------|-----------|
| 000 | 2/17/16 | | Update owner | L Barrett | R SanLuis |
| 000 | 2/17/16 | 1 | Remove outdated battery codes | L Barrett | R SanLuis |
| 000 | 2/17/16 | 3.2 | Specify anticoagulant, add instructions for lipemia | L Barrett A Chini | R SanLuis |
| 000 | 2/17/16 | 4.2 | Add hazard/chemical information | A Chini | R SanLuis |
| 000 | 2/17/16 | 5.2 | Update stability | A Chini | R SanLuis |
| 000 | 2/17/16 | 6.4, 6.6 | Replace LIS with Unity Real Time | L Barrett | R SanLuis |
| 000 | 2/17/16 | 11.3 | Add report comment for lipid panel | A Chini | R SanLuis |
| 000 | 2/17/16 | 16 | Update titles | L Barrett | R SanLuis |
| 000 | 2/17/16 | 17 | Update package insert information | A Chini | R SanLuis |
| 000 | 2/17/16 | Footer | Version # leading zero's dropped due to new EDCS in use as of 10/7/13. | L Barrett | R SanLuis |
| 1 | 5/30/17 | Header | Add WAH | L Barrett | R SanLuis |
| 1 | 5/30/17 | 3.2 | Remove specimen onboard stability | L Barrett | R SanLuis |
| 1 | 5/30/17 | 4,5,6 | Remove individual section labeling instructions and add general one | L Barrett | R SanLuis |
| 1 | 5/30/17 | 6.1, 6.2 | Update QC material and storage | L Barrett | R SanLuis |
| 1 | 5/30/17 | 7.2 | Change freezer upper limit to -50C | L Barrett | R SanLuis |
| 1 | 5/30/17 | 10.5 | Move patient review from section 6 | L Barrett | R SanLuis |
| 1 | 5/30/17 | 10.6 | Add requirement to repeat on different analyzer if below CRR | L Barrett | R SanLuis |
| 1 | 5/30/17 | 15 | Update to new standard wording | L Barrett | R SanLuis |
| 1 | 5/30/17 | 17 | Update QC product and PI dates | L Barrett | R SanLuis |

19. ADDENDA

None