

TRAINING UPDATE

Lab Location: SGMC & WAH
Department: Core Lab

Date Distributed: 8/3/2018
Due Date: 8/31/2018
Implementation: 8/27/2018

DESCRIPTION OF PROCEDURE REVISION

Name of procedure:	
Bilirubin by Dimension Vista® System SGAH.C106 v2 <i>This has been converted to a system SOP</i>	
Description of change(s):	
<i>Revise QC info and add exposure to light comments, Most other changes are format updates</i>	
Section	Reason
Header	Add WAH
3.1	Add protect from light
4,5,6	Remove individual section labeling instructions and add general one
5.2	Add light exposure comment
6.1, 6.2	Update QC material and storage
10.5	Move patient review from section 6
10.6	Remove repeat value below AMR/CRR
15	Update to new standard wording
16	Update policy title
17	Update QC product and PI dates
This revised SOP will be implemented on August 27, 2018	

Document your compliance with this training update by taking the quiz in the MTS system.

Technical SOP

Title	Bilirubin by Dimension Vista® System	
Prepared by	Ashkan Chini	Date: 6/25/2012
Owner	Robert SanLuis	Date: 4/12/2016

Laboratory Approval		Local Effective Date:
Print Name and Title	Signature	Date
<i>Refer to the electronic signature page for approval and approval dates.</i>		

Review		
Print Name	Signature	Date

JUN16C016 (REV) 04/12/2016

TABLE OF CONTENTS

1. Test Information.....	2
2. Analytical Principle.....	3
3. Specimen Requirements.....	3
4. Reagents.....	4
5. Calibrators/Standards.....	5
6. Quality Control.....	7
7. Equipment And Supplies.....	9
8. Procedure.....	10
9. Calculations.....	11
10. Reporting Results And Repeat Criteria.....	11
11. Expected Values.....	13
12. Clinical Significance.....	14
13. Procedure Notes.....	14
14. Limitations Of Method.....	14
15. Safety.....	15
16. Related Documents.....	15
17. References.....	15
18. Revision History.....	16
19. Addenda.....	17

1. TEST INFORMATION

Assay	Method/Instrument	Local Code
Bilirubin, Direct	Dimension Vista® System	DBIL
Bilirubin, Indirect		(Calculated value)
Bilirubin, Total		TBIL
Bilirubin, Neonatal		TBILN, DBILN
Bilirubin, Cord		CBIL

Synonyms/Abbreviations
Bilirubin Direct and Total are included in Batteries/Packages: COMP, L1VP Bilirubin Neonatal is included in Batteries/Packages: NBIL

Department
Chemistry

JUN16C016 (REV) 04/12/2016

2. ANALYTICAL PRINCIPLE

2.1 Total Bilirubin

Diazotized sulfanilic acid is formed by combining sodium nitrite and sulfanilic acid at low pH. Bilirubin (unconjugated) in the sample is solubilized by dilution in a mixture of caffeine/benzoate/acetate/EDTA. Upon addition of the diazotized sulfanilic acid, the solubilized bilirubin including conjugated bilirubins (mono and digluconides) and the delta form (biliprotein-bilirubin covalently bound to albumin) is converted to diazo-bilirubin, a red chromophore representing the total bilirubin which absorbs at 540 nm and is measured using a bichromatic (540, 700 nm) endpoint technique. A sample blank correction is used.

Solubilized bilirubin + Diazotized sulfanilic acid -----> Red chromophore (absorbs at 540 nm)

2.2 Direct Bilirubin

Diazotized sulfanilic acid is formed by combining sodium nitrite and sulfanilic acid at low pH. The sample is diluted in 0.5M HCl. A sample blank reading is taken to eliminate interference from non-bilirubin pigments. Upon addition of the diazotized sulfanilic acid, the conjugated bilirubin is converted to diazo-bilirubin, a red chromophore which absorbs at 540 nm and is measured using a bichromatic (540, 700 nm) endpoint technique.

Conjugated bilirubin + Diazotized sulfanilic acid -----> Red chromophore (absorbs at 540 nm)

3. SPECIMEN REQUIREMENTS

3.1 Patient Preparation

Component	Special Notations
Fasting/Special Diets	N/A
Specimen Collection and/or Timing	Normal procedures for collecting and storing serum and plasma may be used for samples to be analyzed by this method.
Special Collection Procedures	Bilirubin is extremely photosensitive. Care should be taken to protect sample from both daylight and fluorescent light to avoid photodegradation
Other	N/A

3.2 Specimen Type & Handling

Criteria	
Type	Plasma (Lithium Heparin)
-Preferred	Serum, Cord Blood, Plasma (EDTA)
-Other Acceptable	

Criteria	
Collection Container	Plasma: Mint green top tube (PST) Serum: Red top tube, Serum separator tube (SST)
Volume - Optimum	1.0 mL
- Minimum	0.5 mL
Transport Container and Temperature	Plastic vial or spun barrier tube at room temperature, protect from light.
Storage Requirements	Protect from light
Stability Requirements	Room Temperature: 8 hours
	Refrigerated: 5 days
	Frozen: 3 months (Direct Bilirubin) 6 months (Total Bilirubin)
	Instrument on board aliquot stability (T-BIL ONLY) 2 hours (T-BIL ONLY) Not Applicable for D-BIL, for D-BIL use a fresh sample.
Timing Considerations	Serum or plasma should be physically separated from cells as soon as possible with a maximum limit of two hours from the time of collection.
Unacceptable Specimens & Actions to Take	Specimens that are unlabeled, improperly labeled, or those that do not meet the stated criteria are unacceptable. Request a recollection and credit the test with the appropriate LIS English text code for "test not performed" message. Examples: Quantity not sufficient-QNS; Wrong collection-UNAC. Document the request for recollection in the LIS.
Compromising Physical Characteristics	Gross hemolysis. Reject sample and request a recollection. Credit the test with the appropriate LIS English text code explanation of HMT (Specimen markedly hemolyzed)
Other Considerations	Allow Red Top or SST to clot completely prior to centrifugation.

NOTE: Labeling requirements for all reagents, calibrators and controls include: (1) Open date, (2) Substance name, (3) Lot number, (4) Date of preparation, (5) Expiration date, (6) Initials of tech, and (7) Any special storage instructions. Check all for visible signs of degradation. When placed onboard the analyzer, the instrument captures the date / time loaded and calculates and tracks the opened expiration.

4. REAGENTS

The package insert for a new lot of kits must be reviewed for any changes before the kit is used. A current Package Insert is included as a Related Document.

4.1 Reagent Summary

Reagents	Supplier & Catalog Number
Total Bilirubin	Siemens, Flex® reagent cartridge, Cat. No. K1167
Direct Bilirubin	Siemens, Flex® reagent cartridge, Cat. No. K2125

4.2 Reagent Preparation and Storage

Reagent	Total Bilirubin
Container	Reagent cartridge
Storage	Store at 2-8°C
Stability	<ul style="list-style-type: none"> Stable until expiration date stamped on reagent cartridges. Sealed wells on the instrument are stable for 30 days. Open wells: 5 days for wells 1 - 8 Diazotized sulfanilic acid in wells 9 - 10 prepared by the automatic addition of sodium nitrite from wells 11 - 12 is stable for 5 days.
Preparation	All reagents are liquid and ready to use.

Reagent	Direct Bilirubin
Container	Reagent cartridge
Storage	Store at 2-8°C
Stability	<ul style="list-style-type: none"> Stable until expiration date stamped on reagent cartridges. Sealed wells on the instrument are stable for 30 days. Open well stability: <ul style="list-style-type: none"> 2 days for wells 1 - 4 30 days for wells 5 - 6 8 days for wells 7 - 8 Sulfanilic acid in wells 9 – 12 is used immediately to prepare the diazo reagent in wells 1 – 4.
Preparation	All reagents are liquid and ready to use. Diazotized sulfanilic acid is prepared automatically by the instrument in wells 1 - 4 with the addition of sodium nitrite from wells 7 - 8 and sulfanilic acid/HCL from wells 9 - 12.

5. CALIBRATORS/STANDARDS

5.1 Calibrators/Standards Used

Calibrator	Supplier and Catalog Number
BILI CAL	Siemens Dimension Vista®, Cat. No. KC210

5.2 Calibrator Preparation and Storage

Calibrator	BILI CAL
Preparation	Calibrator is ready for use. No preparation is required.

Storage/Stability	<ul style="list-style-type: none"> Store at 2 - 8°C Unopened Calibrator: until expiration date on the box. Opened Calibrator: once the stopper is punctured, assigned values are stable for 14 days when stored on board the Dimension Vista System. Exposure to light may cause significant decrease in bilirubin concentration
--------------------------	---

5.3 Calibration Parameter

Criteria	Special Notations
Reference Material	BILI CAL
Assay Range	Total Bilirubin: 0.1 – 25.0 mg/dL Direct Bilirubin: 0.1 – 16.0 mg/dL
Suggested Calibration Level	See Reagent Package Insert for lot specific assigned values in mg/dL
Frequency	<ul style="list-style-type: none"> Every new reagent cartridge lot. Every 90 days for any one lot When major maintenance is performed on the analyzer. When control data indicates a significant shift in assay.
Calibration Scheme	2 levels, n = 5

5.4 Calibration Procedure

Auto Calibration:

- Place the required calibrator vials in a carrier. Make sure the barcode labels are entirely visible through the slots.
- Place the carrier in the loading area.
- Position the carrier with the labels facing away from the user.
- Press the **Load** button.
- Automatic calibration requires that calibrators be on the instrument. As the time for processing approaches, the instrument reviews onboard inventory for the appropriate calibrators.

Manual Calibration:

- Verify that calibrators and reagents are in inventory on the instrument.
- Press **System > Method Summary > Calibration**.
- Select a method from the sidebar menu. Press the **Order Calibration** button on the screen.
- Verify that the information on the screen is correct. Verify that the calibrator lot is correct using the drop-down menu.
 - When calibrating using Vials press **OK**.
 - When calibrating using Cups, check the Use Cups box. This displays the rack and cup position fields. For additional cups use the positions in ascending order. Be sure to use the number of calibration levels and cups as specified in

- the method IFU. Scan the rack barcode and place calibrator cups in an adapter in position 1 on a rack. Press **OK** and load the rack on the instrument.
- The status field in the calibration screen changes sequentially to Awaiting Scheduling, Preparing Calibrators and Processing.

5.5 Tolerance Limits

IF.....	THEN.....
If result fall within assay-specific specification, and QC values are within acceptable limits,	proceed with analysis
If result falls outside assay-specific specification, or QC values are out of Acceptable limits,	troubleshoot the assay and/or instrument and repeat calibration

6. QUALITY CONTROL

6.1 Controls Used

Controls	Supplier and Catalog Number
Liquid Assayed Multiqua® Levels 1 and 3	Bio-Rad Laboratories Cat. No. 337 and 339
Liquichek™ Pediatric Control Level 2	Bio-Rad Laboratories Cat. No. 177

6.2 Control Preparation and Storage

Control	Liquid Assayed Multiqua, Levels 1 and 3
Preparation	Allow the frozen control to stand at room temperature (18-25°C) for 30 minutes or until completely thawed. Swirl the contents gently to ensure homogeneity. (Do not use a mechanical mixer) Use immediately. After each use, promptly replace the stopper and return to 2-8°C storage.
Storage/Stability	Frozen: stable until the expiration date at -20 to -50°C. Thawed and Unopened: When stored at 2-8°C and the stopper is not punctured, it will be stable for 11 days for DBIL This product can be used for 7 days when stored on-board the Siemens Dimension Vista at 2-8°C. Thawed and Opened: Once the stopper is punctured, all analytes will be stable for 5 days when stored at 2-8°C. Store away from light.

Control	Liquichek Pediatric Control, Level 2
Preparation	Allow the frozen control to stand at room temperature (18-25°C) until completely thawed. Swirl the contents gently to ensure homogeneity. (Do not use a mechanical mixer) Load vials on the instrument immediately.

Liquichek Pediatric Control

Storage/Stability	Frozen: stable until the expiration date at -20 to -50°C. Thawed and Unopened: When stored at 2-8°C and the stopper is not punctured, it will be stable for 45 days. Thawed and Opened: Once the stopper is punctured, all analytes will be stable for 14 days when stored at 2-8°C.
--------------------------	---

6.3 Frequency

Analyze all levels of QC material after every calibration and each day of testing (notated on the QC frequency sheets posted on the instruments).

Refer to the Dimension Vista® QC Schedule in the Laboratory policy Quality Control Program and in the Dimension Vista® Quick Reference Guide.

6.4 Tolerance Limits and Criteria for Acceptable QC

Step	Action
1	Acceptable ranges for QC are programmed into the instrument's Quality Control software system and Unity Real Time, and may be posted near the instrument for use during computer downtime.
2	Run Rejection Criteria <ul style="list-style-type: none"> Anytime the established parameters are exceeded (if one QC result exceeds 2 SD), the run is considered out of control (failed) and patient results must not be reported. The technologist must follow the procedure in the Laboratory QC Program to resolve the problem.
3	Corrective Action: <ul style="list-style-type: none"> All rejected runs must be effectively addressed through corrective action. Steps taken in response to QC failures must be documented. Patient samples in failed analytical runs must be <u>reanalyzed according to the Laboratory QC Program</u>. Supervisors may override rejection of partial or complete runs only with detailed documentation and criteria for overrides that are approved by the Medical Director. Consult corrective action guidelines in Laboratory QC Program. Follow corrective action guidelines in the Laboratory QC Program. Corrective action documentation must follow the Laboratory Quality Control Program.

Liquichek Pediatric Control

Step	Action
4	<p>Review of QC</p> <ul style="list-style-type: none"> QC must be reviewed weekly by the Group Lead or designee and monthly by the Supervisor/Manager or designee. If the SD and/or CV are greater than established ranges, investigate the cause for the imprecision and document implementation of corrective actions.

6.5 Documentation

- QC tolerance limits are programmed into the instrument and Unity Real Time; it calculates cumulative mean, SD and CV and stores all information for easy retrieval.
- Quality control records are reviewed daily at the bench, weekly by the Group Lead or designee, and monthly by the Supervisor/Manager or designee.
- Refer to complete policies and procedures for QC documentation and for record retention requirements in the Laboratory QC Program.

6.6 Quality Assurance Program

- Each new lot number of reagent or new shipment of the same lot of reagent must be tested with external control materials and previously analyzed samples. Performance of the new lot must be equivalent to the previous lot; utilize published TEA for acceptability criteria.
- Training must be successfully completed and documented prior to performing this test. This procedure must be incorporated into the departmental competency assessment program.
- The laboratory participates in CAP proficiency testing. All proficiency testing materials must be treated in the same manner as patient samples.
- Monthly QC must be presented to the Medical Director or designee for review and signature.
- Monthly QC mean and SD are sent to Bio-Rad Laboratories for peer group comparison.
- Consult the Laboratory QC Program for complete details.

7. EQUIPMENT and SUPPLIES

7.1 Assay Platform

Dimension Vista® System

7.2 Equipment

- Refrigerator capable of sustaining 2–8°C.

- Freezer capable of sustaining range not to exceed -20 to -50°C.
- Centrifuge

7.3 Supplies

- Aliquot Plates
- System Fluids
- Assorted calibrated pipettes (MLA or equivalent) and disposable tips

8. PROCEDURE

TBIL Flex® reagent cartridge Cat. No. K1167 and DBIL Flex® reagent cartridge Cat. No. K2125 is required to perform this test.

Bilirubin is performed on the Dimension Vista® System after the method is calibrated (see Reference Material in Calibration section) and Quality Controls are acceptable.

NOTE: For all procedures involving specimens, buttoned lab coats, gloves, and face protection are required minimum personal protective equipment. Report all accidents to your supervisor.

8.1	Sample Processing
1.	A sample rack holding tubes or cups is placed on the rack input lane.
2.	The sample shuttle moves the rack to the barcode reader which identifies the rack and samples to the system.
3.	The rack moves into the sample server and to the rack positioner.
4.	At the same time, aliquot plates move from the aliquot loader into position.
5.	The aliquot probe aspirates the sample from the tubes or cups and dispenses it into the wells of the aliquot plates.
6.	After each aspirate-dispense action, the probe is thoroughly rinsed inside and out to prevent sample carryover.
7.	When sample aspiration is completed, the sample server moves the rack back to the sample shuttle, where it is placed on the output lane and can be removed by the operator.

8.2	Specimen Testing
1.	For QC placement and frequency, refer to the Dimension Vista® QC Schedule in the Laboratory QC Program.
2.	Follow the instructions, outlined in the Dimension Vista® Operator’s Manual
3.	The instrument reporting system contains error messages to warn the user of specific malfunctions. Results followed by such error messages should be held for follow-up. Refer to the Dimension Vista® system manual “Error messages” section for troubleshooting.

8.2	Specimen Testing
4.	Follow protocol in Section 10.5 "Repeat criteria and resulting" for samples with results above or below the Analytical Measurement Range (AMR). Investigate any failed delta result and repeat, if necessary.
5.	Append the appropriate English text code qualifier messages to any samples requiring a comment regarding sample quality and/or any other pertinent factors.

Test Conditions		
	TBIL	DBIL
Sample Volume:	5 µL	5 µL
Reagent 1 Volume:	125 µL	12.5 µL
Reagent 2 Volume:	23.5 µL	25 µL
Reaction Time:	7.5 minutes	4.9 minutes
Test Temperature:	37°C	
Wavelength:	540 & 700 nm	
Type of measurement:	Bichromatic endpoint	

NOTE: In the event that the test system becomes inoperable, notify supervision or designee for further direction. Patient specimens must be stored in a manner that maintains the integrity of the specimen.

9. CALCULATIONS

The instrument automatically calculates the concentration of Bilirubin in mg/dL.

10. REPORTING RESULTS AND REPEAT CRITERIA

10.1 Interpretation of Data

None required

10.2 Rounding

No rounding is necessary. Instrument reports results up to one decimal point.

10.3 Units of Measure

mg/dL

10.4 Clinically Reportable Range (CRR)

Direct Bilirubin: 0.1 – 80.0 mg/dL

Total Bilirubin: 0.1 – 125.0 mg/dL

JUN101012 09:53:53 AM

10.5 Review Patient Data

Each result is reviewed for error messages. Refer to the Dimension Vista system manual "Error messages" section for troubleshooting. Resolve any problems noted before issuing patient reports.

10.6 Repeat Criteria and Resulting

All repeats must replicate the original result within the total allowable error (TEa) of the assay. Refer to TEa listing for specific information.

Values that fall below or within the AMR or CRR may be reported without repeat. Values that exceed the upper ranges must be repeated.

Direct Bilirubin:

IF the result is ...	THEN...
< 0.1 mg/dL	Assure there is sufficient sample devoid of bubbles, cellular debris, and/or fibrin clots. Report as: < 0.1 mg/dL
≥ 16.0 mg/dL	On Board Automated Dilution: Results ≥ 16.0 mg/dL will automatically have repeat testing performed into the instrument using dilution factor of 4. No multiplication is necessary.
> 64.0 mg/dL	Manual Dilution: Using the primary tube, make the smallest dilution possible to bring the raw data within the AMR. Maximum allowable dilution: x 5 DILUENT: Reagent Grade Water Enter dilution factor as a whole number. Re-assay. Readout is corrected for dilution.
> 80.0 mg/dL	If the recommended dilution does not give results within the clinically reportable range, report as: "> 80.0 mg/dL-REP" Bring to the attention of your supervisor prior to releasing result.

Total Bilirubin:

IF the result is ...	THEN...
< 0.1 mg/dL	Assure there is sufficient sample devoid of bubbles, cellular debris, and/or fibrin clots. Report as: < 0.1 mg/dL
≥ 25.0 mg/dL	On Board Automated Dilution: Results ≥ 25.0 mg/dL will automatically have repeat testing performed into the instrument using dilution factor of 4. No multiplication is necessary.

JUN101012 09:53:53 AM

IF the result is ...	THEN...
> 100.0 mg/dL	Manual Dilution: Using the primary tube, make the smallest dilution possible to bring the raw data within the AMR. Maximum allowable dilution: x 5 DILUENT: Reagent Grade Water Enter dilution factor as a whole number. Re-assay. Readout is corrected for dilution.
> 125.0 mg/dL	If the recommended dilution does not give results within the clinically reportable range, report as: "> 125.0 mg/dL-REP" Bring to the attention of your supervisor prior to releasing result.

Message	Code
Verified by repeat analysis	Append -REP to the result.

11. EXPECTED VALUES

11.1 Reference Ranges

Total Bilirubin

Age	Male/Female
Adult (> 17 years):	< 1.0 mg/dL
Pediatric:	
1 month – 17 years	< 0.8
3 days – 30 days	< 10.3
1 – 2 days	< 7.2
0 – 24 hours	< 5.1

Direct Bilirubin, all ages 0.0 – 0.2 mg/dL

Cord Blood Bilirubin < 2.0 mg/dL

11.2 Critical Values

Total Bilirubin, all ages > 17.9 mg/dL

Cord Blood Bilirubin > 17.9 mg/dL

11.3 Standard Required Messages

None established

12. CLINICAL SIGNIFICANCE

Measurements of bilirubin are used in the diagnosis and treatment of liver, hemolytic hematological and metabolic disorders, including hepatitis and gall bladder disease. There are at least four distinct bilirubin species that make up the total bilirubin in serum. The direct reacting species are mono- and diconjugated bilirubin (β - and γ -bilirubin) and the delta fraction (δ -bilirubin), which is tightly bound to albumin. Unconjugated bilirubin (α -bilirubin) is water-insoluble and reacts only after addition of an accelerator such as caffeine.

13. PROCEDURE NOTES

- **FDA Status:** FDA Approved/cleared
- **Validated Test Modifications:** None

The instrument reporting system contains error messages to warn the operator of specific malfunctions. Any report slip containing such error messages should be held for follow-up. Refer to your Dimension Vista Operator's Guide.

The expected maximum observed standard deviations for repeatability using n = 5 replicates at the following bilirubin concentrations are:

TBIL Concentration	Acceptable S.D. Maximum
1.4 mg/dL	0.2 mg/dL
18.3 mg/dL	1.3 mg/dL
BIL Concentration	Acceptable S.D. Maximum
1.0 mg/dL	0.1 mg/dL
14.8 mg/dL	0.9 mg/dL

14. LIMITATIONS OF METHOD

14.1 Analytical Measurement Range (AMR)

Total Bilirubin: 0.1 – 25.0 mg/dL
Direct Bilirubin: 0.1 – 16.0 mg/dL (lower value adjusted to one decimal)

14.2 Precision

Material	Mean mg/dL	Standard Deviation (%CV)	
		Repeatability	Within-Lab
TBIL, Serum Pool	1.0	0.03 (2.7)	0.03 (3.3)
TBIL, MAS bilirubin 3	19.3	0.31 (1.6)	0.56 (2.9)
DBIL, Serum Pool 1	0.4	0.02 (5.7)	0.02 (5.7)
DBIL, Serum Pool 2	13.8	0.22 (1.6)	0.37 (2.7)

14.3 Interfering Substances

HIL Interference:

The TBIL and DBIL method was evaluated for interference according to CLSI/NCCLS EP7-A2. Bias, defined as the difference between the control sample (does not contain interferent) and the test sample (contains interferent), is shown in the table below. Bias exceeding 10% is considered “interference”.

Substance tested	Substance Concentration	TBIL mg/dL	Bias %
Hemoglobin (hemolysate)	1000 mg/dL	1.1, 25	<10
	800 mg/dL	1.1	10
Lipemia (Intralipid®)	1000 mg/dL	1.2	20
	3000 mg/dL	1.4, 21	40, <10

Substance tested	Substance Concentration	DBIL mg/dL	Bias %
Hemoglobin (hemolysate)	≥ 50 mg/dL	≤ 16	
Lipemia (Intralipid®)	1000 mg/dL	0.4	<10
	3000 mg/dL	0.4, 5.1, 14.1	50, <10

14.4 Clinical Sensitivity/Specificity/Predictive Values

Not available

15. SAFETY

Refer to your local and corporate safety manuals and Safety Data Sheet (SDS) for detailed information on safety practices and procedures and a complete description of hazards.

16. RELATED DOCUMENTS

1. Dimension Vista® Clinical Chemistry System Operator’s Manual
2. Dimension Vista® Calibration/Verification Procedure
3. Dimension Vista® Cal Accept Guidelines
4. Dimension Vista® Calibration summary
5. Dimension Vista® Sample Processing, Startup and Maintenance procedure
6. Laboratory Quality Control Program
7. QC Schedule for Siemens Dimension Vista®
8. Laboratory Safety Manual
9. Safety Data Sheets (SDS)
10. Dimension Vista® Limits Chart (AG.F200)
11. Quest Diagnostics Records Management Procedure
12. Dimension Vista® System Error Messages Chart
13. Centrifuge Use, Maintenance and Functions Checks (Lab policy)
14. **Specimen Acceptability Requirements** (Lab policy)
15. Repeat Testing Requirement (Lab policy)
16. Current Allowable Total Error Specifications at http://questnet1.qdx.com/Business_Groups/Medical/qc/docs/qc_bpt_tea.xls
17. Current package insert TBIL Flex® Reagent Cartridge K1167

18. Current package insert DBIL Flex® Reagent Cartridge K2125

17. REFERENCES

1. Ghoshal, Amit K. and Soldin, Steven J., Evaluation of the Dade Behring Dimension® RxL: Integrated chemistry system-pediatric reference ranges. Clinica Chimica Acta 2003; 331:144.
2. Package Insert, TBIL Flex® Reagent Cartridge K1167, Siemens Healthcare Diagnostics Inc., 03/01/2018.
3. Package Insert, DBIL Flex® Reagent Cartridge K2125, Siemens Healthcare Diagnostics Inc., 03/01/2018.
4. Package Insert, BILI CAL, Siemens Healthcare Diagnostics Inc., 03/2015.
5. Package Insert, **Liquid Assayed Multiqua®** Chemistry Controls, Bio-Rad Laboratories, 05/2017.
6. Package Insert, Liquichek Pediatric Control, Bio-Rad Laboratories, 8/2017.

18. REVISION HISTORY

Version	Date	Section	Reason	Reviser	Approval
000	4/12/16		Update title page	L Barrett	R SanLuis
000	4/12/16	3.2	Specify anticoagulant	L Barrett	R SanLuis
000	4/12/16	4.2	Update stability of total bili reagent wells 9-10 to 5 days, add safety statement	A Chini	R SanLuis
000	4/12/16	5.2	Removed uncapped calibrator stability	A Chini	R SanLuis
000	4/12/16	6.1	Update pediatric QC catalog number	A Chini	R SanLuis
000	4/12/16	6.2	Update pediatric QC preparation and stability	A Chini	R SanLuis
000	4/12/16	6.4, 6.6	Replace LIS with Unity Real Time	A Chini	R SanLuis
000	4/12/16	7.2	Change freezer requirements	L Barrett	R SanLuis
000	4/12/16	11.1,11.2	Add Cord Blood	A Chini	R SanLuis
000	4/12/16	11.2	Reformat value to eliminate ≥ sign	L Barrett	R SanLuis
000	4/12/16	17	Update package insert dates	A Chini	R SanLuis
000	4/12/16	Footer	Version # leading zero’s dropped due to new EDCS in use as of 10/7/13	L Barrett	R SanLuis
1	7/25/18	Header	Add WAH	L Barrett	R SanLuis
1	7/25/18	3.1	Add protect from light	L Barrett	R SanLuis
1	7/25/18	4,5,6	Remove individual section labeling instructions and add general one	L Barrett	R SanLuis
1	7/25/18	5.2	Add light exposure comment	L Barrett	R SanLuis
1	7/25/18	6.1, 6.2	Update QC material and storage	L Barrett	R SanLuis

Version	Date	Section	Reason	Reviser	Approval
1	7/25/18	10.5	Move patient review from section 6	L Barrett	R SanLuis
1	7/25/18	10.6	Remove repeat value below AMR/CRR	L Barrett	R SanLuis
1	7/25/18	15	Update to new standard wording	L Barrett	R SanLuis
1	7/25/18	16	Update policy title	L Barrett	R SanLuis
1	7/25/18	17	Update QC insert, update PI dates	L Barrett	R SanLuis

19. ADDENDA

None

Form revised 2/02/2007