

TRAINING UPDATE

| | | | |
|----------------------|----------------------------------|--------------------------|------------------|
| Lab Location: | GEC, SGMC & WAH | Date Distributed: | 8/19/2019 |
| Department: | Technical Specialists, Mgmt & QA | Due Date: | 9/15/2019 |
| | | Implementation: | 9/15/2019 |

DESCRIPTION OF PROCEDURE REVISION

| |
|--|
| Name of procedure: |
| Comparison of Intra-laboratory Test Results SGMC.QA978 v1 Instrument to Instrument Comparison Study for 2 Equivalent Instruments (Same or equivalent instrument, same reference range, 2 instruments) AG.F458 Instrument to Instrument Comparison Study: Intra-Lab (Different method, same reference range, within laboratory) AG.F459 BNP Comparison Study AG.F460 Method to Method: Known Bias, Different Method Comparison AG.F461 Qualitative/Semi-Quantitative Comparison Study AG.F462 |
| Description of change(s): |
| Header: Updated parent facility, changed WAH to WOMC Section 2: Replaced Iris with AUWi PRO Section 5: Updated references to worksheets Section 6: Added worksheets Section 9: Moved worksheets to 'forms'; Added POCT comparison studies Note: The forms are NOT attached. They are excel worksheets that perform the necessary calculations & are nearly identical to the QD versions. They will be placed on the G drive and emailed. This SOP and Forms will be implemented on September 15, 2019 |

Document your compliance with this training update by taking the quiz in the MTS system.

Non-Technical SOP

| | | |
|--------------------|--|-----------------|
| Title | Comparison of Intra-laboratory Test Results | |
| Prepared by | Leslie Barrett | Date: 9/28/2017 |
| Owner | Cynthia Bowman-Gholston, Robert SanLuis | Date: 9/28/2017 |

| Laboratory Approval | | |
|--|------------------|-----------------------|
| Print Name and Title | Signature | Date |
| <i>Refer to the electronic signature page for approval and approval dates.</i> | | |
| | | |
| | | |
| Local Issue Date: | | Local Effective Date: |

| Review: | | |
|-------------------|------------------|-------------|
| Print Name | Signature | Date |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |

TABLE OF CONTENTS

| | |
|--------------------------------|---|
| 1. PURPOSE..... | 2 |
| 2. SCOPE | 2 |
| 3. RESPONSIBILITY..... | 2 |
| 4. DEFINITIONS..... | 2 |
| 5. PROCEDURE..... | 3 |
| 6. RELATED DOCUMENTS | 5 |
| 7. REFERENCES | 5 |
| 8. REVISION HISTORY..... | 6 |
| 9. ADDENDA AND APPENDICES..... | 6 |

1. PURPOSE

This procedure describes the process for periodic instrument and/or method comparison by providing the steps for verifying that an acceptable relationship exists between test results using the same or different methodologies or instruments within a laboratory.

2. SCOPE

This procedure applies to test procedures that are performed

- On multiple instruments within the same laboratory
- Using more than one method within the same laboratory
- Using the same or a different method to test an analyte in point of care and the main laboratory

Note: Examples of systems that require method comparison are:

- Automated vs. manual ABO, Rh, and antibody screening
- BNP on the Centaur vs. BNP by Triage Meter
- Multiple chemistry, hematology, coagulation, etc. analyzers
- Specific gravity by refractometer, dipstick, Clinitek, AUWi PRO

3. RESPONSIBILITY

The Technical Supervisor/Technical Consultant is responsible for implementing this process, ensuring it is performed at the defined frequency and for reviewing all comparison data and initiating corrective action, as necessary.

4. DEFINITIONS

Allowable Total Error (TEa): The amount of error that can be tolerated without invalidating the medical usefulness of the analytical result or the maximum amount of error defined for successful performance in proficiency testing.

Analytical Measurement Range (AMR): The AMR is the range of analyte values that a method can directly measure on the specimen without any dilution, concentration, or other pretreatment not part of the usual assay process.

Estimate of Bias: The difference in results obtained by two different methods. It is calculated as the difference in the mean values from multiple analyses of each method.

Instrument/method mean: The average value of multiple samples run on a single instrument or by a single method.

Sample mean: The average value of the same sample when analyzed on multiple instruments or by multiple methods.

5. PROCEDURE

A. Intra-laboratory Process

1. The department selects a minimum of five (5) specimens that are appropriate for the test method.
2. For quantitative methods obtain specimens with results that span the assay's AMR (low, medium, high).
3. For qualitative methods obtain specimens with positive and negative results.
4. For semi-quantitative methods obtain specimens with positive, negative and equivocal (when applicable) results.
5. If insufficient samples that span the assay range are available, then 'simulated' (analyte-spike) specimen samples may be used. For example, positive specimens for urine drugs of abuse.
6. If possible, analyze the same aliquots on all instruments and by all methods on which the test is performed.
7. Data is submitted to the Technical Supervisor, or designee, for evaluation and review. Refer to appendices for worksheets with calculations.

B. Data Evaluation Criteria

1. Same Analyte and Same / Equivalent Instrument Model with Same Reference Range, Two Instrument Comparison

| |
|--|
| Quantitative (use AG.F458 worksheet) |
| Individual Result Evaluation: <ol style="list-style-type: none">1) Select one instrument as the reference for purposes of comparison.2) The difference between individual sample results should be $< TEa$. |
| Estimate of Bias: The difference between the instrument / method means should be $< TEa/4$. For example: If the from two (2) instruments are 100 and 106 and the $TEa = 24$: <ol style="list-style-type: none">1) $TEa/4 = 6$ and2) The difference in means = 6.3) The result comparison passes. |

Form revised 3/31/00

| |
|---|
| Qualitative (use AG.F462 worksheet) |
| <p>Qualitative: Results are expected to achieve 100% concordance.</p> <ol style="list-style-type: none"> 1) An equivocal specimen is acceptable if it remains equivocal or reads “high” negative or “low” positive. <ol style="list-style-type: none"> a) A high negative is defined as a result that is not $< 70\%$ of the cutoff signal b) A low positive is defined as a result that is not $> 130\%$ of the cutoff signal 2) Semi-Quantitative results that are converted from an OD or Index (specimen signal \div cutoff signal) to a qualitative result are evaluated as qualitative results. 3) Results with a titered or graded result should duplicate within one (+/-1) dilution or grade. |

2. Same Analyte by Different Instrument / Method with Same Reference Range

| |
|---|
| Quantitative (use AG.F459 or BNP-specific worksheet [AG.F460]) |
| <p>Individual Result Evaluation:</p> <ol style="list-style-type: none"> 1) Select one instrument / method as the reference for purposes of comparison. 2) The difference between individual sample results should be $< TEa$. |
| <p>Estimate of Bias: Select the data from the instrument/method having the higher test volume as the reference method. The alternate instrument/method mean must be within $TEa/3$ of the reference instrument/ method mean.</p> |
| Qualitative (use AG.F462 worksheet) |
| <p>Qualitative: Results are expected to achieve 100% concordance.</p> <ol style="list-style-type: none"> 1) Semi-Quantitative results that are converted from an OD or Index (specimen signal \div cutoff signal) to a qualitative result are evaluated as qualitative results. 2) Results with a titered or graded result should duplicate within one (+/-1) dilution or grade. |

3. Same Analyte Different Method with Different Reference Range

| |
|---|
| Quantitative (use AG.F461 worksheet) |
| <p>Individual Result Evaluation: Individual results from alternate platform must be within main platform results +/- TEa, after adjustment for the known bias</p> |
| <p>Estimate of Bias: The observed bias for the alternate platform should be within the instrument/method mean +/- $TEa/2$, after adjustment for the known bias.</p> |

C. Frequency

The minimum frequency for result comparison is every six (6) months.

D. Corrective Action

1. Same Analyte, Same/Equivalent Instrument Model, Same Reference Range:
Service the instrument as needed to bring the comparison data into specifications.
If troubleshooting does not bring data into specification, estimate bias at TEA/3
and obtain approval from the medical director.
2. Same Analyte, Different Instrument/Method, Same Reference Range:
Initiate appropriate corrective action that may include instrument/method
replacement.
3. Same Analyte, Different Method, Different Reference Range:
Corrective action is not needed if the known relationship remains as expected.
If the relationship varies from the expected, initiate an investigation to determine
which method is at fault. Implement corrective action to bring the methods into
specifications.
4. Patient testing will not be performed on any analyte using any test system that
does not provide acceptable comparison data.

E. Documentation

1. Documentation will be maintained of the result comparison studies as well as any
corrective action that is required should the comparison study not meet the
acceptability requirements.
2. The QA Recurring Calendar is utilized as a tool to facilitate this process.

6. RELATED DOCUMENTS

Current Allowable Total Error Specifications at

http://questnet1.qdx.com/Business_Groups/Medical/qc/docs/qc_bpt_tea.xls

Instrument to Instrument Comparison Study for 2 Equivalent Instruments (Same or
equivalent instrument, same reference range, 2 instruments) (AG.F458)

Instrument to Instrument Comparison Study: Intra-Lab (Different method, same
reference range, within laboratory) (AG.F459)

BNP Comparison Study (AG.F460)

Method to Method: Known Bias, Different Method Comparison (AG.F461)

Qualitative/Semi-Quantitative Comparison Study (AG.F462)

7. REFERENCES

Process for Comparison of Intra/Interlaboratory Test Results, Quality Assurance Best
Practice Team, 6/30/2005

8. REVISION HISTORY

| Version | Date | Reason for Revision | Revised By | Approved By |
|---------|--------|---|-----------------------|--------------------------------|
| | | Supersedes SGAH/WAH/GEC.QA16.1 | | |
| 0 | 8/8/19 | Header: Updated parent facility, changed WAH to WOMC Section 2: Replaced Iris with AUWi PRO Section 5: Updated references to worksheets Section 6: Added worksheets Section 9: Moved worksheets to 'forms'; Added POCT comparison studies | L Barrett S Codina | R SanLuis C Bowman-Gholston |
| | | | | |

9. ADDENDA AND APPENDICES

| Appendix | Title |
|----------|----------------------------------|
| A | Point of Care Comparison Studies |
| B | Test and Analyzer List |

Appendix A

Point of Care Comparison Studies

| Step | Action |
|------|--|
| 1 | Point of care meters will be compared to laboratory instruments at least once every six months if a particular analyte is tested in both locations. |
| 2 | A single meter will be compared to a single laboratory instrument that performs equivalent testing using this procedure. Additional meters will be compared to the test meter. |
| 3 | When performing comparisons using the Roche AccuChek II glucose and Abbott iStat meters, a specimen will be collected in a vacutainer tube containing lithium heparin. The specimen will be tested on the point of care meter using whole blood first. Then, the specimen will be centrifuged and the plasma will be tested using the laboratory method. |
| 4 | Point of care comparisons will be reported using the POCT Instrument Comparison Approval Form. |

Appendix B

TEST and ANALYZER LIST

Chemistry - Dimension Analyzers

| | |
|-------|----------------------------|
| ACTM | Acetaminophen |
| ALTI | Alanine Aminotransferase |
| ALB | Albumin |
| ETOH | Alcohol (Ethyl) |
| ALPI | Alkaline Phosphatase |
| AMM | Ammonia |
| AMY | Amylase |
| AST | Aspartate Aminotransferase |
| DBIL | Bilirubin, Direct |
| TBIL | Bilirubin, Total |
| CA | Calcium |
| CRBM | Carbamazepine |
| CTNI | Cardiac Troponin-I |
| CL | Chloride |
| HDLC | Cholesterol, HDL |
| CHOL | Cholesterol, Total |
| CRP | C-Reactive Protein |
| CKI | Creatine Kinase |
| CREA | Creatinine |
| DGNA | Digoxin |
| CO2 | Enzymatic Carbonate |
| FERR | Ferritin |
| FOLAC | Folate |
| FT4 | Free T4 |
| GGT | Gamma Glutamyl Transferase |
| GENT | Gentamicin |
| GLUC | Glucose |
| HA1C | Hemoglobin A1C |
| | |

| | |
|-------|-----------------------------------|
| HCG | Human Chorionic Gonadotropin |
| IRON | Iron |
| TIBC | Iron Binding Capacity, Total |
| LA | Lactic Acid (Lactate) |
| LDI | Lactic Dehydrogenase |
| LIPA | Lipase |
| LITH | Lithium |
| MG | Magnesium |
| MMB | Mass Creatine Kinase MB Isoenzyme |
| MYO | Myoglobin |
| PHNO | Phenobarbital |
| PTN | Phenytoin |
| PHOS | Phosphorus |
| K | Potassium |
| PRALB | Prealbumin |
| PSAT | PSA Total |
| SAL | Salicylate |
| NA | Sodium |
| THEO | Theophylline |
| TSH | Thyroid Stimulating Hormone |
| TOBR | Tobramycin |
| TP | Total Protein |
| TGL | Triglycerides |
| BUN | Urea Nitrogen |
| URCA | Uric Acid |
| VALP | Valproic acid |
| VANC | Vancomycin |
| VB12 | Vitamin B12 |

| | | |
|-------------|------|--|
| QUAL | | Urine Amphetamine/Methamphetamine Screen |
| | AMPH | Urine Amphetamine/Methamphetamine Screen |
| | BARB | Urine Barbiturate Screen |
| | BENZ | Urine Benzodiazepines Screen |
| | COC | Urine Cocaine Metabolite Screen |
| | OPI | Urine Opiates |
| | PCP | Urine Phencyclidine Screen (PCP) |
| | THC | Urine Cannabinoids Screen (THC) |

| | | |
|-----------|------|------------------------|
| UR | CREA | Creatinine, Urine |
| | K | Potassium, Urine |
| | NA | Sodium, Urine |
| | UCFP | Protein, Urine and CSF |
| | GLUC | Glucose, Urine |

| | | |
|-------------|------------|---------------------|
| Calc | % Iron Sat | % Iron Saturation |
| | A/G Ratio | A/G Ratio |
| | IBIL | Bilirubin, Indirect |
| | AGAP | Anion Gap |
| | ALDL | Cholesterol, LDL |

Other Chemistry

| | |
|------|----------------------------|
| BNP | Triage vs. Centaur |
| CTNI | Xpand vs. iSTAT (GEC only) |
| | |

Hematology

Sysmex

| | |
|-------|---|
| WBC | White Blood Cell |
| RBC | Red Blood Cell |
| HGB | Hemoglobin |
| HCT | Hematocrit |
| MCV | Mean Cell Volume |
| MCH | Mean Corpuscular Hemoglobin |
| MCHC | Mean Corpuscular Hemoglobin Concentration |
| RDW | Red Cell distribution Width |
| DIFF | Differential Count |
| PLT | Platelet |
| MPV | Mean Platelet Volume |
| RETIC | Reticulocyte Count |

| | |
|--|-------------------|
| Manual vs Automated Instruments | |
| | Sysmex body Fluid |

| | |
|-----------------|-----------------------------|
| GEC only | |
| | Sysmex vs. back up analyzer |

Urinalysis

Manual vs Automated Instruments

| |
|--|
| Glucose |
| Bilirubin |
| Ketone |
| Blood |
| Protein |
| Nitrite |
| Leukocytes |
| Specific Gravity (Refractometer, Iris, Dipstick) |
| pH |
| Urobilinogen |
| UA Microscopic |

Coagulation

Stagos

| | |
|---------|---------------------------------------|
| PT | Prothrombin Time and INR |
| APPT | Activated Partial Thromboplastin Time |
| Fibro | Fibrinogen |
| D-Dimer | D Dimer |
| | |
| | |
| | |