

TRAINING UPDATE

Lab Location: SGMC & WOMC
Department: Core Lab

Date Distributed: 10/1/2020
Due Date: 10/31/2020
Implementation: 10/21/2020

DESCRIPTION OF PROCEDURE REVISION

Name of procedure:	
Lipase by Dimension Vista® System SGAH.C117 v4	
Description of change(s):	
<p>Note: The information added about the calibration process is already in practice.</p>	
Section	Reason
Header	Changed WAH to WOMC
5.3	Added lot specific correlation factors
8.2	Corrected section reference in step 4
10.6	Updated instruction if above CRR
17	Added customer notification, updated PI
<p>This revised SOP will be implemented on October 21, 2020</p>	

Document your compliance with this training update by taking the quiz in the MTS system.

Technical SOP

Title	Lipase by Dimension Vista® System	
Prepared by	Ashkan Chini	Date: 6/25/2012
Owner	Robert SanLuis	Date: 6/25/2012

Laboratory Approval	Local Effective Date:	
Print Name and Title	Signature	Date
<i>Refer to the electronic signature page for approval and approval dates.</i>		

TABLE OF CONTENTS

1.	Test Information.....	2
2.	Analytical Principle	2
3.	Specimen Requirements.....	2
4.	Reagents	3
5.	Calibrators/Standards	4
6.	Quality Control	6
7.	Equipment And Supplies	8
8.	Procedure	8
9.	Calculations.....	9
10.	Reporting Results And Repeat Criteria.....	10
11.	Expected Values.....	11
12.	Clinical Significance	11
13.	Procedure Notes	11
14.	Limitations Of Method	12
15.	Safety	12
16.	Related Documents	12
17.	References.....	12
18.	Revision History	13
19.	Addenda	14

1. TEST INFORMATION

Assay	Method/Instrument	Test Code
Lipase	Dimension Vista® System	LIPA


Synonyms/Abbreviations
LIPL

Department
Chemistry

2. ANALYTICAL PRINCIPLE

The LIPL method uses a chromogenic ester of methylresorufin as a substrate. Colipase and alkaline pH in the reaction specifically activates pancreatic lipase, the bile salts emulsify the substrate, and cholates suppress other esterase activities in the sample.

Lipase hydrolyzes the substrate 1,2-O-dilauryl-rac-glycero-3-glutaric acid-(6'-methylresorufin) ester to an unstable intermediate: glutaric acid-6-methylresorufin ester in the presence of colipase, bile salt, and CaCl₂. This intermediate is then hydrolyzed by H₂O to yield free methylresorufin which absorbs at 577 nm. Lipase activity is measured as a bichromatic rate reaction at 577/700 nm. The rate of the reaction is proportional to the amount of lipase in the sample.

**3. SPECIMEN REQUIREMENTS****3.1 Patient Preparation**

Component	Special Notations
Fasting/Special Diets	N/A
Specimen Collection and/or Timing	Normal procedures for collecting and storing serum and plasma may be used for samples to be analyzed by this method.
Special Collection Procedures	N/A
Other	N/A

3.2 Specimen Type & Handling

Criteria	
Type -Preferred -Other Acceptable	Plasma (Lithium Heparin) Serum
Collection Container	Plasma: Mint green top tube (PST) Serum: Red top tube, Serum separator tube (SST)
Volume - Optimum - Minimum	1.0 mL 0.5 mL
Transport Container and Temperature	Collection container or Plastic vial at room temperature
Stability & Storage Requirements	Room Temperature: 24 hours
	Refrigerated: 7 days
	Frozen: 1 year
Timing Considerations	Serum or plasma should be physically separated from cells as soon as possible with a maximum limit of two hours from the time of collection.
Unacceptable Specimens & Actions to Take	Specimens that are unlabeled, improperly labeled, or those that do not meet the stated criteria are unacceptable. Request a recollection and credit the test with the appropriate LIS English text code for “test not performed” message. Examples: Quantity not sufficient-QNS; Wrong collection-UNAC. Document the request for recollection in the LIS.
Compromising Physical Characteristics	Gross hemolysis. Reject sample and request a recollection. Credit the test with the appropriate LIS English text code explanation of HMT (Specimen markedly hemolyzed)
Other Considerations	Allow Red Top or SST to clot completely prior to centrifugation.

NOTE: Labeling requirements for all reagents, calibrators and controls include: (1) Open date, (2) Substance name, (3) Lot number, (4) Date of preparation, (5) Expiration date, (6) Initials of tech, and (7) Any special storage instructions. Check all for visible signs of degradation. When placed onboard the analyzer, the instrument captures the date / time loaded and calculates and tracks the opened expiration.

4. REAGENTS

The package insert for a new lot of kits must be reviewed for any changes before the kit is used. A current Package Insert is included as a Related Document.

4.1 Reagent Summary

Reagents	Supplier & Catalog Number
Lipase	Siemens, Flex® reagent cartridge, Cat. No. K3056

4.2 Reagent Preparation and Storage

Reagent	Lipase
Container	Reagent cartridge
Storage	Store at 2-8°C. Protect from light after opening.
Stability	<ul style="list-style-type: none"> • Stable until expiration date stamped on package. • Sealed wells on the instrument are stable for 30 days. • Open well stability: 7 days for wells 1 - 12
Preparation	All reagents are liquid and ready to use.

5. CALIBRATORS/STANDARDS

5.1 Calibrators/Standards Used

Calibrator	Supplier and Catalog Number
ENZ 1 CAL	Siemens Dimension Vista®, Cat. No. KC310

5.2 Calibrator Preparation and Storage

Calibrator	ENZ 1 CAL
Preparation	Calibrator is ready for use. No preparation is required.
Storage/Stability	<ul style="list-style-type: none"> • Store at 2-8°C • Unopened Calibrator: until expiration date on the box. • Opened Calibrator: once the stopper is punctured, assigned values are stable for 7 days when stored on board the Dimension Vista System. • Opened Calibrator: once cap is removed, assigned values are stable for 30 days when recapped immediately after use and stored at 2-8°C. Do not use this vial on board the instrument.

5.3 Calibration Parameter

Criteria	Special Notations
Reference Material	ENZ 1 CAL
Assay Range	10 – 1500 U/L
Suggested Calibration Level	<p>See Reagent Package Insert for lot specific assigned values in U/L.</p> <p>Note: Each reagent carton contains an alert card with lot specific correlation factors. These factors must be entered into the Vista system prior to processing QC and patient samples. To apply the correlation factors correctly:</p> <ul style="list-style-type: none"> • Only one lot of LIPL reagents can be in use at any given time.

	<ul style="list-style-type: none"> • Auto-calibration of new lots must be de-selected in the Calibration Triggers on Method Configuration Screen.
Additional information	<ul style="list-style-type: none"> • The new lot of Lipase reagent is calibrated on one analyzer. QC is performed at the time of the calibration. Coefficient values are not changed at this time. Correlation / patient sample comparison is performed between the new lot and the old lot. • Once all the data is collected for the new lot studies (calibration, QC and correlation), the new lot of Lipase is removed off of the instrument. • When the old lot is used up and the new lot is loaded, the coefficients are changed at this time.
Frequency	<ul style="list-style-type: none"> • Every new reagent cartridge lot. • Every 45 days for any one lot • When major maintenance is performed on the analyzer. • When control data indicates a significant shift in assay.
Calibration Scheme	2 levels, n = 5

5.4 Calibration Procedure

Auto Calibration:

1. Place the required calibrator vials in a carrier. Make sure the barcode labels are entirely visible through the slots.
2. Place the carrier in the loading area.
3. Position the carrier with the labels facing away from the user.
4. Press the **Load** button.
5. Automatic calibration requires that calibrators be on the instrument. As the time for processing approaches, the instrument reviews onboard inventory for the appropriate calibrators.

Manual Calibration:

1. Verify that calibrators and reagents are in inventory on the instrument.
2. Press **System > Method Summary > Calibration**.
3. Select a method from the sidebar menu. Press the **Order Calibration** button on the screen.
4. Verify that the information on the screen is correct. Verify that the calibrator lot is correct using the drop-down menu.
 - a. When calibrating using Vials press **OK**.
 - b. When calibrating using Cups, check the Use Cups box. This displays the rack and cup position fields. For additional cups use the positions in ascending order. Be sure to use the number of calibration levels and cups as specified in the method IFU. Scan the rack barcode and place calibrator cups in an adapter in position 1 on a rack. Press **OK** and load the rack on the instrument.
5. The status field in the calibration screen changes sequentially to Awaiting Scheduling, Preparing Calibrators and Processing.

5.5 Tolerance Limits

IF.....	THEN.....
If result fall within assay-specific specification, and QC values are within acceptable limits,	proceed with analysis
If result falls outside assay-specific specification, or QC values are out of Acceptable limits,	troubleshoot the assay and/or instrument and repeat calibration

6. QUALITY CONTROL

6.1 Controls Used

Controls	Supplier and Catalog Number
Liquid Assayed Multiquel® Levels 1 and 3	Bio-Rad Laboratories Cat. No. 337 and 339

6.2 Control Preparation and Storage

Control	Liquid Assayed Multiquel® Levels 1 and 3
Preparation	Allow the frozen control to stand at room temperature (18-25°C) for 30 minutes or until completely thawed. A precipitate may be present that dissolves upon mixing. Before loading vials, gently swirl the contents to ensure homogeneous with no visible sign of precipitate. (Do not use a mechanical mixer)
Storage/Stability	Frozen: stable until the expiration date at -20 to -50°C. Thawed and Unopened: When stored at 2-8°C and the stopper is not punctured, it will be stable for 30 days for Lipase. This product can be used for 7 days when stored on-board the Siemens Dimension Vista at 2- 8°C. Thawed and Opened: Once the stopper is punctured, all analytes will be stable for 5 days when stored at 2- 8°C. Store away from light.

6.3 Frequency

Analyze all levels of QC material after every calibration and each day of testing (notated on the QC frequency sheets posted on the instruments).

Refer to the Dimension Vista® QC Schedule in the Laboratory policy Quality Control Program and in the Dimension Vista® Quick Reference Guide.

6.4 Tolerance Limits and Criteria for Acceptable QC

Step	Action
1	Acceptable ranges for QC are programmed into the instrument's Quality Control software system and Unity Real Time, and may be posted near the instrument for use during computer downtime.
2	<p>Run Rejection Criteria</p> <ul style="list-style-type: none"> • Anytime the established parameters are exceeded (if one QC result exceeds 2 SD), the run is considered out of control (failed) and patient results must not be reported. • The technologist must follow the procedure in the Laboratory QC Program to resolve the problem.
3	<p>Corrective Action:</p> <ul style="list-style-type: none"> • All rejected runs must be effectively addressed through corrective action. Steps taken in response to QC failures must be documented. Patient samples in failed analytical runs must be <u>reanalyzed according to the Laboratory QC Program</u>. Supervisors may override rejection of partial or complete runs only with detailed documentation and criteria for overrides that are approved by the Medical Director. Consult corrective action guidelines in Laboratory QC Program. Follow corrective action guidelines in the Laboratory QC Program. • Corrective action documentation must follow the Laboratory Quality Control Program.
4	<p>Review of QC</p> <ul style="list-style-type: none"> • QC must be reviewed weekly by the Group Lead or designee and monthly by the Supervisor/Manager or designee. • If the SD and/or CV are greater than established ranges, investigate the cause for the imprecision and document implementation of corrective actions.

6.5 Documentation

- QC tolerance limits are programmed into the instrument and Unity Real Time; it calculates cumulative mean, SD and CV and stores all information for easy retrieval.
- Quality control records are reviewed daily at the bench, weekly by the Group Lead or designee, and monthly by the Supervisor/Manager or designee.
- Refer to complete policies and procedures for QC documentation and for record retention requirements in the Laboratory QC Program.

6.6 Quality Assurance Program

- Each new lot number of reagent or new shipment of the same lot of reagent must be tested with external control materials and previously analyzed samples.

Performance of the new lot must be equivalent to the previous lot; utilize published TEA for acceptability criteria.

- Training must be successfully completed and documented prior to performing this test. This procedure must be incorporated into the departmental competency assessment program.
- The laboratory participates in CAP proficiency testing. All proficiency testing materials must be treated in the same manner as patient samples.
- Monthly QC must be presented to the Medical Director or designee for review and signature.
- Monthly QC mean and SD are sent to Bio-Rad Laboratories for peer group comparison.
- Consult the Laboratory QC Program for complete details.

7. EQUIPMENT and SUPPLIES

7.1 Assay Platform

Dimension Vista® System

7.2 Equipment

- Refrigerator capable of sustaining 2–8°C.
- Freezer capable of sustaining range not to exceed -20 to -50°C.
- Centrifuge

7.3 Supplies

- Aliquot Plates
- System Fluids
- Assorted calibrated pipettes (MLA or equivalent) and disposable tips

8. PROCEDURE

LIPL Flex® reagent cartridge Cat. No. K3056 is required to perform this test.

Lipase is performed on the Dimension Vista® System after the method is calibrated (see Reference Material in Calibration section) and Quality Controls are acceptable.

NOTE: For all procedures involving specimens, buttoned lab coats, gloves, and face protection are required minimum personal protective equipment. Report all accidents to your supervisor.

8.1	Sample Processing
1.	A sample rack holding tubes or cups is placed on the rack input lane.

8.1	Sample Processing
2.	The sample shuttle moves the rack to the barcode reader which identifies the rack and samples to the system.
3.	The rack moves into the sample server and to the rack positioner.
4.	At the same time, aliquot plates move from the aliquot loader into position.
5.	The aliquot probe aspirates the sample from the tubes or cups and dispenses it into the wells of the aliquot plates.
6.	After each aspirate-dispense action, the probe is thoroughly rinsed inside and out to prevent sample carryover.
7.	When sample aspiration is completed, the sample server moves the rack back to the sample shuttle, where it is placed on the output lane and can be removed by the operator.

8.2	Specimen Testing
1.	For QC placement and frequency, refer to the Dimension Vista® QC Schedule in the Laboratory QC Program.
2.	Follow the instructions, outlined in the Dimension Vista® Operator's Manual
3.	The instrument reporting system contains error messages to warn the user of specific malfunctions. Results followed by such error messages should be held for follow-up. Refer to the Dimension Vista® system manual "Error messages" section for troubleshooting.
4.	Follow protocol in Section 10.6 "Repeat criteria and resulting" for samples with results above or below the Analytical Measurement Range (AMR). Investigate any failed delta result and repeat, if necessary.
5.	Append the appropriate English text code qualifier messages to any samples requiring a comment regarding sample quality and/or any other pertinent factors.

Test Conditions	
Sample Volume:	1.62 µL
Reagent 1 Volume:	100 µL
Reagent 2 Volume:	61 µL
Reaction Time:	2.3 minutes
Test Temperature:	37°C
Wavelength:	577 & 700 nm
Type of measurement:	Bichromatic rate

NOTE: In the event that the test system becomes inoperable, notify supervision or designee for further direction. Patient specimens must be stored in a manner that maintains the integrity of the specimen.

9. CALCULATIONS

The instrument automatically calculates the concentration of Lipase in U/L.

10. REPORTING RESULTS AND REPEAT CRITERIA

10.1 Interpretation of Data

None required

10.2 Rounding

No rounding is necessary. Instrument reports results as a whole number.

10.3 Units of Measure

U/L

10.4 Clinically Reportable Range (CRR)

10 – 30,000 U/L

10.5 Review Patient Data

Each result is reviewed for error messages. Refer to the Dimension Vista system manual “Error messages” section for troubleshooting. Resolve any problems noted before issuing patient reports.

10.6 Repeat Criteria and Resulting

All repeats must replicate the original result within the total allowable error (TEa) of the assay. Refer to TEa listing for specific information.

Values that fall below or within the AMR or CRR may be reported without repeat. Values that exceed the upper ranges must be repeated.

IF the result is ...	THEN...
< 10 U/L	Assure there is sufficient sample devoid of bubbles, cellular debris, and/or fibrin clots. Report as: < 10 U/L
≥ 1500 U/L	<p>On Board Automated Dilution: Results ≥ 1500 U/L will automatically have repeat testing performed into the instrument using dilution factor of 20. No multiplication is necessary.</p>
> 30,000 U/L	<p>If the recommended dilution does not give results within the clinically reportable range, report as: “> 30,000 U/L-REP” Bring to the attention of Tech in Charge (TIC) or Group Lead to check for integrity issues prior to release of results.</p>

Message	Code
Verified by repeat analysis	Append –REP to the result.

11. EXPECTED VALUES

11.1 Reference Ranges

Lipase	Male and Female
Adult (>19 years):	63 - 286 U/L
Pediatric:	
17 - 18 years	58 - 260 U/L
11 - 16 years	46 - 211 U/L
7 - 10 years	44 - 199 U/L
2 - 6 years	44 - 199 U/L
1 year	44 - 199 U/L
3 - 11 months	43 - 190 U/L
0 - 2 months	44 - 174 U/L

11.2 Critical Values

None established

11.3 Standard Required Messages

None established

12. CLINICAL SIGNIFICANCE

Pancreatic lipase degrades dietary triglycerides to glycerol and free fatty acids in the duodenum in the presence of bile salts.

Lipase measurements are used to diagnose and monitor treatment of diseases of the pancreas, such as acute and chronic pancreatitis and obstruction of the pancreatic duct.

13. PROCEDURE NOTES

- **FDA Status:** FDA Approved/cleared
- **Validated Test Modifications:** None

The instrument reporting system contains error messages to warn the operator of specific malfunctions. Any report slip containing such error messages should be held for follow-up. Refer to your Dimension Vista Operator's Guide.

The expected maximum observed standard deviations for repeatability using n = 5 replicates at the following lipase concentrations are:

LIPL Concentration	Acceptable S.D. Maximum
72 U/L	8 U/L
150 U/L	11 U/L
519 U/L	32 U/L

14. LIMITATIONS OF METHOD

14.1 Analytical Measurement Range (AMR)

10 – 1500 U/L

14.2 Precision

Material	Mean U/L	Standard Deviation (%CV)	
		Repeatability	Within-Lab
Multiquel Unassayed Control			
Level 1	72	2 (3)	2 (3)
Level 2	150	3 (2)	4 (3)
Level 3	519	8 (2)	12 (2)

14.3 Interfering Substances

HIL Interference:

The LIPL method was evaluated for interference according to CLSI/NCCLS EP7-A2.

Substance tested	Substance Concentration	LIPL U/L	Bias %
Hemoglobin (hemolysate)	1000 mg/dL	200, 1200	<10
Bilirubin (unconjugated)	80 mg/dL	200, 1200	<10
Bilirubin (conjugated)	80 mg/dL	200, 1200	<10
Lipemia Intralipid®	3000 mg/dL	200, 1200	<10

14.4 Clinical Sensitivity/Specificity/Predictive Values

Not available

15. SAFETY

Refer to your local and corporate safety manuals and Safety Data Sheet (SDS) for detailed information on safety practices and procedures and a complete description of hazards.

16. RELATED DOCUMENTS

1. Dimension Vista® Clinical Chemistry System Operator's Manual
2. Dimension Vista® Calibration/Verification Procedure
3. Dimension Vista® Cal Accept Guidelines

4. Dimension Vista® Calibration summary
5. Dimension Vista® Sample Processing, Startup and Maintenance procedure
6. Laboratory Quality Control Program
7. QC Schedule for Siemens Dimension Vista®
8. Laboratory Safety Manual
9. Safety Data Sheets (SDS)
10. Siemens Dimension Vista® Limits Chart (AG.F200)
11. Quest Diagnostics Records Management Procedure
12. Dimension Vista® System Error Messages Chart
13. Centrifuge Use, Maintenance and Functions Checks (Lab policy)
14. Specimen Acceptability Requirements (Lab policy)
15. Repeat Testing Requirement (Lab policy)
16. Current Allowable Total Error Specifications at
http://questnet1.qdx.com/Business_Groups/Medical/qc/docs/qc_bpt_tea.xls
17. Current package insert LIPL Flex® Reagent Cartridge K3056

17. REFERENCES

1. Ghoshal, Amit K. and Soldin, Steven J., Evaluation of the Dade Behring Dimension® RxL: Integrated chemistry system-pediatric reference ranges. Clinica Chimica Acta 2003; 331:144.
2. Package Insert, LIPL Flex® Reagent Cartridge K3056, Siemens Healthcare Diagnostics Inc., 01/24/2020.
3. Package Insert, ENZ 1 CAL, Siemens Healthcare Diagnostics Inc., 6/2019.
4. Package Insert, Liquid Assayed Multiquel® Chemistry Controls, Bio-Rad Laboratories, 08/2019.
5. Customer Notification, *Follow-up Information Regarding Lot-to-Lot Variability When Using Dimension Vista Lipase (LIPL) Flex reagent cartridges*, DC-18-05.B.US-DV, Siemens Healthcare Diagnostics Inc., 11/5/2018

18. REVISION HISTORY

Version	Date	Section	Reason	Reviser	Approval
000	3/9/2013	10.4	Revised CRR upper limit	A Chini	R SanLuis
000	3/9/2013	10.5	Removed manual dilution, revised dilution factor for on board automatic dilution	A Chini	R SanLuis
001	3/2/2017	Header	Add WAH	L Barrett	R SanLuis
001	3/2/2017	3.2	Specify anticoagulant, remove specimen onboard stability	L Barrett	R SanLuis
001	3/2/2017	4,5,6	Remove individual section labeling instructions and add general one	L Barrett	R SanLuis
001	3/2/2017	6.1, 6.2	Update QC material and storage	L Barrett	R SanLuis
001	3/2/2017	6.4, 6.6	Replace LIS with Unity Real Time	L Barrett	R SanLuis

Version	Date	Section	Reason	Reviser	Approval
001	3/2/2017	7.2	Change freezer upper limit to -50C	L Barrett	R SanLuis
001	3/2/2017	10.5	Move patient review from section 6	L Barrett	R SanLuis
001	3/2/2017	15	Update to new standard wording	L Barrett	R SanLuis
001	3/2/2017	17	Update QC product and PI dates	L Barrett	R SanLuis
001	3/2/2017	Footer	Version # leading zero's dropped due to new EDCS in use as of 10/7/13	L Barrett	R SanLuis
2	3/22/19	Header	Update parent facility	L Barrett	R SanLuis
2	3/22/19	10.6	Remove repeat value below AMR/CRR	L Barrett	R SanLuis
2	3/22/19	16	Update policy title	L Barrett	R SanLuis
2	3/22/19	17	Update PI dates	L Barrett	R SanLuis
3	9/28/20	Header	Changed WAH to WOMC	L Barrett	R SanLuis
3	9/28/20	5.3	Added lot specific correlation factors	L Barrett	R SanLuis
3	9/28/20	8.2	Corrected section reference in step 4	L Barrett	R SanLuis
3	9/28/20	10.6	Updated instruction if above CRR	L Barrett	R SanLuis
3	9/28/20	17	Added customer notification, updated PI	L Barrett	R SanLuis

19. ADDENDA

None