

x Germantown Emergency Centerx Shady Grove Medical Center

x White Oak Medical Center

## **Planned Deviation Documentation**

Procedure:	D-Dimer			
Reason for Deviation:	<ul> <li>Shortage of Quality Control - D-Dimer testing volumes have exponentially increased due to COVID-19. The vendor has reported a shortage of D-Dimer Quality Control material.</li> <li>To decrease the amount of QC testing for D-Dimer, the QC setting on each Stago instrument has changed from QC being triggered to run every 4 hours to every 8 hours. (CAP requirement is every 8 hours with patient testing.)</li> <li>When QC is reconstituted, 500 uL aliquot will be removed, frozen, and labeled. 500 uL will be loaded on the instrument in microcups.</li> </ul>			
Medical Director or Designee:	Dr. Cacciabeve / Rob SanLuis – <i>refer to electronic</i> <i>signatures</i>		Date:	
Supervisor:	Hollie Genser		Date: 2/18/21	
Deviation Start Date:	2/18/21		Deviation End Date: Revise SOP to match change	
Testing Personnel Signatures	Start Deviation Signatures		End Deviation Signatures	
Printed Name / Tech Code	Signature	Date	Signature	Date
	To be captured electronically via MTS			