TRAINING UPDATE

Lab Location: Department: WOMC and FWMC Blood Bank

Date Implemented:
Due Date:

3/10/2022 3/16/2022

DESCRIPTION OF PROCEDURE REVISION

Name of procedure:

Referral of Blood Bank Workups from FWMC to WOMC

Description of change(s):

Effective Wednesday, March 16, the FWMC T&S order will generate 2 labels.

- The TS label will be used to collect the T&S specimen using the BB/TS armbanding system
- The BBWO label will be used to collect an extra pink tube to be sent to WOMC for storage. The extra pink tube will be used to complete AbID and ABO discrepancy testing if needed.

This change was made to help decrease the time delays related to sample recollection and transport when a workup must be referred to WOMC. We will send all samples with the next available courier. If blood is ordered and additional workup is needed, a STAT courier will be called to transport the specimen.

FWMC BB staff members will complete the "Blood Bank Workup Referral Form" and fax to WOMC BB when antibody ID or ABO discrepancy testing is required.

- 1. The FWMC BB staff member will notify the patient care area that the sample is being referred and compatible blood products may be delayed.
- 2. The FWMC BB staff member will also obtain the patient's pregnancy, RhIG, and transfusion history.
- 3. The completed form will be faxed to WOMC BB when workup is required.



Blood Bank Workup Referral Form				
Patient Information				
Name			Birthdate	
MRN			Gender	
T&S Accession			Collection Date	
Diagnosis			Transfuse Orders? Is blood needed STAT?	
	FV	NMC Ho	spital Notification	
Nam	ne of person notified			
Date and	time of notification			
	Tech code			
		Obtain	Patient History	
What is the	patient's pregnancy history?	□ Curre □ Has b □ Neve	ently pregnant been pregnant in the past r pregnant nancy history unknown/unavaila	able
Has the patient r	received RhIG within the past 90 days?	□ No	ate given: history unknown/unavailable	
Has the patient been transfused?		 □ Yes, within the previous 90 days □ Yes, but not within the past 90 days □ Never transfused □ Transfusion history unknown/unavailable 		
Has the patient been hospitalized within the past 90 days?		 □ Yes Hospital Name: Approximate Date: □ No □ Hospitalization history unknown/unavailable 		
Name and Phon	e Number of FWMC Hospital Contact			
F	ax This Form AND Belo	ow Infor	mation to WOMC BB at 301-388	3-7506
	tibody Identification V		Antibody screen result	

Fax This Form AND Below Information to WOMC BB at 301-388-7506		
Antibody Identification Workup	Antibody screen results on antigram	
	Original ABO/Rh typing results	
ABO Discrepancy Workup	Repeat ABO/Rh typing results	

Non-Technical SOP

Title	Referral of Blood Bank Workups from FWMC to WOMC		
Prepared by	Stephanie Codina	Date: 3/9/22	
Owner	Stephanie Codina	Date: 3/9/22	

Laboratory Approval		
Print Name and Title Signature Date		
Refer to the electronic signature page for approval and approval dates.		
Local Issue Date:	Local Effective Date:	

TABLE OF CONTENTS

1.	PURPOSE	1
	SCOPE	
3.	RESPONSIBILITY	. 1
4.	DEFINITIONS	. 2
5.	PROCEDURE	. 2
	RELATED DOCUMENTS	
	REFERENCES	
8.	REVISION HISTORY	. 3
	ADDENDA AND APPENDICES	

1. PURPOSE

To define the process for referring blood bank workups from Fort Washington Medical Center to White Oak Medical Center.

2. SCOPE

This procedure applies to any Fort Washington Medical Center Blood Bank specimen that needs to be referred to White Oak Medical Center for additional testing.

3. RESPONSIBILITY

All blood bank staff members must understand and adhere to this procedure when referring blood bank specimens for testing.

4. **DEFINITIONS**

NA

5. PROCEDURE

Step	Action
1	 Each T&S order placed at FWMC will generate a separate "BBWO" order. A. The BBWO order requires the collection of a 6mL pink tube of blood. B. The BBWO tube will be referred to WOMC with the next available courier using routine tracking. C. Blood bank staff will result the BBWO test using the code "ATWO" which translates to "Additional Tube at White Oak Medical Center." The BBWO tube will be used to complete additional blood testing as indicated based on the T&S results. The original T&S tube will remain at FWMC for crossmatch testing after the workup is complete. The specimen should be sent via STAT courier if the patient is awaiting
	transfusion.
2	When additional testing is necessary, FWMC staff will complete the "Blood Bank Workup Referral Form" and fax the completed form to WOMC.
3	 Complete the patient information section of the form. A. Document the patient's full name, medical record number, birthdate, gender, and diagnosis. B. Document the T&S accession number and collection date. C. Indicate if transfuse orders have been placed.
4	 Immediately notify the patient care area when testing indicates the need to refer a specimen for additional testing. This includes, but is not limited to, a positive antibody screen or ABO discrepancy. A. Document the name of the person notified and time of notification on the form. B. Inform the patient care area that crossmatch compatible blood products will be delayed. C. Offer emergency release blood products as needed.
5	Ask the patient's primary nurse or provider the following questions and document the responses on the form. A. What is the patient's pregnancy status? B. Has the patient received RhIG in the past 90 days? C. Has the patient been transfused within the previous 90 days? D. Has the patient been hospitalized within the previous 90 days? If yes, determine the hospital name/location and approximate admission date. E. Does the workup need to be completed STAT?

Step	Action		
6	Fax the completed form and the supporting documents to the WOMC laboratory.A. Send the completed antigram for positive antibody screens.B. Send downtime documentation of the original and repeat ABO testing for ABO discrepancy workups.		
7	 WOMC staff members will complete the workup (including resulting and billing) per procedure. A. Blood products may be provided by WOMC blood bank or ordered from Inova Blood Donor Services. Avoid sending ARC units to FWMC due to contract and billing issues. B. Crossmatching may be performed to ensure blood products are compatible, but IS and gel crossmatched must be performed by FWMC staff using the original T&S specimen prior to issue. 		

6. RELATED DOCUMENTS

Blood Bank Referral Form (AG.FW31)

7. REFERENCES

NA

8. REVISION HISTORY

Version	Date	Reason for Revision	Revised By	Approved By
	· · · · · · · · · · · · · · · · · · ·			

9. ADDENDA AND APPENDICES

NA