

TRAINING UPDATE

Lab Location: Core Lab
Department: Microbiology

Date Distributed: 8/8/2022
Due Date: 8/22/2022

DESCRIPTION OF PROCEDURE REVISION

| |
|---|
| Name of procedure: |
| Technical SOP: Resulting Microbiology Direct Exams (AHC.M28 v.4) and Wet Prep (AHC.M23 v.8). |
| Description of change(s): |
| <p>The wet prep keyboard has been updated to allow reporting to match SOP.</p> <p>The Wet Prep SOP, section 10.1, states; “Report the <u>presence or absence</u> of Trichomonas, yeast and clue cells”. Up until the update, we were not able to report the absence of those organism individually. The keyboard update now allows for reporting “Yeast not seen”, “Clue cells not seen” and “Trichomonas not seen”. We will now report the absence, in addition to presence, of those organisms, as necessary per the SOP.</p> <p>For details on the updated keyboard, refer to the SOP “Resulting Direct Microbiology Exams”. See highlighted sections on attached copy of SOP for changes.</p> <p>This revised SOP will be implemented on August 23, 2022</p> |

Document your compliance with this training update by taking the quiz in the MTS system.

AHC.M23 Wet Prep

Copy of version 8.0 (in review)

Effective Date 8/23/2022

Uncontrolled Copy printed on 8/8/2022 1:06 PM

Organization Adventist HealthCare

Approval and Periodic Review Signatures

| Type | Description | Date | Version | Performed By | Notes |
|---|---------------------------|-----------|---------|--------------------|--|
| Approval | Lab Director | 7/8/2022 | 7.0 | Nicolas Cacciabeve | |
| Approval | Microbiology Lab Director | 7/7/2022 | 7.0 | Ronald Master | |
| Approval | Lab Director | 7/8/2020 | 6.0 | Nicolas Cacciabeve | |
| Approval | Micro Director approval | 7/7/2020 | 6.0 | Ronald Master | |
| Approval | QA approval | 7/6/2020 | 6.0 | Leslie Barrett | |
| Approval Captured outside MediaLab | Lab Director | 7/16/2018 | 5.0 | Nicolas Cacciabeve | Recorded on 11/21/2018 by Leslie Barrett when document added to MediaLab |
| Periodic review Captured outside MediaLab | Designated Reviewer | 7/16/2018 | 5.0 | Nicolas Cacciabeve | Recorded on 11/21/2018 by Leslie Barrett when document added to MediaLab |

Approvals and periodic reviews that occurred before this document was added to the MediaLab Document Control system may not be listed.

Version History

| Version | Status | Type | Date Added | Date Effective | Date Retired |
|---------|----------------------|-----------------------------------|------------|----------------|--------------|
| 7.0 | Approved and Current | Major revision | 7/7/2022 | 7/8/2022 | Indefinite |
| 6.0 | Retired | Major revision | 7/6/2020 | 7/8/2020 | 7/8/2022 |
| 5.0 | Retired | First version in Document Control | 11/21/2018 | 8/7/2018 | 7/8/2020 |

Technical SOP

| | | |
|--------------------|-----------------|-----------------|
| Title | Wet Prep | |
| Prepared by | Ron Master | Date: 8/25/2009 |
| Owner | Ron Master | Date: 8/25/2009 |

| Laboratory Approval | Local Effective Date: | |
|--|------------------------------|------|
| Print Name and Title | Signature | Date |
| <i>Refer to the electronic signature page for approval and approval dates.</i> | | |

TABLE OF CONTENTS

| | | |
|-----|--|---|
| 1. | Test Information..... | 2 |
| 2. | Analytical Principle | 2 |
| 3. | Specimen Requirements..... | 2 |
| 4. | Reagents..... | 3 |
| 5. | Calibrators/Standards..... | 3 |
| 6. | Quality Control | 3 |
| 7. | Equipment And Supplies | 3 |
| 8. | Procedure | 4 |
| 9. | Calculations..... | 4 |
| 10. | Reporting Results And Repeat Criteria..... | 4 |
| 11. | Expected Values..... | 5 |
| 12. | Clinical Significance..... | 6 |
| 13. | Procedure Notes | 6 |
| 14. | Limitations Of Method | 6 |
| 15. | Safety | 7 |
| 16. | Related Documents..... | 7 |
| 17. | References..... | 7 |
| 18. | Revision History | 7 |
| 19. | Addenda | 8 |

1. TEST INFORMATION

| Assay | Method/Instrument | Test Code |
|----------|-------------------|-----------|
| Wet Prep | N/A | WETP |

| Synonyms/Abbreviations |
|------------------------|
| N/A |

| Department |
|--------------|
| Microbiology |

2. ANALYTICAL PRINCIPLE

Yeast and *Trichomonas* sp. can be found in urine or vaginal discharges. *Trichomonas* move by a progressive undulating, whipping of flagella and pseudopodial movement. Outside the body *Trichomonas* rapidly succumbs at temperatures higher than 40°C; therefore it is imperative to examine specimens requested for *Trichomonas* immediately upon arrival in the Laboratory. Yeast and clue cells can also be detected from vaginal discharge.

3. SPECIMEN REQUIREMENTS**3.1 Patient Preparation**

| Component | Special Notations |
|-----------------------------------|--|
| Fasting/Special Diets | N/A |
| Specimen Collection and/or Timing | Specimen should be obtained on a swab and submitted in a tube containing 0.5 mL 0.85-0.9% saline (0.85-0.9% sodium chloride). A red top vacutainer tube (without additives) may be used. Deliver to Laboratory immediately. |
| Special Collection Procedures | N/A |
| Other | N/A |

3.2 Specimen Type & Handling

| Criteria | |
|-------------------|---------------------|
| Type | Vaginal discharges. |
| -Preferred | |
| -Other Acceptable | None |

| Criteria | |
|---|---|
| Collection Container | Swab in tube or a red-top vacutainer (without additives) containing 0.5 mL 0.85-0.9% saline |
| Volume - Optimum - Minimum | Swab in 0.5mL 0.85-0.9% saline Swab in 0.5mL 0.85-0.9% saline |
| Transport Container and Temperature | Collection container at room temperature |
| Stability & Storage Requirements | Room Temperature: Do not let stand, test immediately. |
| | Refrigerated: Unacceptable |
| | Frozen: Unacceptable |
| Timing Considerations | Process immediately. |
| Unacceptable Specimens & Actions to Take | Dry swab. Call and request another sample |
| Compromising Physical Characteristics | N/A |
| Other Considerations | N/A |

4. REAGENTS

N/A

5. CALIBRATORS/STANDARDS

N/A

6. QUALITY CONTROL

N/A

7. EQUIPMENT and SUPPLIES**7.1 Assay Platform**

N/A

7.2 Equipment

Microscope

7.3 Supplies

Sterile transfer pipette

Paper towels

Glass slide

Coverslip

Gloves

8. PROCEDURE

NOTE: For all procedures involving specimens, buttoned lab coats, gloves, and face protection are required minimum personal protective equipment. Report all accidents to your supervisor.

| 8.1 | Action |
|-----|--|
| 1. | Check order and verify patient name in the LIS matches name on specimen. |
| 2. | Place a drop or two of the specimen on a slide using a plastic transfer pipette. |
| 3. | Gently place a coverslip over the surface of the material on the slide. |
| 4. | Examine with low (10x objective) and high dry (40x objective) power objectives. |
| 5. | Enter results in computer utilizing Microbiology Result Entry. |

9. CALCULATIONS

N/A


10. REPORTING RESULTS AND REPEAT CRITERIA

** See procedure AHC.M28- Resulting Microbiology Direct Exams on how to result.



10.1 Interpretation of Reporting of Results

Report the presence or absence of *Trichomonas*, yeast and clue cells.


- In wet mounts the trophozoites of *T. vaginalis* move with a rapid, jerky motion and possess an undulating membrane that extends half the length of the organism. They are pear-shaped, 7-23 μm long and 5-15 μm wide. There are 3-5 anterior flagellae and 1 posterior flagella.



- Clue cells are squamous epithelial cells from the vagina with many small rod shaped bacteria adherent to their surfaces. **NOTE: Bacteria must extend beyond the cytoplasmic margins to be considered a clue cell.**



- Yeast are small, ovoid cells which often have small buds or pseudohyphae attached.



DO NOT REPORT SPERMATOZOA.

11. EXPECTED VALUES

11.1 Reference Ranges

No Trichomonas, yeast, or clue cells seen

11.2 Critical Values

None established

11.3 Standard Required Messages

None established

12. CLINICAL SIGNIFICANCE

Trichomonas vaginalis is a parasite that is distinguished by its rapid movement and flagella. Motile *Trichomonas* trophozoites may be identified in a vaginal sample by its characteristic structure.

Bacterial vaginosis is the most common type of vaginal infection and can sometimes be detected by the presence of “clue cells”. Clue cells are epithelial cells entirely covered with bacteria giving them a “furlike” appearance. If the organisms sticking to the edges or on top of the cell, without extending past the cytoplasmic margins, a diagnosis of clue cells cannot be made. Note: Certain anaerobic, non-pathogenic, species tend to adhere to the epithelial surface.

Yeast vaginitis is primarily caused by *Candida albicans*, although other *Candida* species are becoming increasingly important as disease agents. *Candida albicans*, in low numbers, is considered part of the normal vaginal flora, but may proliferate to cause an infection.

13. PROCEDURE NOTES

- **FDA Status:** LDT without message
 - **Validated Test Modifications:** None
1. Examine wet preparation immediately as motility disappears rather rapidly (within 35-40 minutes) making it impossible to detect the parasite in wet preparations.
 2. *Trichomonas vaginalis* can also be observed on Gram stain of the specific discharge or sedimented urine. Here the *Trichomonas* appear larger than the polymorphonucleated white cells, but smaller than epithelial cells. *Trichomonas* cytoplasm is typically foamy in appearance; the parasite stains slightly pink, demonstrating its typical pear to oval or any shape. The elliptical nucleus is clearly visible as it stains darker than the cytoplasm of the organism. Flagella are sometimes visible but not always.

14. LIMITATIONS OF METHOD**14.1 Analytical Measurement Range (AMR)**

N/A

14.2 Precision

N/A

14.3 Interfering Substances

N/A

14.4 Clinical Sensitivity/Specificity/Predictive Values

N/A

15. SAFETY

Refer to your local and corporate safety manuals and Safety Data Sheet (SDS) for detailed information on safety practices and procedures and a complete description of hazards.

16. RELATED DOCUMENTS

AHC.M28 Resulting Microbiology Direct Exams, Microbiology procedure

17. REFERENCES

- Baron, Ellen Jo, Sydney Finegold, Bailey and Scott's Diagnostic Microbiology, C.V. Mosby Co., St. Louis, .2002.
- P.C. Beaver, R.C. Jung, E.W. Cupp. *Clinical Parasitology*, 9th Edition, p. 49-51, Lee and Febiger Publishers, Philadelphia, PA, 1984.
- Henry, J.B, *Clinical Diagnosis and Management by Laboratory methods*, 19th ed., W. B. Saunders Company, Philadelphia, 1996
- Quest Diagnostics Incorporated procedure *Wet Mount (Vaginal) for Clue Cells, Yeast and Trichomonas vaginalis*, QDRRL820v2

18. REVISION HISTORY

| Version | Date | Section | Reason | Reviser | Approval |
|---------|-----------|---------|--|------------|-----------|
| | | | Supersedes SOP M010.002 | | |
| 000 | 10/12/09 | 8.1 | LIS update to GUI system | A. Sears | R. Master |
| 000 | 10/12/09 | 16 | Added procedure for resulting | L. Barrett | R. Master |
| 001 | 10/4/2011 | 3.2 | Deleted sources other than vaginal | R. Master | R. Master |
| 002 | 5/16/2012 | 10.1 | Deleted hyphae, deleted redundant report comments | R. Master | R. Master |
| 003 | 7/6/16 | Header | Added other sites | L. Barrett | R. Master |
| 003 | 7/6/16 | 3.2 | Specify saline concentration | R. Master | R. Master |
| 003 | 7/6/16 | 11.1 | Corrected reference range | R. Master | R. Master |
| 003 | 7/6/16 | Footer | Version # leading zero's dropped due to new EDCS in use as of 10/7/13. | L. Barrett | R. Master |
| 4 | 7/10/18 | 10.1 | Added photos | L. Barrett | R. Master |
| 4 | 7/10/18 | 15 | Updated to new standard wording | L. Barrett | R. Master |
| 4 | 7/10/18 | 17 | Added RRL SOP | L. Barrett | R. Master |
| 5 | 7/6/20 | Header | Changed WAH to WOMC | L. Barrett | R. Master |
| 6 | 7/7/22 | Header | Changed site to All Laboratories | D Collier | R. Master |
| 6 | 7/7/22 | Footer | Changes prefix to AHC | D Collier | R. Master |

| | | | | | |
|---|---------|-------------------|---|------------------------|-----------|
| 7 | 7/12/22 | Section 10 | Added SOP reference on how to result Added more info defining Clue Cells | M Sabonis D Collier | R. Master |
| 7 | 7/12/22 | Section 10 and 16 | Added related document and changed SOP prefix to AHC | M Sabonis D Collier | R. Master |

19. ADDENDA

None

Retired or Not Yet Effective
Uncontrolled copy

AHC.M28 Resulting Microbiology Direct Exams

Copy of version 4.0 (approved, not yet effective)

Last Approval or Periodic Review Completed 8/2/2022

Uncontrolled Copy printed on 8/8/2022 1:04 PM

Next Periodic Review Needed On or Before 8/2/2024

Organization Adventist HealthCare

Effective Date 8/23/2022

Approval and Periodic Review Signatures

| Type | Description | Date | Version | Performed By | Notes |
|---|--------------------------------|------------|---------|---|--|
| Approval | Lab Director | 8/2/2022 | 4.0 | Nicolas Cacciabeve | |
| Approval | Laboratory Operations Director | 8/2/2022 | 4.0 | <i>Robert SanLuis</i> Robert SanLuis | |
| Approval | Lab Director | 11/9/2020 | 3.0 | Nicolas Cacciabeve | |
| Approval | Micro Director approval | 11/6/2020 | 3.0 | Ronald Master | |
| Approval | QA approval | 11/6/2020 | 3.0 | Leslie Barrett | |
| Periodic review | Designated Reviewer | 12/26/2018 | 2.0 | <i>Robert SanLuis</i> Robert SanLuis | |
| Approval Captured outside MediaLab | Lab Director | 3/24/2017 | 2.0 | Nicolas Cacciabeve | Recorded on 11/21/2018 by Leslie Barrett when document added to MediaLab |
| Periodic review Captured outside MediaLab | Designated Reviewer | 3/24/2017 | 2.0 | Nicolas Cacciabeve | Recorded on 11/21/2018 by Leslie Barrett when document added to MediaLab |

Approvals and periodic reviews that occurred before this document was added to the MediaLab Document Control system may not be listed.

Version History

| Version | Status | Type | Date Added | Date Effective | Date Retired |
|---------|-----------------------------|-----------------------------------|------------|----------------|--------------|
| 4.0 | Approved, Not Yet Effective | Major revision | 8/1/2022 | 8/23/2022 | Indefinite |
| 3.0 | Approved and Current | Major revision | 11/6/2020 | 11/9/2020 | 8/23/2022 |
| 2.0 | Retired | First version in Document Control | 11/21/2018 | 4/10/2017 | 11/9/2020 |

Non-Technical SOP

| | | |
|--------------------|--|-----------------|
| Title | Resulting Microbiology Direct Exams | |
| Prepared by | Marie Sabonis | Date: 10/8/2009 |
| Owner | Marie Sabonis, Ronald Master | Date: 10/8/2009 |

| Laboratory Approval | | |
|--|------------------|-----------------------|
| Print Name and Title | Signature | Date |
| <i>Refer to the electronic signature page for approval and approval dates.</i> | | |
| Local Issue Date: | | Local Effective Date: |

TABLE OF CONTENTS

1. PURPOSE..... 1

2. SCOPE..... 1

3. RESPONSIBILITY 2

4. DEFINITIONS 2

5. PROCEDURE..... 2

6. RELATED DOCUMENTS 7

7. REFERENCES 7

8. REVISION HISTORY 7

9. ADDENDA AND APPENDICES..... 7

1. PURPOSE

Describe the procedure for resulting Microbiology direct exams.

2. SCOPE

Applies to the resulting of gram stains (non-blood cultures), wet preps and malaria smears. To assist in resulting, certain orderables are tied to resulting keyboard.

Example:

| Orderable | Keyboard |
|-------------------------------------|----------|
| Gram Stain(GS) | MDE |
| CSF Culture and gram stain (XCSFC) | MDE |
| Fluid culture and gram stain(XFLC) | MDE |
| Wound culture and gram stain(XWDCG) | MDE |
| Malaria(MAL) | DES2 |
| Wet Prep(WET) | WPDE |

3. RESPONSIBILITY

It is the responsibility of all personnel assigned to Microbiology and Lead technologists to read, understand and perform this procedure.

4. DEFINITIONS


None

5. PROCEDURE


1. Access the GUI Microbiology Result Entry via the Misys Gateway application.
2. The first screen is divided into three sections
 - Data Entry Mode
 - Data Search
 - Accession / Battery search
3. In the Data Entry Mode section, keep the default of specimen mode
4. In the Data Search section, key in the accession number you want to result in the "Value" field and click the search button or press enter. See below

Type in the Accession number in the Value field
and then click the Search button

5. The following will now display on the screen. The center section will provide the patient information and the Order code. The lower section will display information in regards to the accession.



- If this is the accession that you want to result, click on the “Select” button at the bottom of the screen or press “ALT + “S”. You will now enter into the actual result screen for the direct exam.
- Microbiology Result Screen display:



Micro Result Entry Screen is divided into three components.

- a. Section 1: Patient Demographics
Located on top of screen, includes patient current location (Curr Evnt Loc)
 - b. Section 2: Accession number information, includes
 - 1) Collect/Receive date/time
 - 2) Source (Spec Desc, Spec Req)
 - 3) Report: This provides the status of the report, i.e., Prelim or Final
 - 4) Patients location at time of order is noted here (Ord Loc)
 - 5) Attending physician
 - 6) Order Mod: Displays the priority code for the order.
 - c. Section 3: Resulting
This is the result section. Note, that on the top of this section there are tabs **bolded**, denoting that they are associated with the order.
 - 1) **Misc Update** - if you click on this tab, the specimen description and special request display. If you need to change the specimen description you can do it from here.
 - 2) **Direct Exam** - This tab is where you would result a gram stain, wet prep or malaria smear.
8. Resulting
- a. To result, your cursor must be in the result field.
 - b. Make sure your keyboard is on (keyboard icon will say “on” with a green dot). This allows you to use the preprogrammed keyboard codes. If the keyboard is not displayed on the screen you can press “F8” to pull it up. You can move the keyboard on the screen to any location that you wish by moving your mouse to the banner section of the keyboard (“Direct Exam Result/Modifier keys”). Hold left mouse button down and “drag and drop” the keyboard where you want it.
 - c. The drop down menu labeled ‘keyboard’ will default to the correct keyboard code for the test to be resulted.
 - The MDE keyboard is used to result Gram stains.
 - The DES2 keyboard is used to result Malaria.
 - The WPDE keyboard is used to result wet preps.

See addenda for detailed keyboard descriptions.

To access a keyboard other than the default, open the drop down menu and select from the list.

Gram Stain Keyboard (MDE):



Malaria Keyboard (DES2):**Wet Prep(WPDE):**

- d. Enter your first observation either by using the actual keyboard **or** point and click the keyboard displayed on the screen.

Example: If you want to enter “few WBCs”, select the number **2** key and then select the **W** key. Note that both entries display on the same observation line but are linked together.

Enter one observation (**quantity + organism or cell type**) per observation line (always enter the quantity first.)

The translation of what you entered displays on the right hand side of the screen.


- e. To get to the next Observation, use the tab key to tab over to the next result field or use your mouse or click on the next line.
- f. To enter English Text codes or Free Text:
- 1) English Text - Enter “;” (semicolon) in an observation line. This turns off the keyboard codes (keyboard icon says OFF) and brings up the ellipsis box . Click on the ellipsis to bring up the Select Text Code screen and enter the correct English text code.
 - 2) Free Text - Enter “;” (semicolon) to turn off the keyboard codes and then a second “;” to allow Free Texting. Enter result.
 - 3) When you are finished, click on another result field to turn keyboard entry back on and continue result entry.
- g. After entering all your observations, you must save your results as a preliminary or final result.
- 1) To save the direct exam as a preliminary:

Click on **save** at the bottom of the screen or press “**ALT**”+“**S**”. Do NOT select the final key (/) from your keyboard or the keyboard displaying on the screen.

- 2) To finalize the direct exam:
 - Make sure the cursor is in a blank result field, then select the final key “/” from your keyboard or the keyboard displaying on the screen.
 - i. A message will display on the screen stating the direct exam has been finalized and the status will change on the second section of the screen.
 - ii. Select ‘**Save the results**’ or press “**ALT**”+“**S**” Before exiting.
 - iii. If you attempt to exit before saving, a message will prompt: “**You have unsaved data. Do you want to discard it?**” If you select “Yes” then everything that you entered will be lost.

Notes:

1. If you have entered a direct exam observation and you want to change it
 - a. Press the DELETE key while in the result field that contains the observation results that you want to delete. Example, if you entered in the wrong quantity or organism.
 - b. If you want to delete everything on an observation line then click on “Delete Observation [J]” or “ALT”+“J”



2. TAB headers. Note that there is an underscore associated with the Tab headers for Direct Exam and Misc Updates. The underscore lets you know that you can use the ALT key plus the underscored letter to navigate over to that section/area.
3. If you try to enter an observation and the following message displays:



Check the CAPs Lock on your keyboard. It may be off.

4. Once you have verified a result you return to the screen where you can select the next patient or accession number to result. A "Previous" button now appears next to the "Select" button. This will take you back to the last patient or accession number you worked on depending on whether you are under Patient ID or Acc Number. This will reopen the patient in the Misc Updates tab, to return to the Direct Exam screen click the Direct Exam tab (see above).
5. If a Direct Exam is finalized in error you can return to the Result entry screen and enter "/" to unfinalize the results. When corrected results are entered, press "/" again to re-finalize.

6. RELATED DOCUMENTS

None

7. REFERENCES

Sunquest Microbiology Admin Manual

8. REVISION HISTORY

| Version | Date | Reason for Revision | Revised By | Approved By |
|---------|---------|--|------------|-------------|
| | | Supersedes SOP LIS047.001 | | |
| 000 | 6/6/13 | Section 5: In step 8.d - change order for reporting an observation by replacing (organism or cell type + quantity) with (quantity + organism or cell type) Notes 1.b - update Result screen print. | M Sabonis | R Master |
| 001 | 3/20/17 | Header: add other sites Section 5: update keyboard screen shots in step 8.c Section 9: update Appendix A & B Footer: version # leading zero's dropped due to new EDCS in use as of 10/7/13 | M Sabonis | R Master |
| 2 | 11/6/20 | Header: changed WAH to WOMC | L Barrett | R Master |
| 3 | 8/1/22 | Section 4 Definition removed- n/a Section 7 Reference-removed version # of user manual Section 8c Added updated screen shot of Wet Prep keyboard Added appendix C- Keyboard translation for Wet preps | M Sabonis | R SanLuis |
| 3 | 8/1/22 | Header: Changed site to All Laboratories Footer: Changed prefix to AHC | D Collier | R SanLuis |

9. ADDENDA AND APPENDICES

- A: Gram Stain and Wet Prep Keyboard Codes
- B: Malaria Keyboard Codes
- C: Wet Prep Codes

Appendix A Gram Stain and Wet Prep Keyboard Codes

| KEY | TRANSLATION | CODE |
|-----|---|---------|
| / | <finalize> | (fnl) |
| — | <do not report> | HIDE |
| ; | <free text> <English text codes> | (other) |
| F | Branching gram positive rods | BRPR |
| J | chains | CHAIN |
| Q | Clue cells | CLUE |
| K | clusters | CLUST |
| ? | Corrected report, previously reported as | PREVR |
| E | Epithelial cells | EPIT |
| 2 | Few | FEW |
| D | Gram negative diplococci | GNDC |
| C | Gram negative cocci | GNC |
| M | Gram negative coccobacilli | GNCB |
| I | Gram negative rods | GNR |
| H | Gram positive cocci | GPC |
| G | Gram positive coccobacilli | GPCB |
| U | Gram positive rods | GNR |
| S | Gram positive rods resembling diptheroids | GPRD |
| V | Gram variable coccobacilli | GVCB |
| B | Gram variable rod | GVR |
| 7 | Hyphal elements seen | HYPH |
| 4 | Many | MANY |
| 3 | Moderate | MOD |
| N | Negative | NEG |
| 5 | No cell or organisms seen | NOCO |
| O | No organisms seen | NOS |
| 6 | No Trichomonas, yeast or clue cells seen | NTY |
| L | pairs | PAR |
| P | Positive | POSIT |
| 1 | Rare | RARE |
| R | RBC's | RBCP |
| Z | The sputum specimen submitted contains 25 or more squamous epithelial cells per low power field and is unacceptable for culture due to oropharyngeal contamination. Please resubmit another specimen if clinically indicated. | SPUU |
| T | Trichomonas | TRIC |
| W | WBC's | WBC |
| Y | Yeast | YST |
| 8 | Yeast with pseudohyphae | YPSU |

Appendix B Malaria Keyboard Codes

| KEY | MALARIA KEYBOARD TRANSLATION | CODE |
|-----|--|--------|
| M | Parasitemia | INF2 |
| / | <FINAL> | (fnl) |
| N | negative | NEG |
| B | No parasites seen. One set of blood films cannot exclude the diagnosis of malaria. | NMAL1 |
| A | Plasmodium falciparum | PLAF |
| D | Plasmodium malariae | PLAM |
| R | Plasmodium species not P. falciparum | PLSNF |
| E | Plasmodium vivax/ovale | PLAVO |
| F | Plasmodium ovale | PLAO |
| G | Plasmodium species | PLA |
| I | Unable to identify species due to low parasitemia. | UNABLE |
| S | Plasmodium vivax | PLAV |
| L | Blood Flagellate | BFLAG |
| P | Positive | POSIT |
| T | Microfilaria | MFIL |
| K | Plasmodium species, unable to identify species | PSUID |
| H | Presumptive positive, confirmation and identification to follow | PMAL1 |
| ? | Corrected report, previously reported as | PREVR |
| J | Thin smear presumptive negative, thick smear and final report to follow | NMLP1 |

Appendix C**Wet Prep Keyboard Codes**

| KEY | TRANSLATION | CODE |
|-----|--|---------|
| / | <finalize> | (fnl) |
| — | <do not report> | HIDE |
| ; | <free text> <English text codes> | (other) |
| 6 | No Trichomonas, yeast or clue cells seen | NTY |
| E | Trichomonas present | TRP |
| D | Yeast present | YSTP |
| C | Clue cells present | CSP |
| U | Trichomonas not seen | TRKN |
| K | Yeast not seen | YSTN |
| M | Clue cells not seen | CNS |
| ? | Corrected report, previously reported as | PREVR |