



TRAINING UPDATE

Lab Location: All sites
Department: Urinalysis

Date Distributed: 7/23/24
Due Date: 8/12/24
Implementation: **Immediately**

DESCRIPTION OF PROCEDURE REVISION

Name of procedure:
<p>SOPs:</p> <ol style="list-style-type: none"> 1. Urinalysis, Multistix 10 SG Reagent Strips (AHC.U11) 2. Urinalysis, Clinitek Status Plus (AHC.U948) 3. Urinalysis, Clinitek Advantus (GEC.E102)Attached <p>Forms:</p> <ol style="list-style-type: none"> 1. AG.F531 Clinitek Advantus QC Log Attached 2. AG.F370 Clinitek Status QC Log
Description of change(s):
<p>All 3 SOPs: Section 6.3 QC Frequency: Added: “and when opening a new bottle of test strips”.</p> <p>QC forms: 5 lines were added to the QC logs to provide space to record QC when opening a new canister of the same lot of multistix 10 SG test strips.</p>

Document your compliance with this training update by taking the quiz in the MTS system.



Clinitek Advantus QC Log

- Germantown Emergency Center
- Fort Washington Medical Center

Month: _____
 Control Level: _____ I - High Abnormal
 Reagent Strip: _____ Multistix 10 SG

Year : _____
 Lot #: _____
 Lot #: _____

Exp Date: _____
 Exp Date: _____

Day			Glu	Bili	Ket	SG	Blood	PH	Pro	Urob	Nit	Leuk	Color	Clarity	RBCs	WBCs	TechCode
√ to verify lot	↓	Range:															
1																	
2																	
3																	
4																	
5																	
6																	
7																	
8																	
9																	
10																	
11																	
12																	
13																	
14																	
15																	
16																	
17																	
18																	
19																	
20																	
21																	
22																	
23																	
24																	
25																	
26																	
27																	
28																	
29																	
30																	
31																	

Example

New Multistix Canister (same lot) Opened: (Date)

Note: Reconstitute with 15 mL distilled H2O. Once reconstituted stable for 7 days.

Weekly Review:	Weekly Review:	Weekly Review:
Weekly Review:	Weekly Review:	Monthly Review:

GEC.E 102 Urinalysis, Clinitek Advantus

Copy of version 5.0 (approved and current)

Last Approval or
Periodic Review Completed 7/9/2024

Uncontrolled Copy printed on 7/23/2024 11:45 AM

Next Periodic Review
Needed On or Before 7/9/2026

Printed By Demetra Collier (110199)

Organization Adventist HealthCare

Effective Date 7/9/2024

Approval and Periodic Review Signatures

Type	Description	Date	Version	Performed By	Notes
Approval	Lab Director	7/9/2024	5.0	<i>Nicolas Cacciabeve MD</i> Nicolas Cacciabeve	
Approval	Core lab approvals	7/8/2024	5.0	<i>Robert SanLuis</i> Robert SanLuis	
Approval	Lab Director	6/5/2023	4.0	<i>Nicolas Cacciabeve MD</i> Nicolas Cacciabeve	
Approval	Core lab approvals	6/5/2023	4.0	<i>Robert SanLuis</i> Robert SanLuis	
Approval	Lab Director	4/27/2022	3.0	<i>Nicolas Cacciabeve</i> Nicolas Cacciabeve	
Approval	Core lab approvals	4/27/2022	3.0	<i>Robert SanLuis</i> Robert SanLuis	
Periodic review	Lab Service director	6/6/2022	2.0	<i>Robert SanLuis</i> Robert SanLuis	
Approval	Lab Director	7/13/2020	2.0	<i>Nicolas Cacciabeve</i> Nicolas Cacciabeve	
Approval	Core lab approvals	7/10/2020	2.0	<i>Robert SanLuis</i> Robert SanLuis	
Approval	QA approval	7/9/2020	2.0	Leslie Barrett	
Approval	Lab Director	6/16/2020	1.0	<i>Nicolas Cacciabeve</i> Nicolas Cacciabeve	
Approval	Core lab approvals	6/16/2020	1.0	<i>Robert SanLuis</i> Robert SanLuis	
Approval	QA approval	6/10/2020	1.0	Leslie Barrett	

Version History

Version	Status	Type	Date Added	Date Effective	Date Retired
---------	--------	------	------------	----------------	--------------

5.0	Approved and Current	Major revision	6/28/2024	7/9/2024	Indefinite
4.0	Retired	Major revision	5/31/2023	6/5/2023	7/9/2024
3.0	Retired	Major revision	4/15/2022	6/14/2022	6/5/2023
2.0	Retired	Major revision	7/9/2020	7/15/2020	6/14/2022
1.0	Retired	Initial version	6/10/2020	6/16/2020	7/15/2020

Linked Documents

- AG.F 531 Clinitek Advantus QC Log
- AG.F 532 Clinitek Advantus Daily Maintenance Log

Uncontrolled copy
Current as of 7/23/2024 11:45 AM

Adventist HealthCare

Site: Germantown Emergency Center, Fort Washington Medical Center

Title: **Urinalysis, Clinitek Advantus**

Technical SOP

Title	Urinalysis, Clinitek Advantus	
Prepared by	Demetra Collier	Date: 6/5/2020
Owner	Robert SanLuis	Date: 6/5/2020

Laboratory Approval		Local Effective Date:
Print Name and Title	Signature	Date
<i>Refer to the electronic signature page for approval and approval dates.</i>		

TABLE OF CONTENTS

1.	Test Information.....	2
2.	Analytical Principle	2
3.	Specimen Requirements.....	3
4.	Reagents.....	4
5.	Calibrators/Standards.....	5
6.	Quality Control	5
7.	Equipment And Supplies	6
8.	Procedure	7
9.	Calculations.....	8
10.	Reporting Results And Repeat Criteria.....	8
11.	Expected Values.....	12
12.	Clinical Significance.....	13
13.	Procedure Notes	14
14.	Limitations Of Method	14
15.	Safety	16
16.	Related Documents	16
17.	References.....	16
18.	Revision History	17
19.	Addenda	17

1. TEST INFORMATION

Assay	Method/Instrument	Test Code
Urinalysis	Clinitek Advantus	UAI
Urinalysis with reflex to culture	Clinitek Advantus	UAIRX

Synonyms/Abbreviations
Urine Macroscopic, UA, UA with reflex to culture

Department
Urinalysis

2. ANALYTICAL PRINCIPLE

The Clinitek Advantus is a reflectance spectrophotometer that analyzes color and intensity of light reflected from the reagent areas on the Multistix 10-SG and reports the results in clinically meaningful units.

- a. **Protein:** This test is based on the protein-error-of-indicators principle. At a constant pH, the development of any green color is due to the presence of protein. Colors range from yellow for "Negative" through yellow-green and green to green-blue for "Positive" reaction.
- b. **Occult Blood:** This test is based on the peroxidase-like activity of hemoglobin, which catalyzes the reaction of diisopropylbenzene dihydroperoxide and 3,3',5,5'-tetramethylbenzidine. The resulting color ranges from orange through green; very high levels of blood may cause the color development to continue to blue.
- c. **Leukocytes:** Granulocytic leukocytes contain esterases that catalyze the hydrolysis of the derivatized pyrrole amino acid ester to liberate 3-hydroxy-5-phenyl pyrrole. This pyrrole then reacts with a diazonium salt to produce a purple product.
- d. **Nitrite:** This test depends upon the conversion of nitrate to nitrite to by action of Gram negative bacteria in the urine. At the acid pH of the reagent area, nitrite in the urine reacts with p-arsanilic acid to form a diazonium compound. This diazonium compound in turn couples with 1,2,3,4-tetrahydrbenzo(h)quinolin-3ol to produce a pink color.
- e. **Glucose:** This test is based on a double sequential enzyme reaction. One enzyme, glucose oxidase, catalyzes the formation of gluconic acid and hydrogen peroxide from the oxidation of glucose. A second enzyme, peroxidase, catalyzes the reaction of hydrogen peroxide with a potassium iodide chromogen to oxidize the chromogen to colors ranging from green to brown.
- f. **Ketone:** This test is based on the development of colors ranging from buff-pink, for a negative reading, to maroon when acetoacetic acid reacts with nitroprusside.
- g. **pH:** The test is based on the double indicator principle that gives a broad range of colors covering the entire urinary pH range. Colors range from orange through yellow and green to blue.

- h. Specific Gravity: This test is based on the apparent pKa change of certain pretreated polyelectrolytes in relation to ionic concentration. In the presence of an indicator, colors range from deep blue-green in urine of low ionic concentration through green and yellow-green in urines of increasing ionic concentration.
- i. Bilirubin: This test is based on the coupling of bilirubin with diazotized dichloroaniline in a strongly acid medium. The color ranges through various shades of tan.
- j. Urobilinogen: This test is based on a modified Ehrlich reaction in which p-diethylaminobenzaldehyde in conjunction with a color enhancer reacts with urobilinogen in a strongly acid medium to produce a pink-red color.

3. SPECIMEN REQUIREMENTS

3.1 Patient Preparation

Component	Special Notations
Fasting/Special Diets	N/A
Specimen Collection and/or Timing	Normal procedures for collecting urine may be used for samples to be analyzed by this method. Transfer contents to Urine Collection Kit to better preserve the sample.
Special Collection Procedures	A first-morning specimen is preferred but random collections are acceptable.
Other	If Urine Collection Kit is not used, submit to Laboratory within 2 hours of collection.

3.2 Specimen Type & Handling

Criteria	
Type -Preferred -Other Acceptable	Urine, freshly voided None
Collection Container	Clean or sterile container
Volume - Optimum - Minimum	12 mL 1 mL
Transport Container and Temperature	Urine Collection Kit (Urine Analysis Preservative Tube preferred) or container at room temperature. *If order is UAIRX then specimen must be placed/received in a marble and gray collection tubes
Stability & Storage Requirements	Room Temperature: 24 hours in Urine Analysis Preservative Tube 2 hours for other containers
	Refrigerated: 24 hours
	Frozen: Unacceptable
Timing Considerations	Test the urine within two hours after voiding, sooner if testing for bilirubin or urobilinogen.

Adventist HealthCare

Site: Germantown Emergency Center, Fort Washington Medical Center

Title: **Urinalysis, Clinitek Advantus**

Criteria	
Unacceptable Specimens & Actions to Take	Specimens that are unlabeled, improperly labeled, or those that do not meet the stated criteria are unacceptable. Request a recollection and credit the test with the appropriate LIS English text code for “test not performed” message. Examples: Quantity not sufficient-QNS; Wrong collection-UNAC. Document the request for recollection in the LIS.
Compromising Physical Characteristics	If specimen refrigerated, let it return to room temperature before testing. The container should allow for complete dipping of all reagent strip areas.
Other Considerations	After testing, samples will be held until the next successful QC performance.

NOTE: Labeling requirements for all reagents, calibrators and controls include: (1) Open date, (2) Substance name, (3) Lot number, (4) Date of preparation, (5) Expiration date, (6) Initials of tech, and (7) Any special storage instructions. Check all for visible signs of degradation.

4. REAGENTS

The package insert for a new lot of kits must be reviewed for any changes before the kit is used. A current Package Insert is included as a Related Document.

4.1 Reagent Summary

Primary Reagent	Supplier & Catalog Number
Multistix 10 SG Reagent Strips	Siemens Reagent Strips Cat. No. 2161

4.2 Reagent Preparation and Storage

Reagent	Multistix 10 SG Reagent Strips
Container	Plastic Bottle
Storage & Stability	<p>Store at temperatures between 15-30°C.</p> <ul style="list-style-type: none"> All unused strips must remain in the original bottle. Transfer to any other container may cause reagent strips to deteriorate and become un-reactive. Do not use strips after the expiration date printed on the original bottle. Do not store the bottle in direct sunlight and do not remove the desiccant from the bottle. Never leave the container uncapped.
Preparation	None

Adventist HealthCare

Site: Germantown Emergency Center, Fort Washington Medical Center

Title: **Urinalysis, Clinitek Advantus****5. CALIBRATORS/STANDARDS**

Calibration is performed automatically each time a Reagent Strip is analyzed.

6. QUALITY CONTROL**6.1 Controls Used**

Controls	Supplier and Catalog Number
Human Urinalysis Control level I	KOVA-Trol™ HYCOR® Cat. No. 91017
Human Urinalysis Control level II	KOVA-Trol™ HYCOR® Cat. No. 87128
Human Urinalysis Control Level III	KOVA-Trol™ HYCOR® Cat. No. 87328

6.2 Control Preparation and Storage

Control	Level I Urine control
Preparation	Reconstitute the vial of control with exactly 15 mL of Reagent Grade water. Allow the reconstituted material to stand at room temperature for 15 minutes and gently rotate the bottle intermittently until all of the material has dissolved.
Storage/Stability	Once reconstituted, the controls remain stable for 7 days at 2-8°C in its original capped vial.

Control	Level II and Level III Urine controls
Preparation	Reconstitute each vial of control with exactly 60 mL of Reagent Grade water. Allow the reconstituted material to stand at room temperature for 15 minutes and gently rotate the bottle intermittently until all of the material has dissolved.
Storage/Stability	Once reconstituted, the controls remain stable for 7 days at 2-8°C in its original capped vial.

6.3 Frequency

All three levels of Human Urinalysis Control are tested once per day **and when opening a new bottle of test strips.**

The analyzer will prompt for QC after 24 hours.

Daily QC Procedure:

1. From the HOME page select MENU.
2. Select QC.
3. Enter the QC ID and press ENTER.
4. Dip QC and place on the platform
5. Repeat steps 1-4 for each level.

6.4 Tolerance Limits and Criteria for Acceptable QC

All QC Values must be within acceptable limits listed in manufacture's package insert.

IF the result is ...	THEN...
not acceptable	<ul style="list-style-type: none"> • Verify it is the correct control/reagent. • Verify the control/reagent has not expired. • Check for technical/clerical errors. • Visually inspect the condition of the control/reagent. • Inspect the instrument status, do maintenance and troubleshoot. • Repeat the QC test. • Notify the Supervisor if these results are not acceptable.

6.5 Documentation

- Save the instrument printed paper. Record results on "Clinitek Advantus QC Log", located in Urinalysis Quality Control binder.
- Quality control records are reviewed daily at the bench, weekly by the Group Lead or designee, and monthly by the Supervisor/Manager or designee.
- Refer to complete policies and procedures for QC documentation and for record retention requirements in the Laboratory QC Program.

6.6 Quality Assurance Program

- Each new lot number of reagent or new shipment of the same lot of reagent must be tested with external control materials and previously analyzed samples. Performance of the new lot must be equivalent to the previous lot.
- Training must be successfully completed and documented prior to performing this test. This procedure must be incorporated into the departmental competency assessment program.
- The laboratory participates in CAP proficiency testing. All proficiency testing materials must be treated in the same manner as patient samples.
- Consult the Laboratory QC program for complete details.

7. EQUIPMENT and SUPPLIES**7.1 Assay Platform**

Clinitek Advantus

7.2 Equipment

- Centrifuge, 1600 RPM
- Refractometer

7.3 Supplies

- Disposable pipettes
- Plastic Conical Urinalysis tubes

8. PROCEDURE

NOTE: For all procedures involving specimens, buttoned lab coats, gloves, and face protection are required minimum personal protective equipment. Report all accidents to your supervisor.

8.1	Test Run
1.	Verify that lot number and expiration date stored in the instrument matches the lot number of Multistix SG10 in use. Update the Lot number and expiration date of the Mutistix SG whenever a new lot is started. <ul style="list-style-type: none"> • Select MENU • Select the “ Primary Strip Lot Number” to change the lot number • Enter the lot # and expiration date, ENTER • Verify that the information entered is correct • Return to Home menu by selecting the Home Icon (Advantus will save new lot information).
2.	Select ID and Scan or enter patient’s accession number.
3.	Select the color and clarity description for each specimen. Use “OTHER” for colors not listed. If “OTHER” is selected, the result will hold in DI. Use DI to enter your result using “insert coded entry”. See Addendum B
4.	Completely immerse all reagent areas of a Multistix 10 SG Reagent Strip in fresh, well-mixed, un-centrifuged urine.
5.	Immediately remove the Reagent Strip. While removing, slowly run the edge of the entire length of the Reagent Strip against the side of the urine container to remove excess urine. Do not blot the edge of the strip against a paper towel.
6.	Place the Reagent Strip, with reagent areas facing up , onto the strip supports of the strip loading station.
7.	The presence of the reagent strip is detected as soon as it is placed on the loading station. The push bar moves the strip along the loading station to the read area.

Some medications cause urine to become abnormally colored (GREEN, AMBER, ORANGE or PINK) and the Clinitek Advantus will report false positive results. For urines that are abnormally colored:

8.2	Color Interference
1.	Run the strip through the Clinitek Advantus
2.	Verify the specific gravity by manual refractometer (rounding to the nearest .005). Report the results of the manual refractometer.

Adventist HealthCare

Site: Germantown Emergency Center, Fort Washington Medical Center

Title: **Urinalysis, Clinitek Advantus**

8.2	Color Interference
3.	Tests that are NEGATIVE on the Clinitek Advantus can be reported as negative.
4.	Report the Color and Clarity as you see it.
5.	Enter the comment COLINT, which expands out to “Results not reported due to color interference”, for the remainder of the tests.
6.	Perform a microscopic exam on all abnormally colored urines.
8.3	Bloody Urines
1.	Measure the specific gravity by manual refractometer (rounding to the nearest .005). Report the results of the manual refractometer.
2.	Report the Color as BLOODY and the Clarity as you see it.
3.	Centrifuge the specimen. Pour the supernatant into a plastic conical urinalysis tube
4.	Perform dipstick testing on the supernatant and run through the Clinitek Advantus.
5.	Report the remaining results of the supernatant from the Clinitek Advantus (GLU, BIL, KET, PH, PRO, URO, NIT, and LEU).
6.	Perform a microscopic exam on the sediment.
8.4	RESEND or REPRINT a result
1.	From the HOME screen select MENU
2.	Select MEMORY.
3.	Select result to recall
4.	Select RESEND (a circle with an arrow icon) or REPRINT (a printer icon).
8.5	Instrument Maintenance
1.	Refer to Addendum A for maintenance instructions.
2.	Record maintenance on the appropriate log.

NOTE: In the event that the test system becomes inoperable, notify supervision or designee for further direction. Patient specimens must be stored in a manner that maintains the integrity of the specimen.

9. CALCULATIONS

Not applicable

10. REPORTING RESULTS AND REPEAT CRITERIA

10.1 Interpretation of Data

Macroscopic Analysis

Adventist HealthCare

Site: Germantown Emergency Center, Fort Washington Medical Center

Title: **Urinalysis, Clinitek Advantus**

Test	Report As
Color	Yellow, Orange, Pink, Green, Amber, Brown, Bloody, Dark Yellow, Straw
Appearance	Clear, Cloudy, Slightly Cloudy, Turbid
Specific Gravity	1.005 – 1.030 (in increments of 0.005)
pH	5.0 – 9.0 (in increments of 0.5)
Glucose	Negative, Trace, 1+, 2+, 3+, 4+
Bilirubin	Negative, 1+, 2+, 3+
Urobilinogen	0.2, 1.0, 2.0, 4.0, 8.0
Ketone	Negative, Trace, 1+, 2+, 3+, 4+
Occult Blood	Negative, Trace, 1+, 2+, 3+
Protein	Negative, Trace, 1+, 2+, 3+, 4+
Nitrite	Negative, Positive
Leukocytes Esterase	Negative, Trace, 1+, 2+, 3+

Microscopic Analysis

Power Field Instructions for Microscopy	
High Power Field (HPF)	Low Power Field (LPF)
RBCs and WBCs	Squamous Epithelial Cells
Renal & Transitional Epithelial Cells	All Casts
Bacteria / Yeast / Crystals	Mucus

Test	# seen	LIS translation
WBC (average # / HPF)	0 - 2	O0
	3-5	O3
	6-10	O6
	11-20	O11
	21-100	O21
	>100	TNTC
RBC (average # / HPF)	0 - 2	O0
	3-5	O3
	6-10	O6
	11-20	O11
	21-100	O21
	>100	TNTC

Adventist HealthCare

Site: Germantown Emergency Center, Fort Washington Medical Center

Title: **Urinalysis, Clinitek Advantus**

Test	# seen	LIS translation
Epithelial (average # / LPF)	0 - 2	Rare
	3-5	Occasional
	6-10	1+
	11-20	2+
	21-100	3+
	> 100	4+
Casts (average # / LPF)	0-1	O01
	2-5	O2
	6-10	O6
	11-20	O11
	21-100	O21
	TNTC	TNTC
Bacteria / HPF	None seen	Negative
	Few	1+
	Small	2+
	Moderate	3+
	Large	4+
	Packed	TNTC

Only report these analytes if seen during microscopic review:		
Test	# seen	LIS translation
Transitional Epithelial Cells (average # / HPF)	1-2	Rare
	3-5	Occasional
	6-10	1+
	11-20	2+
	21-100	3+
	> 100	4+
Renal Epithelial Cells (average # / HPF)	1-2	Rare
	3-5	Occasional
	6-10	1+
	11-20	2+
	21-100	3+
	> 100	4+
Crystals (average # / HPF)	1-5	Few
	6-10	1+
	11-20	2+
	>21	3+
Mucus / LPF	Occasional	Occasional
	Small	1+
	Moderate	2+
	Large	3+
	Packed	4+

Only report these analytes if seen during microscopic review:		
Test	Test	Test
Yeast / HPF	Occasional	Occasional
	Small	1+
	Moderate	2+
	Large	3+
	Packed	4+
Trichomonas	No quantitation – report “present” if seen	
Enterobius Vermicularis	No quantitation – report “present” if seen. Consult with pathologist prior to releasing results.	
Schistoma Haematobium	No quantitation – report “present” if seen. Consult with pathologist prior to releasing results.	
Oval Fat Bodies	No quantitation – report “present” if seen	

10.2 Rounding

N/A

10.3 Units of Measure

Urobilinogen EU/dL

10.4 Clinically Reportable Range (CRR)

N/A

10.5 Review Patient Data

Technologist must check for unusual patterns, trends, or distributions in patient results (such as an unusually high percentage of abnormal results). Resolve any problems noted before issuing patient reports. Repeat patient samples with other methodologies if necessary.

10.6 Repeat Criteria and Resulting

Test	If the result is...	Then...
Bilirubin	1+, 2+ and 3+	The ETC (English Text Code) of UPPB will be appended to the result by LIS. The code translates to “Presumptive positive bilirubin. Consider confirmation by serum bilirubin if clinically indicated.”

- Microscopic Exam:
 1. Review the results. The following macroscopic abnormalities trigger a microscopic exam:
 - a. Occult Blood: any positive
 - b. Protein: any positive

Adventist HealthCare

Site: Germantown Emergency Center, Fort Washington Medical Center

Title: **Urinalysis, Clinitek Advantus**

- c. Nitrite: any positive
 - d. Leukocytes: any positive
 - e. Clarity (Character): Slightly Cloudy, Cloudy or Turbid
2. Centrifuge the specimens that require a microscopic exam at 1600 RPM for 5 minutes.
 3. Refer to procedure “Microscopic Examination of Urine” for instructions on performing microscopic examination of urine.
 4. Enter Microscopic results using DI (See Addendum B).
- Test UMM (UR Microscopic Added?) is added to the Urine Chemistry group. If criteria to perform a microscopic is met then UMM is resulted with TADD (Test Added). If not met then it is resulted with NIND (Not Indicated)
 - Urinalysis with reflex to culture (UAIRX): Any of the following macroscopic or microscopic abnormalities will trigger a reflex to Urine culture (XURNC) by Sunquest(LIS):
 - a. Nitrite: positive
 - b. Leukocyte: 2+, 3+
 - c. WBC: >10

11. EXPECTED VALUES

11.1 Reference Ranges

Glucose	Negative
Bilirubin	Negative
Ketones	Negative
Specific gravity (SG)	1.005 – 1.030
Occult Blood	Negative
pH	5.0 – 9.0
Protein	Negative
Urobilinogen (URO)	0.2 – 1.0 EU/dL
Nitrite	Negative
Leukocyte	Negative
Color	Yellow
Clarity	Clear

11.2 Critical Values

None established

11.3 Standard Required Messages

None established

12. CLINICAL SIGNIFICANCE

The strips are intended for use in at-risk patient groups to assist diagnosis in the following areas: Kidney function, urinary tract infections, carbohydrate metabolism and liver function. The strips also measure physical characteristics, including acid-base balance and urine concentration. Test results can be used along with other diagnostic information to rule out certain disease states and to determine if microscopic analysis is needed.

Protein: In normal urine, less than 150 mg of total protein is excreted per day. Clinical proteinuria is indicated at greater than 500 mg of protein per day. Positive results may also indicate tubular or overflow proteinuria in the absence of any glomerular abnormality or proteins of renal origin that may be excreted during infection. Urinary protein excretions can be temporarily elevated in the absence of renal abnormality by strenuous exercise, orthostatic proteinuria, dehydration, urinary tract infections, and acute illness with fever.

Occult Blood: Normally, no hemoglobin is detectable in urine. Occult blood occurs in urine as intact erythrocytes and hemoglobin, which can occur during urological, nephrological and bleeding disorders. Small amounts of blood are sufficiently abnormal to require further investigation. The significance of the Trace reaction may vary among patients, and clinical judgment is required for assessment in an individual case.

Leukocytes: Normal urine specimens generally yield negative results. An increase in leukocytes is an indication of pyuria and is found in nearly all diseases of the kidney and urinary tract; however, pyuria may often be present in non-infection conditions. A strip result of small or greater is a useful indicator of infection. Trace results may be of questionable clinical significance; however, Trace results observed repeatedly may be clinically significant.

Nitrite: Normally no nitrite is detectable in urine. Many enteric gram-negative organisms give positive results when their number is greater than 10^5 /mL.

Glucose: Small amounts of glucose are normally excreted by the kidney. These amounts are usually below the sensitivity of this test but on occasion may produce a color between the Negative and the 100 mg/dL color blocks and that is interpreted by the instrument as a positive result. Results at the first positive level may be significantly abnormal if found consistently.

Ketone: Normally, no ketone is detectable in urine. In ketoacidosis, starvation or with other abnormalities of carbohydrate or lipid metabolism, ketones may appear in urine at levels of 10 mg/dL or higher before serum ketone levels are elevated. Clinical judgment is needed to determine the significance of trace results, which may occur during physiological stress conditions such as fasting, pregnancy and frequent strenuous exercise.

pH: The normal pH of urine can range from 4.6 to 8.0. Certain dietary conditions can produce acid or alkaline urines, which can be useful in the treatment of some calculi.

Specific Gravity: The normal SG of urine ranges from 1.001 – 1.035. If the specific gravity of random urine is 1.023 or greater, the concentrating ability of the kidneys can be considered normal.

Bilirubin: Normal adult urine contains about 0.02 mg/dL of bilirubin, which is not detectable by even the most sensitive methods. Even trace amounts of bilirubin are sufficiently abnormal to require further investigation. Since very small amounts of bilirubin may be found in the earliest phases of liver disease, the user must consider whether the sensitivity of Siemens Reagent Strips to bilirubin is sufficient for the intended use.

Urobilinogen: Urobilinogen is normally present in urine at concentrations up to 1.0 mg/dL. A result of 2.0 mg/dL represents the transition from normal to abnormal, and the patient and/or urine specimen should be evaluated further for hemolytic and hepatic disease.

13. PROCEDURE NOTES

- **FDA Status:** Approved/cleared
- **Validated Test Modifications:** None

14. LIMITATIONS OF METHOD

14.1 Analytical Measurement Range (AMR)

N/A

14.2 Precision

N/A

14.3 Interfering Substances

Bloody urine and color interference explained in sections 8.2 and 8.3

For all tests, false positive results and/or false negative results can occur when substances that cause abnormal urine color are present, such as:

- visible levels of blood or bilirubin
- drugs containing dyes
- nitrofurantoin
- riboflavin

14.4 Clinical Sensitivity/Specificity/Predictive Values

Sensitivities listed in the following table depend upon the presence or absence of inhibitory and matrix factors typically found in urine, such as specific gravity and pH.

<i>Test Name</i>	<i>False Positive or Increased values</i>	<i>False Negative or Decreased values</i>
Glucose	Temperature	<ul style="list-style-type: none"> ▪ Ascorbic acid ($\geq 50\text{mg/dL}$) may affect a 75 to 125 mg/dL glucose level ▪ Ketones ($\geq 40\text{mg/dL}$) may affect a 75 to 125 mg/dL glucose level ▪ High specific gravity ▪ Temperature

Adventist HealthCare

Site: Germantown Emergency Center, Fort Washington Medical Center

Title: **Urinalysis, Clinitek Advantus**

Test Name	False Positive or Increased values	False Negative or Decreased values
Bilirubin	<ul style="list-style-type: none"> ▪ Indican (indoxyl sulfate) may impart a yellow-orange to red color on the pad ▪ Metabolites of Lodine (etodolac) 	<ul style="list-style-type: none"> ▪ Ascorbic acid (≥ 25mg/dL). ▪ Urine specimen was more than one hour old (instability of bilirubin). ▪ Contamination with chlorhexidine (found in some skin cleansers)
Ketone	<ul style="list-style-type: none"> ▪ Highly pigmented urines ▪ Large amounts of levodopa (L-dopa) metabolites ▪ Compounds that contain sulfhydryl groups 	
Specific Gravity	<ul style="list-style-type: none"> ▪ Moderate (100 – 750 mg/dL) quantities of protein ▪ Contamination with chlorhexidine (found in some skin cleansers) 	<ul style="list-style-type: none"> ▪ Highly buffered/alkaline urines
Occult Blood	<ul style="list-style-type: none"> ▪ Oxidizing contaminants (e.g. bleach) ▪ Microbial peroxidase from urinary tract infections 	<ul style="list-style-type: none"> ▪ High specific gravity ▪ Capoten® (Captopril)
pH	<ul style="list-style-type: none"> ▪ Bacterial growth that converts urea to ammonia 	<ul style="list-style-type: none"> ▪ Run-over from the protein reagent pad
Protein	<ul style="list-style-type: none"> ▪ Highly buffered or alkaline urines ▪ Contamination with quaternary ammonium compounds (from some antiseptics and detergents) or Chlorhexidine (found in some skin cleansers) 	
Urobilinogen	<ul style="list-style-type: none"> ▪ Temperature $> 26^{\circ}\text{C}$ (79°F) ▪ p-aminosalicylic acid (PAS) and sulfonamides ▪ p-aminobenzoic acid (PABA) may cause atypical color development 	<ul style="list-style-type: none"> ▪ Temperature $< 22^{\circ}\text{C}$ (72°F) ▪ Formalin
Nitrite		<ul style="list-style-type: none"> ▪ Infections caused by organisms that don't contain reductase ▪ Urine was not in bladder long enough (at least 4 hours) ▪ Absence of dietary nitrate ▪ High specific gravity ▪ Ascorbic acid (≥ 25 mg/dL) may affect a low positive nitrate level (< 0.06 mg/dL nitrate ion)
Leukocytes	<ul style="list-style-type: none"> ▪ Formalin ▪ Temperature $>26^{\circ}\text{C}$ (79°F) 	<ul style="list-style-type: none"> ▪ Elevated glucose ($\geq 3,000$ mg/dL) ▪ High specific gravity ▪ Cephalexin (Keflex®) or Cephalothin (Keflin®) ▪ High concentrations of oxalic acid ▪ Tetracycline ▪ Temperature $<22^{\circ}\text{C}$ (72°F)

Adventist HealthCare

Site: Germantown Emergency Center, Fort Washington Medical Center

Title: **Urinalysis, Clinitek Advantus**

<i>Test Name</i>	<i>False Positive or Increased values</i>	<i>False Negative or Decreased values</i>
Color	<ul style="list-style-type: none"> ▪ Concentration ▪ Food Pigments ▪ Dyes ▪ Blood ▪ Various pathological conditions 	<ul style="list-style-type: none"> ▪ These all can affect negatively as well.

15. SAFETY

Refer to your local and corporate safety manuals and Safety Data Sheet (SDS) for detailed information on safety practices and procedures and a complete description of hazards.

16. RELATED DOCUMENTS

1. Laboratory Quality Control Program
2. Laboratory Safety Manual
3. Safety Data Sheets (SDS)
4. Quest Diagnostics Records Management Procedure
5. Specific Gravity Using the Refractometer, Urinalysis procedure
6. Microscopic Examination of Urine, Urinalysis procedure
7. Clinitek Advantus QC Log (AG.F531)
8. Clinitek Advantus Daily Maintenance Log (AG.F532)
9. Current Allowable Total Error Specifications at http://questnet1.qdx.com/Business_Groups/Medical/qc/docs/qc_bpt_tea.xls
10. Current package insert Multistix 10 SG (manufacturer provides alert when changes are made)

17. REFERENCES

1. Operator's Guide, Siemens Clinitek Advantus, Siemens Healthcare Diagnostics, Inc., revised 8/2013 (*a copy is located on the AHC G drive at LD USERS, GEC, Advantus Operator Guide and Multistix 10 SG package insert*)
2. Package Insert, Siemens Multistix 10 SG, Siemens Healthcare Diagnostics, Inc. revised 7/2017 (*a copy is located on the AHC G drive*)
3. Package Insert, KOVA-Trol™ HYCOR, P/N 91017-09, 10/2016.
4. CLINITEK ADVANTUS Technical Procedure, doc # 035103. National Committee for Clinical Laboratory Standards (NCCLS). Clinical Laboratory Procedure Manuals-3rd Edition (GP2-A3), 1996.

Adventist HealthCare

Site: Germantown Emergency Center, Fort Washington Medical Center

Title: **Urinalysis, Clinitek Advantus****18. REVISION HISTORY**

Version	Date	Section	Reason	Reviser	Approval
1	7/6/20	10.1	Added microscopic tables	L Barrett	R SanLuis
1	7/6/20	10.6	Deleted instruction for pH >8.0	L Barrett	R SanLuis
1	7/6/20	Add. B	Deleted spec. gravity confirmation	L Barrett	R SanLuis
2	3/17/22	1	Added UAIRX(UA with reflex to culture) and updated Synonyms	M Sabonis	R SanLuis
		3.2	Add UAIRX aliquoting to gray tube	M Sabonis	R SanLuis
		10.1	Replaced "Blood" with "Occult Blood"	M Sabonis	R SanLuis
		Addendum B	Replaced Coded Entry with drop down resulting from DI	M Sabonis	R SanLuis
		Addendum B	Added new "REQUIRED ELEMENTS" screen shot for DI Added FMWC UA keyboard	M Sabonis	R SanLuis
		Addendum B	Added auto-release of urine chemistry and updated order of release	M Sabonis	R SanLuis
		Addendum B	Added URTYP description and updated screen shot Added info and screen on process for UAIRX(UA with reflex to culture)	M Sabonis	R SanLuis
		Header	Added site FWMC	D Collier	R SanLuis
3	5/25/23	Addendum C	Added Macroscopic and Microscopic result entry when instrument down	M Sabonis	R SanLuis
		Addendum B	Addend #5 Urine culture message	M Sabonis	R SanLuis
4	6/18/24	Section 6.3	Added: and when opening a new bottle of test strips.	D Collier	R SanLuis

19. ADDENDA

- A. Clinitek Advantus Maintenance
- B. DI (Data Innovations) Information and Actions
- C. Macroscopic and Microscopic result entry when instrument down

Addendum A**Clinitek Advantus Maintenance****A. Daily Maintenance**

1. From the main screen, press the back key until you are at the Ready/Run screen.
2. Turn off the instrument. The on/off switch is located in the lower left, in the rear of the instrument.
3. Remove the push bar by tilting it slightly upwards and pulling straight out.
4. Remove the waste bin and discard the used strips, into the appropriate container. Inspect the liner. If it has any cracks or is extremely dirty it should be replaced.
5. Remove the fixed platform by pulling the entire assembly towards you. Remove the moving table in the same manner.
6. Remove the hold-down plate from the fixed platform by pressing up against the tab at the back of the plate. Then pull the other end from its retaining hole.
7. Clean all parts with water and mild soap. A toothbrush may be used if sediment accumulation is observed.
Note: DO NOT use solvent or alcohol
8. When cleaning the fixed platform, Do NOT wipe across the two white calibration bars. The white calibrator bars should be GENTLY cleaned with water on a cotton-tipped applicator.
9. Rinse and dry all parts with paper towel except the calibrator pads. Use mild soap if necessary. The calibrator pads should be allowed to air dry. Check the white calibration bars for scratches or discoloration. Notify the supervisor/designee if they appear overly scratched.
10. Reinstall the moving table as follows:
 - a) Hold the table with the small rectangle to the back.
 - b) Align the two grooves on the bottom of the table with the edges of the platform on which the table rests.
 - c) Gently push the table in as far as it will go. It must be pushed past a dent in order to be correctly in position.
11. Reinstall the hold-down on the fixed platform.
12. Position the hold-down with the arrow side facing up and the arrow pointing to the back. Place the pin on the front of the hold-down into the hole at the front of the platform. Then align the tab at the back of the hold-down with the slot at the back of the platform and snap the hold-down into place. Make sure the white calibration bars are visible.

13. Reinstall the clean fixed platform by:
 - a) Aligning the two flared grooves on the bottom of the fixed platform with the arms extending from the instrument
 - b) Gently push the platform in as far as it will go. (It must be pushed past a slight dent to be correctly positioned.)
 - c) Ensure the moving platform is correctly positioned.
14. Hold the push bar by its flattened end and, with this end slightly upward; insert the peg on the other end of the bar into the hole in the pusher mechanism. Lower the push bar into place.
15. Clean the display screen with lens paper dampened with water.

Notes:

 - Dry with lens paper. Do NOT use Kimwipes or paper towels as this may scratch the screen.
 - Do NOT put water directly on the screen.
 - Do NOT use bleach
16. Turn the instrument on. The Clinitek will go through a verification check that all parts have been correctly positioned.

Note: If the instrument gives an error (e.g; "table not positioned properly"), refer to the Clinitek Advantus Urine Chemistry Analyzer Operating manual - Troubleshooting section.
17. Run quality control according to section 6 of this procedure
18. Complete the daily Maintenance log sheet to document that maintenance was performed.

Addendum B

DI (Data Innovations) Information and Actions

1. Result Processing

A. Changing Result -Click on the result cell that you want to edit.

- If applicable, a drop down button displays. Click on the drop down button and result options display. Click on result you want to change to.

- (none)	
Color	YELLOW
Clarity	
Protein	AMB - Amber
Bilirubin	BLDY - Bloody
Urobilinogen	BRWN - Brown
pH	DYEL - Dark Yel...
Ketone	GRN - Green
Nitrite	ORNG - Orange
Specific Gravity	RED - Red
Glucose	STRAW - Straw
Leukocytes	YEL - Yellow
Occult Blood	

Example: of drop down displaying result options to select.

B. Positive Bilirubin

If Bilirubin is resulted with 1+, 2+ or 3+, DI will add the English text code **UPPB** to the test comment. UPPB translates to “Presumptive positive bilirubin. Consider confirmation by serum bilirubin if clinically indicated.”

Run Worksheet								
	Test Name ^	Result (1)	Test Statu...	Units (1)	Test Ins...	Error Cod...	Error ...	Test Comm...
- URINE CHEMISTRY								
	Color	YEL	Held for V...		GCA			
	Clarity	CLER	Held for V...		GCA			
	Protein	3+	Held for V...	mg/dL	GCA			
	Bilirubin	1+	Held for V...	mg/dL	GCA			UPPB
	Urobilinogen	4.0	Held for V...	EU/dL	GCA			

2. Performing Manual Microscopy using the Urinalysis Keyboard

A. Before you begin, you must select a “keyboard”:

- Select the “Cell COUNTER” tab
- If at GEC, then, select “GEC UA Keyboard” from the drop-down menu
- If at FWMC, then select "FWMC UA Keyboard" from the drop-down menu

B. Once you have selected your keyboard, right-click on your macroscopic results, and select “Verify with Cell Counter” from the drop down menu.

The Urinalysis keyboard is used to enter the observational results from the manual microscopy. Each test on the keyboard can be resulted by left- clicking on the result field for that test and selecting the drop-down arrow to reveal a list of available results.

Under the 'REQUIRED ELEMENTS' section there are four elements denoted with "???". These **MUST ALWAYS BE REPORTED.**

- REQUIRED ELEMENTS		
RBC_Urine	???	/HPF
WBC_Urine	???	/HPF
Bacteria	???	/HPF
Squamous_Epith...	???	/LPF

- REQUIRED ELEMENTS		
RBC_Urine	???	/HPF
WBC_Urine	???	/HPF
Bacteria	???	/HPF
Squamous_Epith...	???	/LPF
- CAST		
Hyaline_Cast		/LPF
Broad_Cast		/LPF
Cellular_Cast		/LPF
Epithelial_Cast		/LPF
Fatty_Cast		/LPF
Granular_Cast		/LPF
Hemoglobin_Cast		/LPF
RBC_Cast		/LPF
Waxy_Cast		/LPF
WBC_Cast		/LPF
- FORMED ELEMENTS		
Renal_Epithelial		/HPF
Transitional_Epith...		/HPF
Mucus		/LPF
Trichomonas		
Yeast		/HPF
Oval_Fat_Body		/HPF
Enterobius_Vermi...		
Schistosoma_Hae...		
- CRYSTAL		
Ammonium_Biurate		/HPF
Calcium_Carbonat...		/HPF
Calcium_Oxalate_...		/HPF
Calcium_Phospha...		/HPF
Calcium_Sulfate_...		/HPF
Cholesterol_Crystal		/HPF
Cystine_Crystal		/HPF

3. Order of Release

The Urinalysis results consist of 2 to 3 groups in DI. They must be released in DI in a certain order to ensure proper filing into Sunquest. The order is explained below.

Note:

- The Urine Chemistry group is ALWAYS auto-released to Sunquest, unless held on DI
- If you have to manually release the results ALWAYS release the URINE CHEMISTRY first as noted in the header description (see below).

+	Release First - URINE CHEMISTRY
+	AUTOMATED MICROSCOPY
+	MANUAL MICROSCOPY

Results with just Urine Chemistry

- Release the Urine Chemistry group

Results with Urine Chemistry and Manual Microscopy

- Release the Chemistry group
- Release the Manual Microscopy group

To release or reject a group, follow the steps below:

- Right click on the test within the group to be released/rejected and select the appropriate action. Example: If you select "Release URINE CHEMISTRY/Reject Other Runs," DI will release the selected Urine Chemistry group and reject other Urine Chemistry group from a different run

To Reject a test, follow the steps below:

- Right click on the test within the group to be rejected and select the appropriate action. Example: if you select "Reject Result," DI will reject that result. Once rejected, that result can no longer be released from DI.

4. Microscopic Billing

DI will add a billing testcode of Microscopic completed? to the Manual Microscopy group whenever there is a microscopic test done. This test code is resulted with "DONE." This test code must be released together with the rest of the group.

Adventist HealthCare

Site: Germantown Emergency Center, Fort Washington Medical Center

Title: **Urinalysis, Clinitek Advantus**

- AUTOMATED MICROSCOPY	
WBC Urine	06
RBC Urine	00
Squamous Epithelial	2+
Hyaline Cast	001
Bacteria	NEG
Microscopic Completed?	DONE

5. Urine Culture message

The Clinitek Advantus is not able to query the orders that the LIS (Sunquest) has sent to DI. When the results from the Advantus have crossed over to DI, DI is not able to differentiate between UAI or UAIRX orders. As a result, whenever the criteria is met to reflex to a urine culture, DI will always display, “If order UAIRX, XURNC reflexed” to the test that met triggered the rule. If the order is for UAI, this error message can be ignored. Continue to process the urine as UAI. If the order is for UAIRX, once the results are released from DI, the reflex order for XURNC (Urine culture) will be generated and the UR Culture (UCRX) will be resulted with UCADD (Urine culture has been reflexed) by the LIS (Sunquest). If the criteria is not met and the order is for UAIRX, the UR Culture will be resulted with UCNAD (Urine culture reflex not indicated) by the LIS.


Nitrite	NEG	Released	5/20/2023 12:38:00 ...	GCA	
Specific Gravity	1.025	Released	5/20/2023 12:38:00 ...	GCA	
Glucose	NEG	Released	5/20/2023 12:38:00 ...	GCA	
Leukocytes	2+	Released	5/20/2023 12:38:00 ...	GCA	If order=UAIRX, XURNC reflexed...
Occult Blood	2+	Released	5/20/2023 12:38:00 ...	GCA	Perform Microscopic
UR Microscopic Added?	TADD	Released	5/20/2023 12:11:17 ...	GCA	Perform Microscopic,Perform Micro...

Also, during LIS downtime, and the Patient Order Management is in use, there is no need to order UAI or UAIRX to a specimen that gets run on the Clinitek Advantus.

**Collection label for XURNC will print on Sunquest lab printer at GEC and FWMC.


Addendum C**Macroscopic and Microscopic result entry when instrument down**

1. We have two urinalysis orderables.
 - UAI- Urinalysis. UAI may reflex to a UMIC-Microscopic
 - UAIRX- Urinalysis with reflex to culture. UAIRX may reflex to UMIC and/or display message to order a XURNC-Urine culture, if one or the following criteria is met.
 - ULEUK 2+ or greater
 - UNIT Positive
 - UWBC >10
 - Both orderables include test UMM(Urine Microscopy)- this notifies physicians if we are performing a microscopic or not.
 - UAIRX has a test UCRX that will notify physicians if a urine culture has been added or not.
2. Log into Sunquest GUI, select the **Urinalysis Result Entry**. The following information dialog box displays demonstrating the different keyboards options
 - **Choose the keyboard based on the order code, **UAI or UAIRX** then click on OK.




3. **Macroscopic Resulting**, type in the Accession # and press **ENTER**.
 - *****IMPORTANT- MAKE SURE THAT YOUR CAP LOCK IS ALWAYS ON.**

- Tests that **MUST** be resultated are denoted in **YELLOW**.




- Click on test/key you want to result. The result key options display. Click on result key then press ENTER key on PC keyboard OR click on ENTER on Sunquest keyboard. If result keys do not display, checked to see if you Caps Lock is on.


- For Specific gravity and pH a pop-up windows displays for you to enter the numeric values. Enter value then click on OK.




- Once done resulting, select **QA Review** tab to review the results and click on the **SAVE** button to save and file the results.



- When you click on the QA Review Tab, QA checking of the special logic rules behind the scene fire off.
 - If results for mandatory tests are missing, a message will display denoting the test(s) that are missing results. Click on OK then click on the Resulting TAB and enter in results. Once completed then click on the QA tab again.



- If microscopic is indicated then UMIC is automatically ordered and 2 messages displays.




- Test code UMM(Urine Microscopic) is automatically resulted with TADD (has been added).

- If UMIC is NOT indicated then test code UMM (Urine Microscopic) is automatically resulted with NIND (Not indicated) and message displays.




- If UAIRX is ordered and if one or more of the following criteria is met then message displays to order a urine culture.
 - (UNIT) is Positive
 - (ULEUK) is 2+ or greater




- **YOU MUST order urine culture in Sunquest using collect date/time and source of UAIRX.**

- If Urine Culture indicated then UCRX is automatically resulted with UCADD(urine culture has been added)

- Click on OK another message display denoting tests that are **mandatory and must be reported.**



- **Once you click on OK the MICROSCOPIC tests/resulting keys display on the same keyboard.**
****Make sure that your Caps Lock is ON**



- The urine macroscopic results will display on the left side and the microscopic results display on the right side.

4. Resulting urine Microscopic (UMIC)

- There are FIVE (5) tests that are mandatory for each microscopic analysis (these tests are denoted in YELLOW): [they are circled in RED below].
 - (UWBC)White blood cells (UWBC)
 - (URBC)Red blood cells (URBC)
 - (UBACT)Bacteria (UBACT)
 - (USEPI)Squamous Epithelial cells(USEPI)
 - (UCAST)- Must report out Hyaline casts (HYAL)- Note you MUST select UCAST then HYAL then quantity

Adventist HealthCare

Site: Germantown Emergency Center, Fort Washington Medical Center


Title: **Urinalysis, Clinitek Advantus**

- Click on test you want to result. The result key options display. Click on result key then press ENTER key on PC keyboard OR click on ENTER on Sunquest keyboard.
- Resulting casts(UCAST) and/or crystals(UCRYS) consists of two parts
 - Type of cast or crystal then Quantity
****must be resulted in this order to post correctly into Cerner**
 - **NOTE:** You must always report out Hyaline casts
- To append a comment, select the test code, click on the **EDIT/COMMENT** button and enter free text and/or an English text code in the Comment box.
 Note: Use UOMIC to add any applicable observations that are not on the keyboard.
- Click on QA Review tab. To save and file the urine microscopic click on the **SAVE** button. Special rule logic checking will fire off.
- If UAIRX is ordered and if UWBC >10 then message displays to order a urine culture.

Adventist HealthCare

Site: Germantown Emergency Center, Fort Washington Medical Center

Title: **Urinalysis, Clinitek Advantus**



- **YOU MUST order urine culture in Sunquest using collect date/time and source of UAIRX.**
 - If Urine Culture indicated then UCRX is automatically resulted with UCADD(urine culture has been added)
- If Urine Culture is not indicated then UCRX is automatically resulted with NIND(Not indicated)

Uncontrolled copy