



TRAINING UPDATE

Lab Location: SGMC and WOMC
Department: Microbiology

Date Distributed: 4/15/25
Due Date: 4/30/25
Implementation: **Immediately**

DESCRIPTION OF PROCEDURE REVISION

Name of procedure:
Malaria (AHC.M06) and Malaria/Blood Parasite History Form (F289)
Description of change(s):
<p>Malaria specimens are collected in EDTA lavender tubes, instead of by fingerstick. Smears are made by technical staff within 60 minutes. (30 minutes is preferred)</p> <p>A Malaria/Blood Parasite History Form must be completed and accompany the EDTA tube / patient specimen.</p>

Document your compliance with this training update by taking the quiz in the MTS system.

Adventist HealthCare

Site: All Laboratories

Title: **Malaria**

Technical SOP

Title	Malaria	
Prepared by	Ron Master	Date: 5/11/2009
Owner	Vittal Ponraj	Date: 8/11/2023

Laboratory Approval		Local Effective Date:
Print Name and Title	Signature	Date
<i>Refer to the electronic signature page for approval and approval dates.</i>		

TABLE OF CONTENTS

1.	Test Information.....	2
2.	Analytical Principle	2
3.	Specimen Requirements.....	2
4.	Reagents.....	4
5.	Calibrators/Standards	5
6.	Quality Control	5
7.	Equipment And Supplies	6
8.	Procedure	6
9.	Calculations.....	7
10.	Reporting Results And Repeat Criteria.....	7
11.	Expected Values.....	21
12.	Clinical Significance.....	21
13.	Procedure Notes	22
14.	Limitations Of Method	23
15.	Safety	23
16.	Related Documents	23
17.	References.....	23
18.	Revision History	24
19.	Addenda	25

1. TEST INFORMATION

Assay	Method/Instrument	Test Code
Malaria, thick and thin smears	Manual	MAL

Synonyms/Abbreviations
Malaria smear, Malaria ID, Malaria Parasites, <i>Plasmodium</i> species

Department
Microbiology

2. ANALYTICAL PRINCIPLE

Examination of stained peripheral blood smears is used for screening and identifying malarial parasites, *Babesia*, trypanosomes, and microfilaria. Malarial and *Babesia* parasites infect circulating red cells and undergo various stages of development within the red cell. The Wright Giemsa stain highlights morphologic features of these stages.

3. SPECIMEN REQUIREMENTS**3.1 Patient Preparation**

Component	Special Notations
Fasting/Special Diets	None
Specimen Collection and/or Timing	<p>Slides are to be prepared when the patient presents with symptoms of malaria, and every 6 hours for 36 hours. Specimens obtained during the febrile state yield the greatest number of parasites in circulating blood.</p> <p>Malaria specimens are collected in an EDTA lavender tube. Smears must be made within 60 minutes but preferred to be made within 30 minutes. This reduces distortion of the parasites and RBCs.</p> <p>Prepare four thin smears first: Mix the blood thoroughly prior to making the slides. Place one drop of blood near one end of a slide and then spread the blood over the surface with a second slide.</p> <ul style="list-style-type: none"> Strive for a thin smear that is rounded, feathered and progressively thinner toward the center of the slide.

Component	Special Notations
	<ul style="list-style-type: none"> The thin, feathered end should be at least 2 cm long, and the film should occupy the central area of the slide, with free margins on each side. <p>Prepare four thick smears: Mix the blood thoroughly prior to making slides. Place a small drop of blood in the center of the slide. Using an applicator stick, spread the drop in a circular pattern until it is about the size of a dime. Allow for complete air-drying of the smears.</p> <ul style="list-style-type: none"> Care should be taken not to make the thick smear too thick, as the blood will flake off when dried. A good test to determine if the blood is thick enough, is to hold printed material under the film. If the print is readable, the smear has the correct thickness. <p>Label the frosted end of the slides using a pencil. Include the patient name, medical record number, accession number, tech code, and date.</p> <p>Allow all slides to air-dry in a flat position.</p>
Special Collection Procedures	<p>Smears must be made on site within 60 minutes, but 30 minutes is preferred, to reduce the distortion of the parasites and RBCs. Slides are read and tracked as follows:</p> <p>GEC → SGMC (both thick and thin smears) using LIS template GLAB</p> <p>FWMC → SGMC (both thick and thin smears) Using LIS template FMANS</p> <p>WOMC → SGMC (thin is read on site, and thick smear is sent immediately, if it cannot be read within 24 hours) using LIS template WMAN.</p> <p>Smears, the EDTA tube, and the Malaria/Blood Parasite form (AG.F289) will be sent to Shady Grove, via STAT courier, where the slides will be stained and examined-</p> <p>Notify the Microbiology tech at SGMC that slides are on the way.</p>
Other	A Malaria/Blood Parasite History Form MUST be completed for each patient.

3.2 Specimen Type & Handling

Criteria

Criteria	
Type -Preferred -Other Acceptable	Two thin and two thick smears
Collection Container	See section 3.1
Volume - Optimum - Minimum	N/A N/A
Transport Container and Temperature	Slide holder at room temperature
Stability & Storage Requirements	Room Temperature: Slides must be made from the EDTA tube within 60 minutes, but 30 minutes is preferred. Slides are stable for 1 month.
	Refrigerated: Unacceptable
	Frozen: Unacceptable
Timing Considerations	N/A
Unacceptable Specimens & Actions to Take	If specimen is too old test must not be performed. Improperly prepared or improperly labeled slides. Reject specimen and request recollection.
Compromising Physical Characteristics	N/A
Other Considerations	Treatment with anti-malarial or other antiparasitic drugs may reduce the sensitivity of the test.

NOTE: Labeling requirements for all reagents, calibrators and controls include: (1) Open date, (2) Substance name, (3) Lot number, (4) Date of preparation, (5) Expiration date, (6) Initials of tech, and (7) Any special storage instructions. Check all for visible signs of degradation.

4. REAGENTS

The package insert for a new lot of kits must be reviewed for any changes before the kit is used. A current Package Insert is included as a Related Document.

4.1 Reagent Summary

Reagents	Supplier and Catalog Number
Giemsa Stain	Epredia 89002
Buffer	Alphatec Giemsa (Malaria) Stain Buffer – 033-25

4.2 Reagent Preparations and Storage

Reagent	Giemsa Stain
Container	1 L bottle
Storage	15-30°C
Stability	Stable until expiration date
Preparation	None

Reagent	Alphatec Giemsa (Malaria) Stain Buffer
Container	125 mL bottle
Storage	15-30°C
Stability	Stable until expiration date
Preparation	None

5. CALIBRATORS/STANDARDS

N/A

6. QUALITY CONTROL

6.1 Controls

Appearance of blood cells is noted every time a patient's smear for malaria is performed.

Romanowsky Color Range

Chromatin of white blood cells	purple
Nuclei of parasitic protozoa	red
Basophilic cytoplasm of lymphocytes, monocytes, and parasitic protozoa	blue
Eosinophilic granules	pink
Neutrophilic granules	purple
Red blood cells	salmon pink (to bluish)
Bacteria	deep blue

Record QC results on Malaria Stain QC Form.

If controls are unacceptable do not report patient results, notify supervisor.

6.2 Control Preparations and Storage

N/A

6.3 Frequency

Each batch of patient smears is evaluated for proper staining characteristics.

6.4 Tolerance Limits and Criteria for Acceptable QC

A run is rejected if the WBCs, RBCs, and platelets on the thin smear are not stained adequately.

Rejected runs must be effectively addressed by corrective action. Steps taken in response to QC failures must be documented. Patient samples in failed analytical runs must be reanalyzed.

6.5 Documentation

Steps taken in response to QC failures must be documented.

6.6 Quality Assurance Program

N/A

7. EQUIPMENT and SUPPLIES**7.1 Assay Platform**

None

7.2 Equipment

Microscope
Wescor Hematology Slide Stainer

7.3 Supplies

Immersion oil
Glass slides
Geimsa Stain
Buffer

8. PROCEDURE

NOTE: For all procedures involving specimens, buttoned lab coats, gloves, and face protection are required minimum personal protective equipment. Report all accidents to your supervisor.

8.1	Thin Smears
1.	Allow smear to dry thoroughly before staining.
2.	Place patient thin smears on Wescor Hematology Slide Stainer. (see Hematology procedure for stainer instructions)
3.	Examine thin smears under 10X and 100X (oil immersion) to screen for the presence of malarial parasites, <i>Babesia</i> , microfilaria, and trypanosomes.
4.	At least 300 fields must be viewed with a 100X oil immersion lens for adequate assessment.
8.2	Thick Smears
1.	Allow smear to dry thoroughly before staining (at least 2 hours).
2.	Do not fix with alcohol or heat or dry in an incubator. Heat will prevent RBC lysis.

8.2	Thick Smears
3.	Obtain working Giemsa solution. Into a Coplin jar add 49ml of the phosphate buffer, 1 ml of the Giemsa Blood Stain. Mix well before use. The working stain solution is stable for 24 hours.
4.	Place the thick smears directly into the working solution for 45-60 minutes. The water-based Giemsa stain disrupts the red cell membrane (laking) during the staining procedure exposing the parasites.
5.	Wash the smears by rinsing them with buffer (pH 7.0 to 7.2) for 3-5 minutes.
6.	Record pH of buffer on QC sheet.
7.	Air-dry in a vertical position.
8.3	Reading
1.	Scan the smear under low power first to detect presence of microfilaria or trypanosomes.
2.	Next read under oil immersion (100X objective).
3.	At least 300 fields under oil immersion must be examined.
4.	All shifts will stain thin and thick smears and screen thin smears for malaria, <i>Babesia</i> , microfilaria, and trypanosomes.
5.	If slides must be sent to another site for interpretation, keep 1 set of slides at the originating site.
6.	If positive and the species cannot be determined after review by a microbiology lead tech, supervisor or director, thin and thick smears may be sent to Washington Adventist Hospital to Dr. Beltaifa for pathologist review if she is available. If Dr. Beltaifa is not available, refer the slides to Chantilly.

9. CALCULATIONS

Parasitemia: In areas of the slide where the RBCs are evenly spread out over the entire field and not overlapping count the number of infected cells per field of 200 cells. Do this on 10 different areas on each thin smear. Calculate the average and divide by 2. The resulting number is the percentage of RBC's infected.

For parasitemia lower than 0.01%, report the parasitemia as <0.01%.

10. REPORTING RESULTS AND REPEAT CRITERIA

10.1 Interpretation of Smears

Smears are examined utilizing 10X and oil immersion lens (100x).

Read a minimum of 300 fields under oil immersion before determining that the thin smears are negative.

Thick smears are to be read before finalizing the report as negative. Thick smears are also a guide to the intensity of the infection. Thick smears allow a large amount of blood to be examined, increasing the detection of parasites in light infections. If parasites are detected on the thick smears, species determination must be made using the thin smear examination. This is determined by the recovery and identification of life cycle stages observed on the thin smear.

10.2 Reporting

10.2.1 General Information

Call both positive preliminary and final results to the nursing unit or physician. The call back information must be documented in the LIS.

Preliminary Reports:

If thin smears are negative, report: "Thin smear presumptive negative, thick smear and final report to follow". (NMLP1)

If thin smears are positive, report: "Presumptive positive, confirmation and identification to follow." (PMAL1)

If *Plasmodium falciparum* can be ruled out, report *Plasmodium* species, not *P. falciparum*. Send all smears (thin and thick) to Fort Washington Medical Center to Dr. Beltaifa for pathologist review.

If microfilaria or trypanosomes are seen, report their presence and send to Chantilly for confirmation.

*** Do not finalize the thin smear preliminary report in the LIS.**

Final Reports:

If negative, report: "No parasites seen. One set of blood films can not exclude the diagnosis of malaria." (NMAL1)

If positive report:

1. report genus and species for malaria or "*Babesia* species" if *Babesia* is seen.
2. report the level of parasitemia

Parasitemia:

Report the percentage of cells infected on all positive *Plasmodium* species or *Babesia* species. See section 9 for instructions on performing the calculation.

Report: "x.xx % Parasitemia"

Enter the number and % sign, then enter the code INF2 or enter using the Sunquest keyboard as in 10.2.2 below.

The call back information must be documented in the LIS.

*** More than one technologist must review all initial positive malaria smears. Repeated positive smears on the same patient do not require review by a second technologist.**

Document both tech codes on the LIS workcard.

Reporting to Maryland DHMH:

Smears positive for malaria or *Babesia* species must be reported by the technologist who reported the result to the Maryland Department of Health and Mental Hygiene by completing DHMH form 1281. Reports must be submitted within 1 working day (fax or mail).

Fax the form to the Montgomery County Health Department (240-777-4680).

If mailing, mark the sealed envelope “confidential” and send to

Montgomery County Health Department
2000 Dennis Ave
Suite 238
Silver Spring, MD 20992

One set of thick and thin smears must be sent to the Maryland DHMH.

- Place slides in a plastic slide holder and package in a padded shipping envelope or box.
- Form DHMH 4676 must be completed and must accompany the slides.
- Address the package to Maryland Department of Health and Mental Hygiene, 201 W Preston St., Baltimore, MD 21201.
- Place package with the newborn screen samples for courier pickup

10.2.2 LIS resulting

1. Use GUI function **Micro Result Entry**.
2. Key in the accession # and click on **SELECT**.
3. Press on F8 to display the resulting keyboard. Note: to turn off the keyboard press F8 again.

4. Click on the Micro keyboard to enter in your results
5. If you have a positive malaria smear then result as follows:
 - a. **Observation #: Organism** - Click on the organism from the Micro keyboard, then press the tab key until you are at the next observation line in the **result** field.
 - b. **Observation #: Infectivity rate** - Press **;** twice (the first ; will not display on the screen but the second one will) and then free text the infectivity rate (example 2.0 %), then press the tab key. From the Micro keyboard, click on the **M** key. This will add **-Parasitemia** to your infectivity rate.
 - c. **Observation #:** ;CBACK <cr> (expands to 'Called to and read back by:'); ; add free text call documentation <cr> <cr>
6. Press the tab key twice so that your cursor is in the next result field.
7. If this is a prelim, click on **SAVE**
8. If this is a final, click on the **/ (final)** key from the Micro Keyboard and then click on **SAVE**.


* Each call must be documented. Do not delete previous call back information.

10.3 Interpretation of data


Morphological Characteristics of *Plasmodium falciparum* and *Babesia* species:

Appearance of parasite	<i>Plasmodium falciparum</i>	<i>Babesia</i> species
Size	Small (1/4 to 1/3 RBC diameter,	Tiny to small (1/8 to 1/4 RBC

	3-5 μ m)	diameter, 1-5 μ m)
Shape	Consistent oval to round ring	Pleomorphic: pear-shaped to round ring
Appliqué Forms	Common, either marginal or bulging forms	Common, either marginal or bulging forms
Number of Chromatin dots	1 to 2	1 to many ("string of pearls")
Multiple rings/RBC	Common	Common; two adjacent parasites may appear to be split into mirror images
Tetrads	No	Rarely seen
Appearance of RBCs	Normal size and shape	Normal size and shape
Parasite stages present	Ring: trophozoite with pigment (in heavy infections); banana-shaped gametocytes (rarely found)	Ring: late ring or trophozoite with no pigment, may contain a white central vacuole not seen in <i>Plasmodium</i>



A: *Babesia microti* infection, Giemsa-stained thin smear. The organisms resemble *Plasmodium falciparum*; however *Babesia* parasites present several distinguishing features: they vary more in shape and in size; and they do not produce pigment.




B. and C: Infection with *Babesia*. Giemsa-stained thin smears. Note in **B** the tetrad (left side of the image), a dividing form pathognomonic for *Babesia*. Note also the variation in size and shape of the ring stage parasites (compare **B** and **C**), and the absence of pigment.

Morphology of *Plasmodium* species in Wright-Giemsa stained smears:

Characteristics	<i>P. falciparum</i>	<i>P. vivax</i>	<i>P. ovale</i>	<i>P. malariae</i>
-----------------	----------------------	-----------------	-----------------	--------------------

Characteristics	<i>P. falciparum</i>	<i>P. vivax</i>	<i>P. ovale</i>	<i>P. malariae</i>
Size and shape of infected erythrocytes	Normal size and shape	Enlarged up to twofold, may be oval	Normal to enlarged, frequently oval, may be fimbriated	Small to normal size, normal shape
Stippling	Occasional Mauer's dots, less numerous than Schuffner's	Schuffner's dots (stippling) usually present, except in rings	James' stippling, darker than Schuffner's, present in all stages, including rings	Zeiman's dots rarely seen; requires deliberate over staining
Stages seen in peripheral blood	Rings and gametocytes	All	All	All
Multiply infected erythrocytes	Common	Occasional	Occasional	Rare
Early trophozoites	Delicate ring, frequently with two small chromatin dots, often at the edge of the erythrocyte	Ring up to 1/3 diameter of the erythrocyte; larger chromatin dot than <i>P. falciparum</i>	Similar to <i>P. vivax</i>	Smaller than <i>P. vivax</i> ; otherwise similar
Mature trophozoites	Not seen in peripheral blood	Amoeboid shape, fine golden brown pigment	Similar to <i>P. vivax</i> except less amoeboid, pigment darker brown	Compact cytoplasm, oval, round, or band-shaped, dark brown pigment
Schizonts	Not seen in peripheral blood	12-24 merozoites	8-12 merozoites	6-12 merozoites often radially arranged around central pigment in a rosette form
Gametocytes	Crescent or banana-shaped	Round to slightly oval	Round to slightly oval	Round to slightly oval
Most characteristic findings	Absence of mature trophozoites and schizonts; normal size of infected erythrocytes; multiply infected RBCs; appliqué forms; banana-shaped gametocytes	Enlarged infected erythrocytes; Schuffner's dots frequently present; amoeboid trophozoite; 12-14 merozoites in each schizont	Normal to enlarged, oval or fimbriated infected RBCs, James' stippling may be seen in rings; schizonts with 8-12 merozoites	Normal size of infected erythrocytes; no stippling; "band" trophozoite; rosette schizont with 6-12 merozoites


Malaria Identification




Always calculate and report % parasitemia

Ring Stage Parasites

Plasmodium falciparum: Rings



Plasmodium vivax: Rings



Plasmodium ovale: Rings




***Plasmodium malariae*: Rings**




Schizonts


***Plasmodium falciparum*: Schizonts (usually not seen in blood)**




***Plasmodium vivax*: Schizonts**



***Plasmodium ovale*: Schizonts**




***Plasmodium malariae*: Schizonts**




Trophozoites


***Plasmodium falciparum*: Trophozoites** (early forms may be seen but later forms usually not seen in blood)



***Plasmodium vivax*: Trophozoites**




***Plasmodium ovale*: Trophozoites**




Gametocytes


Plasmodium falciparum: Gametocytes



Plasmodium vivax: Gametocytes



Plasmodium ovale: Gametocytes





Plasmodium malariae: Gametocytes




Morphology of trypanosomes and microfilaria.

Trypanosomes



Microfilaria



10.5 Rounding / Units of Measure / Clinically Reportable Range (CRR)

N/A

10.6 Review Patient Data

Review patient results for unusual patterns, trends or distributions in patient results such as an unusually high percentage of abnormal results.

10.7 Repeat Criteria and Resulting

N/A

11. EXPECTED VALUES**11.1 Reference Ranges**

No parasites seen.

11.2 Critical Values

Any positive smear

11.3 Standard Required Messages

None

12. CLINICAL SIGNIFICANCE

Malaria is a disease of worldwide importance characterized by fever, anemia and splenomegaly. Although four species of the genus *Plasmodium* (*P. falciparum*, *P. vivax*, *P. malariae*, and *P. ovale*) infect humans, malaria is clinically two diseases; the benign type due to *P. vivax*, *P. malariae*, and *P. ovale*, and the malignant type due to *P. falciparum*.

Determination of parasitemia becomes important when therapy is initiated. The patient's parasitemia is monitored so that possible cases involving drug-resistant strains of *P. falciparum* may be detected. In those cases where the patient is hospitalized, monitoring of the parasitemia should be performed at 24, 48 and 72 hours after initiating therapy. Generally, if the malarial strain is susceptible to the therapeutic regime, the parasitemia will drop significantly within the first 24 hours (often by 50% or more).

Babesia is a malaria-like disease characterized by fever, chills, headache, lethargy and myalgia. Hemolytic anemia and hemoglobinuria are typical and may be severe. The disease is transmitted by the bite of hard ticks of the family Ixodidae. This disease is suspected when individuals have traveled through tick-infested areas and present with a malaria-like illness. The disease becomes apparent 1-3 weeks after the bite of an infectious tick. In splenectomized and immunocompromised patients this disease may be fatal. Determination of % parasitemia helps direct therapy. In severe parasitemia (>10%), exchange transfusion may be considered.

Parasitemia and Clinical Correlation

Parasitemia	Parasites /uL	Clinical Correlation
0.0001-0.0004%	5 – 20	Number of organism that are required for a positive thick film (sensitivity)
0.002%	100	Patients may be symptomatic below this level
0.2%	10,000	Level above which immune patients will exhibit systems
2%	100,000	Maximum parasitemia of <i>P. vivax</i> and <i>P. ovale</i> (infect young RBSs only)
2-5%	100,000 – 250,000	Hyperparasitemia, severe malaria, increased mortality
10%	500,000	Exchange transfusion may be considered, high mortality

13. PROCEDURE NOTES

- **FDA Status:** LDT without message
- **Validated Test Modifications:** None

Any alcohol left on the skin prior to collection may fix the red cells and then they will not clear in the staining procedure.

Do not dry smears using heat, as this will fix the red cells.

Slides prepared from EDTA blood are not optimal as they may cause distortion in the parasites, making identification difficult. However, the Emergency Center at Germantown is the ONLY location where an EDTA specimen is acceptable rather than fingertip smears.

Organisms are most likely to be detected if the smears are obtained immediately upon the onset of fever, or immediately before the fever is anticipated. In patients with a strong clinical history, but repeatedly negative results, multiple sampling throughout the fever may prove successful.

Platelets sitting on top of red blood cells may have the appearance of a ring form of malaria.

Precipitated stain may obscure malarial forms on the smear.

Identification to species should not be based solely on the examination of the thick smear preparation. Both thick and thin smears are required for a comprehensive blood parasite examination.

The patient's travel history may provide helpful information in the identification of malaria, *Babesia* species, and other blood parasites. Blood parasites are endemic to certain regions of the world; knowing what countries the patient has visited will aid in diagnosis.

The chart below can be used as a guide for diagnosis. It is not to be used as the primary diagnostic factor:

Blood Parasite	Endemic Area(s)
<i>Plasmodium falciparum</i>	Africa, Asia, Indian subcontinent, South America (Tropical areas worldwide)
<i>Plasmodium vivax</i>	Tropical and Temperate areas worldwide
<i>Plasmodium malariae</i>	Africa, Asia, Indian subcontinent, South America (Tropical areas worldwide)
<i>Plasmodium ovale</i>	West Africa, India, South America, some South Pacific Islands

14. LIMITATIONS OF METHOD

14.1 Analytical Measurement Range (AMR)

N/A

14.2 Precision

N/A

14.3 Interfering Substances

N/A

14.4 Clinical Sensitivity/Specificity/Predictive Values

N/A

15. SAFETY

Refer to your local and corporate safety manuals and Safety Data Sheet (SDS) for detailed information on safety practices and procedures and a complete description of hazards.

16. RELATED DOCUMENTS

Hematology Slide Stainer Cytocentrifuge, Wescor Aerospray; Hematology procedure
 Resulting Microbiology Direct Exams, Microbiology procedure
 Malaria Smear Collection- GEC and FWMC
 Tracking Specimens between AHC Lab Sites
 Malaria/Blood Parasite History Form (AG.F289)
 Reportable Results to State and Outside Agencies, Laboratory policy

17. REFERENCES

1. Jacobs DS, et al, Laboratory Test Handbook, 4th edition, Hudson, OH: Lexi-Comp, Inc., 1999, pp. 332-333.
2. Kjeldsberg C, et al, Practical Diagnosis of Hematologic Disorders, 2nd edition, Chicago, IL: ASCP Press, 1995, pp. 172-173.
3. Atlas of Human Parasitology, 3rd edition, Chicago, IL: ASCP Press.

4. Hansheid, T. 1999. Diagnosis of malaria: A review of alternatives to conventional microscopy. Clin. Lab. Haematol. 21:235-245.
5. Wilkinson, R.J., J.L. Brown, G. Pasvol, P.L. Chiodini, and R.N. Davidson. 1994. Severe falciparum malaria: predicting the effect of exchange transfusion. Q.J. Med. 87:553-557.

18. REVISION HISTORY

Version	Date	Section	Reason	Reviser	Approval
			Supersedes SOP M028.005		
000	10/12/09	10.2.2	LIS upgrade to GUI system	A. Sears	R. Master
000	10/12/09	16	Added procedure for resulting	L. Barrett	R. Master
001	9/19/2011	3.1, 13	Added use of EDTA specimen at GEC	C. Reidenauer	R. Master
001	9/19/2011	4.2	Changed storage temperature for buffer	R. Master	R. Master
001	9/19/2011	8	Remove statement regarding pkg insert	L. Barrett	R. Master
001	9/19/2011	8.3	Added trypanosomes and microfilaria	R. Master	R. Master
001	9/19/2011	11.2	Update title to local terminology	L. Barrett	R. Master
002	11/19/12	9 10.2.1 10.2.2	Change report to "Parasitemia followed by the % infectivity". Changed steps on how to report the % infectivity (English Text code first then free text the rate %)	M. Sabonis	R. Master
003	7/17/14	3.1, 3.2	Changed EDTA time to 30 min. Removed sending all thick smears to WAH.	R. Master	R. Master
003	7/17/14	8.2	Add stability of working solution. Removed comment to send all think smears to WAH.	R. Master	R. Master
003	7/17/14	8.3, 10.2	Add to send to Dr. Beltaifa if species could not be determined.	R. Master	R. Master
003	7/17/14	9	Clarified calculation	R. Master	R. Master
003	7/17/14	10.2	Add preliminary report of <i>Plasmodium</i> species, not <i>P. falciparum</i> .	R. Master	R. Master
003	7/17/14	10.2.1	Change order of reporting parasitemia, deleted redundant calculation	R. Master	R. Master
003	7/17/14	10.2.2	Add instructions for entering results in GUI version of LIS	R. Master	R. Master
003	7/17/14	10.3	Added flow chart	R. Master	R. Master
003	7/17/14	16	Update titles, add form number	L. Barrett	R. Master
003	7/17/14	Footer	Version # leading zero's dropped due to new EDCS in use as of 10/7/13	L. Barrett	R. Master
4	9/20/16	Header	Add WAH	L. Barrett	R. Master
4	9/20/16	3.1	Add date and MR# to slide label	R. Master	R. Master

Adventist HealthCare

Site: All Laboratories

Title: **Malaria**

4	9/20/16	4	Update to new standard labeling instruction	L Barrett	R. Master
4	9/20/16	8.3	Add stain thin and thick smears. Add instructions to keep 1 set of slides at originating site. Clarified options if species cannot be determined	R. Master	R. Master
4	9/20/16	10.2	Add reporting of identification and parasitemia of Babesia	R. Master	R. Master
4	9/20/16	10.2	Clarified reporting of patients with repeat positive results. Added reporting and sending slides to Maryland DHMH	R. Master	R. Master
4	9/20/16	10.6	Move patient review from section 6	R. Master	R. Master
4	9/20/16	15	Update to new standard wording	L Barrett	R. Master
4	9/20/16	16	Added SOP Reportable Results to State and Outside Agencies	R. Master	R. Master
4	9/20/16	19	Added Maryland DHMH forms	R. Master	R. Master
5	2/18/20	Header	Changed WAH to WOMC	L Barrett	R. Master
5	2/18/20	9	Added minimum parasitemia level of <0.01%	R. Master	R. Master
6	2/3/22	Header	Changed site to All Laboratories	D. Collier	R. Master
6	2/3/22	Footer	Changed prefix to AHC	D. Collier	R. Master
6	2/3/22	3.1	Added FWMC, changed the number of slides required, and added a reference to site specific collection SOP	D. Collier	R. Master
6	2/3/22	16	Added FWMC to related documents	D. Collier	R. Master
7	8/11/23	3.1, 3.2	Updated FWMC tracking to SGMC Updated Owner of SOP	D. Collier	V. Ponraj
8	12/13/24	4.1	Updated brand name of Giemsa Stain to EpreDia	M. Belay	V. Ponraj
9	4/14/25	3.1	Removed fingerstick; changed sample to EDTA tube; updated slide prep instructions for clarity Added tracking information	H. Genser D. Collier	V. Ponraj
9	4/14/25	3.2	Updated stability info	H. Genser	V. Ponraj
9	4/14/25	16	Updated name of form AG.F289 Removed Phlebotomy SOP and GEC and added tracking SOP.	H. Genser D Collier	V. Ponraj

19. ADDENDA

Maryland form DHMH 1281

Maryland form DHMH 4676

Adventist HealthCare

Site: All Laboratories

Title: **Malaria****CONFIDENTIAL REPORT: LABORATORY EVIDENCE OF CERTAIN COMMUNICABLE DISEASES
USE FOR REPORTING TO: MARYLAND STATE DEPARTMENT OF HEALTH AND MENTAL HYGIENE****USE FOR ALL COMMUNICABLE CONDITIONS EXCEPT HIV and CD4. (Use form DHMH 4492 for HIV and CD4.)****(PLEASE TYPE OR PRINT USING BLACK INK.)**

PATIENT LAST NAME		FIRST	MIDDLE INITIAL		HOSPITAL NUMBER		PREGNANT? (FEMALE) YES <input type="checkbox"/> NO <input type="checkbox"/>	
DATE OF BIRTH			AGE	SEX	ETHNICITY HISPANIC <input type="checkbox"/> NON-HISPANIC <input type="checkbox"/>		RACE	
NUMBER STREET		APT	CITY	STATE	ZIP	COUNTY	(AREA CODE) PHONE	
ORDERING PROVIDER NAME								
NUMBER STREET		SUITE	CITY	STATE	ZIP	COUNTY	(AREA CODE) PHONE	
(AREA CODE) FAX								
ORDERING FACILITY NAME								
NUMBER STREET		SUITE	CITY	STATE	ZIP	COUNTY	(AREA CODE) PHONE	
DATE SPECIMEN COLLECTED		DATE SPECIMEN RECEIVED		DATE RESULTED		LAB ACCESSION NUMBER		
TYPE OF SPECIMEN Sputum <input type="checkbox"/> Stool <input type="checkbox"/> Pharyngeal Swab <input type="checkbox"/> Discharge <input type="checkbox"/> Blood <input type="checkbox"/> CSF <input type="checkbox"/> Washing <input type="checkbox"/> Other (Specify) _____								
SITE OF SPECIMEN (CERVIX, EYE, ETC.)								
NAME OF TEST						TEST NUMBER OR CODE		
RESULT WITH REFERENCE RANGE & INTERPRETATION								
(IF AN ORGANISM RESULT: INCLUDE SPECIES, SEROGROUPING, OR OTHER SUBTYPING IF KNOWN)								
IF A HEPATITIS C RESULT:								
Signal to Cut-Off Ratio (SCO)		Critical Value for SCO		Hepatitis A IgM Result		Hepatitis B Core IgM Result		
LAB NAME (LAB PERFORMING THE TEST)						LAB CLIA NUMBER		
LAB ADDRESS								
LAB DIRECTOR			LAB (AREA CODE) PHONE			DATE OF REPORT		

DHMH 1281

SEND TO YOUR LOCAL HEALTH DEPARTMENT


Revised JAN 26, 2012

For more forms or information, go to <http://ideha.dhmh.maryland.gov/SitePages/what-to-report.aspx>

Adventist HealthCare

Site: All Laboratories

Title: **Malaria**

 Laboratories Administration MD DHMH 201 W. Preston St. • Baltimore, MD 21201 P.O. Box 2355 • Baltimore, MD 21203-2355 410-767-6100 www.dhmh.state.md.us/labs Robert A. Myers, Ph.D., Director		STATE LAB Use Only
INFECTIOUS AGENTS: CULTURE/DETECTION		
TYPE OR PRINT REQUIRED INFORMATION OR PLACE LABELS ON ALL FOUR COPIES	<input type="checkbox"/> DEH <input type="checkbox"/> FFP <input type="checkbox"/> MTY/PN <input type="checkbox"/> NOD <input type="checkbox"/> STD <input type="checkbox"/> TB <input type="checkbox"/> CD <input type="checkbox"/> COR	
	Patient SS# (last 4 digits):	
	Health Care Provider: Last Name: <input type="checkbox"/> SR <input type="checkbox"/> JR <input type="checkbox"/> Other	
	Address: First Name: M.I. Maiden:	
	City: County: Date of Birth (mm/dd/yyyy): / /	
	State: Zip Code: Address:	
	Contact Name: City: County:	
	Phone#: Fax#: State: Zip Code:	
	Test Request Authorized by:	
	Sex: <input type="checkbox"/> Male <input type="checkbox"/> Female <input type="checkbox"/> Transgender M to F <input type="checkbox"/> Transgender F to M Ethnicity: Hispanic or Latino Origin? <input type="checkbox"/> yes <input type="checkbox"/> no	
Race: <input type="checkbox"/> American Indian/Alaska Native <input type="checkbox"/> Asian <input type="checkbox"/> Black/African American <input type="checkbox"/> Native Hawaiian/other Pacific Islander <input type="checkbox"/> White		
Case # DOC# Outbreak # Submitter Lab#		
Collect Date: Collect Time: <input type="checkbox"/> am <input type="checkbox"/> pm Onset Date:		
Reason for Test: <input type="checkbox"/> Screening <input type="checkbox"/> Diagnosis <input type="checkbox"/> Contact <input type="checkbox"/> Test of Cure <input type="checkbox"/> 2-3 Months Post Rx <input type="checkbox"/> Suspected Carrier <input type="checkbox"/> Isolate for ID <input type="checkbox"/> Release		
Therapy/Drug Treatment: <input type="checkbox"/> No <input type="checkbox"/> Yes Therapy/Drug Type: Therapy/Drug Date:		
SPECIMEN CODE BACTERIOLOGY Bacterial Culture - Routine Additional specimen codes: <input type="checkbox"/> Bordetella pertussis <input type="checkbox"/> Group A Strep <input type="checkbox"/> Group B Strep Screen <input type="checkbox"/> C. difficile Toxin <input type="checkbox"/> Diphtheria <input type="checkbox"/> Foodborne Pathogens (B. cereus, C. perfringens, S. aureus) <input type="checkbox"/> Gonorrhea Culture: Incubated? <input type="checkbox"/> yes <input type="checkbox"/> no Hrs. incubated: Add'l specimen codes: <input type="checkbox"/> MRSA (rule out) <input type="checkbox"/> VRE (rule out) ENTERIC INFECTIONS <input type="checkbox"/> Campylobacter <input type="checkbox"/> E. coli O157 typing <input type="checkbox"/> Enteric Culture - Routine (Salmonella, Shigella, E. coli O157, Campylobacter) <input type="checkbox"/> Salmonella typing <input type="checkbox"/> Shigella typing <input type="checkbox"/> V. parahaemolyticus <input type="checkbox"/> Yersinia REFERENCE MICROBIOLOGY <input type="checkbox"/> ABC'S (BIDS) # Organism: <input type="checkbox"/> Bacteria Referred Culture for ID Specify:	SPECIMEN CODE SPECIAL BACTERIOLOGY <input type="checkbox"/> Legionella Culture <input type="checkbox"/> Leptospira <input type="checkbox"/> Mycoplasma MYCOBACTERIOLOGY/AFB/TB <input type="checkbox"/> AFB/TB Culture and Smear <input type="checkbox"/> AFB/TB Referred Culture for ID <input type="checkbox"/> M. tuberculosis Referred Culture for Genotyping <input type="checkbox"/> Nucleic Acid Amplification test for M. tuberculosis Complex (MTC) PARASITOLOGY <input type="checkbox"/> Blood Parasites: Country visited outside US: <input type="checkbox"/> Cvs & Parasites: Immigrant? <input type="checkbox"/> yes <input type="checkbox"/> no <input type="checkbox"/> Cryptosporidium <input type="checkbox"/> Cyclospora/Isospora <input type="checkbox"/> Microsporidium <input type="checkbox"/> Pinworm VIRUS ISOLATION/CHLAMYDIA <input type="checkbox"/> Adenovirus* <input type="checkbox"/> Arbovirus Panel (WNV, EEEV, SLEV) <input type="checkbox"/> Chlamydia trachomatis <input type="checkbox"/> Cytomegalovirus (CMV) <input type="checkbox"/> Enterovirus (Inc. Echo & Coxsackie) <input type="checkbox"/> Herpes Simplex Virus (Types 1 & 2) <input type="checkbox"/> Influenza (Types A & B)* <input type="checkbox"/> Parainfluenza (Types 1, 2 & 3)* <input type="checkbox"/> Respiratory Syncytial Virus (RSV)* <input type="checkbox"/> Varicella (VZV) *MAY INCLUDE RESPIRATORY SCREENING PANEL. Comments:	SPECIMEN CODE RESTRICTED TESTS Pre-approved submitters only <input type="checkbox"/> Chlamydia trachomatis GC NAAT <input type="checkbox"/> Chlamydia trachomatis only NAAT <input type="checkbox"/> Norovirus ** (see comment on back) OTHER TESTS FOR INFECTIOUS AGENTS Test name: Prior arrangements have been made with the following DHMH Laboratories Administration employee:
SPECIMEN CODE: PLACE CODE IN BOX NEXT TO TEST B Blood BW Bronchial Washing CSF Cerebrospinal Fluid CX Cervix/Endocervix E Eye F Feces N Nasopharynx/Nasal P Penis R Rectum SP Sputum T Throat URE Urethra UPV Urine (First Void) UCC Urine (Clean Catch) V Vagina W Wound O Other:		

DHMH 4676 Revised 04/12

ORIGINAL

Malaria/Blood Parasite History Form



- ☐ Germantown Emergency Center
- ☐ Shady Grove Medical Center
- ☐ White Oak Medical Center
- ☐ Fort Washington Medical Center

Attach Specimen Label
or complete the following:

Patient Name (last, first): _____

Hospital Record No.: _____

Age (yrs): _____ (mos): _____ Gender: Male
Female

Date of symptom onset of **this** attack (mm/dd/yyyy): ____/ ____/ ____

Physician Name (last, first): _____

Has the patient traveled or lived outside the U.S. during the past 4 years? Yes / No

If yes, specify Country or Countries:

Was malaria chemoprophylaxis taken? Yes / No

If yes, which drugs were taken?
Chloroquine / Mefloquine / Doxycycline / Primaquine / Malarone® / Other: _____

History of Malaria, or any blood parasite in last 12 months (prior to this report) Yes / No

If yes, date of previous illness: ____/ ____/ ____

If yes, species (check all that apply):
☐ *P. falciparum*
☐ *P. vivax*
☐ *P. malariae*
☐ *P. ovale*
☐ Other (eg: Babesia) _____
☐ Not Determined

LABORATORY USE ONLY:

Malarial Antigen result: _____

Smear Preliminary report: _____

Tech code: _____
2nd Tech code (if positive): _____ CBACK: _____

Smear Final report: _____
Identification (if positive): _____
% Parasitemia (if applicable): _____

Tech code: _____
2nd Tech code (if positive): _____ CBACK: _____

Refer to Pathologist? Yes / No

Pathologist result: _____

Pathologist result entered by Tech: _____