

## High-Sensitivity Cardiac Troponin I

**Signature Healthcare** 

August 2021

**Steven Noel, PhD, DABCC, MT(ASCP)** Senior Clinical Consultant Scientific & Clinical Affairs



## High Sensitivity Troponin I (TNIH)

#### SIEMENS ... Healthineers

#### **Circulation**

#### **ESC/ACC/AHA/WHF EXPERT CONSENSUS DOCUMENT**

Fourth Universal Definition of Myocardial Infarction (2018)

## **Definition of High Sensitivity Troponin**

- The %CV at the 99<sup>th</sup> percentile should be  $\leq$ 10%.
- At least 50% of healthy individuals should have measurable concentrations above the Limit of Detection (LoD) and below the 99<sup>th</sup> percentile.
- Sex-specific 99<sup>th</sup> percentiles are established.

#### Why do we need high sensitivity troponin?


- 1. Better precision offers improved discrimination at diagnostic decision levels
- 2. Allows for the detection of smaller infarcts
- 3. Allows for the potential use of accelerated diagnostic protocols (ADP) for the more rapid triage of patients
- 4. The use of ADPs can reduce ED wait times and shorten patient turnaround times in the ED
- 5. More efficient and potentially better patient care

## How does a contemporary troponin assay compare to a high sensitivity troponin?

#### 99<sup>th</sup> percentile

SIEMENS ...

Healthineers



Troponin below the LOQ

**Detectable with CV between 10 and 20%** 

Detectable with  $CV \le 10\%$ 


Source: Atellica IM High-Sensitivity Troponin I & Atellica IM TnI-Ultra IFUs

## How does a contemporary troponin assay compare to a high sensitivity troponin?

#### 99<sup>th</sup> percentile

SIEMENS .

Healthineers



Troponin below the LOQ

Detectable with CV between 10 and 20%

Detectable with  $CV \le 10\%$ 

#### High Sensitivity Troponin Units Change


- From ng/mL to **ng/L**
- 1000 fold change
- Whole numbers

## Information from the Fourth Universal Definition of MI (2018)



- Cardiac troponin is the preferred biomarker for the evaluation of myocardial injury, and high-sensitivity assays are recommended for routine clinical use.
  - The definition of MI includes a significant rise and/or fall pattern of troponin results
  - Used along with **clinical interpretation** (e.g. History, ECG, HEART Score, etc.)
- High Sensitivity Troponin assays support the use accelerated diagnostic protocols of 0 to 1 hour or 0 to 2 hours
- Definition of "significant change" (rise and/or fall) relative to 99<sup>th</sup> percentile value
  - If the patient's result is <99<sup>th</sup> percentile: delta at least 50-60% of the 99<sup>th</sup> percentile
  - If the patient's result is >99<sup>th</sup> percentile: delta at least 20% of the 99<sup>th</sup> percentile
  - Absolute changes are assay dependent but appear superior to relative percent changes

# **Comparison: Guideline Acceptable vs. Clinically Usable Troponin Assays**



SIEMENS

Healthineers


**Options for Using the 99th Percentile** 

- 1. **Combined** male and female
  - Atellica 45 ng/L (ng/L)

## 2. Sex specific

- Atellica male: 53 ng/L
- Atellica female: 34 ng/L
- Assays are not alike ... different analytical sensitivity and 99<sup>th</sup> percentile

## **Pros and Cons of using sex-Specific 99th Percentiles**



### **High-Sensitivity Troponin Using Sex-Specific Cutpoints**




	Patients with Suspected ACS (n = 1126)		
	Men (n = 622; 55%)		Women (n = 504; 45%)
	High Sensitivity Troponin Classification Using a Single Cutpoint		
p < 0.021	——— n = 142 (23%)	Type 1 MI	n = 80 (16%)
	n = 30 (5%)	Type 2 MI	n = 25 (5%)
	n = 15 (2%)	Myocardial Injury	n = 23 (5%)
	n = 20 (3%)	Unstable Angina	n = 8 (2%)
	n = 425 (67%)	Others	n = 368 (73%)
	High Sensitivity Troponin Classification Using Sex-Specific Cutpoints		
	—— n = 131 (21%)	Type 1 MI	n = 111 (22%)
	n = 28 (5%)	Type 2 MI	n = 28 (6%)
	n = 12 (2%)	Myocardial Injury	n = 28 (6%)
	n = 22 (4%)	Unstable Angina	n = 6 (1%)
	n = 429 (69%)	Others	n = 331 (66%)

Diagnosis by consensus between the two cardiologists through in-depth review of source documents.

### The Clinical Significance of High Sensitivity Troponin



Cardiac troponin assays are the preferred biomarkers for the evaluation of myocardial injury, and high-sensitivity assays are recommended for routine clinical use.




4<sup>th</sup> Universal Definition of MI JACC August 2018

### **NSTEMI Evaluation Algorithms using High-Sensitivity cTn**

- Analytical sensitivity and precision are both important considerations when selecting a evaluation protocol.
- If you are dealing with very early presenters (<3 hours CPO), you may be able to make an earlier diagnosis based on the trajectory of the troponin result for those with high clinical suspicion; however, AMI should not be ruled out until sufficient time has passed from onset of symptoms.
- The institution needs to **determine and/or "fine tune" the algorithm** used based on their individual patient population and the clinical goals of the medical staff.

# Clinical Performance of the Centaur hs-cTnI Assay using the 0/2 hr Testing Protocol - APACE



Note: The Centaur TNIH assay is not intended to be used in isolation; results should be interpreted in conjunction with other diagnostic tests and clinical information. Results may vary according to the prevalence of AMI in the population being tested and other factors that may influence assay performance.

\*If chest pain onset >3h before presentation to the ED


12

SIEMENS


Healthinee

#### **Massachusetts General Hospital Pathway**





#### **Signature Healthcare Algorithm**




SIEMENS

Healthineers

#### Higher Troponin Values are More Consistent with AMI



A broader differential diagnosis associated with lower-range elevations of hs-cTn begins to narrow as concentrations are higher.



Diagnosis	Peak hsTnI (ng/L)	
Type 1 MI	855 (104 – 6755)	
Type 2 MI	125 (48 – 604)	
Acute Injury	74 (37 – 307)	
Chronic Injury	55 (34 – 145)	

9,115 patients with an elevated hsTnI results Chapman et al Circulation. 2020;141:161–171


### "False Positive" Troponin Results

- In the absence of overt myocardial ischemia, elevated cTn levels are often labelled as 'false-positive' hs-cTn results.
- This term should be avoided, as most of these unexpected hs-cTn elevations are 'true positive' for myocardial injury (rather than MI) and reflect previously undetected or underestimated cardiac disease including valvular heart disease, heart failure, and chronic coronary artery disease.
- Many cardiac and non-cardiac disorders may lead to substantial amounts of cardiomyocyte injury and thereby hs-cTn elevations.
- It is important to note that <u>cTn elevations universally portend a worse prognosis</u> than otherwise similar patients without a cTn elevation, irrespective of the underlying disease.
- This is true regardless of whether the patient has <u>heart failure, renal dysfunction, gastrointestinal</u> bleeding, sepsis, respiratory disease, pulmonary embolism, subarachnoid hemorrhage, or stroke or whether the patient is asymptomatic without known cardiovascular disease.
- Obviously, the medical consequences of cardiomyocyte injury as quantified by cTn elevations will be highly individualized and different from that in patients with MI.

#### **Procedures after the Implementation of hs-cTn**



## 2,544 patients presenting to the ED with symptoms suggestive of AMI before and after the introduction of hs-cTn



#### **Contact Information**



#### Steven Noel, PhD, DABCC, MT(ASCP)

Senior Clinical Consultant Scientific & Clinical Affairs (267) 810-4621 <u>steven.noel@siemens-healthineers.com</u>

