

DXC (AST) ASPARTATE AMINOTRANSFERASE

St. Joseph Medical Center Tacoma, WA
 St. Clare Hospital Lakewood, WA
 St. Elizabeth Hospital Enumclaw, WA
 St. Francis Hospital Federal Way, WA
 St. Anthony Hospital Gig Harbor, WA
 Highline Medical Center Burien, WA
 PSC

PURPOSE

To provide instructions for the quantitative determination of aspartate aminotransferase on the DXC 600/800.

PRINCIPLE

AST reagent, when used in conjunction with UniCel® DxC 600/800 System(s), is intended for the quantitative determination of Aspartate Aminotransferase activity in human serum or plasma.

BACKGROUND


Clinical Significance

Aspartate aminotransferase measurements are used in the diagnosis and treatment of certain types of liver and heart disease.

Methodology

The AST reagent is used to measure aspartate aminotransferase activity by an enzymatic rate method. In the assay reaction, the AST catalyzes the reversible transamination of L-aspartate and α -ketoglutarate to oxaloacetate and L-glutamate. The oxaloacetate is then reduced to malate in the presence of malate dehydrogenase (MDH) with the concurrent oxidation of β -Nicotinamide Adenine Dinucleotide (reduced form) (NADH) to β -Nicotinamide Adenine Dinucleotide (NAD).

The SYNCHRON® System(s) automatically proportions the appropriate sample and reagent volumes into a cuvette. The ratio used is one part sample to 11 parts reagent. The system monitors the rate of change in absorbance at 340 nanometers over a fixed-time interval. This rate of change in absorbance is directly proportional to the activity of AST in the sample and is used by the SYNCHRON® System(s) to calculate and express the AST activity.



RELATED DOCUMENTS

| | |
|-------------|--|
| R-PO-CH0810 | Quality Control Program General Laboratory |
| R-PO-CH0809 | Quality Control Westgard Rules Statistics |
| R-PR-AD0540 | Specimen Rejection/Cancellation Protocol |
| J-F-CH0820 | DXC 800 Controls |
| M-F-CH0820 | Chemistry Controls |
| J-F-CH0826 | DXC 800 Calibrators |

| | |
|------------|--|
| M-F-CH0826 | Chemistry Calibrators |
| M-F-CH1940 | DXC 600 (AMR) Analytical Measurement Range |
| J-F-CH1940 | DXC 800 (AMR) Analytical Measurement Range |
| R-W-CH0815 | DXC Reagent Lot to Lot Correlations |
| R-F-CH0814 | Lot-to-Lot Correlation |

SPECIMEN

Type of Specimen

Biological fluid samples should be collected in the same manner routinely used for any laboratory test. Freshly drawn serum or plasma is the specimens of choice. Acceptable anticoagulants are listed in PROCEDURAL NOTES section of this chemistry information sheet. Whole blood is not recommended for use as a sample.

Specimen Storage and Stability

1. Tubes of blood are to be kept closed at all times and in a vertical position. It is recommended that the serum or plasma be physically separated from contact with cells within two hours from the time of collection.
2. Separated serum or plasma should not remain at room temperature longer than 8 hours. If assays are not completed within 8 hours, serum or plasma should be stored at +2°C to +8°C. If assays are not completed within 48 hours, or the separated sample is to be stored beyond 48 hours, samples should be frozen at -15°C to -20°C. Frozen samples should be thawed only once. Analyte deterioration may occur in samples that are repeatedly frozen and thawed.

| Sample Type | Volume | Sample Stability |
|--------------|--------|---|
| Plasma/Serum | 0.5mL | <ul style="list-style-type: none"> • Separate serum from cells within 2 hours • Room Temp 8 hours • Refrigerated 48 hours • Frozen 3 months |

Criteria for Unacceptable Specimens

See Specimen Rejection/Cancellation Protocol

SAMPLE VOLUME

A filled 0.5 mL sample cup is the optimum volume. For optimum primary sample tube volumes in primary tube samples and minimum volumes, refer to the Primary Tube Sample Template for your system.

REAGENTS

Contents

Each kit contains the following items:

Two Aspartase Aminotransferase Reagent Cartridges (2 x 200 tests) or (2 x 400 tests and two bottles of AST [A-reagent])

| Volume per Test | |
|----------------------|----------------------------|
| Sample Volume | 23 µL |
| Ordac Sample Volume | 3 µL |
| Total Reagent Volume | 250 µL |
| Cartridge Volumes | A 242 µL B 8 µL C -- |

| Reactive Ingredients | |
|----------------------------|-------------|
| α-Ketoglutarate | 16 mmol/L |
| Malate dehydrogenase (MDH) | >2300 IU/L |
| L-Aspartate | 218 mmol/L |
| NADH | 0.18 mmol/L |

Also non-reactive chemicals necessary for optimal system performance.

Reagent Preparation

For P/N 442665 (200 tests):

Transfer all of the contents of the smallest reagent compartment (C) into the largest reagent compartment (A).

For P/N 476831 (400 tests):

Transfer all the contents of one AST (A-reagent) bottle into the largest reagent compartment (A).

Replace the cartridge caps and gently invert the cartridge several times to ensure adequate mixing.

Acceptable Reagent Performance

The acceptability of a reagent is determined by ensuring that quality control results are within your facility's acceptance criteria.

NOTE: New lots of reagent require lot to lot correlation studies. Refer to Related Documents section for related work instructions/forms.

Reagent Storage and Stability

AST reagent when stored unopened at +2°C to +8°C will obtain the shelf-life indicated on the cartridge label. Once prepared, the reagent is stable for 30 days at +2°C to +8°C. Do not use beyond the manufacturer's expiration date. DO NOT FREEZE.

CALIBRATION

Calibrator Required

Calibration is not required.

Traceability

This measurand (analyte) is traceable to the manufacturer's selected Measurement Procedure as described in the Methodology section.

QUALITY CONTROL

See Related Documents J-F-CH0820 DXC 800 Controls & M-F-CH0820 Chemistry Controls

STEPS

1. If necessary prepare reagent as defined in the Reagent Preparation section of this chemistry information sheet and load the reagent onto the system.
2. Program controls for analysis.
3. After loading controls onto the system, follow the protocols for system operation. To load samples manually refer to the FHS DXC Series Manual Sample Programming procedure. For detailed testing procedures, refer to the UniCel DxC 600/800 System *Instructions For Use* (IFU) manual.

CALCULATIONS

SYNCHRON® System(s) perform all calculations internally to produce the final reported result. The system will calculate the final result for sample dilutions made by the operator when the dilution factor is entered into the system during sample programming.

ANTICOAGULANT TEST RESULTS

If plasma is the sample of choice, the following anticoagulants were found to be compatible with this method based on a study of 20 healthy volunteers:

| Anticoagulant | Level Tested for In Vitro Interference |
|------------------|--|
| Ammonium Heparin | 29 Units/mL |
| Lithium Heparin | 29 Units/mL |
| Sodium Heparin | 29 Units/mL |

The following anticoagulants were found to be incompatible with this method:

| Anticoagulant | Level Tested for In Vitro Interference |
|-------------------|--|
| Potassium Oxalate | 2 Units/mL |
| Sodium Fluoride | 2.5 Units/mL |
| Sodium Citrate | 6.6 Units/mL |

PERFORMANCE CHARACTERISTICS

Reference Range

| Sample Type | Age (years) | Conventional Units |
|-----------------|-------------|--------------------|
| Serum or Plasma | 0-7 | 20-60 U/L |
| Serum or Plasma | 7-11 | 20-40 U/L |
| Serum or Plasma | 11-18 | 14-40 U/L |
| Serum or Plasma | >18 | 10 –45 U/L |

Analytic Range

The SYNCHRON[®] System(s) method for the determination of this analyte provides the following analytical range:

| Sample Type | Conventional Units |
|-------------------------|--------------------|
| Serum or Plasma | 5 – 400 IU/L |
| Serum or Plasma (Ordac) | 350-2600 IU/L |

Samples with activities exceeding the high end of the analytical range should be rerun with ORDAC enabled or diluted with saline and reanalyzed.

Reporting results outside of analytical range

| | | |
|--------------------------|-----------|--|
| Lower limit of detection | 5 IU/L | Results below 5; Report as <5 IU/L (See Limitations below for "OIR LO" results) |
| Upper limit of range | 2600 IU/L | Result >2600 should be diluted with 0.9% saline, reanalyzed and dilution factor applied. The maximum allowable dilution is X5. Results >13,000 are reported as >13,000 IU/L. |

Sensitivity

Sensitivity is defined as the lowest measurable concentration which can be distinguished from zero with 95% confidence. Sensitivity for AST determination is 5 IU/L (0.08 µkat/L).

LIMITATIONS

Samples with extremely high enzyme activity (>12,000 IU/L) may consume all of the NADH substrate before the first absorbance measurement is taken after sample addition. These samples can report either very low enzyme activities or suppress the result as "OIR LO". These samples should be diluted 1:20 with saline and rerun.

Interferences

1. The following substances were tested for interference with this methodology:

| Substance | Source | Level Tested | Observed Effect |
|------------|----------------|-------------------------|---|
| Hemoglobin | RBC hemolysate | INDEX OF 1 | AVOID HEMOLYSIS |
| Bilirubin | Bovine | 30 mg/dL INDEX of 20 | No significant interference (within ± 6 IU/L or 7%) |
| Lipemia | Intralipid | 320 mg/dL INDEX of 8 | +7 @ 67 IU/L |
| | | Airfuge recommended | No significant interference (within ± 6 IU/L or 7% at 307 IU/L) |
| Pyruvate | Pyruvic Acid | 2.4 mg/dL | +8 @ 61 IU/L |
| | | 6.0 mg/dL | No significant interference (within ± 6 IU/L or 7% at 286 IU/L) |

2. Samples showing evidence of hemolysis should not be used. Hemolysis may cause falsely elevated results.

3. Refer to References (9,10,11) for other interferences caused by drugs, disease and preanalytical variables.

ADDITIONAL INFORMATION

For more detailed information on the UniCel DxC Systems, refer to the appropriate system manual.

REFERENCES

1. Kamen, A., Wroblewski, F., LaDue, J. E., *J. Clin. Inv.*, 34 :126 133 (1955).
2. Henry, R. J., et al., *Amer. J. Clin. Path.*, 34:381 (1960).
3. Tietz, N. W., "Specimen Collection and Processing; Sources of Biological Variation", *Textbook of Clinical Chemistry*, 2nd Edition, W. B. Saunders, Philadelphia, PA (1994).
4. National Committee for Clinical Laboratory Standards, *Procedures for the Handling and Processing of Blood Specimens*, Approved Guideline, NCCLS publication H18-A, Villanova, PA (1990).
5. Tietz, N. W., *Clinical Guide to Laboratory Tests*, 3rd Edition, W. B. Saunders, Philadelphia, PA (1995).
6. National Committee for Clinical Laboratory Standards, *How to Define, Determine, and Utilize Reference Intervals in the Clinical Laboratory*, Approved Guideline, NCCLS publication C28-A, Villanova, PA (1995).
7. Tietz, N. W., ed., *Fundamentals of Clinical Chemistry*, 3rd Edition, W. B. Saunders, Philadelphia, PA (1987).
8. Henry, J. B., *Clinical Diagnosis and Management by Laboratory Methods*, 18th Edition, W. B. Saunders Company, Philadelphia, PA (1991).
9. Young, D. S., *Effects of Drugs on Clinical Laboratory Tests*, 4th Edition, AACC Press, Washington, D. C. (1995).
10. Friedman, R. B., Young, D. S., *Effects of Disease on Clinical Laboratory Tests*, 3rd Edition, AACC Press, Washington, D.C. (1997).
11. Young, D. S., *Effects of Preanalytical Variables on Clinical Laboratory Tests*, 2nd Edition, AACC Press, Washington, D. C. (1997).
12. National Committee for Clinical Laboratory Standards, *Method Comparison and Bias Estimation Using Patient Samples*, Approved Guideline, NCCLS publication EP9-A, Villanova, PA (1995).
13. National Committee for Clinical Laboratory Standards, *Precision Performance of Clinical Chemistry Devices*, Tentative Guideline, 2nd Edition, NCCLS publication EP5-T2, Villanova, PA (1992).

DOCUMENT APPROVAL Purpose of Document / Reason for Change:

Updated for current process

| | | | |
|--------------------------------|--|---|---------------------------------------|
| Committee Approval Date | <input checked="" type="checkbox"/> Date: 8/26/2014 <input type="checkbox"/> NA – revision of department-specific document which is used at only one facility | Medical Director Approval <i>(Electronic Signature)</i> | <i>Katie Wilkinson, MD</i> 8/26/14 |
|--------------------------------|--|---|---------------------------------------|