

## DXC (CAR) CARBAMAZEPINE (TEGRETOL)

- St. Joseph Medical Center Tacoma, WA   
  St. Clare Hospital Lakewood, WA   
  St. Elizabeth Hospital Enumclaw, WA  
 St. Francis Hospital Federal Way, WA   
  St. Anthony Hospital Gig Harbor, WA   
  Highline Medical Center Burien, WA   
 PSC

### PURPOSE

To provide instruction for the quantitative determination of carbamazepine (Tegretol) on the DXC 600/800.

### PRINCIPLE

CAR reagent, when used in conjunction with the UniCel® DxC 600/800 System(s) and SYNCHRON® Systems Drug Calibrator 1 set, is intended for quantitative determination of Carbamazepine concentration in human serum or plasma.

### BACKGROUND

#### Clinical Significance

Carbamazepine is indicated for the treatment of psychomotor and grand mal seizures as well as trigeminal neuralgia. Carbamazepine therapy is monitored for suspected inadequate dose or toxicity.

#### Methodology

CAR reagent is used to measure the carbamazepine concentration by a particle enhanced turbidimetric inhibition immunoassay method. Particle-bound drug (PBD) binds carbamazepine specific antibody (Ab) resulting in the formation of insoluble aggregates causing light scatter. Non-particle-bound carbamazepine in the patient sample competes with the PBD for the antibody binding sites, inhibiting the formation of insoluble aggregates. The rate and amount of particle aggregation is inversely proportional to the concentration of carbamazepine in the sample.

The SYNCHRON® System(s) automatically proportions the appropriate sample and reagent volumes into a cuvette. The ratio used is one part sample to 97 parts reagent. The system monitors aggregate formation by measuring the change in absorbance at 340 nanometers. This change in absorbance is inversely proportional to the concentration of CAR in the sample and is used by the System to calculate and express the CAR concentration based upon a multi-point calibration curve.



### RELATED DOCUMENTS

R-PO-CH-0810	Quality Control Program General Laboratory
R-PO-CH-0809	Quality Control Westgard Rules Statistics
R-PR-AD-0540	Specimen Rejection/Cancellation Protocol
J-F-CH-0820	DXC 800 Controls
J-F-CH-0826	DXC 800 Calibrators
J-F-CH-1940	DXC Analytical Measurement Range
M-F-CH-0820	Chemistry Controls
M-F-CH-0826	Chemistry Calibrators

**Type of Specimen**

Biological fluid samples should be collected in the same manner routinely used for any laboratory test. Freshly drawn serum or plasma are the preferred specimens. Acceptable anticoagulants are listed in the PROCEDURAL NOTES section of this chemistry information sheet. Whole blood or urine are not recommended for use as a sample.

**Specimen Storage and Stability**

1. Tubes of blood are to be kept closed at all times and in a vertical position. It is recommended that the serum or plasma be physically separated from contact with cells within two hours from the time of collection.
2. Separated serum or plasma should not remain at room temperature longer than 8 hours. If assays are not completed within 8 hours, serum or plasma should be stored at +2°C to +8°C. If assays are not completed within 48 hours, or the separated sample is to be stored beyond 48 hours, samples should be frozen at -15°C to -20°C. Frozen samples should be thawed only once. Analyte deterioration may occur in samples that are repeatedly frozen and thawed.

Sample Type	Volume	Sample Stability
Plasma/Serum	0.5mL	<ul style="list-style-type: none"> <li>• Separate serum from cells within 2 hours</li> <li>• Room Temp 8 hours</li> <li>• Refrigerated 48 hours</li> <li>• Frozen 3 months</li> </ul>

**Criteria for Unacceptable Specimens**

See Specimen Rejection/Cancellation Protocol

**Sample Volume**

A filled 0.5 mL sample cup is the optimum volume. For optimum primary sample tube volumes in primary tube samples and minimum volumes, refer to the Primary Tube Sample Template for your system.

**REAGENTS**

**Contents**

Each kit contains the following items:  
 Two CAR Reagent Cartridges (2 x 100 tests) Kit reorder #469112

Volume per Test	
Sample Volume	3 uL
Total Reagent Volume	292 uL
Cartridge Volumes	A 230uL B 30uL C 32uL

Reactive Ingredients	
Carbamazepine Particle Reagent	4.8 mL
Monoclonal anti-Carbamazepine Antibodies (Mouse)	5.4 mL
Carbamazepine Reaction Buffer	34.8 mL

Also non-reactive chemicals necessary for optimal system performance

**Reagent Preparation**

No preparation is required. Do not mix.

### **Acceptable Reagent Performance**

The acceptability of a reagent is determined by ensuring that quality control results are within your facility's acceptance criteria.

### **Reagent Storage and Stability**

CAR Reagent when stored unopened at +2°C to +8°C, will remain stable until the expiration date printed on the cartridge label. Once opened, the reagent is stable for 42 days at +2°C to +8°C. Do not use beyond the manufacturer's expiration date. DO NOT FREEZE. Do not expose reagent to temperatures above +35°C or to direct sunlight.

## **CALIBRATION**

### **Calibrator Required**

SYNCHRON<sup>®</sup> Systems Drug Calibrator 1 set

### **Calibrator Preparation**

No preparation is required.

### **Calibrator Storage and Stability**

SYNCHRON<sup>®</sup> Systems Drug Calibrator 1 set is stable until the expiration date printed on the calibrator bottle if stored capped in the original container at +2°C to +8°C.

### **Calibration Information**

1. The system must have a valid calibration curve in memory before control or patient samples can be run.
2. Under typical operating conditions the CAR reagent cartridge must be calibrated every 14 days and also with certain parts replacements or maintenance procedures, as defined in the UniCel DxC 600/800 System *Instructions For Use* (IFU) manual. This assay has within-lot calibration available. Refer to the UniCel DxC 600/800 System *Instructions For Use* (IFU) manual for information on this feature.
3. For detailed calibration instructions, refer to the UniCel DxC 600/800 System *Instructions For Use* (IFU) manual.
4. The system will automatically perform checks on the calibration and produce data at the end of calibration. In the event of a failed calibration, the data will be printed with error codes and the system will alert the operator of the failure. For information on error codes, refer to the UniCel DxC 600/800 System *Instructions For Use* (IFU) manual.

### **Traceability**

For Traceability information refer to the Calibrator instructions for use.

## **QUALITY CONTROL**

**STEPS**

1. If necessary, load the reagent onto the system.
2. After reagent load is completed, calibration is required.
3. Program quality control for analysis.
4. After loading controls onto the system, follow the protocols for system operation. To load samples manually refer to the FHS DXC Series Manual Sample Programming procedure. For detailed testing procedures, refer to the UniCel DxC 600/800 System *Instructions For Use* (IFU) manual.

**CALCULATIONS**

SYNCHRON<sup>®</sup> System(s) perform all calculations internally to produce the final reported result. The system will calculate the final result for sample dilutions made by the operator when the dilution factor is entered into the system during sample programming.

**ANTICOAGULANT TEST RESULTS**

1. If plasma is the sample of choice, the following anticoagulants were found to be compatible with this method based on a study of 20 healthy volunteers:

Anticoagulant	Level Tested for In Vitro Interference	Average Plasma-Serum Bias (µg/mL)
Lithium Heparin	14 Units/mL	No Significant Interference (within ± 1.2 µg/mL or 10%)
Sodium Heparin	14 Units/mL	No Significant Interference (within ± 1.2 µg/mL or 10%)

2. The following anticoagulants were found to be incompatible based on the same study:

Anticoagulant	Level Tested for In Vitro Interference	Plasma-Serum Bias (µg/mL)
EDTA	1.5 mg/mL	+1.8

**PERFORMANCE CHARACTERISTICS****Reference Range**

Therapeutic	4 – 12 ug/mL
Critical	>15 ug/mL

For Critical Value reporting protocol, refer to the FHS Critical Value Policy.

**Analytic Range**

The SYNCHRON<sup>®</sup> System(s) method for the determination of this analyte provides the following analytical ranges:

Sample Type	Conventional Units
Serum or Plasma	2.0 – 20.0 µg/mL

Samples with concentrations outside of the analytical range will be reported as "<2.0 µg/mL" Samples reported out as greater than the analytical range may be confirmed by diluting with saline and reanalyzing. The appropriate dilution factor should be applied to the reported result.

**Reporting results outside of analytical range**

Lower limit of detection: serum / plasma	2.0 µg/mL	Results below 2.0, report as <2.0
Upper limit of range: serum / plasma	20.0 µg/mL	Results >20.0 should be diluted using 0.9% saline, reanalyzed and dilution factor applied. The maximum allowable dilution is X 2. Results >40 are reported as >40 µg/mL

Samples reported out as "**SUPPRESSED**" due to **RXN ERROR** should be reanalyzed.

**Sensitivity**

Sensitivity is defined as the lowest measurable concentration which can be distinguished from zero with 95% confidence. Sensitivity for CAR determination is 2.0 µg/mL

**LIMITATIONS**

None identified.

**Interferences**

1. The following substances were tested for interference with this methodology:

Substance	Source	Level tested	Observed effect
Hemoglobin	RBC hemolysate	500 mg/dL INDEX of 10	No Significant Interference (within ± 1.2 µg/mL or 10%)
Bilirubin	Porcine	30 mg/dL INDEX of 20	No Significant Interference (within ± 1.2 µg/mL or 10%)
Rheumatoid Factor	Human	300 IU/mL	No Significant Interference (within ± 1.2 µg/mL or 10%)
Lipemia	Human	320 mg/dL INDEX of 8 Airfuge Recommended	No Significant Interference (within ± 1.2 µg/mL or 10%)
Paraprotein (IgM)	Human	500 mg/dL	No Significant Interference (within ± 1.2 µg/mL or 10%)

2. Refer to References (9,10,11) for other interferences caused by drugs, disease and preanalytical variables
3. For assays employing mouse antibodies, the possibility exists for interference by human anti-mouse antibodies (HAMA) in the sample. Human anti-mouse antibodies may be present in samples from patients who have received immunotherapy or diagnostic procedures utilizing monoclonal antibodies or in individuals who have been regularly exposed to animals.<sup>12,13</sup> Additionally, other heterophile antibodies, such as human anti-goat antibodies may be present in patient samples. Interpretation of results should be done in the context of the overall clinical presentation of the patient, including symptoms, clinical history, data from additional tests and other appropriate information.

**ADDITIONAL INFORMATION**

For more detailed information on UniCel DxC Systems, refer to the appropriate system manual.

## REFERENCES

1. Newman, D. J., Henneberry, H., Price, C. P., "Particle Enhanced Light Scattering Immunoassay", *Ann. Clin. Biochem.*, 29:22-42 (1992).
2. Tietz, N. W., "Specimen Collection and Processing; Sources of Biological Variation", *Textbook of Clinical Chemistry*, 2nd Edition, W. B. Saunders, Philadelphia, PA (1994).
3. National Committee for Clinical Laboratory Standards, *Procedures for the Handling and Processing of Blood Specimens*, Approved Guideline, NCCLS publication H18-A, Villanova, PA (1990).
4. CDC-NIH manual, *Biosafety in Microbiological and Biomedical Laboratories*, U.S. Government Printing Office, Washington, D.C. (1984).
5. Burtis, C. A., Ashwood, E. R., *Tietz Textbook of Clinical Chemistry*, 2nd Edition, W. B. Saunders, Philadelphia, PA (1994).
6. National Committee for Clinical Laboratory Standards, *How to Define, Determine, and Utilize Reference Intervals in the Clinical Laboratory*, Approved Guideline, NCCLS publication C28-A, Villanova, PA (1995).
7. Tietz, N. W., ed., *Fundamentals of Clinical Chemistry*, 3rd Edition, W. B. Saunders, Philadelphia, PA (1987).
8. Henry, J. B., *Clinical Diagnosis and Management by Laboratory Methods*, 18th Edition, W. B. Saunders Company, Philadelphia, PA (1991).
9. Young, D. S., *Effects of Drugs on Clinical Laboratory Tests*, 4th Edition, AACC Press, Washington, D. C. (1995).
10. Friedman, R. B., Young, D. S., *Effects of Disease on Clinical Laboratory Tests*, 3rd Edition, AACC Press, Washington, D.C. (1997).
11. Young, D. S., *Effects of Preanalytical Variables on Clinical Laboratory Tests*, 2nd Edition, AACC Press, Washington, D. C. (1997).
12. Bjerner, J., et al., "Immunometric Assay Interference: Incidence and Prevention", *Clin. Chem.* 48:613-621 (2002).
13. Kricka, L. J., "Interferences in Immunoassays-Still a Threat", *Clin. Chem.*, 46:1037-1038 (2000).
14. National Committee for Clinical Laboratory Standards, *Method Comparison and Bias Estimation Using Patient Samples*, Approved Guideline, NCCLS publication EP9-A, Villanova, PA (1995).
15. National Committee for Clinical Laboratory Standards, *Precision Performance of Clinical Chemistry Devices*, Tentative Guideline, 2nd Edition, NCCLS publication EP5-T2, Villanova, PA (1992).

**DOCUMENT APPROVAL Purpose of Document / Reason for Change:**

Updated for current process

<b>Committee Approval Date</b>	<input checked="" type="checkbox"/> Date: 8/26/14 <input type="checkbox"/> NA – revision of department-specific document which is used at only one facility	<b>Medical Director Approval</b> <i>(Electronic Signature)</i>	<i>Katie Wilkinson, MD</i> 8/26/14
--------------------------------	--	---	---------------------------------------