

DXC (BARB) BARBITURATES

St. Joseph Medical Center Tacoma, WA
 St. Clare Hospital Lakewood, WA
 St. Elizabeth Hospital Enumclaw, WA
 St. Francis Hospital Federal Way, WA
 St. Anthony Hospital Gig Harbor, WA
 Highline Medical Center Burien, WA
 PSC

PURPOSE

To provide instructions for the qualitative determination of barbiturates on the DXC 600/800.

PRINCIPLE

BARB reagent, in conjunction with UniCel® DxC 600/800 System(s) and SYNCHRON® Systems Drugs of Abuse Testing (DAT) Urine Calibrators, is intended for the qualitative determination of Barbiturates in human urine at a cutoff value of 200 ng/mL.

The BARB assay provides a rapid screening procedure for determining the presence of Barbiturates (BARB) and its metabolites in urine. This test provides only a preliminary analytical result; a positive result by this assay should be confirmed by another generally accepted non-immunological method such as thin layer chromatography (TLC), gas chromatography (GC), or gas chromatography/mass spectrometry (GC/MS). GC/MS is the preferred confirmatory method.

Clinical consideration and professional judgement should be applied to any drug of abuse test result, particularly when preliminary positive results are used.

BACKGROUND

Clinical Significance


Barbiturates are a class of central nervous system depressants that are used as sedatives and hypnotics. Measurements of barbiturates are used in the diagnosis and treatment of barbiturate use or overdose, and in monitoring the presence of barbiturates to ensure appropriate therapy.

Barbiturate-derivatives such as phenobarbital are useful in the management of epileptic seizures. When barbiturate is ingested, it is rapidly metabolized and excreted into urine. The ratio of unchanged drug to metabolites varies depending upon the duration of action. Short-acting barbiturates will generally be excreted in urine as metabolites, while the long-acting barbiturates will primarily appear unchanged. The most commonly abused barbiturates are short- to intermediate-acting agents (24 hours or less), such as secobarbital, pentobarbital and amobarbital. Barbiturate overdose produces shock syndrome and can result in death from respiratory depression.

Methodology

This assay utilizes a homogenous enzyme immunoassay method. The BARB reagent is comprised of a specific antibody which can detect most analyte in urine. A drug-labeled glucose-6-phosphate dehydrogenase (G6PDH) conjugate competes with any free drug from the urine sample for a fixed amount of antibody binding sites. In the absence of free drug from the sample, the drug-labeled G6PDH conjugate is bound by the specific antibody and the enzyme activity is inhibited. This reaction creates a direct relationship between the presence of drug and enzyme activity. The G6PDH enzyme activity is determined spectrophotometrically by measuring its ability to convert nicotinamide adenine dinucleotide (NAD) to NADH (reduced form).

The SYNCHRON® System(s) automatically proportions the appropriate sample and reagent volumes into a cuvette. The ratio for BARB is one part sample to 25 parts reagent. The System monitors the change in absorbance at 340 nanometers to calculate and express a reaction rate. A qualitative result is reported based on a comparison of the sample rate to the calibrated cutoff rate.



RELATED DOCUMENTS

| | |
|-------------|--|
| R-PO-CH0810 | Quality Control Program General Laboratory |
| R-PO-CH0809 | Quality Control Westgard Rules Statistics |
| R-PR-AD0540 | Specimen Rejection/Cancellation Protocol |
| J-F-CH0820 | DXC 800 Controls |
| J-F-CH0826 | DXC 800 Calibrators |
| M-F-CH0820 | Chemistry Controls |
| M-F-CH0826 | Chemistry Calibrators |

SPECIMEN

Type of Specimen

Freshly collected urine samples should be used for testing. Collect urine samples in glass or plastic (i.e., polypropylene, polycarbonate, polyethylene) containers. Urine samples should be collected in the manner routinely used for drug screening analysis. Samples should be at room temperature for testing.

Specimen Storage and Stability

If the sample cannot be analyzed immediately, it may be stored at +2°C to +8°C for up to 7 days.² If longer storage is required or when a split sample collection method is used, samples should be stored frozen at -20°C or less.⁵

| Sample Type | Volume | Sample Stability |
|-------------|--------|--|
| Urine | 0.5mL | <ul style="list-style-type: none"> • Analyze at Room Temp • Refrigerated 7 days • Frozen (-20°C) >7 days |

Criteria for Unacceptable Specimens

See Specimen Rejection/Cancellation Protocol

Sample Volume

A filled 0.5 mL sample cup is the optimum volume. For optimum primary sample tube volumes in primary tube samples and minimum volumes, refer to the Primary Tube Sample Template for your system.

REAGENTS

Contents

Each kit contains the following items:
One BARB Reagent Cartridge (1 x 250 tests)

| Volume per Test | |
|----------------------|-----------------------------|
| Sample Volume | 10 uL |
| Total Reagent Volume | 250 uL |
| Cartridge Volumes | A 200 uL B 50 uL C -- |

| Reactive Ingredients | |
|---|-------|
| Antibody/Substrate Reagent | 69 ml |
| Monoclonal anti-Barbiturates Antibody (mouse) | |
| Glucose-6-phosphate | |
| Nicotinamide adenine dinucleotide (NAD) | |
| Tris buffer | |
| Enzyme Conjugate Reagent | 18 ml |
| G6PDH labeled with barbituric acid | |
| Tris buffer | |

Also non-reactive chemicals necessary for optimal system performance

Reagent Preparation

No preparation is required.

Acceptable Reagent Performance

The acceptability of a reagent is determined by successful calibration and by ensuring that quality control results are within acceptance criteria.

Reagent Storage and Stability

BARB reagent when stored unopened at +2°C to +8°C, will remain stable until the expiration date printed on the cartridge label. Once opened, the reagent is stable for 90 days at +2°C to +8°C. Do not use beyond manufacturer's expiration date. DO NOT FREEZE.

CALIBRATION

Calibrator Required

SYNCHRON Systems DAT Negative Urine Calibrator (0 ng/mL secobarbital)
SYNCHRON Systems DAT Multi-Drug Low (cutoff) Urine Calibrator (200 ng/mL secobarbital)
SYNCHRON Systems DAT Multi-Drug High Urine Calibrator (1000 ng/mL secobarbital)

Calibrator Preparation

No preparation is required.

Calibrator Storage and Stability

SYNCHRON[®] Systems Drugs of Abuse Testing (DAT) Urine Calibrators are stable until the expiration date printed on the calibrator bottles if stored capped in the original containers at +2°C to +8°C.

Calibration Information

1. The DAT assays require three levels of calibrators. The calibration measures the separation between calibrators to ensure reagent integrity. NOTICE: The calibration factor generated is non-functional for sample result calculation.
2. The system must have a valid calibrator cutoff value in memory before controls or patient samples can be run. The cutoff value for each DAT chemistry represents the mean reaction rate of the Low Calibrator, and is reported in mA/min units on patient and control reports. Cutoff values are stored in memory until the next successful calibration.
3. Under typical operating conditions the BARB reagent cartridge must be calibrated every 14 days and also with certain parts replacements or maintenance procedures, as defined in the UniCel DxC 600/800 System *Instructions For Use* (IFU) manual. This assay has within-lot calibration available. Refer to the UniCel DxC 600/800 System *Instructions For Use* (IFU) manual for information on this feature.
4. For detailed calibration instructions, refer to the UniCel DxC 600/800 System *Instructions For Use* (IFU) manual.
5. The system will automatically perform checks on the calibration and produce data at the end of calibration. In the event of a failed calibration, the data will be printed with error codes and the system will alert the operator of the failure. For information on error codes, refer to the UniCel DxC 600/800 System *Instructions For Use* (IFU) manual.

TRACEABILITY

For Traceability information refer to the Calibrator instructions for use.

QUALITY CONTROL

RELATED DOCUMENTS

J-F-CH-0820 DXC 800 Controls
M-F-CH-0820 Chemistry Controls

STEPS

1. If necessary, load the reagent onto the system.
2. After reagent load is completed, calibration may be required.

3. Program quality control for analysis
4. After loading controls onto the system, follow the protocols for system operation. To load samples manually refer to the FHS DXC Series Manual Sample Programming procedure. For detailed testing procedures, refer to the UniCel DxC 600/800 System *Instructions For Use* (IFU) manual.

CALCULATIONS

The system performs all calculations internally to produce the final qualitative result, reported as POSITIVE or NEGATIVE. The qualitative result is based on a comparison of the sample rate to the calibrated cutoff rate; a sample rate greater than or equal to the cutoff rate is reported as POSITIVE. A POSITIVE result (≥ 200 ng/mL) from this assay indicates only the presence of Barbiturates and does not necessarily correlate with the extent of physiological and psychological effects. A NEGATIVE test result indicates that Barbiturates are either not present, or are present at levels below the cutoff threshold of the test.

PERFORMANCE CHARACTERISTICS

Reference Range

| Sample Type | Normal |
|-------------|----------------|
| Urine | "Not Detected" |

| Analyte | Normal Results | Suspect Results (Adulterated) |
|------------|----------------|-------------------------------|
| Creatinine | >20 mg/mL | <20 mg/mL |

If results for Creatinine are outside normal parameters for urine drug testing, **DIL1** comment will be appended to screen results.

Reporting results

| Result | Reported As |
|----------|----------------|
| Negative | "Not Detected" |
| Positive | "Detected" |

Relative Sensitivity and Specificity

The cutoff value of the SYNCHRON Systems Barbiturates assay is 200 ng/mL.

LIMITATIONS

1. The test is designed for use with human urine only.
2. Do not dilute the urine samples since this is a qualitative assay. Dilution of samples may produce erroneous results.
3. Interference has been demonstrated from mefenamic acid, a nonopioid analgesic.

4. Adulteration of the urine sample may cause erroneous results. Alteration of a urine specimen may be detected by checking the appearance, temperature, pH specific gravity, and creatinine levels of a sample. If adulteration is suspected, obtain another sample and forward both specimens to the laboratory for testing.
5. An effort should be made to keep pipetted samples free from gross debris. It is recommended that highly turbid specimens be centrifuged before analysis.

Cross Reactivity

Barbiturates and various potential interfering substances in a human urine matrix were tested for cross-reactivity with the SYNCHRON Systems BARB assay. The following table summarizes the results obtained at the concentrations tested for each potential cross-reactant.

| Compound | Concentration (µg/mL) | Effect |
|---------------------------------------|-----------------------|----------|
| Secobarbital (cutoff) | 0.2 | Positive |
| Alphenal | 0.25 | Positive |
| Amobarbital | 0.3 | Positive |
| Aprobarbital | 0.2 | Positive |
| Barbital | 1.5 | Positive |
| Butobarbital | 0.25 | Positive |
| Butalbital | 0.4 | Positive |
| Butethal | 0.3 | Positive |
| Diallylbarbital | 0.6 | Positive |
| Pentobarbital | 0.5 | Positive |
| Phenobarbital | 0.8 | Positive |
| Talbutal | 0.08 | Positive |
| Thiopental | 0.8 | Positive |
| Acetaminophen | 1000 | Negative |
| Acetylsalicylic Acid | 1000 | Negative |
| Albuterol | 1000 | Negative |
| d-Amphetamine | 1000 | Negative |
| Caffeine | 100 | Negative |
| Codeine | 1000 | Negative |
| Diphenhydramine | 500 | Negative |
| Glutethimide | 80 | Negative |
| 5-OH-Phenyl-5-phenyl-hydantoin (HPPH) | 500 | Negative |
| Meperidine | 1000 | Negative |
| Methadone | 1000 | Negative |
| Methsuximide | 100 | Negative |
| Morphine | 1000 | Negative |
| Normethsuximide | 100 | Negative |
| Oxazepam | 500 | Negative |
| Phencyclidine | 1000 | Negative |
| Phenytoin (DPH) | 500 | Negative |
| Propoxyphene | 1000 | Negative |

ADDITIONAL INFORMATION

For more detailed information on UniCel DxC Systems, refer to the appropriate system manual.

REFERENCES

1. National Institute on Drug Abuse Research, "Urine Testing for Drugs of Abuse", Monograph 73 (1986).
2. National Institute on Drug Abuse, "Mandatory Guidelines for Federal Workplace Drug Testing Programs", *Federal Register*, Vol. 53, No. 69 (1988).
3. Rubenstein, K. E., Schneider, R. S., Ullman, E. F., "Homogenous Enzyme Immunoassay: A New Immunochemical Technique", *Biochem. Biophys. Res. Commun.*, Vol. 47, 846 851 (1972).
4. McEnvoy, G. K., Litvak, K., *AHFS Drug Information*, American Society of Hospital Pharmacists, Inc., Bethesda, MD (1992).
5. National Committee for Clinical Laboratory Standards, *Urine Drug Testing in the Clinical Laboratory*, Proposed Guideline, Vol.13, No. 7, NCCLS publication T/DM8-P, Villanova, PA (1993).
6. "USP XXII, NF XVII", United States Pharmacopeial Convention, Inc., Rockville, MD (1990).
7. CDC-NIH manual, *Biosafety in Microbiological and Biomedical Laboratories*, U.S. Government Printing Office, Washington, D.C. (1984).
8. Substance Abuse and Mental Health Service Administration, "Mandatory Guidelines for Federal Workplace Drug Testing Programs", *Federal Register*, Vol. 58, No. 14, (1993).
9. Crane, T., et al., "Mefenamic Acid Prevents Assessment of Drug Abuse with EMIT™ Assays", *Clin. Chem.*, Vol. 39, No. 3, 549 (1993).
10. Tietz, N. W., *Fundamentals of Clinical Chemistry*, 4th Edition, W. B. Saunders Company, Philadelphia, PA (1996).

DOCUMENT APPROVAL Purpose of Document / Reason for Change:

Updated for current process

| | | | |
|--------------------------------|--|---|--|
| Committee Approval Date | <input checked="" type="checkbox"/> Date: 8/26/2014 <input type="checkbox"/> NA – revision of department-specific document which is used at only one facility | Medical Director Approval <i>(Electronic Signature)</i> |  8/26/14 |
|--------------------------------|--|---|--|