

Quality Control in the Coagulation Laboratory

Our Passion.
Your Results.



Objectives


- Define Quality
- Understand how to Establish a QC range
- Review examples of QC errors
- Learn how to use the Levey-Jennings chart to guide troubleshooting of the system
- Review commonly used Westgard rules
- Review how to determine data to be submitted to Accutrak

Responsibilities

- It is the responsibility of the Medical Director of the Laboratory or their designee to determine each laboratory's quality control policy.

What is a Quality Result?

- **Accurate** - shows the “true” answer; on target
- Reflects the *in vivo* state of the patient
- Provided in a **timeframe** that is useful to the clinician




Goal of Quality Control

Catch **all** significant errors without repeating tests unnecessarily

Significant Error: A wrong answer that causes a change in the diagnosis or treatment of a patient; or a proficiency testing failure

Levey-Jennings Plot


- 1950 - Levey and Jennings proposed that quality could be monitored in a quantitative, statistical way
- Recommended that labs use the mean of a set of repetitive values, with acceptable limits set at $\pm 3SD$



Statistical Definitions

- Mean - average of a series of measurements
- SD - standard deviation; A measure of the imprecision; variability of results over time
- $\%CV = SD/mean \times 100$
- Bias - a measure of inaccuracy; How far a series of measurement is from the “true” value

Precision (%CV) and Accuracy (bias)



Precision: How closely multiple values match each other

Accuracy: How close values are to the “true” value

Establishing QC Ranges

- **New Method** : Determine method precision and bias.
- A minimum of 20 values over 20 days optimal, run at various times of the day, by various operators, with various bottles of reagent and controls.
- Screen for outliers and exclude from statistics.
- Calculate mean and SD

New Method QC Range

| QC Level 1 | PT | | | | | | | | |
|------------|----|------|--|------------------|--|--|-------------|--|--|
| | 1 | 11.8 | | | | | | | |
| | 2 | 11.2 | | | | | | | |
| | 3 | 11.9 | | | | | | | |
| | 4 | 11.7 | | | | | | | |
| | 5 | 11.9 | | | | | | | |
| | 6 | 11.5 | Mean | =Average(B2:B21) | | | 11.7 | | |
| | 7 | 11.9 | | | | | | | |
| | 8 | 11.7 | SD | =stdev(B2:B21) | | | 0.20 | | |
| | 9 | 11.6 | 2SD | | | | 0.40 | | |
| | 10 | 11.4 | | | | | | | |
| | 11 | 11.6 | Calculate mean +/- 2SD | | | | 11.3 - 12.1 | | |
| | 12 | 11.9 | | | | | | | |
| | 13 | 11.8 | Exclude outliers >3SD and recalculate if necessary | | | | | | |
| | 14 | 12 | | | | | | | |
| | 15 | 11.7 | %CV | =(H9/H7)*100 | | | 1.7 | | |
| | 16 | 11.8 | Compare %CV to package insert value and Accutrak | | | | | | |
| | 17 | 11.6 | | | | | | | |
| | 18 | 11.8 | | | | | | | |
| | 19 | 11.5 | | | | | | | |
| | 20 | 11.6 | | | | | | | |

Establishing QC Ranges

New lot numbers for existing method:

- Minimum of 20 values over 20 days optimal
- If not possible, 4 values over 5 days
- Tested with various vials of reagent, times of day, operators
- Calculate mean
- SD – if a well-established method, and mean is similar to the previous, can use the historical SD/CV

Historical SD or %CV

Advantage: takes into account slight variation in instrumentation and reagents over time.

- Cumulative data from at least 6 months
- Requires stable performance
 - No drifts or shifts due to reagent deterioration

Example:

- Old lot PT Level 1 mean = 11.5, SD 0.3, %CV 2.6
- New lot PT Level 1 mean = 12.0
 - Calculate historical SD based on CV
 - $0.026 \times 12.0 = 0.3$

Historical CV Example


| | Normal control | | | |
|----------------|----------------|--|----------------------|--|
| | %CV | | | |
| Jan | 5.2 | | | |
| Feb | 3.6 | | Control's | |
| March | 4.8 | | mean value = 95% | |
| Apr | 5.0 | | | |
| May | 5.5 | | SD = 95 x (4.6/100) | |
| June | 4.8 | | | |
| July | 4.2 | | SD = 4.37 | |
| Aug | 4.0 | | | |
| Sept | 5.3 | | 2SD range = 86 - 104 | |
| Oct | 4.5 | | | |
| Nov | 3.7 | | | |
| Dec | 5.0 | | | |
| Average | 4.6 | | | |

Assayed QC - Expected Values

“The reported ranges were determined over multiple runs on IL Coagulation ...Systems using specific lots of reagents. **The mean of the control range determined in your laboratory may vary due to the lot of reagent used.**”

“Due to differences in reagents and instrumentation, each laboratory should establish its own Target Value and Acceptance Range (mean and standard deviation). However, any properly functioning coagulation system should yield **mean** values within the Acceptance Range on the package insert.”

Is This Excellent QC Performance?




Is This Acceptable QC performance?


| | | | |
|---------------------|---------------------|-----------------------|---------|
| Test Code: | APTT-SS | Data Interval: | 60 Days |
| QC Material: | Normal C. Unassayed | | |

| | |
|----------------|------|
| Target: | |
| Mean: | 30.0 |
| SD: | 1.4 |
| Unit: | s |

| | | | |
|-------------------|------|--------------------|-----|
| Statistic: | | | |
| Mean: | 29.7 | N: | 133 |
| SD: | 1.3 | Omitted: | 0 |
| CV: | 4.4 | Total Data: | 133 |




Levey-Jennings Graph Showing Abrupt Shift



Levey-Jennings Showing a Shift

QC Material: High Abn C Unassayed Data Interval: 30 Points


| | | | | | |
|----------------|------|-------------------|------|--------------------|----|
| Target: | | Statistic: | | | |
| Mean: | 54.0 | Mean: | 53.5 | N: | 30 |
| SD: | 3.0 | SD: | 4.1 | Omitted: | 0 |
| Unit: | s | CV: | 7.6 | Total Data: | 30 |



Common Causes of QC shifts

- Change in reagent preparation, e.g. incorrect reconstitution volume
- New calibration curve
- Change in analyzer, e.g. new syringes, PM done
- Incorrect lot number on board

Levey-Jennings with a Trend




Common Causes of Trends

- Changes in a reagent or control, e.g. deterioration
- Contamination of the analyzer
- Minor trends in control or reagent are expected over the stability of the materials.
- Gradual wear of an instrument part, e.g. syringes seals.

Trends alert us to possible future failures

Random Error - Imprecision


| | | | |
|----------------------------------|-----------------------|------------------------|-----------------------|
| Instrument Coag | File ID PT PTT 2 | Level 2 | Lot Opened 09/13/2009 |
| Lot Number NO987358 | Location PMH Main Lab | Lot Expires 09/30/2010 | Lot Closed 00/00/00 |
| Protime QC (10/01/09 - 10/30/09) | | | 102 points |



File Test Comment: [09/13/09 11:22 1007T] Exp Mean modified. [31.1]->[25.9], Exp Stdev modified. [1.137]->[1.3], Reason: incorrect mean and sd entered

| File/Dates | M(E) | SD(E) | CV(E) | Range(E) | M(O) | SD(O) | CV(O) | Range(O) |
|------------------|-------|-------|-------|-----------|-------|-------|-------|-------------|
| PT PTT 2 Current | 25.90 | 1.30 | 5.01 | 23.3-28.5 | 26.06 | 1.38 | 5.32 | 23.29-28.84 |

Levey-Jennings Graph Showing a Statistical Outlier



Causes of QC Failures

Our Passion.
Your Results.




Mistakes in Reagent or QC preparation

- Storage or shipping temperature failures
- Reconstitution errors:
 - Pipetting
 - Water
 - Mixing
 - TIMING




Incorrect Mean




Incorrect Calibration Curve

Regression Type: Linear regression
X - Transformation: $\ln(x)$
Y - Transformation: $\ln(y)$
Y-Intercept: 8.118E+00
Slope: -1.130E+00
r2: 0.991
Adjusted Y:
NPP:
NPP Value:

● Normal points ● Edited points — Calibration curve
✕ Deleted points ● Omitted points — Adjusted calibrated curve



Mechanical Failures



Improperly Performed Maintenance

- Or...maintenance not performed
- Incorrect solution used for cleaning or rinsing
- Syringes not correctly seated after replacement
- Incorrect probe/needle alignment

Bubbles!



QC Failure Classification and Management

Our Passion.
Your Results.



Corrective Actions

- Rerun – if using 2SD rule only
- Make fresh reagent – both levels out?
- Make fresh control material – one level out?
- Perform maintenance – clean, avoid contamination
- Look for trends/shifts
- Other considerations (Calibration, H₂O, pipettes)
- Consider instrument malfunction
- Consider instrument failure

How to Handle an Outlier (within the +/- 3 SD Range)

If using Westgard multi-rules – not an automatic failure

If using only a 2SD range:


- If the numeric results for the outlier is between 2 SD and 3 SD, this is very likely a statistical outlier. If the graphical display does not show a 3 SD range, the operator must run the control again.
- If the result of the rerun is now “in control” (falls within 2 SD) the analytic run is accepted and **both QC runs (the statistical outlier and the acceptable rerun) are to be included in the control data base.**

How to Handle an Outlier (within the +/- 3 SD Range)...cont'd

- If the rerun of the control reflects the same parameter as still being out of control between the 2 SD and 3 SD limits, (fails 2-2s) this reflects systematic error and the operator should troubleshoot
- Look for shifts or trends, are the mean and SD correct?


Error Classification

- Random – spuriously occurs, not consistent. Indicates imprecision. – high CV



Error Classification

- Systematic – an ongoing error such as a shift or drift in results. Indicates inaccuracy



Westgard Rules Overview

Our Passion.
Your Results.




1980's Westgard Rules

- N = number of control observations
- Subscript = rule that is violated
 - e.g. 1_{2s} - one control level out by more than 2SD
- Multiple rules written with a slash between, e.g.
 - $1_{2s}/1_{3s}/2_{2s}/R_{4s}/4_{1s}/10_x$
- Rules may be designated as rejection rules or as warning rules

Westgard JO, Barry PL, Hunt MR, Groth T. A multi-rule Shewhart chart for quality control in clinical chemistry. *Clin Chem* 1981;27:493-501.

Commonly Used Westgard Rules


1_{3s} - Reject Run Random error



Commonly Used Westgard Rules

22S -Reject Run

Systematic error




2 consecutive values outside the same 2SD (across runs) OR 2 consecutive values outside the same 2SD (within run)

Commonly Used Westgard Rules

4_{1s} – Warning (reject at lab's discretion)


Systematic error



4 consecutive values on one side of the mean and more than 1SD from the mean; one control over 4 runs, or 2 controls across 2 runs

Commonly Used Westgard Rules

10_x – Warning rule Systematic error




10 consecutive values on ONE side of the mean;
can be one control across 10 runs, or 2 controls
over 5 runs

Commonly Used Westgard Rules


R_{4S} - Reject Run (Warning rule at lab's discretion)

Random Error



The range between 2 controls **within run** exceeds 4 SD

Multi-rule Flow Chart 2 Controls per Run



This scheme will lead to <1% false rejects

QC Example 1

Technologist is performing an AT activity assay using 2 levels of controls

- Control 1 = 97% (2SD limit is 95 to 115%)
- Control 2 = 29% (2SD limit is 30 - 40%)

Should the run be accepted?

- Must look at the previous run's result for that control.
- Control 2 read 31% on the previous run. Passes 2_{2S} .
- Looking back, the previous 9 runs have read below the mean. Fails 10x; Accept run, but troubleshoot shift.

QC Example 2

Technologist is performing an APTT assay using 2 levels of controls.

- Control 1 = 35 sec (2SD range = 27 - 33 secs)
- Control 2 = 56 sec (2SD range = 51 - 57 secs)

Should the run be accepted?

- Must look at the previous run's result for that control.
- Previous result for Control 1 was 36 secs.
- Fails 2_{2S} rule. Reject Run.

Important points when using Westgard Rules

- The mean and SD must be current and accurate for the lot numbers in use.
- The mean and SD must have been based on an acceptable number of data points; minimum of 20 values obtained under normal testing conditions
- Westgard rules should not be used when using a manufacturer's package insert range on an assayed control

QC File Management

- The frequency of outlier occurrences should be carefully monitored. A high outlier frequency may indicate excessive or systematic error.
- All control results less than 3 SD should be used in statistical calculations.
- In some labs it is still common practice to omit the outliers; however this practice may artificially narrow the control range if recalculated.

Accutrak Reporting


Our Passion.
Your Results.



Submitting Data to Accutrak Monthly

- Filter QC data to show only the last month
- Review points for outliers ($\pm >3SD$) and omit from statistics
- Enter the ACTUAL mean and SD for the month, not the TARGET.

Actual Mean and SD on an ACL TOP® System



Actual Mean and SD on an ACL ELITE® System

Analysis QC Calibration Diagnostic Setup Utility

QC Review

| LIQUID ID | CONFIGURED TESTS | UNIT | |
|------------------------------------|------------------|-----------------------|---|
| Abnormal Control Plasma - Level II | PT HS | S | |
| Assess High Heparin Control | APCR-V | | ← |
| Abnormal Control Plasma - Level I | PCX | ACTUAL MEAN | ← |
| Normal Control Plasma | TT | TARGET MEAN | ← |
| Assess Low Heparin Control | AT-III | ACTUAL SD | ← |
| | F9SP-P | TARGET SD | ← |
| | | ACTUAL CV | ← |
| | | RESULTS IN STATISTICS | ← |
| | | RESULTS IN DB | |

Plot and Statistics Cumulative Results Setup Show Enabled

✓

Accutrak Data

Using Accutrak to assess accuracy

- Data is now reagent lot number specific in the full report
- A minimum of 20 labs is statistically significant
- Confidence limits of the peer group is a range of means

Accutrak Data

Using Accutrak to assess precision

- Average peer group CV may not reflect all reagent lot numbers in use

Reasons Control Data May Flag in Accutrak

- An instrument problem
- A reagent problem
- A control storage and/or handling problem
- A change in laboratory environment, including pipettes/water
- Has this problem been fixed already?

Accutrak Report Message:

“Your Lab’s 2 SD range is wider than the Group’s 3 SD range.”

- Verify lot numbers
- Instrument maintenance
- Instrument problem
- Isolate cause of Imprecision

Conclusions

- QC ranges must be correctly established and maintained to catch errors before they occur
- One can guide troubleshooting based on the Levey-Jennings graph and type of error – random vs. systematic
- Comparison of mean and CV with peer group labs via Accutrak is a valuable tool

References

- CLSI C24-A3. Statistical quality control for quantitative measurements: Principles and definitions; Approved guideline - second edition. NCCLS, Wayne, PA, 2006.
- CLSI EP23A. Laboratory Quality Control Based on Risk Management. CLSI, Wayne, PA, 2011.
- Westgard JO, Barry PL, Hunt MR, Groth T. A multi-rule Shewhart chart for quality control in clinical chemistry. *Clin Chem* 1981;27:493-501.
- Levey S, Jennings E. The Use of Control Charts in the Clinical Laboratory. *Am J Clin Pathol* 1950;20:1059-66..
- www.westgard.com/mltirule.htm

Additional Reading

- Quality in Laboratory Hemostasis and Thrombosis edited by Kitchen, Steve, et al. Blackwell Publishing Ltd. 2009