



STANDARDIZED TESTING / OPERATING PROTOCOL REQUEST/ANNOUNCEMENT

Mini-Vidas Interface Activation

Description:	Alamance will begin using the interface for reporting from the mini-Vidas
Implementation Date:	2/27/2018
Staff Update:	Staff will use interface ARVID to view results for release in the LIS. Qualifying results will auto-release.
Performing Locations:	Click on the boxes that apply: □ Alamance Cancer Center ☑ Alamance Regional □ Annie Penn Hospital □ Moses Cone Hospital □ Med Center at High Point □ Med Center at Mebane □ Wesley Long Hospital □ Women's Hospital
Affected Locations:	Click on the boxes that apply: □Alamance Cancer Center □Alamance Regional □Annie Penn Hospital □Moses Cone Hospital □Med Center at High Point □Med Center at Mebane □Wesley Long Hospital □Women's Hospital

	Click on the boxes that apply:
Affected Departments:	□Blood Bank □Cytology □Flow Cytometry □Histology □Microbiology □Phlebotomy □Point of Care ⊠Rapid Response Lab □Respiratory Therapy □Specimen Processing
Specimen Type:	Appropriate specimen for each test type
Updated Clinical Lab Procedures:	COAG-715 DDimer COAG-716 PCT LIS-113-CH Result Reporting and Report Format
Retired Clinical Lab Procedures:	None
Training/Competency:	Is training required? ☐ Yes ☒ No If yes, enter date complete: Do training checklists need to be created or updated? ☐ Yes ☒ No If yes, enter date complete: Do training quizzes need to be created or updated? ☐ Yes ☒ No If yes, enter date complete: Do competency/direct observation forms need to be created or updated? ☐ Yes ☒ No If yes, enter date complete:

	Click on the boxes that apply: ⊠Section Not Applicable □Memo Needed
Notification to Client:	Distribution of Memo: Medical Staff Allied Health Professionals (PA, Nurse Practitioners) Anesthesia Annie Penn (Primary Source Physicians) Dentist Emergency Department/Urgent Care Centers Family Practice Infectious Docs #ID Docs (John Campbell, Robert Comer, Jeffrey Hatcher, Cynthia Snider, Kees Van Dam) OB/GYN Pathology Pediatricians Psych Radiology Surgery #Nursing Leadership (Directors, Asst. Directors, Clinical Nurse Manager) Pharmacy - Send to DeAnne Brooks & Jim Hasspacher #IM Residents Kim Helsabeck Phlebotomy Managers and Supervisors Point of Care: Sheila, Kim & Marty
Accreditation Section:	Click on the boxes that apply: Section Not Applicable □CAP Test menu change needed □CMS Analyte form change needed □Proficiency Testing surveys changes needed or ordered
Laboratory IT section:	Click box and type needed changes/additions: □Section Not Applicable □LIS changes □Reference range change/addition Reporting units for D dimer will change from whole numbers to 2 decimal places.

	
	□Technical Failure change/addition
	□Critical Value change/add
	□Text comments needed
	□Specimen collection instructions
	□Need to monitor TAT
	□CPT code for tests(s)
Performing Location Post Implementation Tasks	Monitor interface to ensure proper function. Complete autovalidation verifications annually.
Laboratory IT Post Implementation Tasks	Verify patient from miniVIDAS to Sunquest to Cone HealthLink post go live for both PCT and DDIMR.
STOP Initiator:	Wendy Turner
Alamance Medical	Quality Department will obtain signature:
Director	
Signature/Date:	Julian Rulas 1.22.18
MedCenter Mebane	
Medical Signature/Date:	N/A
Oignature/Date.	
Greensboro/Reidsville	
Medical Director	N/A
Signature/Date:	
	Performing Location Post Implementation Tasks Complete and Acceptable
Post Implementation Complete	Date:
	QA Department Signature:(Attach evidence of completion)
(Completed by QA Department before filing)	Laboratory IT Post Implementation Tasks Complete and Acceptable
	Date:
	QA Department Signature:
	(Attach evidence of completion)

¥392000839 92/07/2018 1348 212973 U1484 TEST . KATE

Alamance Regional Medical Center 1240 Huffman Mill Rd. Burlington, NC 27215

INTERIM REPORT

NAME: TEST, KATE

MRN : 212973 ACCT: 100034160 LOC: LDA ROOM: LDR6-1

W1404 COLL: 02/07/2018 13:48 REC: 02/07/2018 14:18 PHYS: LESCHBER, VALERI

AGE: 29Y SEX: P DOB: 02/21/1988

MD: JACKSON, STEPHEN D

PROCALCITONIN

PROCALCITONIN

0.58

STAT

ng/mL {AH}

Interpretation: PCT > 0.5 ng/mL and <= 2 ng/mL: Systemic infection (sepsis) is possible, but other conditions are known to elevate PCT as well.

(NOTE)

ICU PCT Algorithm

Non ICU PCT Algorithm

PCT < 0.25 ng/mL

PCT < 0.1 ng/mL

Stopping of antibiotics strongly encouraged.

Stopping of antibiotics strongly encouraged.

PCT level decrease by >= 80% from peak PCT OR PCT 0.25 - 0.5 ng/mL PCT < 0.25 ng/mL

Stopping of antibiotics encouraged.

Stopping of antibiotics encouraged.

PCT level decrease by < 80% from peak PCT AND PCT >= 0.5 ng/mL

PCT >= 0.25 ng/mL Continuing antibiotics encouraged.

Continuing antibiotics ' encouraged.

PCT level increase compared

PCT > 0.5 ng/mL

with peak PCT AND PCT >= 0.5 ng/mL

Escalation of antibiotics strongly encouraged.

Escalation of antibiotics

strongly encouraged.

{AH} - Performed at Alamance Hospital Lab, 1240 Huffman Mill Rd., Burlington, NC 27215

TEST, KATE PRINT DATE: 02/07/2018 END OF REPORT

PRINT TIME: 14:24

Mylus Billon 2.22.17

Alamance Regional Medical Center 1240 Huffman Mill Rd. Burlington, NC 27215

INTERIM REPORT

NAME: TESTING, LAB

MRN : 123456789

W1391

LOC: EDA

AGE: 31Y SEX: M

DOB: 07/15/1986 ACCT: 111 MD: DEPAULT. PROVIDER

FIBRIN DERIVATIVES D-DIMER (AR

{AH}

H 547.21

[0.00-499.00]

Justulus 2.22-18

(NOTE)

<> Exclusion of Venous Thromboembolism (VTE) - QUIPATIENT ONLY (Emergency Department or Mebane)

0-499 ng/ml (FEU): With a low to intermediate pretest probability

COLL: 02/07/2018 08:46 REC: 02/07/2018 08:47 PHYS: DEFAULT, PROVIDE

for VTE this test result excludes the diagnosis

of VTE.

>499 ng/ml (FEU) : VTE not excluded; additional work up for VTE is

required.

<> Testing on Impatients and Evaluation of Disseminated Intravascular Coagulation (DIC) Reference Range: 0-499 ng/ml (FEU)

{AH} = Performed at Alamance Hospital Lab, 1240 Huffman Mill Rd., Burlington, NC 27215

> 123456789 EDA TEST ING. LAB 97/15/1986 M ARV1D

62/07/2018 mine

TESTING, LAB PRINT DATE: 02/07/2018

END OF REPORT PRINT TIME: 09:11 PAGE ٦.