

## Performing a Wet Mount Examination

---

**Principle** This procedure provides instructions for performing a microscopic direct examination using phase contrast for identification of clue cells, yeast, and Trichomonas in a saline suspension of vaginal secretions, along with identification of yeast and Trichomonas in a variety of other specimens.

---

**Policy** Wet mounts are to be performed STAT and read within 1 hour of collection.

- **Order Code: WET (Wet Mount Smear)**

---

**Specimen**

- Vaginal specimens, urethral discharges, prostatic secretions and sediments of fresh urine may be examined for the presence of yeast and Trichomonas. Clue cell presence is only reported for vaginal specimens.
- Flocked collection swabs are the preferred swab. Rayon swabs are also acceptable.
  - Dry or cotton swabs will be rejected.
  - Flocked swabs received in Eswab containers or Universal Viral Transport (UVT) media will be rejected.
- Swabs must be received in a properly labeled sterile container with at least 0.5 ml of sterile saline.
- Specimens must be transported to the lab at room temp within 1 hr of collection.

---

**Equipment/Supplies**

- Phase contrast microscope
- Centrifuge (for fluid specimens only)
- Glass microscope slides and coverslips
- Plastic transfer pipettes

---

**Specimen Preparation** Prepare the specimen based upon specimen type:

<b>If</b>	<b>Then</b>
Collection swab in 0.5-1.0 ml of saline	<ol style="list-style-type: none"><li>1. <b>Vortex</b> tube for 3-5 seconds</li><li>2. Place 1-2 drops of liquid onto labeled glass slide (Do not use swab to roll onto slide)</li><li>3. Apply cover slip</li></ol>
Collection swab in >1.0 ml of saline or Body Fluid >1.0 ml	<ol style="list-style-type: none"><li>1. <b>Vortex</b> tube for 3-5 seconds (if swab sample)</li><li>2. Centrifuge specimen for 5 minutes</li><li>3. Remove supernatant to bring sample to at least 1.0 ml</li><li>4. Resuspend sediment and place 1-2 drops of liquid onto labeled glass slide</li><li>5. Apply cover slip</li></ol>
Body Fluid <1.0 ml	<ol style="list-style-type: none"><li>1. Mix sample well</li><li>2. Place 1-2 drops of liquid onto labeled glass slide</li><li>3. Apply cover slip</li></ol>

---

*Continued on next page*

## Performing a Wet Mount Examination

**Procedure** Follow the steps to perform Wet Mount examination.

Step	Action												
1.	Using phase contrast microscopy, examine the <i>entire slide</i> on low power (10X) for the presence of motile Trichomonas and pseudohyphae. Using low power is also helpful to evaluate the distribution of cellular material on the slide.												
2.	Re-examine the <i>entire slide</i> on high power (40X) to evaluate the presence or absence of clue cells, weakly motile Trichomonas, and yeast (budding/non-budding, and pseudohyphae).												
3.	Record findings on the “Wet Mount Report” form, using the following grading criteria: <table border="1" data-bbox="578 753 1162 984" style="margin-left: 40px;"> <thead> <tr> <th>Amount Seen</th> <th>Report</th> </tr> </thead> <tbody> <tr> <td>&gt;2 per High Power Field</td> <td>Many</td> </tr> <tr> <td>1 to 2 per High Power Field</td> <td>Moderate</td> </tr> <tr> <td>4 to 10 per entire slide</td> <td>Few</td> </tr> <tr> <td>1 to 3 per entire slide</td> <td>Rare</td> </tr> <tr> <td>0 per entire slide</td> <td>None Seen</td> </tr> </tbody> </table> <p>Also include the following information on the form:</p> <ul style="list-style-type: none"> <li>• Patient Name/MRN/Accession number</li> <li>• Source</li> <li>• Tech Code</li> <li>• Date/Time testing performed</li> </ul>	Amount Seen	Report	>2 per High Power Field	Many	1 to 2 per High Power Field	Moderate	4 to 10 per entire slide	Few	1 to 3 per entire slide	Rare	0 per entire slide	None Seen
Amount Seen	Report												
>2 per High Power Field	Many												
1 to 2 per High Power Field	Moderate												
4 to 10 per entire slide	Few												
1 to 3 per entire slide	Rare												
0 per entire slide	None Seen												
4.	Discard specimen in the appropriate biohazard waste container.												

*Continued on next page*

## Performing a Wet Mount Examination

**Reporting** Follow the steps to result patient results in LIS.

Step	Action														
1.	In Sunquest Microbiology Result Entry function, enter the patient accession number to be resultted.														
2.	<p>Using designated result keys &lt;&gt;, enter Wet Mount results under the Direct Exam tab. Press &lt;F8&gt; to display the keyboard on screen. Enter all findings that are present before entering findings that are absent using the approved corresponding keys below:</p> <ul style="list-style-type: none"> <li>• The presence or absence of Yeast or Trichomonas must always be reported.</li> <li>• For Vaginal specimens only: also report the presence or absence of clue cells.</li> <li>• Do not report the presence or absence of any other cell type/findings.</li> </ul> <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th style="width: 50%;">If</th> <th style="width: 50%;">Then enter</th> </tr> </thead> <tbody> <tr> <td>Yeast, Trichomonas or Clue Cells are seen</td> <td> <ul style="list-style-type: none"> <li>• Quantity followed by finding description:                             <ul style="list-style-type: none"> <li>○ Rare &lt;5&gt;</li> <li>○ Few &lt;6&gt;</li> <li>○ Moderate &lt;7&gt;</li> <li>○ Many &lt;8&gt;</li> <li>○ Trichomonas &lt;T&gt;</li> <li>○ Yeast &lt;Y&gt;</li> <li>○ Clue Cells &lt;U&gt;</li> </ul> </li> <li>• Repeat for each finding seen</li> </ul> <p><i>Example:</i></p> <ul style="list-style-type: none"> <li>• <i>Many &lt;8&gt; Yeast &lt;Y&gt;</i></li> <li>• <i>Few &lt;6&gt; Trichomonas &lt;T&gt;</i></li> </ul> </td> </tr> <tr> <td>No Yeast or Trichomonas are seen</td> <td>No Yeast or Trichomonas seen &lt;B&gt;</td> </tr> <tr> <td>No Trichomonas, clue cells, or yeast are seen (Vaginal specimens)</td> <td>No Trichomonas, clue cells, or yeast seen &lt;N&gt;</td> </tr> <tr> <td>No Trichomonas are seen</td> <td>No Trichomonas seen &lt;F&gt;</td> </tr> <tr> <td>No Yeast are seen</td> <td>No Yeast seen &lt;H&gt;</td> </tr> <tr> <td>No Clue Cells are seen</td> <td>No Clue Cells seen &lt;G&gt;</td> </tr> </tbody> </table>	If	Then enter	Yeast, Trichomonas or Clue Cells are seen	<ul style="list-style-type: none"> <li>• Quantity followed by finding description:                             <ul style="list-style-type: none"> <li>○ Rare &lt;5&gt;</li> <li>○ Few &lt;6&gt;</li> <li>○ Moderate &lt;7&gt;</li> <li>○ Many &lt;8&gt;</li> <li>○ Trichomonas &lt;T&gt;</li> <li>○ Yeast &lt;Y&gt;</li> <li>○ Clue Cells &lt;U&gt;</li> </ul> </li> <li>• Repeat for each finding seen</li> </ul> <p><i>Example:</i></p> <ul style="list-style-type: none"> <li>• <i>Many &lt;8&gt; Yeast &lt;Y&gt;</i></li> <li>• <i>Few &lt;6&gt; Trichomonas &lt;T&gt;</i></li> </ul>	No Yeast or Trichomonas are seen	No Yeast or Trichomonas seen <B>	No Trichomonas, clue cells, or yeast are seen (Vaginal specimens)	No Trichomonas, clue cells, or yeast seen <N>	No Trichomonas are seen	No Trichomonas seen <F>	No Yeast are seen	No Yeast seen <H>	No Clue Cells are seen	No Clue Cells seen <G>
If	Then enter														
Yeast, Trichomonas or Clue Cells are seen	<ul style="list-style-type: none"> <li>• Quantity followed by finding description:                             <ul style="list-style-type: none"> <li>○ Rare &lt;5&gt;</li> <li>○ Few &lt;6&gt;</li> <li>○ Moderate &lt;7&gt;</li> <li>○ Many &lt;8&gt;</li> <li>○ Trichomonas &lt;T&gt;</li> <li>○ Yeast &lt;Y&gt;</li> <li>○ Clue Cells &lt;U&gt;</li> </ul> </li> <li>• Repeat for each finding seen</li> </ul> <p><i>Example:</i></p> <ul style="list-style-type: none"> <li>• <i>Many &lt;8&gt; Yeast &lt;Y&gt;</i></li> <li>• <i>Few &lt;6&gt; Trichomonas &lt;T&gt;</i></li> </ul>														
No Yeast or Trichomonas are seen	No Yeast or Trichomonas seen <B>														
No Trichomonas, clue cells, or yeast are seen (Vaginal specimens)	No Trichomonas, clue cells, or yeast seen <N>														
No Trichomonas are seen	No Trichomonas seen <F>														
No Yeast are seen	No Yeast seen <H>														
No Clue Cells are seen	No Clue Cells seen <G>														
3.	Enter each observation on a separate line. Press Tab to move from one line to another.														
4.	On a separate line, enter the location of the performing hospital by pressing “,” to turn OFF keyboard. Then enter the appropriate ETC (RVPER = <i>Performed at Sutter Roseville Medical Center</i> ).														

*Continued on next page*

## Performing a Wet Mount Examination

---

### Reporting (continued)

Step	Action
5.	Final the report by clicking the period key <.>.
6.	Before selecting Save, ensure that findings for Trichomonas and yeast, along with clue cells (vaginal specimen), have been reported. If a correction needs to be made, press “.” in the result field to remove final status and make corrections. The repeat Step 5 to finalize result.
7.	Select Save to finish resulting of Wet Mount.
8.	Verify correct manual entry of results.
9.	Retain Wet Mount Report in the designated area of the laboratory.

---

**Reference Range**

- No Yeast or Trichomonas
- No Clue Cells (vaginal specimens)

---

**Procedure Notes**

- **Trichomonas vaginalis:**  
Pear shaped organisms that are larger than WBCs, slightly refractile with a granular appearance and have long flagella. They move with a rapid, jerky, rotating and non-directional leaf-like motion. Their numbers vary in clinical specimens, so it is essential to scan the entire slide. Trichomonas has limited survival outside the host and examining the specimen immediately is critical. As viability decreases they become spherical and non-motile. Identification should be made upon visualization of flagella movement.

- **Yeast/Pseudohyphae:** Yeast will appear refractile and can be seen as a single cell with or without a bud and is about the same size as an RBC. Hyphae elements appear as thick walled tube like structures that may be single or branching. Do not confuse hyphae with natural fibers from a swab. Candida albicans is the most frequently encountered Candida sp. Associated with vaginal infections but the organism does occur in low numbers as part of the normal vaginal flora.

- **Clue Cells:**  
Epithelial cells of vaginal origin covered with bacteria. The bacterium coating the cells gives the cytoplasm a characteristic refractile, stippled, or granular appearance with shaggy or bearded cell borders. Most of the cell surface should be covered by bacteria for it to be identified as a clue cell. Epithelial cells may have occasional irregular kerato-hyaline granules in the cytoplasm that should not be mistaken for adherent bacteria. Bacterial vaginosis (BV) is the most common type of vaginal infection and is associated with the presence of clue cells. A key microscopic feature of BV is the lack of inflammatory cells and the absence of or reduction in normal flora.

---

*Continued on next page*

## Performing a Wet Mount Examination

---

**Attachments** ● Form A: Wet Mount Report

---

- References**
- Bailey and Scott's *Diagnostic Microbiology*, 12<sup>th</sup> Edition, Mosby Elsevier, 2007.
  - Braude, Abraham L. *Medical Microbiology and Infectious Diseases, Volume II*. WB Saunders Company.
  - Henry, JB. *Clinical Diagnosis and Management by Laboratory Methods*, WB Saunders Company, Philadelphia, 2001.
  - Garcia, LS. Buickner, DA. *Diagnostic Medical Parasitology*, 2<sup>nd</sup> Edition, American Society for Microbiology, Washington, D.C., 1993.
  - Kern, M. *Medical Mycology*, FA Davis Company, Philadelphia, 1990.
  - Egan, ME. Lipsky, MS. *Diagnosis of Vaginitis*. *Amer Fam Physician* 2000; 62: 1995.
- 

*End*