




Current Status: Pending		PolicyStat ID: 9061592	
 Sutter Health Sutter Roseville Medical Center	Origination:	N/A	
	Effective:	Upon Approval	
	Final Approved:	N/A	
	Last Revised:	N/A	
	Next Review:	2 years after approval	
	Owner:	Nadera Poirier: Spvr, Transfusion Services	
	Policy Area:	Lab - Transfusion Service	
References:			
Applicability:	Sutter Roseville Medical Center		

Performing a Rh Blood Type by Tube Method TS.ANA 10.17-RV

Performing a Rh Blood Type by Tube Method

Purpose	This procedure describes how to perform an Rh blood type.		
Policy	An Rh Type and ABO group must be performed on each crossmatch specimen, on each admission when FFP is requested and on all patients without a historical blood type when requests for other components are requested. Weak D (Du) testing must be performed on all patients with negative D results in the following situations: <ul style="list-style-type: none"> • Cord Blood/Neonate <7 days old • Fetal Bleed evaluation • Allogeneic BMT (Bone Marrow Transplant) recipient/donor 		
Equipment/ Reagents/ Supplies			
Equipment	Reagents	Supplies	
<ul style="list-style-type: none"> • Centrifuge • Agglutination viewer 	<ul style="list-style-type: none"> • Anti D • ABO+Rh Control 	<ul style="list-style-type: none"> • Test tubes • Test tube rack • Disposable pipettes 	
Specimen Requirements	<ul style="list-style-type: none"> • Adults: One (1) 10 ml EDTA stoppered tube (preferred) or red stoppered tube without separator gel- minimum of 3ml. • Infants: Two EDTA microtainers or 1.0 ml EDTA or red stoppered tube without separator gel • Donor red blood cells 		
Quality Control	Daily Reagent Quality Control		
Procedure: Rh Typing	Follow the steps below to perform Rh typing:		
Step	Action		
	Label one test tube for Anti-D with patient and reagent identifiers.		
	Add one (1) drop of Anti-D in the appropriately labeled test tube.		
	Prepare a 3-5% patient cell suspension.		

Add one (1) drop of the cell suspension to the tube containing the Anti-D and mix gently.		
Centrifuge.		
Gently resuspend the cell button and examine macroscopically for agglutination.		
Grade and immediately record the result.		
Save tube for use in Weak D (Du) testing, if indicated. (See policy)		
If ABO type is type AB, perform an Rh control. Repeat steps 3 and 4 above substituting Anti-D with Monoclonal Control.		
Interpreting Results		Follow the steps in the table to interpret the Rh test.
If		Then
Cells agglutinate and/or are hemolyzed		Positive reaction
Cells do not agglutinate/and or no hemolysis is present		Negative reaction
If:		Then
Anti D	D Control	
> 2+	0 or ND	Report
1-2+	0 or ND	Check microscopically for mixed field agglutination <i>If present, result as "mf" and free text a BBCNC comment with grade of reactivity</i>
If original testing in tube <2+	0 or ND	Result as Rh weak D
If original testing by ECHO <2+		Proceed to <i>Resulting Tests Using the ECHO Interface</i>
+ or 0	+	Invalid test
Procedure: Current results vs. historical or confirmation		Follow the steps below to compare current Rh type to historical/confirmation results:
Step	Action	
1.	Compare current Rh type result to previous/confirmation.	
2.	If: Current Rh result does match previous/confirmation Rh type	Then: Report results
	Current Rh result does not match previous/confirmation Rh type	<i>Refer to Resolving Patient History Discrepancies</i>
Reporting Results		Enter results and Interpretation directly in LIS testing grid. mf=Mixed Field NA= Not applicable ND= Not done NR= Not required INCL= Inconclusive
D	DC	Report Interp
2-4+	0 or ND	P
<2+	0 or ND	

mf	0 or ND	See table below for Interpretation
0	ND	N
1-4+	1-4+	INVAL
If:		Then:
Rh resulted as INVAL		Transfuse Rh negative
Mixed field agglutination is present		<ul style="list-style-type: none"> • Consideration should be given to discrepancies resulting from circulating mixed cell types • If no history of Rh Type prior to non-identical Rh type RBC transfusion, Rh Type cannot be reported. • Transfuse Rh Negative units until Rh type can be determined • If previous history of Rh type prior to transfusion of non-identical Rh type RBC, result Rh same as previously resulted
Related Procedures		<ul style="list-style-type: none"> • Performing ABO Blood Group • Performing Weak D Rh Type • Daily Reagent Quality Control Procedure • Resolving Patient History Discrepancies
Reference	Manufacturer's product insert Blood Grouping Reagent Anti-D and ABORh Control	

All revision dates:

Attachments

No Attachments

Approval Signatures

Step Description	Approver	Date
Laboratory Director	Lindsey Westerbeck: Dir, Lab	pending