

Status **Active** PolicyStat ID 12203659



Sutter Roseville Medical Center

Origination 1/4/2021
Final 9/12/2022
Approved
Effective 9/12/2022
Last Revised 9/12/2022
Next Review 9/11/2024

Owner Alex Alba:
Supervisor,
Laboratory
Analytic
Policy Area Lab - Coag
Applicability Sutter Roseville
Medical Center

Running Quality Control on the TEG 6s Hemostasis Analyzer

PURPOSE

This procedure describes how to run quality control on the TEG 6s hemostasis analyzer

POLICY

- The external QC is run upon receipt of a new shipment or new lot of test cartridges
- Citrated global hemostasis test cartridge uses the following external QC material: Level 1 normal and level 2 abnormal QC vial kits
- Citrated global hemostasis with lysis (Trauma) test cartridge uses the following external QC material: Level 1 normal and level 2 abnormal QC vial kits
- Platelet Mapping test cartridge uses the following external QC material: normal donor and abnormal donor sample
- The appropriate corrective action is performed on any QC outlier and documented on the TEG 6s QC log
- The internal built-in QC validates that all components and functions of the analyzer-test cartridge combination are operating satisfactorily:
 - Power on self test- when analyzer is turned on, a system QC check is performed on all analyzer functions prior to loading of the test cartridge. Any failure will generate an error message and test will be invalidated.
 - Pre Test - when the test cartridge is inserted, a system QC check is performed to validate the ability of the analyzer to correctly identify the test cartridge and validate the integrity of the analyzer-cartridge interface. Any failure will generate an error message and prevent further system operation.
 - In Use - when the sample is added and test started, a system QC check is performed

to monitor critical operational parameters throughout the duration of the test. If any parameter is out of range then test will be invalidated.

SCOPE

All CLS and MLT staff assigned to the Hematology/Coagulation department

PREPARATION AND STABILITY OF EXTERNAL QC

| External QC | Test Cartridge | Storage and Stability | Preparation |
|----------------------------|--|--|---|
| Level 1 Normal QC vial kit | <ul style="list-style-type: none"> • Citrated Global Hemostasis Test Cartridge • Citrated Global hemostasis with Lysis (Trauma) Test Cartridge | <ul style="list-style-type: none"> • Sealed vials are stored at 2-8°C • Use within 2 hrs of reconstitution | <ul style="list-style-type: none"> • Allow QC vial and diluent water vial to equilibrate at RT for 10 min • Tap the top of QC vial a few times to ensure lyophilized material is on the bottom of the QC vial • Remove seal and stopper of the QC vial, avoid the sharp metal edges • Slowly pour the contents of the diluent water into the QC vial and make sure no water drips out • Re-insert stopper in the QC vial • Hold stopper in place, vigorously shake the QC vial until fully reconstituted then let stand at RT for 5 min • Shake QC vial vigorously again and let stand at RT for 5 more min • Repeat until there is |

Level 2
Abnormal
QC vial
kit

- Citrated Global Hemostasis Test Cartridge
- Citrated Global hemostasis with lysis (Trauma) Test Cartridge

- Sealed vials are stored at 2-8°C
- Use within 2 hrs of reconstitution

no undissolved material remaining in the QC vial

- Allow QC vial and diluent water vial to equilibrate at RT for 10 min
- Tap the top of QC vial a few times to ensure lyophilized material is on the bottom of QC vial
- Remove seal and stopper of the QC vial, avoid the sharp metal edges
- Slowly pour the contents of the diluent water into the QC vial and make sure no water drips out
- Re-insert stopper in the abnormal QC vial
- Hold stopper in place, vigorously shake the abnormal QC vial until fully reconstituted then let stand at RT for 5 min
- Shake abnormal QC vial vigorously again and let stand at RT for 5 more min
- Repeat until there is no undissolved material remaining in the abnormal QC vial

Normal
Donor
Sample

- Platelet Mapping Test

- Normal donor sample stored at RT in a horizontal position

- Refer to Haemonetics donor screening criteria for identification and

COPY

Cartridge

- Use after 30 minutes of collection at RT and within 2 hrs of collection for Platelet Mapping Test cartridge

selection of normal donor samples

- Donors taking medications such as oral contraceptives, hormone replacement of any type, aspirin, ibuprofen, naproxen should be excluded
- Refer to the TEG 6s specimen collection procedure to collect the appropriate blood tubes

Abnormal Donor sample

- Platelet Mapping Test Cartridge

- Abnormal donor sample stored at RT in a horizontal position

- Use after 30 minutes of collection at RT and within 2 hrs of collection for Platelet Mapping Test cartridge

- Refer to Haemonetics donor screening criteria for identification and selection of normal donor samples
- Refer to the TEG 6s specimen collection procedure to collect the appropriate blood tubes
- To prepare aspirin spiked normal donor blood:

- Dissolve (1) 325 mg aspirin tablet in 200 ml saline
- add 160 ul of the ASA solution to the 4 ml heparinized blood tube.
- Recap vial. Thoroughly

and gently
mix

- Incubate
solution at
RT for 30
min prior to
testing

- To prepare
Functional
Fibrinogen rgt spiked
normal donor blood:

- Remove 3
vials of FF
rgt from
the
refrigerator
and allow
to reach RT

- tap the
vials to
ensure all
contents
are at the
bottom of
the vial

- Remove
seal and
stopper on
each vial,
avoid
sharp
metal
edges

- Pipette 500
ul of whole
blood into
each of the
FF vials

- On each
vial,
replace the
stopper
then gently
swirl and
invert 5

COPY

times

- Obtain a large sterile plastic container and pool the contents of the 3 vials
- Gently swirl and invert 5 times prior to testing

PROCEDURE

Follow the procedure below to run QC on the TEG 6s

| Step | Action |
|------|--|
| 1. | Remove the test cartridge to use from refrigerated storage |
| 2. | From the <i>Home</i> screen on the TEG analyzer, select new qc |
| 3. | Tear open the cartridge pouch and when prompt, insert the cartridge into the slot as indicated with the bar code on the left side |
| 4. | On the <i>Confirm Test</i> screen, touch continue |
| 5. | After the cartridge pretest has completed and you have verified that the assay is what you intended to run, touch next |
| 6. | On the Test Information screen: <ul style="list-style-type: none">• Using the external QC barcodes located above the TEG 6s analyzers, scan the appropriate barcode of the QC that will be run (Level 1 normal QC, Level 2 abnormal QC, normal donor QC, abnormal donor QC)• Then touch next |
| 7. | Pipette the prepared QC sample into the cartridge sample port, filling up to or above the line marked on the cartridge. Touch next . <ul style="list-style-type: none">• The TEG analyzer starts the test• The results are displayed as they become available |
| 8. | Touch tracings to view a graphic representation of the results <ul style="list-style-type: none">• To cycle through superimposed, offset, and single-tracing views, touch next tracing until the desired view is displayed |

9. When the analyzer displays the "Remove cartridge" prompt, removed the used cartridge from the slot and immediately dispose of it in a biohazard container
10. Print out and compare the QC results against the established QC ranges
 - Reagent performance is verified if the QC test results fall within the established QC ranges
 - Refer to section "Verifying External QC Results"
11. Previously ran QC may be viewed from the *Home* screen by touching **stored qc**
 - Select the desired test, then touch results
 - The status of each test is shown on the right side of the screen
 - A green check mark indicates that the test completed (all parameters were finalized), a red X indicates that the test timed out, and an orange triangle indicates that the test was stopped early

VERIFYING EXTERNAL QC RESULTS

| Step | Action | | | | | | | | | | | | |
|--|---|------|------|---|---|--|--|----|------|---|---|--|---|
| 1. | Print the external QC results for the test cartridge in use | | | | | | | | | | | | |
| 2. | Refer to the TEG 6s external QC reference ranges to review the external QC results for the test cartridge in use and ensure values are within acceptable range | | | | | | | | | | | | |
| | <table border="1"> <thead> <tr> <th>If</th> <th>Then</th> </tr> </thead> <tbody> <tr> <td>QC results within acceptable range</td> <td> <ul style="list-style-type: none"> • Initial the printouts • Place printout in the TEG 6s QC binder • Proceed to step 3 </td> </tr> <tr> <td rowspan="2">QC results outside acceptable range</td> <td> Repeat QC sample on test cartridge in use <table border="1"> <thead> <tr> <th>If</th> <th>Then</th> </tr> </thead> <tbody> <tr> <td>Repeat QC results within acceptable range</td> <td> <ul style="list-style-type: none"> ◦ Initial the printouts ◦ Place printouts in the TEG 6s QC binder ◦ Proceed to step 3 </td> </tr> <tr> <td>Repeat QC results outside acceptable range</td> <td> <ul style="list-style-type: none"> ◦ Notify Supervisor and contact Haemonetics tech support ◦ Do not use affected </td> </tr> </tbody> </table> </td> </tr> </tbody> </table> | If | Then | QC results within acceptable range | <ul style="list-style-type: none"> • Initial the printouts • Place printout in the TEG 6s QC binder • Proceed to step 3 | QC results outside acceptable range | Repeat QC sample on test cartridge in use <table border="1"> <thead> <tr> <th>If</th> <th>Then</th> </tr> </thead> <tbody> <tr> <td>Repeat QC results within acceptable range</td> <td> <ul style="list-style-type: none"> ◦ Initial the printouts ◦ Place printouts in the TEG 6s QC binder ◦ Proceed to step 3 </td> </tr> <tr> <td>Repeat QC results outside acceptable range</td> <td> <ul style="list-style-type: none"> ◦ Notify Supervisor and contact Haemonetics tech support ◦ Do not use affected </td> </tr> </tbody> </table> | If | Then | Repeat QC results within acceptable range | <ul style="list-style-type: none"> ◦ Initial the printouts ◦ Place printouts in the TEG 6s QC binder ◦ Proceed to step 3 | Repeat QC results outside acceptable range | <ul style="list-style-type: none"> ◦ Notify Supervisor and contact Haemonetics tech support ◦ Do not use affected |
| If | Then | | | | | | | | | | | | |
| QC results within acceptable range | <ul style="list-style-type: none"> • Initial the printouts • Place printout in the TEG 6s QC binder • Proceed to step 3 | | | | | | | | | | | | |
| QC results outside acceptable range | Repeat QC sample on test cartridge in use <table border="1"> <thead> <tr> <th>If</th> <th>Then</th> </tr> </thead> <tbody> <tr> <td>Repeat QC results within acceptable range</td> <td> <ul style="list-style-type: none"> ◦ Initial the printouts ◦ Place printouts in the TEG 6s QC binder ◦ Proceed to step 3 </td> </tr> <tr> <td>Repeat QC results outside acceptable range</td> <td> <ul style="list-style-type: none"> ◦ Notify Supervisor and contact Haemonetics tech support ◦ Do not use affected </td> </tr> </tbody> </table> | If | Then | Repeat QC results within acceptable range | <ul style="list-style-type: none"> ◦ Initial the printouts ◦ Place printouts in the TEG 6s QC binder ◦ Proceed to step 3 | Repeat QC results outside acceptable range | <ul style="list-style-type: none"> ◦ Notify Supervisor and contact Haemonetics tech support ◦ Do not use affected | | | | | | |
| | If | Then | | | | | | | | | | | |
| Repeat QC results within acceptable range | <ul style="list-style-type: none"> ◦ Initial the printouts ◦ Place printouts in the TEG 6s QC binder ◦ Proceed to step 3 | | | | | | | | | | | | |
| Repeat QC results outside acceptable range | <ul style="list-style-type: none"> ◦ Notify Supervisor and contact Haemonetics tech support ◦ Do not use affected | | | | | | | | | | | | |

TEG 6s analyzer for patient testing

- Initial the printouts
- Place printouts in the TEG 6s QC binder

3. Proceed to patient testing when external QC results have been verified to be within acceptable ranges for test cartridges in use

EXTERNAL QC REFERENCE RANGES

| Level 1 Normal QC Citratd Global Hemostasis Test Cartridge | R (min) | K (min) | Angle (deg) | MA (mm) | FLEV (mg/dl) |
|--|----------------|----------------|--------------------|----------------|---------------------|
| CK | 4.6 - 13.1 | 0.8- 4.1 | 58 - 81 | 55 - 73 | NA |
| CRT | NA | NA | NA | 55 - 73 | NA |
| CKH | 3.6 - 8.6 | NA | NA | NA | NA |
| CFF | NA | NA | NA | 55 - 73 | 1003 - 1333 |
| Level 2 Abnormal QC Citratd Global Hemostasis Test Cartridge | R (min) | K (min) | Angle (deg) | MA (mm) | FLEV (mg/dl) |
| CK | 1.0 - 1.7 | 0.7 - 1.4 | 61 - 81 | 22 - 34 | NA |
| CRT | NA | NA | NA | 22 - 34 | NA |
| CKH | 1.0 - 1.7 | NA | NA | NA | NA |
| CFF | NA | NA | NA | 22 - 34 | 401 - 621 |
| Level 1 Normal QC Citratd Global Hemostasis Test Cartridge with Lysis | | | R (min) | MA (mm) | LY30 (%) |
| CK | | | 4.6 - 13.1 | NA | 0.0 - 0.0 |
| CRT | | | NA | 55 - 73 | NA |
| CFF | | | NA | 55 - 73 | NA |
| Level 2 Abnormal QC Citratd Global Hemostasis Test Cartridge with Lysis | | | R (min) | MA (mm) | LY30 (%) |
| CK | | | 1.0 - 1.7 | NA | 79 - 95 |
| CRT | | | NA | 22 - 34 | NA |

| CFF | | | | | NA | 22 - 34 | NA | |
|---------------------------------|-------------|--------------|-------------|--------------|-----------|------------|--------------|------------|
| Platelet Mapping Test Cartridge | HKH-MA (mm) | ActF-MA (mm) | ADP-MA (mm) | ADP (%Inhib) | ADP(%agg) | AA-MA (mm) | AA (% inhib) | AA (% Agg) |
| Normal donor QC | 53 - 68 | 2 - 19 | 45 -69 | 0 - 17 | 83 - 100 | 51 -71 | 0 - 11 | 89 - 100 |
| Abnormal donor QC | <53 | N/A | <45 | >17% | <83% | <51 | >11% | <89% |

RELATED DOCUMENTS

- Running Specimens on the TEG 6s Hemostasis Analyzer
- Donor Screening Questions for Reference Range Study

REFERENCES

- Haemonetics TEG 6s User Manual

All Revision Dates

9/12/2022, 4/7/2022, 1/4/2021

Attachments

[Normal donor screening criteria.pdf](#)

Approval Signatures

| Step Description | Approver | Date |
|---------------------|---|-----------|
| Medical Director | Lindsey Westerbeck: Director, Laboratory Services | 9/12/2022 |
| Laboratory Director | Lindsey Westerbeck: Director, Laboratory Services | 8/15/2022 |