 <b>YALE-NEW HAVEN HOSPITAL</b>	<b>TITLE:</b>  <p style="text-align: center;"><b>SUDAN BLACK B</b></p>		<b>DEPT OF LAB MEDICINE CLINICAL HEMATOLOGY Policy and Procedure Manual</b>
			<b>DOCUMENT #</b> H-05-002
			Page 1 of 3
<b>WRITTEN BY:</b> Paula Morris, MT (ASCP)	<b>EFFECTIVE DATE:</b> 08-05-97	<b>REVISION:</b> H-2 5/24/12	<b>SUPERCEDES:</b> H-1 8/5/97

### I. PRINCIPLE:

Sudan Black B stains lipids in cells. Neutrophils stain the most strikingly, followed by monocytes with a sprinkling of black granules. Eosinophilic granules stain as black circles as only the membrane of the granules stain.

### II. SPECIMEN:

Coverslips or 3x1 slides using patient peripheral blood (EDTA) or bone marrow samples.

### III. REAGENTS:


- A. 40% formalin
- B. Stock buffer solution
  - 1. 91 ml of liquid phenol -88%- (A.R.)
  - 2. 150 ml absolute ethanol (200 proof)
  - 3. 1.5 gm Na<sub>2</sub>HPO<sub>4</sub> 12 H<sub>2</sub>O dissolved first in 500 ml distilled water  
Mix phenol and ethanol then add phosphate and water mixture.
- C. Stock Sudan solution
  - 1. 1.50 gm Sudan Black B (Fisher CI #26150)
  - 2. Dissolve in 500 ml absolute ethanol
  - 3. Dissolve completely by shaking continually two hours.
- D. Working incubation mixture-**make fresh**

#### **For 3x1 slides:**

- 1. 20 ml stock buffer
- 2. 30 ml stock Sudan solution
- 3. Filter into coplin jar over slides or onto slides on staining rack.

#### **For coverslips:**

- 1. 8 ml stock buffer
- 2. 12 ml stock Sudan solution
- 3. Filter over coverslips in coplin jar.

 <b>YALE-NEW HAVEN HOSPITAL</b>	<b>TITLE:</b>  <p style="text-align: center;"><b>SUDAN BLACK B</b></p>		<b>DEPT OF LAB MEDICINE CLINICAL HEMATOLOGY Policy and Procedure Manual</b>
			<b>DOCUMENT #</b> H-05-002
			Page 2 of 3
<b>WRITTEN BY:</b> Paula Morris, MT (ASCP)	<b>EFFECTIVE DATE:</b> 08-05-97	<b>REVISION:</b> H-2 5/24/12	<b>SUPERCEDES:</b> H-1 8/5/97

#### IV. QUALITY CONTROL:

- A. Normal patient EDTA coverslip sample with Monocytes >10 and Lymphocytes > 18 with a WBC count >6,000.

#### V. PROCEDURE:


- A. Fix in Formalin Vapor for 10 minutes.
- B. Stain fixed smear in Sudan (working incubation mixture) for 10 minutes.
- C. Wash slide well with 95 % ethanol (95 ml absolute ethanol with 5 ml DI water).
- D. Wash with tap water.
- E. Counterstain with 1:10 Giemsa using Hematek buffer for 20 minutes.  
**(Mix in a red top tube 1ml bottled Giemsa stain with 9 ml of Hematek buffer)**  
Pour over coverslips in coplin jar or onto slides in coplin jar or on staining rack:  
If using large slides in coplin jar adjust volume as follows: **5ml Giemsa with 45ml Hematek buffer.**
- F. Wash with tap water and air dry.
- G. Mount coverslips with permount, coverslip 3x1 slides with permount.

#### VII. RESULTS:

Results are evaluated by an attending physician and resident.

#### VIII. REFERENCES:

1. Sheehan, H.L., Storey, GW: An improved Method of Staining Leukocyte Granule with Sudan Black B, Journal of Pathology 59:336-337. 1947

 <b>YALE-NEW HAVEN HOSPITAL</b>	<b>TITLE:</b>  <b>SUDAN BLACK B</b>		<b>DEPT OF LAB MEDICINE CLINICAL HEMATOLOGY Policy and Procedure Manual</b>
			<b>DOCUMENT #</b> H-05-002
			Page 3 of 3
<b>WRITTEN BY:</b> Paula Morris, MT (ASCP)	<b>EFFECTIVE DATE:</b> 08-05-97	<b>REVISION:</b> H-2 5/24/12	<b>SUPERCEDES:</b> H-1 8/5/97

2. Baillif, RN, Kinbrough, R.T.: Studies on Leukocyte Granules after Staining with Sudan Black B Nag-Grimwald Giemsa, J. Lab. Clin. Med. 32, p. 155-166.

## **IX. HISTORY:**

- H-1 This procedure was written by P. Morris on 08-05-97.
- H-2 This procedure was revised by Susan Richardson on 5/24/12.