

## **Bone Marrow Transplant Spec Handling Technical Procedure #1549.t**

### **Purpose**

A standardized protocol for performing cell counts and differentials on donor bone marrow is required to ensure accurate and prompt testing and reporting of results.

### **Specimen**

Specimens submitted to the hematology laboratory for cell counts and differential consist of liquid cell suspensions. For Progenitor Cell Specimens (HPC), see procedure 1548t.

### **Reagents**

- Bovine Albumin- 22% Protein Concentration ph 7.2-Ortho Diagnostic Systems.

### **Equipment**

- Coulter LH780 with reagents
- Wrights stainer
- Vortex
- Glass microscope slides
- Microscope

### **Quality Control**

- Coulter LH780-refer to LH780 Procedure #1510.
- Stain quality is evaluated on a slide by slide basis by the CLS.

### **Specimen Requirements**

- Send approximately 1 ml of specimen to the laboratory
  - For bone marrow harvests (usually from UCDCM operating room), the tube and requisition slip are labeled as:
    - Interim Bone Marrow Harvest #1, #2, #3. etc. or
    - Final filtered Bone Marrow

**Bone Marrow Transplant Spec Handling Technical Procedure #1549.t**

- The specimen and the accompanying requisition must be delivered directly to the Hematology Laboratory
- Specify callback name and phone or pager for STAT results.

**Requisitioning**

- Bone marrow harvest specimens are handled as Priority One specimens and are requisitioned in the LIS system.
  - Bone Marrow Harvests
    - Specimens labeled as Bone Marrow Harvest #1, #2, etc. are requisitioned as a nucleated cell count. The mnemonic is WBCO.
  - Final filtered Bone Marrow
    - Specimens labeled as Final Filtered Marrow are requisitioned as WBC, HGB, HCT, PLT, WBCO, and BMD.

**Procedure**

- Cell Counts
  - Mix specimen well – both by inversion and gentle vortexing.
  - Check specimen for clot using applicator sticks.
  - Run specimen on Coulter LH780 in the secondary mode, carefully observing the aspiration to ensure that it does not short sample.
- Reporting Results
  - Interim Bone Marrow Harvest
    - Report only the Total Nucleated Cell Count (uncorrected WBC count), using the canned text mnemonic HWBCO and identify the specimen as “Bone Marrow Harvest #\_”
    - Call the automated cell count results to the callback number.
    - This is a preliminary count and must be verified by slide estimate or manual methods, if needed. If result is determined

University of California, Davis Medical Center  
Department of Pathology, Hematology

**Bone Marrow Transplant Spec Handling Technical Procedure #1549.t**

to be lower than originally reported, re-notify via the callback number.

- Final Filtered Marrow
  - Report the WBC count (use the uncorrected WBC) with the canned text mnemonic BMC, to append the following comment: “THIS SPECIMEN IS NOT PERIPHERAL BLOOD, Specimen is Final Filtered Bone Marrow collected for transplant.” Call the automated cell count results to the callback number.
  - This is a preliminary count and must be verified by slide estimate or manual methods, if needed. If result is determined to be lower than originally reported, re-notify via the callback number.
  - If the platelet count exceeds linearity, the instrument count should be reported with the comment “Result exceeds instrument linearity and had not been confirmed due to insufficient specimen.”
- Slide Preparation and Staining
  - Prepare two wedge smears from Final Filtered Bone Marrow.
  - Stain slides using the Bone Marrow protocol
  - Check smears for stain quality. Remake unacceptable smears, using one drop of albumin to 4-5 drops of specimen, if needed, due to cell fragility.
- Differential
  - A qualified CLS will perform the differential. The Hemepath Fellow or Pathologist may be consulted if needed.
  - Perform a 200 WBC cell differential; NRBC’s are included in the differential. Report as percent %.
  - Enter the results of the differential using the canned text BMD. Report the number of NRBC’s at the bottom of the report (as they are not in the canned text).
  - Completed smears are saved in the progenitor box.

**Bone Marrow Transplant Spec Handling Technical Procedure #1549.t**

**Reference Range**

Reference ranges are included in the BMD canned text.

**References**

- “Bone Marrow Procedure”. Clinical Cryobiology Laboratory Fred Hutchinson Cancer Research Center Seattle, Washington.
- “Procedure For Processing Bone Marrow Harvest Specimens”. Clinical Laboratory Services Rush Presbyterian St. Luke’s Medical Center.
- Sacher, R et al; Marrow Transplantation: Practical and Technical Aspects of Stem Cell Reconstitution, pp133-135; AABB, 1992
- Sacher, R et al; Processing of Bone Marrow for Transplantation, pp.42-43; AABB, 1990

University of California, Davis Medical Center  
Department of Pathology, Hematology

**Bone Marrow Transplant Spec Handling Technical Procedure #1549.t**

**Procedure History**

Date	Written/Revised By	Revision	Approved Date	Approved By
6/18/93	C. Sunset		9/17/93	Ed Larkin, MD
			10/12/93	Ed Larkin, MD
			12/2/94	Ed Larkin, MD
3/20/96	J. Cannon	Minor	11/95	Ed Larkin, MD
		Annual Review	11/96	Ed Larkin, MD
		Annual Review	11/97	Ed Larkin, MD
		Annual Review	10;98	Ed Larkin, MD
		Annual Review	11/99	Ed Larkin, MD
		Annual Review	11/00	Ed Larkin, MD
		Annual Review	10/17/01	Ed Larkin, MD
11/02	J Cannon	MS Word	11/5/02	Ed Larkin, MD
		Annual Review	10/17/03	Ed Larkin, MD
		Annual Review	11/27/04	Ed Larkin, MD
		New Director	02/24/05	Kim Janatpour, MD
11/06	J Cannon	Annual Review	11/18/06	Denis Dwyre MD
04/07	J Cannon	Minor updates		
		Annual Review	05/02/07	Denis Dwyre, MD
		Annual Review	11/05/07	Denis Dwyre MD
		Annual Review	07/03/08	Denis Dwyre MD
		Annual Review	10/27/09	D Dwyre MD
		Annual Review	10/15/10	D Dwyre MD
		Biannual Review	08/24/12	D Dwyre MD
		Biannual Review	10/1/14	D Dwyre MD
10/2015	L Gandy	Minor updates		