

Purpose

To define minimum volumes of sample and diluent for the first manual dilution for each analyte, subsequent manual dilution ratios and maximum dilution ratios that can be used to obtain analyte results.

Procedure

When manual dilutions are required to obtain a result, use the following guidelines:

1. Determine which dilution factor to use by looking at the error codes for the analyte in Remisol. Refer to the following manual dilution tables for the dilution factor, minimum sample volume and diluent volume to use for the analyte(s) being diluted. Do not over-dilute or under-dilute. The "raw" diluted value should be in the middle 50% of the AMR (25% - 75% of the AMR). Dilutions should not be run with the ORDAC function. Serial dilutions should be avoided, with the exception of Methotrexate.
2. Program the dilution for the analyte needing a dilution in Remisol. The Remisol barcode printer will print a label after the dilution is programmed.
3. Label a container with the barcode label generated for the programmed dilution.
4. Always confirm the name and sample ID number of the specimen.
5. Take sample for dilution(s) from the primary tube whenever possible.
6. Make appropriate dilution with recommended diluent stated in the analyte procedure.
7. Only use the following calibrated MLA pipettors for making dilutions:

25 µL
50 µL
100 µL
200 µL
250 µL
300 µL
500 µL
1000 µL

8. For the maximum dilution (X201) a 5 mL volumetric pipet may be used for the diluent

A few analytes may require an off-line dilution (ex. low phosphorus, low suppressed urine Na, some TDMs). Refer to the specific analyte procedure for recommended dilutions and required calculations.

Manual Dilution Tables

Dxl Immunoassays

Analyte	Dilution Ratio	Sample Volume	Diluent Volume
Digoxin	X2	200 µL	200 µL
E2 (Estradiol)	X2	100 µL	100 µL
PTH, PTHIO	X10	50 µL	450 µL
TSH	X5	100 µL	400 µL

University of California, Davis Health System
Department of Pathology and Laboratory Medicine
Chemistry and Urinalysis

Guidelines for Specimen Dilutions

Technical Procedure 3045

Centaur Immunoassays

Analyte	Dilution Ratio	Sample Volume	Diluent Volume
Ferritin	X21	25 µL	500 µL
Ferritin	X41	50 µL	2000 µL
Ferritin	X61	50 µL	3000 µL

DxC Chemistries

Analyte	Dilution Ratio	Sample Volume	Diluent Volume
MC Chemistries			
CO2	X2	50 µL	50 µL
CO2	X3	50 µL	100 µL
Calcium	X2	50 µL	50 µL
BUN, Glucose	X3	50 µL	100 µL
BUN, Glucose	X6	50 µL	250 µL
Crea, Phos, ALB, TP	X2	50 µL	50 µL
Crea, Phos, ALB, TP	X3	50 µL	100 µL
Phos	X4	100 µL	300 µL
Enzymes			
ALP, ALT, AMY, AST, CK, GGT, LD, LIP	X2	50 µL	50 µL
ALP, ALT, AMY, AST, CK, GGT, LD, LIP	X3	50 µL	100 µL
ALP, ALT, AMY, AST, CK, GGT, LD, LIP	X6	50 µL	250 µL
ALP, ALT, AMY, AST, CK, GGT, LD, LIP	X11	50 µL	500 µL
ALP, ALT, AMY, AST, CK, GGT, LD, LIP	X21	25 µL	500 µL
ALP, ALT, AMY, AST, CK, GGT, LD, LIP	X41	50 µL	2000 µL
ALP, ALT, AMY, AST, CK, GGT, LD, LIP	X61	50 µL	3000 µL
ALP, ALT, AMY, AST, CK, GGT, LD, LIP	X81	50 µL	4000 µL
ALP, ALT, AMY, AST, CK, GGT, LD, LIP	X101	50 µL	5000 µL
AMY, CK	X201	25 µL	5000 µL
Lipid Tests			
CHOL, TRIG	X2	50 µL	50 µL
CHOL, TRIG	X3	50 µL	100 µL
CHOL, TRIG	X6	50 µL	250 µL
CHOL, TRIG	X11	50 µL	500 µL
CHOL, TRIG	X21	25 µL	500 µL
HDL, LDL	X2	50 µL	50 µL
HDL, LDL	X3	50 µL	100 µL

University of California, Davis Health System
Department of Pathology and Laboratory Medicine
Chemistry and Urinalysis

Guidelines for Specimen Dilutions

Technical Procedure 3045

DxC Chemistries - continued

Analyte	Dilution Ratio	Sample Volume	Diluent Volume
Miscellaneous CC Tests			
DBIL, FE, ETOH, LACT, MG, TBIL, SALY	X2	50 µL	50 µL
DBIL, FE, ETOH, LACT, MG, TBIL, SALY	X3	50 µL	100 µL
DBIL, ETOH, LACT, MG, TBIL, SALY	X6	50 µL	250 µL
AMMONIA, TRANSF	X2	50 µL	50 µL
URIC	X2	50 µL	50 µL
URIC	X3	50 µL	100 µL
URIC	X6	50 µL	250 µL
Therapeutic Drugs			
ACTM, CAR, PHE, TOB, VPA, GEN	X2	50 µL	50 µL
CAR, PHE, TOB, VPA, GEN	X3	50 µL	100 µL
CAR, PHE, TOB, VPA, GEN	X6	50 µL	250 µL
CAR, PHE, VPA, GEN	X11	50 µL	500 µL
Immunology Tests			
C3, C4, HPT, IgA, IgG, IgM, PAB	X2	50 µL	50 µL
C3, C4, PAB	X3	50 µL	100 µL
C3, C4, PAB	X6	50 µL	250 µL
C3, C4, PAB	X11	50 µL	500 µL
Urine Chemistries			
Urine - Ca, Phos, M-TP, MA, UN, Uric	X2	50 µL	50 µL
Urine - Ca, Phos, M-TP, MA, UN, Uric	X3	50 µL	100 µL
Urine - Ca, Phos, M-TP, MA, UN, Uric	X6	50 µL	250 µL
Urine - Ca, Phos, M-TP, MA, UN, Uric	X11	50 µL	500 µL
Urine - Ca, Phos, M-TP, MA, UN, Uric	X21	25 µL	500 µL
Urine - M-TP	X41	50 µL	2000 µL
Urine - M-TP	X61	50 µL	3000 µL
Urine - Crea	X2	50 µL	50 µL
Urine - Crea	X3	50 µL	100 µL
Urine - Crea	X6	50 µL	250 µL
CSF Chemistries			
CSF Glucose	X2	50 µL	50 µL
CSF Glucose	X3	50 µL	100 µL
CSF Glucose	X6	50 µL	250 µL
CSF Protein	X2	50 µL	50 µL
CSF Protein	X3	50 µL	100 µL
CSF Protein	X6	50 µL	250 µL
CSF Protein	X11	50 µL	500 µL
CSF Protein	X21	25 µL	500 µL

Methotrexate Serial Dilutions

Dilution Ratio	Volume 1	Diluent Volume
X10	50 μ L sample	450 μ L saline
X100	50 μ L x10 dilution*	450 μ L saline
X1000	50 μ L x100 dilution*	450 μ L saline

* Dilutions should be well-mixed with a transfer pipette prior to using the diluted sample for the next dilution.
Refer to the Methotrexate Technical Procedure 3146 for further information.

