


QUALITY CONTROL

How to Evaluate QC and Troubleshoot Out of Control Results

Kathy Dagang, CLSS
Automated Chemistry/Urinalysis Department
UC Davis Health

OBJECTIVES

1. Define Quality Control and its basic elements.
2. Define the criteria we use to evaluate Quality Control: Mean, Standard Deviation, and Coefficient of Variation.
3. Review basic calibration theory.
4. Provide a consistent method to troubleshoot out of control results.
5. Describe Lot-to-Lot evaluation.



QUALITY CONTROL



What should I do?



WHAT IS QUALITY CONTROL?

Statistical process that tracks the ability of reagents, instrumentation (and techs) to obtain an accurate result...




WHY DO WE RUN QC?

- Report out good results
- Detect SIGNIFICANT errors in a timely manner
- If there is an error, identify the source of the error



QC STATISTICS



GLUCOSE RESULTS



QC STATISTICS


The "MEAN", or AVERAGE is a good estimate of a test's ACCURACY:

$$\bar{X} = \frac{\sum X}{n}$$




QC STATISTICS

Normal Distribution




↑
MEAN

Skewed Distribution



↑
MEAN

Skewed Distribution




↑
MEAN




QC STATISTICS

$$SD = \sqrt{\frac{\sum_{i=1}^n (x_i - \bar{x})^2}{n - 1}}$$

The "STANDARD DEVIATION", or SD, is a measure of PRECISION.




Statistically, 95% of QC points will fall within $\pm 2SD$ of the mean.



QC CHARTS

QC CHARTS usually appear similar to this, with the line in the middle representing the MEAN, and lines above and below the mean representing the SDs.



QC STATISTICS


The "COEFFICIENT OF VARIATION", or CV, is a ratio of the SD to the Mean, and is expressed as a %.

$$CV (\%) = \left(\frac{\text{Standard deviation}}{\text{Mean}} \right) \times 100$$



CALIBRATION

A calibration establishes the relationship between



ANALYTE
CONCENTRATION
and
INSTRUMENT
RESPONSE.



CALIBRATORS


A **calibrator** is a sample with a **known analyte concentration**.



The "**set point**" on the calibration report is the **calibrator value**.



CALIBRATION



CALIBRATION CRITERIA

Three criteria are used to evaluate a successful calibration:

ACCURACY
PRECISION
SENSITIVITY



ACCURACY

The definition of ACCURACY is
recovering expected results.

The values obtained must fall between a pre-determined minimum and maximum value.



PRECISION

The definition of PRECISION is
reproducibility.

The values obtained must differ by less than a predetermined amount.



SENSITIVITY

SENSITIVITY allows us to see a **difference in system response** with different analyte concentrations.

Without an acceptable span, the instrument may not "resolve" two close, yet different, analyte concentrations.



CALIBRATION

Successful calibration
requires acceptable results
for ACCURACY, PRECISION,
and SPAN.

QC that's run after calibration can help you
confirm that the calibration was successful.



SO, QC IS "OUT"



First, make sure your control is really "out".
Compare acceptable range in BEAKER
with the ranges in the instrument or Remisol.





CHECK THE CONTROL

- Are you using the correct control/level of control?
- Is it outdated or expired, or the end of the bottle?
- Has it been stored properly? Has it been sitting at room temperature for an extended time?
- Is it the correct lot number?
- Are there any signs of contamination or degradation?
- Is it a different shipment of control?



IF NO APPARENT ISSUES WITH THE CONTROL:

Rerun the test on a New Aliquot of Same
Bottle of Control

IF THE BOTTLE OF CONTROL IS QUESTIONABLE, OR THE RERUN IS STILL OUT OF CONTROL:

Open a NEW BOTTLE of control and
rerun the test.





CHECK THE REAGENT

- Was the reagent prepared properly?
- Is it close to its outdate?
- Is it due to be calibrated?
- Are there just a few tests left (say, less than 10%)?
- Is there anything going on in the instrument that could affect the reagents?



IF NO APPARENT ISSUES WITH THE REAGENT

Recalibrate the test and rerun the QC on all controls for the test.

IF THE REAGENT IS QUESTIONABLE

Prepare and/or Load new reagent before you recalibrate the test and rerun the QC on all controls for the test.






CHECK THE CALIBRATOR

Check the expiration date and open stability of the calibrator.

Verify correct calibrator lot number is loaded/used.

(Compare the Lot# on the calibration printout with the Lot# on the bottle of calibrator.)


Load new calibrator data if required.







CHECK THE CALIBRATOR

Review QC charts for a "shift".



CHECK THE CALIBRATOR




RECALIBRATE AND
RERUN CONTROLS



STILL OUT??



Now's the time to talk to a specialist or supervisor...



LOT CHANGES



Changing lots of reagent can often cause a shift in QC THAT IS NOT SEEN WITH PATIENT SAMPLES.



LOT-TO-LOT VERIFICATION

CONFIRMING THAT CHANGING LOTS OF REAGENT DOESN'T AFFECT PATIENT RESULTS IS REQUIRED BY CAP (COM.30450) AND CLIA REGULATIONS.

PATIENT SPECIMENS ARE RUN ON BOTH THE OLD AND NEW LOTS.



LOT-TO-LOT VERIFICATION

DO NOT USE THIS LOT NUMBER
NEW LOT Record # in Log before using
THIS LOT IS READY FOR USE

THAT MEANS WE
NEED TO KEEP
CLOSE TRACK OF
WHICH LOTS
ARE IN USE.
ALWAYS.



LOT-TO-LOT VERIFICATION

WHEN A NEW LOT IS RECEIVED, **ORANGE** STICKERS WILL BE PLACED ON EACH BOX.

WHEN THE NEW LOT HAS BEEN VERIFIED, **YELLOW** STICKERS WILL BE PLACED OVER THE ORANGE STICKERS.

WHEN THE NEW LOT IS PUT IN USE, **GREEN** STICKERS WILL BE PLACED OVER THE YELLOW STICKERS.

NEW SHIPMENTS OF THE SAME LOT WILL GET THE SAME COLOR STICKERS.



LOT-TO-LOT VERIFICATION

TO VERIFY A NEW LOT, SELECT A NORMAL AND ABNORMAL PATIENT THAT WAS RUN ON THE OLD LOT, KEEPING SPECIMEN STABILITY IN MIND.

- ▶ TROPONIN LOT CHANGES REQUIRE AT LEAST THREE PATIENT SAMPLES: A NORMAL, AN ALERT VALUE, AND A CRITICAL VALUE.

RUN APPROPRIATE QC ON THE NEW LOT.

- ▶ EACH REAGENT WILL HAVE ACCEPTABILITY CRITERIA FOR THE PATIENT COMPARISONS.

IF THE QC IS ABOUT THE SAME AS THE PREVIOUS LOT'S QC, AND THE PATIENT RESULTS AGREE WITHIN ACCEPTABLE LIMITS, PLACE A YELLOW STICKER OVER THE ORANGE STICKER.



LOT-TO-LOT VERIFICATION

IF THE QC IS OUT, OR IS SIGNIFICANTLY DIFFERENT (>1 SD), **BUT THE PATIENT RESULTS ARE ACCEPTABLE**, NOTIFY A SPECIALIST OR SUPERVISOR. PLACE YELLOW STICKERS OVER THE ORANGE STICKERS ON THAT LOT.

IF THE PATIENT RESULTS EXCEED THE ALLOWABLE DIFFERENCE, NOTIFY A SPECIALIST OR SUPERVISOR!



HELPFUL HINTS

- **DON'T RUN CONTROLS OVER AND OVER HOPING IT WILL COME IN**, especially if you're getting the same result.
- Limit the time controls and calibrators are kept at room temperature to a minimum. Always.
- Check expiration dates. Always.
- Be aware of LOT#s of controls and calibrators. Always.
- Review the QC chart at the right of the Beaker results. Look for results close to the mean, and scattered fairly evenly on both sides. *Just because your control is "in", doesn't mean that you're "in control".*
- Ask for help if you need it.



**RAN QC AFTER 6
MONTH MAINTENANCE
AND CALIBRATION**



IN ON FIRST TRY

WWW.FACEBOOK.COMLABHUMOR

It can't always be like this, but with good lab practices and everyone's help monitoring QC, running QC won't be torture...



#QC123

#DONT PANIC

#YOU CAN DO IT

REFERENCES

3000 Quality Control Procedure

3001 Lot to Lot Change Procedure

Basic Lessons in Laboratory Quality Control, BioRad QC Workbook, 2008

Basic QC Practices, 3rd Edition, James Westgard, Ph.D., 2010

Tietz Testbook of Clinical Chemistry, 2nd ed., Burtis, Ashwood, editors, 1998